



FIMA INSTRUMENTS

S U R G I C A L



**Manufacturer and Exporter of Surgical / Dental /
Manicure and Veterinary Instruments.**

FIMA DENTAL INSTRUMENTS

We feel immense pleasure to present our latest catalogue in which we have featured a most comprehensive range of all types of quality Dental & Surgical instruments required for day to day use for Dental Clinics. We will be pleased to quote you a competitive Price and Delivery time.

All instruments are manufactured according to the recognized standard patterns currently in use in Dental Clinics throughout the world and are being produced under restrict control of experts to avoid any flow and to refrain any complaint at any stage.

A group of very vast experienced persons in manufacturing handling our quality Control System.

Any instruments proven defective in workmanship or material will be replaced with out charges.

We express our thanks to all our customers for the confidence they have placed us in the past. We further assure them that we will continue our best efforts to gain their confidence in the future.

Certificates



QUALITY SERVICES & SATISFACTION

- * To provide "Excellent services and economically prices" is the motto of our business.*
- * All instruments are made of stainless steel other wise specified.*
- * All instrument are stamped with words "stainless" & "Pakistan any other country to whom manufactured buyers own brand may also be engraved on demand, "Die-Sunk" procedure is generally used.*
- * Special items can also be supplied if sketch, drawing or samples are provided to us or refer any reference from reputed catalogue for our guidance.*
- * All kind of Dental instruments of stainless steel made the specific material according to its use.*

Specification of material used in production.

Cutting Instruments		
For (scissors, Cutter)	Stainless steel AISI-410/420 (Rock-well hardness-48 to 58)	Hot drop forged
For (Forceps, Holding Clamp)	Stainless steel AISI-410/420 (Rock-well hardness-42 to 48)	Hot drop forged
For (Scalpel, Handles)	Stainless steel AISI-420 (Rock-well hardness-48 to 58)	Hot drop forged
For (Speculums, Retractors)	Stainless steel Sheet-300 Series 18/8 (18% Chromium & 8% Nickel)	

Please Note:

The delicate, hand-crafted dental instrument use only for its intended purpose to insure effectiveness and maximum life. Proper care must be given during cleaning, sharpening and sterilization to preserve design and integrity and to prolong usefulness. Do not expose instruments to sterilization temperatures above 350 F/177 C° or to higher than recommended chemical concentrations. "Longer Life for Dental Instrument" for proper cleaning, sterilization and other instrument care.

Inspected the instrument after each use for: pitting; corrosion; discoloration; broken, loose or bent components; cracked handles; dull cutting edges; over sharpened or worn out blades; chipped or peeled plating; damaged hinges, beaks or blades. Instruments found to Have these or other compromising conditions should be removed from service.

DIAGNOSTIC INSTRUMENTS



STETHOSCOPES & MAGILL FORCEPS


Stethoskope und Magillzangen | Estetoscopios y Pinzas Magill
Stéthoscopes et pinces Magill | Stetoscopuri si forcepsul Magill



**Duplex
Stethoscope
01-001**
super
lightweight



**Duplex-Baby
Stethoscope
01-002**
super
lightweight



**Anaestophon
Stethoscope
01-003**
Super
Lightweight



**Clar Mead Mirror
komplett
6 Volt
01-004**



01-007
Ø 50mm



**Pinard Wooden
Stethoscope
01-008** 18 cm



**Pinard Aluminium
Stethoscope
01-009** 15 cm



**Weber
Aesthesiometer
01-010** 9 cm



**Ziegler
01-005**
Ø 90mm
with head band without
protection cap



**Ziegler
01-006**
Ø 90mm
with head band with
protection cap



01-011
16cm
01-012
20cm
01-013
25cm

**Magill
Catheter
Intubating
Forceps**

Magill Catheter Intubating Forceps Gold
01-014
6.5" (2.5" shaft) , Infant , hole (dia 2.00 mm)
01-015
7.5" (3" shaft) , Infant , hole (dia 2.30 mm)
01-016
7.5" (4" shaft) , Infant , hole (dia 2.60 mm)



Magill Catheter Intubating Forceps GOLD
01-017
8" (4" shaft) , Child , hole (dia 3.00 mm)
01-018
8" (4" shaft) , Child , oval tip (5.38 mm X 4.15 mm)
01-019
9" (5" shaft) , Adult , oval (6.60 mm x 5.00 mm)
01-020
10" (5" shaft) , Adult , oval (9.60 mm x 7.70 mm)



PERCUSSION HAMMERS

Perkussionshämmer | Martillos de percusión
 Marteaux-percuteurs | Martelli percussori



Neuro-Duoflex Hammer
 01-021 22cm



Dejerine Hammer
 01-022 20cm



Dejerine Hammer
 01-023 20cm



Troemner Hammer
 01-024 145 g
 01-025 175 g
 01-026 235 g
 24.5 cm
 spare rubber



Troemner-Simplex Hammer
 01-027 21cm



Buck Hammer
 01-028 18.5cm



Buck Hammer
 01-029 18cm



Taylor Hammer
 01-030 18cm
 01-031 20cm



Taylor Hammer
 01-032 18cm



Babinsky Hammer
 01-033 22 cm



Babinsky Hammer
 01-034 22 cm



Rabiner Hammer
 01-035 25.5 cm
 with needle & brush



Berliner Hammer
 01-036 20 cm




Traube Hammer
 01-037 16 cm



Wartenberg Pinwheel
 01-038 19 cm

MC-INTOSH FIBER OPTIC LARYNGOSCOPE BLADE



BLADE

Size	A	B	C	Working Length
0	77mm	53mm	14mm	55mm
1	95mm	72mm	43mm	70mm
2	111mm	86mm	44mm	78mm
3	133mm	109mm	44mm	108mm
4	151mm	129mm	44mm	133mm

Mc-Intosh Laryngoscope Fiber Optic

Halogen illumination. Manufactured from AISI 304 Stainless-Steel. Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Green plastic case.

Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list

Mc-INTOSH FIBER OPTIC LARYNGOSCOPE SET HALOGEN ILLUMINATION

01-039 PEDIATRIC SET

Set of 0,1,2 Blades
Supplied in Green plastic case

01-040 STANDARD SET

Set of 1,2,3 Blades
Supplied in Green plastic case

01-041 ADULT SET

Set of 1,2,3,4 Blades
Supplied in Green plastic case



Laryngoscope Battery
Handle for Adult
C type 2 batteries
01-047



Mc-Intosh Fiber Optic
Laryngoscope Set
Halogen Illumination
Manufactured by AISI:
304 Stainless Steel Material

Fig. 0
01-042


Fig. 1
01-043

Fig. 2
01-044

Fig. 3
01-045

Fig. 4
01-046

MC-INTOSH LARYNGOSCOPE



BLADE

Size	A	B	C	Working Length
0	77mm	53mm	14mm	55mm
1	95mm	72mm	43mm	70mm
2	111mm	86mm	44mm	78mm
3	133mm	109mm	44mm	108mm
4	151mm	129mm	44mm	133mm

Fiber Optic Laryngoscope

Plastic case, Halogen illumination. Manufactured from AISI 304 Stainless Steel Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Green plastic case.

Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list

Mc-INTOSH FIBER OPTIC LARYNGOSCOPE SET HALOGEN ILLUMINATION

AVAILABLE SETS WITH BATTERY HANDLE

01-048 PEDIATRIC SET

Set of 0,1,2 Blades
Supplied in Green plastic case

01-049 STANDARD SET

Set of 1,2,3 Blades
Supplied in Green plastic case

01-050 ADULT SET

Set of 1,2,3,4 Blades
Supplied in Green plastic case



01-056
Laryngoscope Battery
Handle for Adult
C type 2 Batteries



Fig. 0
01-051



Fig. 1
01-052



Fig. 2
01-053



Fig. 3
01-054



Fig. 4
01-055



01-057
Laryngoscope
Bulb Small



01-058
Laryngoscope
Bulb Large


Mc-Intosh Laryngoscope Set
of Fiber Optic
Halogen Illumination
Manufactured From: AISI 304
Stainless Steel Material

LARYNGOSCOPES

Laryngoskope | Laringoscopios
Laryngoscopes | Laringoscopi



FIBER OPTIC LARYNGOSCOPE MILLER



BLADE

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm



01-067
Laryngoscope Battery
Handle for Adult
C type 2 Batteries

FIBER OPTIC LARYNGOSCOPE MILLER AVAILABLE SETS WITH BATTERY HANDLE

**01-059
PEDIATRIC SET**
Set of 0,1,2
Blades
Supplied in
Green plastic case

**01-060
STANDARD SET**
Set of 1,2,3
Blades
Supplied in
Green plastic case

**01-061
ADULT SET**
Set of 1,2,3,4
Blades
Supplied in
Green plastic case



Fig. 0
01-062



Fig. 1
01-063



Fig. 2
01-064




Fig. 3
01-065



Fig. 4
01-066

MC-INTOSH LARYNGOSCOPE



BLADE

Size	A	B	C	Working Length
0	80mm	58mm	30mm	55mm
1	97mm	76mm	36mm	70mm
2	111mm	88mm	36mm	78mm
3	133mm	112mm	36mm	108mm
4	154mm	131mm	38mm	133mm

Mc-Intosh Laryngoscope

Plastic case, Standard illumination. Manufactured from AISI 304 Stainless Steel Matt Type, Standard Illumination 2.5 V. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Gray plastic case.

Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list



01-077
With one Spare
Bulb



01-076
Laryngoscope Battery
Handle for Adult
C type 2 Batteries



Mc-Intosh Laryngoscope Set
Standard Illumination
Manufactured From
Stainless Steel Material

CONVENTIONAL LARYNGOSCOPE SET

Mc-INTOSH LARYNGOSCOPE SET STANDARD ILLUMINATION

AVAILABLE SETS WITH BATTERY HANDLE

01-068 PEDIATRIC SET
Set of 0,1,2 Blades
Supplied in Gray plastic case

01-069 STANDARD SET
Set of 1,2,3 Blades
Supplied in Gray plastic case

01-070 ADULT SET
Set of 1,2,3,4 Blades
Supplied in Gray plastic case



Fig. 0
01-071



Fig. 1
01-072



Fig. 2
01-073



Fig. 3
01-074




Fig. 4
01-075

LARYNGOSCOPES

Laryngoskope | Laringoscopios
Laryngoscopes | Laringoscopi

MC-INTOSH LARYNGOSCOPE



Size	BLADE			Working Length
	A	B	C	
0	80mm	58mm	30mm	55mm
1	97mm	76mm	36mm	70mm
2	111mm	88mm	36mm	78mm
3	133mm	112mm	36mm	108mm
4	154mm	131mm	38mm	133mm

Mc-Intosh Laryngoscope

Plastic case, Standard illumination. Manufactured from AISI 304 Stainless Steel. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Blue plastic case.

Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list



01-087
With two Spare Bulbs



01-086
Laryngoscope Battery Handle for Adult C type 2 batteries



CONVENTIONAL LARYNGOSCOPE SET

Mc-INTOSH LARYNGOSCOPE CONVENTIONAL STANDARD ILLUMINATION

AVAILABLE SETS WITH BATTERY HANDLE

01-078 PEDIATRIC SET

Set of 0,1,2 Blades.
Supplied in Blue plastic case

01-079 STANDARD SET

Set of 1,2,3 Blades.
Supplied in Blue plastic case


01-080 ADULT SET

Set of 1,2,3,4 Blades.
Supplied in Blue plastic case



Mc-Intosh Laryngoscope Conventional Standard Illumination
Manufactured From Stainless Steel material

LARYNGOSCOPE MILLER



Size	BLADE			Working Length
	A	B	C	
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm

Mc-Intosh Laryngoscope

Plastic case, Standard illumination. Manufactured from AISI 304 Stainless Steel. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Blue plastic case.

Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list



01-097
With two Spare Bulbs



01-096
Laryngoscope Battery Handle for Adult C type 2 batteries

CONVENTIONAL LARYNGOSCOPE SET

LARYNGOSCOPE SET OF CONVENTIONAL MILLER STANDARD ILLUMINATION

AVAILABLE SETS WITH BATTERY HANDLE

01-088 PEDIATRIC SET

Set of 0,1,2 Blades
Supplied in Blue plastic case

01-089 STANDARD SET

Set of 1,2,3 Blades
Supplied in Blue plastic case

01-090 ADULT SET

Set of 1,2,3,4 Blades
Supplied in Blue plastic case



Laryngoscope Set of conventional Miller Standard Illumination
Manufactured From Stainless Steel Material

LARYNGOSCOPES

Laryngoskope | Laringoscopios
Laryngoscopes | Laringoscopi



Laryngoscope Flexible (Mc-Coy)



Available blade sizes 2,3 & 4 with
'AA' & 'C' type battery handles
01-100
01-101
01-102

CONVENTIONAL SET
Flexible Laryngoscope
Conventional
Halogen Illumination
Manufactured From
Stainless Steel High
grade 304 material.

01-098

01-099
'C' type handle



Available Blade sizes 2,3 & 4 with
'AA' & 'C' type battery handles
01-105
01-106
01-107
01-103 FIBER OPTIC SET
Flexible Laryngoscope
Set of Fiber Optic
Halogen Illumination Manufactured
From Stainless Steel High grade 304
material packed in green plastic box.

Laryngoscope Battery
Handle for Adult
AA Batteries
01-104



DIAGNOSTIC PENTOSCOPE SET 01-108

2.5 V Standard illumination for use with
2 batteries type AA.

Consisting:

Universal Handle thread type fit for all
Pentoscope Heads like.

- Otoscope Head.
- Ophthalmoscope Head.
- Throat Branch Head.
- Mouth Mirrors Head.

It have two exploring mirrors No. 4 & 5 mm.

Its spare 06 disposable Plastic ear Canulas are
with sizes 2 x 2.50 mm, 2 x 3.50 mm & 2 x 4.5 mm.

On special order it can be as 2 x 1.5 mm for new
born babies 2 x 2.5 mm for pediatric and 2 x 4 mm
for adults, so it means total 06 canula also available
in a plastic tube packing.

Sturdy adjustment mount bracket

2 spare 2.5 V bulb.

Pentoscope set battery handle is operating with
two AA size batteries.

This Pentoscope set Packed in our new beautiful
Blue Colour plastic case in same blue inner packing material.
Spare parts also available of any item.



With 2 spare
Bulbs



Mouth Mirrors Head



Otoscope Head



Sturdy Adjustment Mount Bracket



Standard
illumination



Plastic ear Canulas

DIAGNOSTIC SETS

Diagnostik-Bestestecke
Trousse à diagnose

Juegos para diagnóstico
Astucci diagnostici



Two laryngeal mirrors: n. 4 and 5



Bent-arm illuminator: pin contact fitting.



Specula



With 2 spare Bulbs

Professional Otoscope & Ophthalmology Set ENT Medical Diagnostic Set 01-109

2.5 V Standard illumination for use with 2 batteries C type.

Including:

Otoscope Head: Luir Lock fitting.
Chrome plated.

Fix lens 4 x in plastic mount can be removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included,

Ophthalmoscope head: Luir Lock fitting.

Lens selection wheel permits choice of corrective lenses from 0 to - 20 and 0 to +20 diopters.

Single large circle aperture.

Bent-arm illuminator: pin contact fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

Diagnostic set is supplied in our new beautiful look and quality Blue Case, which protect our quality Products from moisturizing and damaging.



Otoscope Head



Nasal speculum: Chrome plated.



Ophthalmoscope head



ENT Specula

DIAGNOSTIC SETS

Diagnostik-Bestestecke
Trousses à diagnose

Juegos para diagnóstico
Astucci diagnostici



Two laryngeal mirrors: n. 4 and 5



Bent-arm illuminator: pin contact fitting.



Nasal speculum: Chrome plated.



Handle



ENT Specula



Otoscope Head



Ophthalmoscope Head

Universal Oscope & Ophthalmology Set ENT Medical Diagnostic Set 01-110

2.5 V Standard illumination for use with 2 batteries type C.

Including:

Otoscope Head: Luir Lock fitting. Chrome plated.

Fix lens 4 x in plastic mount can be removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included.

Ophthalmoscope head: pin contact fitting.

Lens selection wheel permits choice of corrective lenses from 0 to - 20 and 0 to +20 diopters.

Single large circle operture.

Bent-arm illuminator: Luir Lock fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

The Diagnostic set is supplied in new beautiful look and quality Blue Case, which protect our quality Products from moisturizing and damaging.

DIAGNOSTIC SETS

Diagnostik-Bestestecke
Trousse à diagnose

Juegos para diagnóstico
Astucci diagnostici



Standard Otoscope + Ophthalmoscope Set
01-111

This Set is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Resuable speculas with sizes 2.5, 3.5 & 4.5 mm. A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image A Full Functional Ophthalmoscope Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Plastic Box Packing.



Otoscope
Head



ENT Specula



Ophthalmoscope
head



ENT Specula



'C' type
Battery handle



Standard Ophthalmoscope Set
01-112

A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image Packed in Black Plastic Case.



Ophthalmoscope
head



'C' type Battery Handle



Standard Otoscope Set
01-113

this otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Resuable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing



DIAGNOSTIC SETS

Diagnostik-Bestestecke
Trousses à diagnose

Juegos para diagnóstico
Astucci diagnostici



Standard Otoscope Set 01-114

this otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing



Otoscope + Ophthalmoscope Set 01-115

This otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image A Full Functional Ophthalmoscope Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing.




Standard Ophthalmoscope 01-116

A Full Functional Ophthalmoscope. compact durable, Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image Handle made of Stainless Steel Material.



Mini Ophthalmoscope 01-117

A Full Functional Ophthalmoscope. compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image Made of plastic Material.



Mini Ophthalmoscope 01-118

A Full Functional Ophthalmoscope. compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image made of plastic Material.

DISSECTING KITS

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



01-119

- Dissecting Set of 6 pcs.
- Dropper Oil
 - Probe Plastic Handle
 - Dressing Forceps
 - Dressing Scissors S/s
 - Hook
 - Scalpel Handle 4



01-120

- Dissecting Set of 12 pcs.
- Iris Scissors
 - Scalpel Knives 4
 - Probe Single Ended
 - Scalpel Handle 3
 - Dressing Scissors S/s
 - Dressing Scissors S/b
 - Dressing Forceps
 - College Tweezer
 - Probe Straight Plastic Handle
 - Probe Angular Plastic Handle
 - Dropper Owl
 - Hook



01-121

New Instructor's Dissecting set

Contain quality stainless steel Instrument Enclosed in an attractive zippered case Scissor surgical 4 1/2" scissors fine point 4 1/2" Forceps fine point Str . Forcep fine point Cvd . Mosquito forceps Cvd 5" box lock . One piece scalpel 1 1/2" blade Surgeon handle # 4 . Mall probe and seeker 2 plastic handle str needle also 2 blade # 22 for each set.



01-122

Dissecting Set- Advanced

Stainless steel, consisting of 12 instruments for accurate dissection. Supplied in zipper leather case.

Kit consists of:

- Scissors - 2 pcs. Forceps - 2 pcs. Needle - 2 pcs.
Artery Forceps - 1 pc. Tissue Forceps - 1 pc.
Seeker - 1 pc. Scalpel Handle w/ Blade - 1 pc.
Bone Cutter - 1 pc. Dental Probe & Seeker - 1 pc



01-123

Dissecting Kit

Kit includes:

- 1 4 1/2" Dissecting forceps with curved and serrated tips
1 Curved dissecting needle 1 Glass pipette with rubber bulb
1 4 1/2" Dissecting forceps straight and serrated tips
1 Dissecting scalpel, 6" long, 1 1/2" cutting edge.
1 4 1/2" Dissecting scissors, sharp, straight points
1 Plastic ruler 6", 1/16" div., millimeter and centimeter



01-124

Dissecting Set - Pro

Stainless steel, consisting of 20 instruments for accurate dissection. Supplied in zipper leather case.

Kit consists of:

- Scissors - 2 pcs. Forceps - 2 pcs. Needles - 4 pcs. Scalpel - 1 pc.
B.P. Handle w/ Blade - 1 pc. Razor - 1 pc. Ruler 6" - 1 pc.
Magnifying Glass - 1 pc. Glass Slide - 1 pc. Chain and Hook - 1 pc.
Cover Slips - 1 pc. Blow Pipe - 1 pc. Needle Wooden Handle - 1 pc.
Dropper Glass - 1 pc. Brush - 1 pc.



01-125

Dissecting Set - Popular

Stainless steel, consisting of 14 instruments for accurate dissection. Complete in zipper leather case.

Kit consists of:

- Scissors - 2 pcs. Forceps - 2 pcs. Needles - 4 pcs. Scalpel - 1 pc.
B.P. Handle w/ Blade - 1 pc. Razor - 1 pc. Magnifying Glass - 1 pc.
Chain and Hook - 1 pc. Blow Pipe - 1 pc.



01-126

Dissecting Set - Senior

Stainless steel, consisting of 9 instruments for accurate dissection. Supplied in zipper leather case.

Kit consists of:

- Scissors - 2 pcs. Forceps - 2 pcs. Needle - 2 pcs. Scalpel - 1 pc.
Scalpel Handle w/ Blade - 1 pc. Ruler 6" - 1 pc.

DISSECTING KITS



Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



01-127
Dissecting Set



01-128
Dissecting Set



01-129
Dissecting Set



01-130
Dissecting Set



01-131
Dissecting Set



01-132
Dissecting Set

MEASURING INSTRUMENTS

Messinstrumente

Instruments pour mesurer

Instrumentos de medidas

Strumenti di misura



Jameson Caliper
01-133 8 cm



Townley Caliper
01-134 10 cm



Vernier Caliper
01-135 15 cm



Vernier Caliper
01-136 15 cm



**Carroll
Pocket Goniometer**
01-137 14.5cm



Möeltgen Flexometer
01-138 24cm



**Bertillon Measuring
Device**
01-139 32 cm



Castroviejo-Epker Caliper
01-140 18 cm
Measuring range 40 mm, straight



Castroviejo-Epker Caliper
01-141 18 cm Angled
Measuring range 40 mm



**Castroviejo
Bone Caliper Straight**
01-142 18cm



Castroviejo
01-143 18cm
Bone Caliper Curved



Castroviejo-Epker Caliper
01-144 9 cm Str.
Measuring range, 40 mm
01-145 9 cm Cvd.
Measuring range, 40 mm



Castroviejo
01-146 8 cm, Str.
Eye View Caliper
Measuring range, 20 mm
01-147 8 cm, Cvd.
Measuring range, 20 mm



Castroviejo
01-148 8 cm
Eye View Caliper
Measuring range
20 mm



Tonometer
01-149
Set



TROCARS, SUCTION TUBES, NEEDLE AND PROBES



TROCARS

Trokare | Trócares
Trocarts | Trequarti



Cap Only



Fig. 1 02-03 1.7 mm Ø
Fig. 2 02-04 3 mm Ø
Fig. 3 02-05 4.5 mm Ø
Fig. 4 02-06 6 mm Ø

Universal Abdominal Trocar Set
02-01
(Handle with 4 needles) 14 cm
Fig. 1.7 mm Ø to 6 mm Ø
02-02
Only Handle



Universal Abdominal Trocar Set
02-07
(Handle with 4 needles) 14 cm
Fig. 1.7 mm Ø to 6 mm Ø
02-08
Only Handle




Fig. 1 02-09 1.7 mm Ø
Fig. 2 02-10 3 mm Ø
Fig. 3 02-11 4.5 mm Ø
Fig. 4 02-12 6 mm Ø



Uckermann Trocar
02-14 13cm



Sinuscopy Trocar
02-15 17cm



3.5mm Ø 02-15
5.0mm Ø 02-16
5.0mm Ø 02-17



Nelson Trocar Thorax
02-18 21cm



02-18
6.0mmØ



02-19
8.3mmØ



02-20
9.5mmØ



02-21
11.6mmØ



Luer, Trocar
02-13
Luer Lock Adapter

TROCARS



Trokare Trócares
Trocarts Trequarti



TROCARS

Trokare | Trócares
Trocarts | Trequarti



Buelau, Trocar
02-45 8 mm
24 cm



Fleurant Trocar
02-46 23cm
5 mm Ø



Fleurant Trocar
02-47 23cm
5mm Ø




Fleurant Trocar
02-48 23cm
5 mm Ø



Pool
02-49 22.5 cm
Suction cannula with
outer tube to be screwed

SUCTION TUBES




Pool
02-50 22.5 cm
Suction cannula with
outer tube to be screwed



02-51 18 cm
for tonsillectomy,
oral and laryngeal
operations



02-52 22 cm
Suction Cannula for
tonsillectomy, oral
and laryngeal operations



02-53 18 cm
suction cannula for
tonsillectomy, oral
and laryngeal operations



02-54 17 cm
S-shaped blunt cannula
for aspirating fluid from
the excretory ducts of
large glands and for
keeping the operations
area clear

SUCTION TUBES



Saugrohre

Tubos de aspiracion

Canules d'aspiration

Tubi per aspirazione



Cooley, Suction Tube
02-55 18 cm
6mm Ø



Cooley, Suction Tube
02-56 18 cm
7mm Ø



Cooley, Suction Tube
02-57 18 cm
9mm Ø



Cooley, Suction Tube
02-58 32.7 cm
9mm Ø Curved
Ø 1.5 mm



02-60

02-61

Andrews-Pyncheon
Suction Tube
02-59 23.5 cm
Ø 2.0 mm



Yankauer
Suction Tube
02-62
31 cm, 15 cm
2.5 cm, 10 cm



Yankauer
Suction Tube
02-63
31 cm, 15 cm
2.5 cm, 10 cm



De Bakey
suction tube
02-64
27.5 cm



Yankauer Suction Tube
02-65 27 cm
Baby Yankauer
02-66 21 cm



Coupland
Suction Tube
02-67



Byrd, Suction Tube
02-68 1.5 mm
02-69 2.5mm
02-70 3.5 mm



Adson
suction tubes
02-71
4 mm Ø
15 cm



Adson
suction tubes
02-72
21 cm
4 mm Ø 12 Charr./Fr
02-73
5 mm Ø 15 Charr./Fr



02-74
Wullstein Suction
Tube Complete



Gorney
Suction Tubes
02-75
18 cm
suction dissector



Landolt Dissecting
Suction Tubes
02-76 4 mm wide
02-77 5 mm wide
20 cm
suction raspatories, sharp,



Lopez-Reinke
Suction Tubes
02-78
23.5 cm
tonsil dissector with
suction device

SUCTION TUBES

Saugrohre | Tubos de aspiracion
Canules d'aspiration | Tubi per aspirazione



Plester Daumenplatt

02-79 6 Charr./Fr.
02-80 7 Charr./Fr.
02-81 8 Charr./Fr.
02-82 9 Charr./Fr.
02-83 10 Charr./Fr.
02-84 11 Charr./Fr.
02-85 12 Charr./Fr.
02-86 15 Charr./Fr.

Baron suction tubes 7.5 cm

02-87 1 mm Ø 3 Charr./Fr.
02-88 1 2/3 mm Ø 5 Charr./Fr.
02-89 2 1/3 mm Ø 7 Charr./Fr.

Baron suction tubes 10 cm

02-90 1 mm Ø 3 Charr./Fr.
02-91 1 2/3 mm Ø 5 Charr./Fr.
02-92 2 1/3 mm Ø 7 Charr./Fr.

House suction tubes 7.5 cm

02-93 1 mm Ø 3 Charr./Fr.
02-94 1 2/3 mm Ø 5 Charr./Fr.
02-95 2 1/3 mm Ø 7 Charr./Fr.
02-96 2 2/3 mm Ø 8 Charr./Fr.
02-97 3 1/3 mm Ø 10 Charr./Fr.
02-98 4 mm Ø 12 Charr./Fr.
02-99 4 2/3 mm Ø 14 Charr./Fr.

House suction tubes 10 cm

02-100 1 mm Ø 3 Charr./Fr.
02-101 1 2/3 mm Ø 5 Charr./Fr.
02-102 2 1/3 mm Ø 7 Charr./Fr.
02-103 2 2/3 mm Ø 8 Charr./Fr.
02-104 3 1/3 mm Ø 10 Charr./Fr.
02-105 4 mm Ø 12 Charr./Fr.
02-106 4 2/3 mm Ø 14 Charr./Fr.

Plester Suction Tubes 17 cm

02-107 2.0 mm Ø 6 Charr./Fr.
02-108 2.3 mm Ø 7 Charr./Fr.
02-109 2.5 mm Ø 8 Charr./Fr.
02-110 3.0 mm Ø 9 Charr./Fr.
02-111 3.3 mm Ø 10 Charr./Fr.
02-112 3.5 mm Ø 11 Charr./Fr.
02-113 4.5 mm Ø 12 Charr./Fr.
02-114 5.0 mm Ø 15 Charr./Fr.


Fergusson Suction Tubes 17 cm

02-115 2.0 mm Ø 6 Charr./Fr.
02-116 2.3 mm Ø 7 Charr./Fr.
02-117 2.5 mm Ø 8 Charr./Fr.
02-118 3.0 mm Ø 9 Charr./Fr.
02-119 3.3 mm Ø 10 Charr./Fr.
02-120 3.5 mm Ø 11 Charr./Fr.
02-121 4.5 mm Ø 12 Charr./Fr.
02-122 5.0 mm Ø 15 Charr./Fr.

FIMA INSTRUMENTS

FRAZIER suction tubes

02-123 6 Fr.
02-124 7 Fr.
02-125 8 Fr.
02-126 9 Fr.
02-127 10 Fr.
02-128 11 Fr.
02-129 12 Fr.
02-130 13 Fr.
02-131 15 Fr.
19 cm



FRAZIER 30° suction tubes

02-132 6 Fr.
02-133 7 Fr.
02-134 8 Fr.
02-135 9 Fr.
02-136 10 Fr.
02-137 11 Fr.
02-138 12 Fr.
02-139 13 Fr.
02-140 15 Fr.
18 cm



FRAZIER 45° suction tubes

02-141 6 Fr.
02-142 7 Fr.
02-143 8 Fr.
02-144 9 Fr.
02-145 10 Fr.
02-146 11 Fr.
02-147 12 Fr.
18 cm



Cannula

02-148 6fr
02-149 7fr
02-150 8fr
02-151 9fr
02-152 10fr
02-153 11fr
02-154 12fr
02-155 13fr
02-156 14fr



SUCTION AND IRRIGATION CANNULAS



FIMA INSTRUMENTS

Saug-und Spülkanülen

Canulas para aspiracion y lavado

Canules pour aspiration

Cannule per aspirazione e lavaggio



Luer
02-157
Suction Tube working
length 110 mm



Nasal Suction Tubes
02-158 Ø 4.5 mm
02-159 Ø 5.5 mm
17.5 cm



Eicken
Suction &
Irrigation Cannula
02-160 2.5 mm Ø
02-161 3.0 mm Ø
Slightly curved



Eicken
Suction &
Irrigation Cannula
02-162 2.5 mm Ø
02-163 3.0 mm Ø
Strongly curved



Tonsil Needle
Suction &
Irrigation Cannula
02-164 0.8 mm Ø
Straight



Tonsil Needle
Suction &
Irrigation Cannula
02-165 0.8 mm Ø
Curved



Septum Needle
Suction &
Irrigation Cannula
02-166 0.8 mm Ø
bayonet




Septum Needle
Suction &
Irrigation Cannula
02-167 0.8 mm Ø
angled on flat



Luer-Lock Suction &
Irrigation Cannula
02-168 3.0 mm Ø



Luer-Lock Suction &
Irrigation Cannula
02-169 2.5 mm Ø



Fisch Suction &
Irrigation Cannula
02-170 4.0 mm Ø

Probe, Double Ended
1.0 mm Ø
02-171 11.5cm
02-172 13.0cm
02-173 14.5cm
02-174 16.0cm
02-175 18.0cm
02-176 20.0cm
02-177 25.0cm

Probe, Double Ended 1.5 mm Ø
02-178 11.5cm
02-179 13.0cm
02-180 14.5cm
02-181 16.0cm
02-182 18.0cm
02-183 20.0cm
02-184 25.0cm

Probe, Double Ended 2.0 mm Ø
02-185 11.5cm
02-186 13.0cm
02-187 14.5cm
02-188 16.0cm
02-189 18.0cm
02-190 20.0cm
02-191 25.0cm



Probe With Eye
2.0mm Ø
02-192 11.5cm
02-193 13.0cm
02-194 14.5cm
02-195 16.0cm
02-196 18.0cm
02-197 20.0cm
02-198 25.0cm
with eyelet



Probe
2.0mm Ø
02-199 11.5cm
02-200 13.0cm
02-201 14.5cm
02-202 16.0cm
myrtle-shaped

EUSTACHIAN CATHETERS

Ohrkatheter

Sondes auriculaires

Sondas auriculares

Sonde auricolari



Hartmann, Eustachian Catheters
Fig. 1

02-203 2.0 mm Ø
02-204 2.5 mm Ø
02-205 3.0 mm Ø



Hartmann Eustachian Catheters
Fig. 2

02-206 2.0 mm Ø
02-207 2.5 mm Ø
02-208 3.0 mm Ø



Killian, Eustachian Catheters
02-209 16 cm 2.5 mm Ø
02-210 16 cm 3.0 mm Ø
oval opening



Hartmann Eustachian Catheters
15cm



Fig. 0
02-211



Fig. 1
02-212



Fig. 2
02-213



Fig. 3
02-214



Fig. 4
02-215



Fig. 5
02-216



Tröltsch Eustachian Catheters
15cm



02-217
Tube end piece
for eustachian
catheter



Fig. 1
02-218



Fig. 2
02-219



Fig. 3
02-220



Fig. 4
02-221



Fig. 5
02-222



Zöllner, Suction Tube
02-223 0.7mm Ø
complete with 1 slip-on
cannula, 0.7mm Ø



Iowa, Suction Tube
02-224 0.7mm Ø
16cm



Rosen Suction Tube
X=80 mm
02-225 0.7mm Ø
02-226 0.9mm Ø
02-227 1.2mm Ø




ROSEN, Suction Tube
X=65 mm
02-228 0.7mm Ø
02-229 0.9mm Ø
02-230 1.2mm Ø



Magill Suction Tube
24.5cm
02-231 3.0mm Ø
02-232 3.5mm Ø
02-233 4.0mm Ø



Suction Tube, Guide Needle
19.5cm
02-234 2.2/3mmØ 8 Charr./Fr.
02-235 3.1/3mmØ 10 Charr./Fr.
02-236 4mmØ 12 Charr./Fr.
02-237 4 2/3mmØ 14 Charr./Fr.
02-238 5 1/3mmØ 16 Charr./Fr.
02-239 6mmØ 18 Charr./Fr.



Suction Tube, Guide Needle
19.5cm
02-240 2mmØ 6 Charr./Fr.
02-241 2.2/3mmØ 8 Charr./Fr.
02-242 3 1/3mmØ 10 Charr./Fr.
02-243 4mmØ 12 Charr./Fr.
02-244 4 2/3mmØ 14 Charr./Fr.
02-245 5 1/3mmØ 16 Charr./Fr.
02-246 6 1/3mmØ 18 Charr./Fr.



Suction Tube, Guide Needle
19.5cm
02-247 2.2/3mmØ 8 Charr./Fr.
02-248 3.1/3mmØ 10 Charr./Fr.
02-249 4mmØ 12 Charr./Fr.
02-250 4 2/3mmØ 14 Charr./Fr.
02-251 5 1/3mmØ 16 Charr./Fr.
02-252 6mmØ 18 Charr./Fr.

PROBES AND DIRECTORS



Sonden
Sondes et stylets

Sondas y estiletes
Specilli e sonde

Butterfly Grooved Director With Tongue Tie Plain End

02-253 11.5cm
02-254 13.0cm
02-255 14.5cm
02-256 16.0cm
02-257 18.0cm
02-258 20.0cm



Butterfly Grooved Director With Tongue Tie Plain End

02-259 11.5cm
02-260 13.0cm
02-261 14.5cm
02-262 16.0cm
02-263 18.0cm
02-264 20.0cm



Doyen Probes & Directors 02-265 14cm



Nélaton Probes & Directors 02-266 16cm



Directors 02-267 13cm 02-268 14.5cm 02-269 16cm Larry Probe and Grooved Director, 17.8 cm



Doyen 02-270 14.5cm



Nélaton 02-271 16cm



Director Probe End 02-272 14cm Director with tip Diam. of tip 3.3 mm



Brodie 02-273 20cm Director with tip Diam. of tip 3.5 mm




Lockhart- Mummery 02-274 16.5cm



Lockhart- Mummery 02-275 16.5cm



Lockhart- Mummery 02-276 16.5cm



Lockhart- Mummery 02-277 16.5cm



Director with wire loop, malleable 02-278 19cm





FIMA INSTRUMENTS

PROBES AND DIRECTORS

Sonden

Sondas y estiletes

Sondes et stylets

Specilli e sonde


COTTON APPLICATORS**Jopson-Horne
Cotton Applicator**
02-279 18cm**Fein Cotton
Applicator**
02-280 28cm**Pelterson Cotton
Applicator**
02-281 23cm**Barr Fistula Probe**
02-282
with 20.3cm sterling shaft,
eye at tip, 29.2cm overall**Barr Crypt Hook**
02-283
25.4cm short and
long hook ends**Buie Suction Tube**
02-284
45 cm with finger
valve, 15 Fr. (5 mm)**Bowman Lacrimal Probes Double Ended, Sterling Silver****02-285**
Bowman Probe,
000-0000**02-286**
Bowman Probe,
00-0**02-287**
Bowman Probe,
1-2**02-288**
Bowman Probe,
3-4**02-289**
Bowman Probe,
5-6**02-290**
Bowman Probe,
7-8**William Lacrimal Probes Double Ended, Sterling Silver****02-291**
William Probe,
000-0000**02-292**
William Probe,
00-0**02-293**
William Probe,
1-2**02-294**
William Probe,
3-4**02-295**
William Probe,
5-6**02-296**
William Probe,
7-8

PROBES AND DIRECTORS



Sonden
Sondes et stylets

Sondas y estiletes
Specilli e sonde



SCALPEL HANDLE, SCALPEL BLADES, DISSECTING KNIVES & SAWS



SCALPEL HANDLES



Skalpellsgriffe

Mangos de bisturís

Manches pour bistouris

Manici per bisturi



03-01
#3 Scalpel Handle
12.50 cm
Blade Size 10, 11,
12, 13, 15, 17



03-02
#3 Scalpel Handle/
Graduated 12.75 cm
Blade Size 10, 11,
12, 13, 15, 17



03-03
Scalpel Handle
#4 Graduated
with scale



03-04
#4 Scalpel Handle
12.75 cm Blade
Size 18, 20, 21, 22,
23, 24, 25, 60,70



03-05
#4 Scalpel Handle/
Graduated 12.75 cm
Blade Size 18, 20, 21,
22, 23, 24, 25, 60,70



03-06
4LA Scalpel Handle
21 cm Blade
Size 18, 20, 21, 22,
23, 24, 25, 60,70



03-07
#3, #4 Scalpel
Handle
Double Ended
16 cm



03-08
Scalpel Handle
3L Long 21 cm
Blade Size
10, 11, 12, 15, 17



03-09
#7 Scalpel
Handle
16.5 cm
Blade 10, 11,
12, 13, 15, 17



03-10
Fig. 10 - 17
12.5cm
Solid



03-11
Scalpel Handle
4L Long 21 cm
Blade
Size 18, 20, 21,
22, 23, 24, 25,
60,70



03-12
#3LA Scalpel
Handle / Long
Angled 20.25 cm
Blade Size
10, 11, 12, 15, 17



03-13
Knife Handle
(Smooth) for
Blades Nos.
60 and 70



SCALPEL HANDLES

Skalpеллgriffe
Mangos de bisturíes

Manches pour bistouris
Manici per bisturi



03-14
Siegel Round
Scalpel Handle
#3 With Adjustable
Tip 16 cm
Blade Size 10,
11, 12, 13, 15,
17



03-15
Siegel Round
Scalpel Handle
#4 With Adjustable
Tip 16 cm
Blade Size 18, 20,
21, 22, 23, 24,
25, 60, 70



03-16
Round Knurled
Scalpel Handle
11.5 cm
Blade Size
50, 60 & 70



03-17
Round Knurled
Scalpel Handle
12.7 cm
Blade Size
50, 60 & 70



03-18
Round Knurled
Scalpel Handle
12.7 cm
Blade Size
50, 60 & 70



03-19
Round Knurled
Scalpel Handle
11 cm
Blade Size
50, 60 & 70



03-20
Round Knurled
Scalpel Handle
13.5 cm
Blade Size
50, 60 & 70



03-21
Beaver
Handle



03-22
Round Knurled
Scalpel Handle
12 cm



03-23
#3RK Kaye
Round Scalpel
Handle 15 cm,
Blade
Size 10, 11,
12, 13, 15, 17



03-24
#3R Corwin
Round Scalpel
Handle
14.5 cm
Blade Size 10,
11, 12, 15, 17



03-25
#4R Corwin Round
Scalpel Handle,
14.5 cm
Blade Size 18, 20,
21, 22, 23, 24,
25, 60, 70



03-26
#3RC Scalpel
Handle 14 cm,
Blade
Size 10, 11, 12,
13, 15, 17



03-27
#3P Scalpel
Handle 15.5 cm,
Blade Size
10, 11, 12, 13,
15, 17



03-28
15.5cm
Blade Size
10 - 17
With Round
Hollow Handle



03-29
Scalpel
Handle



03-30
Scalpel Grip

SCALPEL HANDLES, BLADE BREAKERS



Skalpelligriffe, Klingenbrecher

Mangos de bisturí, rompe cuchillas

Poignées de scalpel, disjoncteurs de lame

Manici per bisturi, rompitrama



03-31
Round Scalpel
Handle,
Double Head #3,
Blade Size 10, 11, 12,
13, 15, 17



**Bever Handle
03-32**



03-33
Scalpel Plastic
Handle #8 Blade
Size 18, 20, 21,
22, 23, 24, 25,
60,70



03-34
Collin Scalpel
Handle # 9 Blade
Size 10, 11, 12,
13, 15, 17



Fig. 10



Fig. 11

Blade Sizes



Fig. 12



Fig. 15



03-35
Scalpel Plastic
Handle



03-36
Scalpel Plastic
Handle # 5
Blade Size
10, 11, 12,
13, 15, 17



03-37
Scalpel Plastic
Handle # 6
Blade Size 18,
20, 21, 22, 23,
24, 25, 60,70



Fig. 20



Fig. 22

RAZOR BLADE HOLDER & BREAKERS



03-38
Cole-Parmer Dissecting
Blade Remover For
Blade Handle



03-39
Troutman-Chris
Blade Breaker &
Holder with Push-
Pull Lock 9 cm



03-40
Troutman-Chris
Blade Breaker &
Holder angled, with
Push-Pull Lock
9 cm



03-41
Blade Holder Blade
Breaker & Holder,
with Ratchet Lock
11 cm



03-42
Castroviejo Blade
Breaker & Holder,
with Push-Pull Lock
13 cm



03-43
Razor Blade
Breaker & Holder,
Straight
18.5 cm

SCALPEL BLADES

Skalpellsklngen,

Lames de bistouris, Mangos

Hojas de bisturi, Manches

Lame per bisturi, Manici



03-44
Sterile Scalpel
Blade Size # 6
Stainless Steel



03-45
Sterile Scalpel
Blade Size # 10
Stainless Steel



03-46
Sterile Scalpel
Blade Size # 10A
Stainless Steel



03-47
Sterile Scalpel
Blade Size # 11
Stainless Steel



03-48
Sterile Scalpel
Blade Size # 11P
Stainless Steel



03-49
Sterile Scalpel
Blade Size # E11
Stainless Steel



03-50
Sterile Scalpel
Blade Size # 12
Stainless Steel



03-51
Sterile Scalpel
Blade Size # 13
Stainless Steel



03-52
Sterile Scalpel
Blade Size # 15
Stainless Steel



03-53
Sterile Scalpel
Blade Size # 15A
Stainless Steel



03-54
Sterile Scalpel
Blade Size # 17
Stainless Steel



03-55
Sterile Scalpel
Blade Size # 18
Stainless Steel



03-56
Sterile Scalpel
Blade Size # 20
Stainless Steel



03-57
Sterile Scalpel
Blade Size # 21
Stainless Steel



03-58
Sterile Scalpel
Blade Size # 22
Stainless Steel



03-59
Sterile Scalpel
Blade Size # 23
Stainless Steel



03-60
Sterile Scalpel
Blade Size # 24
Stainless Steel



03-61
Sterile Scalpel
Blade Size # 25
Stainless Steel



03-62
Sterile Scalpel
Blade Size # 36
Stainless Steel



03-63
Sterile Scalpel
Blade Size # 40
Stainless Steel



03-64
Sterile Scalpel
Blade Size # 60
Stainless Steel



03-65
Sterile Scalpel
Blade Size # 70
Stainless Steel

SCALPEL KNIVES



Operations-Skalpelle

Bisturios para operaciones

Bistouris pour opérations

Bisturi per operazioni



DIEFFENBACH
16 cm



03-66
Fig. 1
16 cm

03-67
Fig. 2
16 cm

03-68
Fig. 3
16 cm

03-69
Fig. 4
16 cm

03-70
Fig. 5
16 cm

03-71
Fig. 6
16 cm

03-72
Fig. 7
16 cm

03-73
Fig. 8
16 cm

03-74
Fig. 9
16 cm

03-75
Fig. 10
16 cm



03-76
Fig. 15
16 cm



03-77
Fig. 16
16 cm



**Collin Dissecting
Knives**
14.5 cm



03-78
Fig. 1
14.5 cm

03-79
Fig. 2
14.5 cm

03-80
Fig. 3
14.5 cm

03-84
Fig. 1
13.5 cm

03-85
Fig. 2
13.5 cm

03-86
Fig. 3
13.5 cm

03-87
Fig. 4
13.5 cm

03-88
Fig. 5
13.5 cm

03-89
Fig. 6
13.5 cm

03-81
Fig. 4
14.5 cm

03-82
Fig. 5
14.5 cm

03-83
Fig. 6
14.5 cm

**Dieffenbach
Dissecting Knives**
13.5 cm

03-90
Fig. 7
13.5 cm

03-91
Fig. 8
13.5 cm

03-92
Fig. 9
13.5 cm

03-93
Fig. 10
13.5 cm

SCALPEL KNIVES

Operations-Skalpelle


Bistouris pour operations

Bisturios para operaciones

Bisturi per operazioni



Dissecting Knives
17 cm



03-94
Fig. 1
17 cm

03-95
Fig. 3
17 cm

03-96
Fig. 5
17 cm

03-97
Fig. 7
17 cm

03-98
Fig. 9
17 cm

03-99
Fig. 10
17 cm

03-100
Fig. 12
17 cm

03-101
Fig. 14
17 cm

03-102
Fig. 16
17 cm

03-103
Fig. 18
17 cm

DISSECTING KNIVES



DISSECTING KNIFE
15 cm
with metal handle



03-104
15 cm
30 mm

03-105
15 cm
36 mm

03-106
15.5 cm
40 mm

03-107
16 cm
45 mm

03-108
16 cm
50 mm

03-109
15 cm
27 mm

03-110
15 cm
35 mm

03-111
16 cm
40 mm

03-112
16 cm
47 mm

03-113
17 cm
50 mm

DISSECTING KNIVES



Anatomische skalpelle
Bisturíes para disección

Bistouris à dissection
Bisturi para dissezione



**Collin Dissecting
Knives Blade Size**
03-114 13 cm
03-115 16 cm
03-116 19 cm
03-117 22 cm
03-118 28 cm



**Catlin Dissecting
Knives Blade Size**
03-119 16 cm
03-120 19 cm
03-121 22 cm



**Liston Dissecting
Knives Blade Size**
03-122 16 cm
03-123 19 cm
03-124 22 cm



**Langenbeck Dissecting
Knives Blade**
03-125 12 cm



**Amputation
Knife**
03-126 29 cm



**Amputation
Knife**
03-127 32.5 cm



**Catlin Amputation
Knife**
03-128 29 cm



**Resection
Knife**
03-129 18 cm



**Langenbeck
Amputation**
03-130 18 cm



**Langenbeck
Amputation**
03-131 18 cm



**Dieffenbach
Dissecting Knives**
03-132 13 cm



**Dieffenbach
Dissecting Knives**
03-133 13 cm



**Simon Dissecting
Knives**
23.5 cm



03-134
Fig. 1



03-135
Fig. 2



03-136
Fig. 3



03-137
Fig. 4



FIMA INSTRUMENTS

DISSECTING KNIVES

Anatomische skalpelle
Bisturíes para disección

Bistouris à dissection
Bisturi para dissezione



Autopsy Knife
03-138
25.5 cm



Walb Dissecting
Knives
03-139
11 cm



Walb Dissecting
Knives
03-140
14 cm



Walb Dissecting
Knives
03-141
17 cm



Virchow Dissecting
Knives
03-142
16 cm



Virchow Dissecting
Knives
03-143
20 cm



Virchow Dissecting
Knives
03-144
24 cm



Ayre Cone Knife
03-145
23.5 cm



Cone Knife
03-146
23 cm



Segond
Myomatome
03-147
27 cm



Olivecrona
Trigeminal knife
03-148
20.5 cm



Trigeminal knife
03-149
23 cm



Fowler
Tonsil Knife
03-150
20 cm



Abraham
Tonsil Knife
03-151
21.5 cm



Brunings
Tonsil knife
03-152
21.5 cm

DERMATOMES - TRANSPLANTATION KNIVES

Dermatome - Transplantationsmesser


Dermatomos - Cuchillos para transplantaciones

Dermatomes - Bistouries pour transplantation

Dermatomi - coltello per trapianto



AMPUTATION KNIVES & RESECTION SAWS



RAZOR HOLLOW
Blades On
Both Sides
Blade Flat / Hollow
03-153



**Charriere Amputation
& Resection**
03-155 27cm
03-156 30cm



**Charriere Amputation
& Resection**
03-157 35cm



**Satterlee Amputation
& Resection**
03-158 31cm



Thiersch
03-159
20 cm



Schink
03-160 30cm
Complete with 1
Blade & 2 plates



03-161
Spare Blade



03-162
Skin fixing plate



Humby Dermatomes
03-163 31cm
Complete with 1 Blade



03-164
spare blade



Watson Dermatomes
03-165 30cm
Complete with 1 Blade,
fixed bar



Silver Dermatomes
03-166 19cm

GENERAL OPERATING SCISSORS



OPERATING SCISSORS

Chirurgische scheren | Ciseaux chirurgicaux
Tijeras quirúrgicas | Forbici per chirurgia



Operating Scissors (Flat Pattern)



Operating Scissors

04-001	11.5 cm	S/S Str.
04-002	12.5 cm	S/S Str.
04-003	13 cm	S/S Str.
04-004	14.5 cm	S/S Str.
04-005	16 cm	S/S Str.
04-006	18 cm	S/S Str.
04-007	20 cm	S/S Str.
04-008	23 cm	S/S Str.
04-009	25 cm	S/S Str.
04-010	28 cm	S/S Str.
04-011	30 cm	S/S Str.



Operating Scissors

04-012	11.5 cm	S/S Cvd
04-013	12.5 cm	S/S Cvd
04-014	13 cm	S/S Cvd
04-015	14.5 cm	S/S Cvd
04-016	16 cm	S/S Cvd
04-017	18 cm	S/S Cvd
04-018	20 cm	S/S Cvd
04-019	23 cm	S/S Cvd
04-020	25 cm	S/S Cvd
04-021	28 cm	S/S Cvd
04-022	30 cm	S/S Cvd



Operating Scissors

04-023	11.5 cm	B/B Str.
04-024	12.5 cm	B/B Str.
04-025	13 cm	B/B Str.
04-026	14.5 cm	B/B Str.
04-027	16 cm	B/B Str.
04-028	18 cm	B/B Str.
04-029	20 cm	B/B Str.
04-030	23 cm	B/B Str.
04-031	25 cm	B/B Str.
04-032	28 cm	B/B Str.
04-033	30 cm	B/B Str.



Operating Scissors

04-034	11.5 cm	B/B Cvd
04-035	12.5 cm	B/B Cvd
04-036	13 cm	B/B Cvd
04-037	14.5 cm	B/B Cvd
04-038	16 cm	B/B Cvd
04-039	18 cm	B/B Cvd
04-040	20 cm	B/B Cvd
04-041	23 cm	B/B Cvd
04-042	25 cm	B/B Cvd
04-043	28 cm	B/B Cvd
04-044	30 cm	B/B Cvd



Operating Scissors

04-045	11.5 cm	S/B Str.
04-046	12.5 cm	S/B Str.
04-047	13 cm	S/B Str.
04-048	14.5 cm	S/B Str.
04-049	16 cm	S/B Str.
04-050	18 cm	S/B Str.
04-051	20 cm	S/B Str.
04-052	23 cm	S/B Str.
04-053	25 cm	S/B Str.
04-054	28 cm	S/B Str.
04-055	30 cm	S/B Str.



Operating Scissors

04-056	11.5 cm	S/B Cvd
04-057	12.5 cm	S/B Cvd
04-058	13 cm	S/B Cvd
04-059	14.5 cm	S/B Cvd
04-060	16 cm	S/B Cvd
04-061	18 cm	S/B Cvd
04-062	20 cm	S/B Cvd
04-063	23 cm	S/B Cvd
04-064	25 cm	S/B Cvd
04-065	28 cm	S/B Cvd
04-066	30 cm	S/B Cvd

OPERATING SCISSORS

Chirurgische scheren
Tijeras Para Cirugia

Ciseaux à chirurgie
Forbici per chirurgia



Operating Scissors

04-067	11.5 cm	S/S Str.
04-068	12.5 cm	S/S Str.
04-069	13 cm	S/S Str.
04-070	14.5 cm	S/S Str.
04-071	16 cm	S/S Str.
04-072	18 cm	S/S Str.
04-073	20 cm	S/S Str.
04-074	23 cm	S/S Str.
04-075	25 cm	S/S Str.
04-076	28 cm	S/S Str.
04-077	30 cm	S/S Str.



Operating Scissors

04-078	11.5 cm	S/S Cvd
04-079	12.5 cm	S/S Cvd
04-080	13 cm	S/S Cvd
04-081	14.5 cm	S/S Cvd
04-082	16 cm	S/S Cvd
04-083	18 cm	S/S Cvd
04-084	20 cm	S/S Cvd
04-085	23 cm	S/S Cvd
04-086	25 cm	S/S Cvd
04-087	28 cm	S/S Cvd
04-088	30 cm	S/S Cvd



Operating Scissors

04-089	11.5 cm	B/B Str.
04-090	12.5 cm	B/B Str.
04-091	13 cm	B/B Str.
04-092	14.5 cm	B/B Str.
04-093	16 cm	B/B Str.
04-094	18 cm	B/B Str.
04-095	20 cm	B/B Str.
04-096	23 cm	B/B Str.
04-097	25 cm	B/B Str.
04-098	28 cm	B/B Str.
04-099	30 cm	B/B Str.



Operating Scissors

04-100	11.5 cm	B/B Cvd
04-101	12.5 cm	B/B Cvd
04-102	13 cm	B/B Cvd
04-103	14.5 cm	B/B Cvd
04-104	16 cm	B/B Cvd
04-105	18 cm	B/B Cvd
04-106	20 cm	B/B Cvd
04-107	23 cm	B/B Cvd
04-108	25 cm	B/B Cvd
04-109	28 cm	B/B Cvd
04-110	30 cm	B/B Cvd



Operating Scissors

04-111	11.5 cm	S/B Str.
04-112	12.5 cm	S/B Str.
04-113	13 cm	S/B Str.
04-114	14.5 cm	S/B Str.
04-115	16 cm	S/B Str.
04-116	18 cm	S/B Str.
04-117	20 cm	S/B Str.
04-118	23 cm	S/B Str.
04-119	25 cm	S/B Str.
04-120	28 cm	S/B Str.
04-121	30 cm	S/B Str.



Operating Scissors

04-122	11.5 cm	S/B Cvd
04-123	12.5 cm	S/B Cvd
04-124	13 cm	S/B Cvd
04-125	14.5 cm	S/B Cvd
04-126	16 cm	S/B Cvd
04-127	18 cm	S/B Cvd
04-128	20 cm	S/B Cvd
04-129	23 cm	S/B Cvd
04-130	25 cm	S/B Cvd
04-131	28 cm	S/B Cvd
04-132	30 cm	S/B Cvd

OPERATING & DISSECTING SCISSORS



FIMR INSTRUMENTS

Bedienen & Präparierscheren
Funcionamiento y disección Tijeras

D'exploitation & Ciseaux de dissection
Forbici di funzionamento e di dissezione



04-133
12.5 cm
B/B Str.

**LANDOLT OPERATING
SCISSORS**



04-134
12.5 cm
B/B Cvd.



04-135
14.5 cm
S/B Str.

**DEAVER OPERATING,
DISSECTING SCISSORS**



04-136
14.5 cm
S/B Cvd.



04-137
14.5 cm
S/B Str.
Saw Edge.



04-138
14.5 cm
S/B Cvd.
Saw Edge.



04-139
14.5 cm
S/B Str.

**LONGE OPERATING,
DISSECTING SCISSORS**



04-140
14.5 cm
S/B Cvd.



04-141
14.5 cm
S/B Str.
Saw Edge.



04-142
14.5 cm
S/S Cvd.
Saw Edge.



04-143 13.0 cm
04-144 14.5 cm
S/S Str.

INCISION OPERATING SCISSORS



**MAYO DISSECTING &
GYNECOLOGICAL
SCISSORS**



04-145
14.5 cm
04-146
15.0 cm
04-147
17.0 cm
04-148
19.5 cm
04-149
21.0 cm
04-150
23.0 cm
B/B, Str.



04-151
14.5 cm
04-152
15.0 cm
04-153
17.0 cm
04-154
19.5 cm
04-155
21.0 cm
04-156
23.0 cm
B/B, Cvd.



**MAYO-STILLE DISSECTING &
GYNECOLOGICAL SCISSORS**



04-157
15.0 cm
04-158
17.0 cm
04-159
19.5 cm
04-160
21.0 cm
B/B, Str.



04-161
15.0 cm
04-162
17.0 cm
04-163
19.5 cm
04-164
21.0 cm
B/B, Cvd.



**MAYO-HARRINGTON DISSECTING
& GYNECOLOGICAL SCISSORS**



04-165
22.5 cm
B/B, Str.



04-166
22.5 cm
B/B, Cvd.



**MAYO-NOBLE DISSECTING
& GYNECOLOGICAL SCISSORS**



04-167
17.0 cm
B/B, Cvd.



04-168
17.0 cm
B/B, Str.

OPERATING & DISSECTING SCISSORS



Bedienen & Präparierscheren

D'exploitation & Ciseaux de dissection

Funcionamiento y disección Tijeras

Forbici di funzionamento e di dissezione



**METZENBAUM OPERATING,
DISSECTING SCISSORS
BLUNT/BLUNT STR.**



04-169
14.5 cm
B/B Str.



04-170
18 cm
B/B Str.



04-171
20 cm
B/B Str.



04-172
23 cm
B/B Str.



04-173
25 cm
B/B Str.



04-174
28 cm
B/B Str.



**METZENBAUM OPERATING,
DISSECTING SCISSORS
BLUNT/BLUNT CVD.**



04-175
14.5 cm
B/B Cvd.



04-176
18 cm
B/B Cvd.



04-177
20 cm
B/B Cvd.



04-178
23 cm
B/B Cvd.



04-179
25 cm
B/B Cvd.



04-180
28 cm
B/B Cvd.


OPERATING & DISSECTING SCISSORS

Bedienen & Präparierscheren

D'exploitation & Ciseaux de dissection

Funcionamiento y disección Tijeras

Forbici di funzionamento e di dissezione




**METZENBAUM-NELSON
SCISSORS**

04-181 18cm
04-182 20cm
04-183 23cm
04-184 25cm
04-185 28cm
04-186 30cm
B/B Str.

04-187 18cm
04-188 20cm
04-189 23cm
04-190 25cm
04-191 28cm
04-192 30cm
S/B Str.

04-193 18cm
04-194 20cm
04-195 23cm
04-196 25cm
04-197 28cm
04-198 30cm
S/S Str.



**METZENBAUM-NELSON
SCISSORS**

04-199 18cm
04-200 20cm
04-201 23cm
04-202 25cm
04-203 28cm
04-204 30cm
B/B Cvd.

04-205 18cm
04-206 20cm
04-207 23cm
04-208 25cm
04-209 28cm
04-210 30cm
S/B Cvd.

04-211 18cm
04-212 20cm
04-213 23cm
04-214 25cm
04-215 28cm
04-216 30cm
S/S Cvd.



**METZENBAUM-FINO
SCISSORS**

04-217
14.5 cm
B/B Str.

04-218
14.5 cm
S/S Str.

04-219
14.5 cm
B/B Cvd.

04-220
14.5 cm
S/S Cvd.



**METZENBAUM-FINO
SCISSORS**

04-221
18 cm
04-222
220 cm
04-223
230 cm
B/B Str.

04-224
18 cm
04-225
220 cm
04-226
230 cm
B/B Cvd.

04-227
18 cm
04-228
220 cm
04-229
230 cm
B/B Slightly
Cvd.

OPERATING & DISSECTING SCISSORS

Bedienen & Präparierscheren

Funcionamiento y disección Tijeras

D'exploitation & Ciseaux de dissection

Forbici di funzionamento e di dissezione




**RICHTER OPERATING,
DISSECTING SCISSORS**
04-230 14.5cm
Angled



**JOSEPH OPERATING,
DISSECTING SCISSORS**



04-231
14cm
S/S Str.



04-232
14cm
S/S Cvd



**KELLY OPERATING,
DISSECTING SCISSORS**



04-233
16cm
04-234
18cm
S/S Str.



04-235
16cm
04-236
18cm
S/S Cvd.



**DOYEN OPERATING,
DISSECTING SCISSORS**
04-237 17.5 cm
B/B, Str.



**DOYEN OPERATING,
DISSECTING SCISSORS**
04-238 17.5 cm
B/B, Cvd.





SIMS DISSECTING & GYNECOLOGICAL SCISSORS

04-239 20cm
04-240 23cm
Str, B/B



04-241 20cm
04-242 23cm
Str, S/B



04-243 20cm
04-244 23cm
Str, S/S



04-245 20cm
04-246 23cm
Cvd, B/B



04-247 20cm
04-248 23cm
Cvd, S/B



04-249 20cm
04-250 23cm
Cvd, S/S



Iris Scissors



04-251 09 cm
04-252 10.5 cm
04-253 11.5 cm
04-254 12.0 cm
S/S Str.



04-255 09 cm
04-256 10.5 cm
04-257 11.5 cm
04-258 12.0 cm
S/S Cvd.



04-259
11.5 cm
Angular Left



04-262 9 cm
04-263 10cm
04-264 11.5cm
04-265 12cm
S/S Str.



04-260
11.5 cm
Laterally Cvd.



04-261
11.5 cm
Angular Right



Iris Scissors



04-266 9 cm
04-267 10cm
04-268 11.5cm
04-269 12cm
S/S Cvd.

OPERATING SCISSORS

Chirurgische scheren
Tijeras Para Cirugia

Ciseaux à chirurgie
Forbici per chirurgia



VASCULAR SCISSORS
04-270 19cm
S/S Angled 125°



COOLEY REVERSE CUT SCISSORS
04-271 17cm
S/S angled 125°



POTTS-SMITH VASCULAR SCISSORS
S/S Angled



04-272
19cm
25°



04-273
19cm
45°



04-274
19cm
60°



04-275
19cm
125°



HEGEMANN-DIETRICH VASCULAR SCISSORS
S/S Angled



25°
04-276 16cm
04-277 19cm



45°
04-278 16cm
04-279 19cm




60°
04-280 16cm
04-281 19cm



90°
04-282 16cm
04-283 19cm




125°
04-284 16cm
04-285 19cm




DE BAKHEY VASCULAR SCISSORS
S/S Angled



25°
04-286 16cm
04-287 19cm
04-288 23cm
04-289 28cm



45°
04-290 16cm
04-291 19cm
04-292 23cm
04-293 28cm



60°
04-294 16cm
04-295 19cm
04-296 23cm
04-297 28cm


VASCULAR SCISSORS

Gefäßschere
Tijeras vasculares

Ciseaux vasculaires
Forbici vascolari



O'BRIEN VASCULAR SCISSORS
04-298 9cm
S/S Angled



HEATH VASCULAR SCISSORS
04-299 15cm
S/S Double Cvd.




COOLEY VASCULAR SCISSORS
04-300 18.4 cm, 1 mm
04-301 18.4 cm, 1.5 mm
probe point on tip.
S/S Angled 45°




04-302
17cm
Serrated
Blade S/S
Double Cvd.

DISSECTING & NASAL SCISSORS



Locklin Dissecting Scissors
04-303 16.5cm
S/S



FAVALORO CORONARY
04-304 10.5 cm
S/S Angled 90°



ENT-DEAN SCISSORS
04-305 17cm
S/S Double Cvd.
Serrated



04-306
15.5cm
S/S

04-307
15.5cm
S/S

DE BAKY VASCULAR SCISSORS


STITCH SCISSORS

stitch Schere | point de ciseaux
Coser Tijeras | Stitch Forbici




Saw Edge


STITCHES SCISSORS
04-308 12.0 cm
Str.



NORTHBENT STITCHES SCISSORS
04-309 12.5 cm
Cvd.



LITTAUER STITCH SCISSORS
04-310 14 cm
Str.




04-311
14.5 cm
Str.

BUCK STITCH SCISSORS



04-312
18.0 cm
Str.

SPENCER STITCH SCISSORS



04-313
9 cm
Str.



04-314
11.5 cm
Str.



SHORTBENT STITCH SCISSORS
04-315 9 cm
Cvd.

TONSIL, NASAL SCISSORS

Tonsillensteinen Schere, Nasenschere
Tijeras de amígdalas, tijeras nasales

Ciseaux d'amygdale, Ciseaux nasaux
Tonsille forbici, forbici nasali



COTTLE SCISSORS
04-316 16 cm
 B/B, Angled



FOMON SCISSORS
04-317 14.5 cm
 B/B, Angled



SAILER SCISSORS
04-318 17.0 cm
 B/B, Angled



HEYMANN SCISSORS
04-319 18 cm
04-320 19.5 cm
 B/B, Angled



HEYMANN SCISSORS
04-321 18.0 cm
 B/B, Angled



OLIVECRONA VASCULAR SCISSORS
04-322 23.0 cm
 B/B

GUM & OPERATING SCISSORS



Chirurgische scheren, Gum Scheren
Goma Tijeras y tijeras de funcionamiento

Gum Ciseaux et Ciseaux d'exploitation
Gum Forbici & Forbici Operativi



GOOD SCISSORS
04-323 19.0 cm
Probe End, Cvd.



PRINCE SCISSORS
04-324 17.5 cm
B/B, Cvd.




HALLE SCISSORS
04-325 19.0 cm
B/B, Cvd.



04-326
14.5 cm
Pointed Tip
S/S, Str.



04-327
14.5 cm
Pointed Tip
S/S, Cvd.



MOD. BONN SCISSORS
04-328 8.0 cm
S/S, Str.



MOD. BONN SCISSORS
04-329 8.0 cm
S/S, Cvd.

WIRE CUTTING SCISSORS

Draht-Ausschnitt-Scheren
Tijeras de corte de alambre

Ciseaux de coupe de câbles
Filo da taglio Forbici



04-330
10.5 cm
04-331
12.0 cm
B/B Str.
04-332
10.5 cm
04-333
12.0 cm
B/B Str.
Saw Edge



04-334
10.5 cm
04-335
12.0 cm
S/S Str.
04-336
10.5 cm
04-337
12.0 cm
S/S Str.
Saw Edge



04-338
10.0 cm
04-339
12.0 cm
B/B Cvd.
04-340
10.0 cm
04-341
12.0 cm
B/B Cvd.
Saw Edge



04-342
10.0 cm
04-343
12.0 cm
S/S Cvd.
04-344
10.0 cm
04-345
12.0 cm
S/S Cvd.
Saw Edge

BEEBEE WIRE CUTTING SCISSORS

BEEBEE WIRE CUTTING SCISSORS



SYSTRUNK
04-346
13.0 cm
S/S, Str.



SYSTRUNK
04-347
13.0 cm
S/S, Cvd.



UNIVERSAL WIRE CUTTING SCISSORS
04-348 12 cm

SYSTRUNK SCISSORS

NAIL NIPPERS, NAIL SCISSORS



Nagelzangen, Nagelscheren
Pince à ongles, ciseaux à ongles

Pinzas para uñas, tijeras de uñas
Tronchesi per unghie, forbici per unghie



04-349
10 cm
04-350
11.5 cm
S/S Str.

STEVENS SCISSORS




04-351
10 cm
04-352
11.5 cm
S/S Cvd.



04-353
10 cm
04-354
11.5 cm
B/B Str.

STEVENS SCISSORS



04-355
10 cm
04-356
11.5 cm
B/B Cvd.



NAIL SCISSORS
04-357 9 cm S/S, Cvd.



NAIL SCISSORS
04-358 9 cm
S/S, Cvd.



NAIL SCISSORS
04-359 9 cm
Arrow Point Curved



NAIL NIPPERS
04-360 11 cm
04-361 13 cm
04-362 15 cm



NAIL NIPPERS
04-363 14 cm



NAIL NIPPERS
04-364 14 cm



TURNBULL NIPPERS
04-365 11.5 cm
04-366 13 cm



NAIL NIPPERS
04-367 10.5 cm



COHEN NIPPERS
04-368 10.5 cm

MICRO SCISSORS

Mikroscheren | micro Ciseaux
Micro Tijeras | Micro Forbici



**ARTERY
SCISSORS**
Probe end
Angled




04-369
18.0 cm
25°



04-370
18.0 cm
45°



04-371
18.0 cm
60°



04-372
18.0 cm
90°



04-373
18.0 cm
125°



04-374
18.0 cm
25°



04-375
18.0 cm
45°



04-376
18.0 cm
60°



04-377
18.0 cm
90°



04-378
18.0 cm
125°



CASTROVIEJO SCISSORS
Spring Type



04-379
9.0 cm
S/B Str.



04-380
9.0 cm
S/B Cvd.



04-381
9.0 cm
S/S Str.



04-382
9.0 cm
S/S Cvd.



04-383
11.4 cm
B/B Str.



04-384
11.4 cm
B/B Cvd.



04-385
11.4 cm
S/B Str.



04-386
11.4 cm
S/B Cvd.



04-387
11.4 cm
S/S Str.



04-388
11.4 cm
S/S Cvd.



04-389
14 cm
B/B Str.



04-390
14 cm
B/B Cvd.



04-391
14 cm
S/B Str.



04-392
14 cm
S/B Cvd.



04-393
14 cm
S/S Str.



04-394
14 cm
S/S Cvd.



MICRO SCISSORS

Mikroscheren | micro Ciseaux
Micro Tijeras | Micro Forbici



NOYES SCISSORS
Spring Handle



04-395
12 cm
S/S



04-396
12 cm
S/B



04-397
12 cm
B/B



NOYES SCISSORS
Spring Handle
04-398 12cm
Angled,



WECKER SCISSORS
Angled



04-399
11.0 cm
S/B



04-400
11.0 cm
B/B



04-401
11.0 cm
S/S



BARRAQUER SCISSORS
04-402 Sharp Points
04-403 Blunt Points
5.5 cm



CASTROVIEJO CORNEAL MICRO SCISSORS
04-404 11cm
Left Angled,
Spring Handle



CASTROVIEJO CORNEAL MICRO SCISSORS
04-405 11cm
Right Angled,
Spring Handle



TROUTMAN MICRO SCISSORS
04-406 10cm
Left Angled,
Spring Handle



TROUTMAN MICRO SCISSORS
04-407 10cm
Right Angled,
Spring Handle



WESTCOTT MICRO SCISSORS



04-408
11.0 cm
S/S Str.



04-410
11.0 cm
S/S Cvd.



04-409
13.3 cm
S/S Str.



04-411
13.3 cm
S/S Cvd.

MICRO SCISSORS

Mikroscheren
Micro Tijeras

micro Ciseaux
Micro Forbici



MICRO SCISSORS
04-412 12cm
S/S, Cvd.



MICRO SCISSORS
S/S, Angled



04-413
18.0 cm
25°



04-414
18.0 cm
45°



04-415
18.0 cm
60°



04-416
18.0 cm
90°



04-417
18.0 cm
125°



MICRO SCISSORS
S/S, Angled pointed probe



04-418
18.0 cm
25°



04-419
18.0 cm
45°



04-420
18.0 cm
60°



04-421
18.0 cm
90°



04-422
18.0 cm
125°



YASARGILL SCISSORS



04-423 16 cm
04-424 20 cm
04-425 22.5cm
S/S, Str.



04-426 16cm
04-427 20cm
04-428 22cm
S/S, Cvd.



04-429 16cm
04-430 20cm
04-431 22cm
B/B, Str.



04-432 16cm
04-433 20cm
04-434 22cm
S/S, Str.



04-435 16cm
04-436 20cm
04-437 22cm
S/S, Cvd.



04-438 16cm
04-439 20cm
04-440 22cm
B/B, Cvd.



NOYES MICRO SCISSORS
04-441 15cm
S/S



NOYES MICRO SCISSORS
04-442 14.5cm
S/S

PLASTER SCISSORS



Gips Schere
Tijeras Yeso

plâtre Ciseaux
intonaco Forbici



ESMARCH KNIFE
04-443 18cm



BRUNS SCISSORS



04-444
24 cm
Cvd, Probe End



04-445
24cm
Cvd, Serrated
Probe End



**LISTER-EXCENTRIC
PLASTER SCISSORS**
04-446 16 cm
04-447 19 cm
Angled, Probe End



Saw Edge Upper
Teeth is Serrated



**LISTER BANDAGE PLASTER
Scissors, one large ring**
04-448 9 cm
04-449 11 cm
04-450 14 cm
04-451 15 cm
04-452 18 cm
04-453 20 cm
Angled, Probe End



BERGMANN SCISSORS
04-454 23 cm
Angled, Probe End



ESMARCH SCISSORS
04-455 20 cm
04-456 23 cm
Angled, Probe End



**LISTER BANDAGE
PLASTER SCISSORS**
04-457 9 cm
04-458 11 cm
04-459 14 cm
04-460 15 cm
04-461 18 cm
04-462 20 cm
Angled, Probe End



PLASTER SCISSORS
04-463 17 cm
Angled



**BURNHAM BANDAGE
PLASTER SCISSORS**
04-464 14 cm
Angled

PLASTER SCISSORS

Gips Schere
Tijeras Yeso

plâtre Ciseaux
intonaco Forbici



BANDAGE SCISSORS
04-465 16.0 cm



**SMITH (Mod. Usa)
PLASTER SCISSORS**
04-466 16 cm
04-467 18 cm
04-468 20 cm
04-469 23 cm



**LORENZ PLASTER
SCISSORS**
04-470 23 cm
B/B, Cvd



**LANGE BANDAGE
SCISSORS**
04-471 18.0 cm
B/B



**BRAUN-STADLER
PLASTER SCISSORS**
04-472 14cm
04-473 22cm
B/B, Angled



**BUSCH PLASTER
SCISSORS**
04-474 13 cm
04-475 16 cm
B/B, Cvd.



**KNOWLES PLASTER
SCISSORS**
04-476 14 cm
Str. Probe End



**KNOWLES PLASTER
SCISSORS**
04-477 14 cm
Angled, Probe End



**WALDMANN PLASTER
SCISSORS**
04-478 18 cm
Probe End

PLASTER SCISSORS

Gips Schere
Tijeras Yeso

plâtre Ciseaux
intonaco Forbici



WALDMANN SCISSORS
04-479 18 cm



SCHUMACHER PLASTER SCISSORS
04-480 16 cm



MODIFIED-USA PLASTER SCISSORS
04-481 10.5 cm



MARBACH PLASTER SCISSORS
04-482 15.5 cm
Probe End, Serrated blade



SEUTIN PLASTER SCISSORS
04-483 20cm
Probe End



SEUTIN PLASTER SCISSORS
04-484 23cm
Probe End

PLASTER SHEARS

GipsScheren
Tijeras Yeso

plâtres Ciseaux
Cesoie in gesso



Stille Plaster Shears
04-485 23 cm
04-486 26 cm
04-487 37 cm
04-488 Spare Blade




STILLE PLASTER SHEARS
04-489 24 cm
04-490 Spare Blade



**HENNIG PLASTER
CAST SPREADER**
04-491 28 cm



MODEL USA CAST SPREADER PLASTER
04-492 23 cm 3 Prong



SPREADER PLASTER
04-493 29 cm

With removable jaws



WOLFF SPREADER PLASTER
04-494 18 cm
04-495 24 cm

SUPER CUT SCISSORS



Super Cut Schere
Súper Cut Tijeras

Ciseaux SuperCut
Forbici SuperCut



**Standard Operating
Super Cut Scissors**



04-496 14.5 cm
04-497 16.5 cm
Str. B/B



04-498 14.5 cm
04-499 16.5 cm
Str. S/B



04-500 14.5 cm
04-501 16.5 cm
Str. S/S



**Standard Operating
Super Cut Scissors**



04-502 14.5 cm
04-503 16.5 cm
Cvd. B/B



04-504 14.5 cm
04-505 16.5 cm
Cvd. S/B



04-506 14.5 cm
04-507 16.5 cm
Cvd. S/S



**Lexer Operating
Super Cut Scissors**



04-508 16 cm
04-509 21 cm
B/B, Str.



04-510 16.5 cm
04-511 21 cm
B/B, Cvd.



**Safety Super
Cut Scissors**
04-512 Str.
04-513 Cvd.
9.0 cm, Probe
Ends



**Goldman-fox Super
Cut Scissors**
04-514 Str.
04-515 Cvd.
9.0 cm, S/S



**Stevens Tenotomy
Super Cut Scissors**
04-516 Str.
04-517 Cvd.
9.0 cm, B/B



**Stevens Tenotomy
Super Cut Scissors**
04-518 Str.
04-519 Cvd.
11.5 cm, S/S



**BABY METZENBAUM
SUPER CUT SCISSORS**
04-520 11.5 cm
B/B, Str.



**BABY METZENBAUM
SUPER CUT SCISSORS**
04-521 11.5 cm
B/B, Cvd.



IRIS SUPER CUT SCISSORS
04-522 9.0 cm
04-523 11.5 cm
S/S, Str.



IRIS SUPER CUT SCISSORS
04-524 9.0 cm
04-525 11.5 cm
S/S, Cvd.

OPERATING SCISSORS T.C.

Chirurgische scheren
Tijeras quirúrgicas

Ciseaux chirurgicaux
Forbici per chirurgia



Operating Scissors, S/S Str.
T/C Insert

04-526	11.5 cm
04-527	13.0 cm
04-528	14.0 cm
04-529	15.0 cm
04-530	16.0 cm
04-531	17.5 cm
04-532	18.5 cm
04-533	20.0 cm



Operating Scissors, S/S Cvd.
T/C Insert

04-534	11.5 cm
04-535	13.0 cm
04-536	14.0 cm
04-537	15.0 cm
04-538	16.0 cm
04-539	17.5 cm
04-540	18.5 cm
04-541	20.0 cm



Operating Scissors, B/B Str.
T/C Insert

04-542	11.5 cm
04-543	13.0 cm
04-544	14.0 cm
04-545	15.0 cm
04-546	16.0 cm
04-547	17.5 cm
04-548	18.5 cm
04-549	20.0 cm



Operating Scissors, B/B Cvd.
T/C Insert

04-550	11.5 cm
04-551	13.0 cm
04-552	14.0 cm
04-553	15.0 cm
04-554	16.0 cm
04-555	17.5 cm
04-556	18.5 cm
04-557	20.0 cm



Operating Scissors, S/B
Standard Str. T/C Insert

04-558	11.5 cm
04-559	13.0 cm
04-560	14.0 cm
04-561	15.0 cm
04-562	16.0 cm
04-563	17.5 cm
04-564	18.5 cm
04-565	20.0 cm



Operating Scissors, S/B
Standard Cvd. T/C Insert

04-566	11.5 cm
04-567	13.0 cm
04-568	14.0 cm
04-569	15.0 cm
04-570	16.0 cm
04-571	17.5 cm
04-572	18.5 cm
04-573	20.0 cm

OPERATING SCISSORS T.C.

Chirurgische scheren
Tijeras quirúrgicas

Ciseaux chirurgicaux
Forbici per chirurgia



Metzenbaum Operating
Scissors T/C Insert

		
04-574 11.5 cm 04-575 14.5 cm B/B Str.	04-576 11.5 cm 04-577 14.5 cm S/B Str.	04-578 11.5 cm 04-579 14.5 cm S/S Str.



Metzenbaum Operating
Scissors T/C Insert

		
04-580 11.5 cm 04-581 14.5 cm B/B Cvd.	04-582 11.5 cm 04-583 14.5 cm S/B Cvd	04-584 11.5 cm 04-585 14.5 cm S/S Cvd



Metzenbaum-Fino
Operating Scissors
T/C Insert

	
04-586 14.5 cm 04-587 18 cm 04-588 20 cm 04-589 23 cm B/B Str.	04-590 14.5 cm 04-591 18 cm 04-592 20 cm 04-593 23 cm B/B Cvd.



Metzenbaum-Fino
Operating Scissors
T/C Insert

	
04-594 14.5 cm 04-595 18 cm 04-596 20 cm 04-597 23 cm S/S Str.	04-598 14.5 cm 04-599 18 cm 04-600 20 cm 04-601 23 cm S/S Cvd.



Metzenbaum-Nelson
Operating Scissors
T/C Insert

		
04-602 18 cm 04-603 20 cm 04-604 23 cm 04-605 25 cm 04-606 28 cm 04-607 30 cm B/B Str.	04-608 18 cm 04-609 20 cm 04-610 23 cm 04-611 25 cm 04-612 28 cm 04-613 30 cm S/B Str.	04-614 18 cm 04-615 20 cm 04-616 23 cm 04-617 25 cm 04-618 28 cm 04-619 30 cm S/S Str.



Metzenbaum-Nelson
Operating Scissors
T/C Insert

		
04-620 18 cm 04-621 20 cm 04-622 23 cm 04-623 25 cm 04-624 28 cm 04-625 30 cm B/B Cvd.	04-626 18 cm 04-627 20 cm 04-628 23 cm 04-629 25 cm 04-630 28 cm 04-631 30 cm S/B Cvd.	04-632 18 cm 04-633 20 cm 04-634 23 cm 04-635 25 cm 04-636 28 cm 04-637 30 cm S/S Cvd.

OPERATING SCISSORS T.C.

Chirurgische scheren
Tijeras quirúrgicas

Ciseaux chirurgicaux
Forbici per chirurgia



**Mayo Dissecting &
Gynecological Scissors
T/C Insert**



04-638 14.5 cm	04-642 14.5 cm
04-639 15.0 cm	04-643 15.0 cm
04-640 17.0 cm	04-644 17.0 cm
04-641 23.0 cm B/B, Str.	04-645 23.0 cm B/B, Cvd.



**Mayo-stille Dissecting &
Gynecological Scissors
T/C Insert**



04-646 15.0 cm	04-650 15.0 cm
04-647 17.0 cm	04-651 17.0 cm
04-648 19.5 cm	04-652 19.5 cm
04-649 21.0 cm B/B, Str.	04-653 21.0 cm B/B, Cvd.



ENT-DEAN SCISSORS
04-654 17 cm
S/S Double Cvd.
Serrated



**KELLY OPERATING,
DISSECTING SCISSORS**



04-655 16.0 cm	04-657 16.0 cm
04-656 18.0 cm S/S Str.	04-658 18.0 cm S/S Cvd.



**Goldmann-Fox
Operating Scissors
T/C Insert**



04-659 13 cm S/S Str.	04-660 13 cm S/S Cvd.
------------------------------------	------------------------------------



WIRE CUTTING SCISSORS
04-661 12 cm

OPERATING SCISSORS T.C.

Chirurgische scheren
Tijeras quirúrgicas

Ciseaux chirurgicaux
Forbici per chirurgia



Lister Bandage Plaster Scissors T/C Insert, one large ring
04-662 8.9cm
04-663 11.43 cm
04-664 14 cm
04-665 16 cm
04-666 18 cm
04-667 20 cm




Lister Bandage Plaster Scissors, T/C Insert
04-668 8.9cm
04-669 11.43 cm
04-670 14 cm
04-671 16 cm
04-672 18 cm
04-673 20 cm



Castroviejo Scissors T/C Insert



04-674
14 cm
B/B Str.



04-675
14 cm
B/B Cvd.



04-676
14 cm
S/B Str.



04-677
14 cm
S/B Cvd.



04-678
14 cm
S/S Str.



04-679
14 cm
S/S Cvd.



Joseph Operating Scissors T/C Insert



04-680
14 cm
S/S Str.



04-681
14 cm
S/S Cvd.



Reynolds Operating Scissors T/C Insert



04-682
16 cm
04-683 18 cm
S/S Str.



04-684
16 cm
04-685 18 cm
S/S Cvd.



Foman Operating Scissors T/C Insert
04-686
13 cm Full Bent



Stevens Operating Scissors T/C Insert



04-687
10 cm
04-688 11 cm
Str.



04-689
10 cm
04-690 11 cm
Cvd.

Universal Paramedic Shear



DRESSING & TISSUE FORCEPS



TISSUE & DRESSING FORCEPS

Pinzas para tejidos, Pincas à dissection
Pinze per tessuti, Pinze per dissezione

Pincas a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección



Standard Dressing Forceps
serrated
05-01 10.5cm
05-02 11.5cm
05-03 13cm
05-04 14.5cm
05-05 16cm
05-06 18cm
05-07 20cm
05-08 25cm
05-09 30cm



Standard Tissue Forceps
1x2 Teeth
05-10 10.5cm
05-11 11.5cm
05-12 13cm
05-13 14.5cm
05-14 16cm
05-15 18cm
05-16 20cm
05-17 25cm
05-18 30cm



Mod. USA, Dressing Forceps
serrated
05-19 10.2cm
05-20 12cm
05-21 13cm
05-22 14.5cm
05-23 15.5cm
05-24 18cm
05-25 21cm
05-26 25cm
05-27 31cm



Mod. USA, Tissue Forceps
1x2 teeth
05-28 10.2cm
05-29 11.4cm
05-30 12.7cm
05-31 14.0cm
05-32 15.2cm
05-33 17.8cm
05-34 20.3cm
05-35 25.4cm
05-36 30.5cm



Mod. Martin Dressing Forceps
05-37 12.7cm
05-38 14.0cm



Standard Dressing Forceps
Curved, serrated
05-39 11.4cm
05-40 12.7cm
05-41 14.0cm
05-42 15.2cm
05-43 17.8cm
05-44 20.3cm
05-45 25.4cm
05-46 30.5cm



Standard Dressing Forceps
wide pattern cross serration
05-47 11.4cm
05-48 12.7cm
05-49 14.0cm
05-50 15.2cm
05-51 17.8cm
05-52 20.3cm
05-53 25.4cm
05-54 30.5cm



Standard Dressing Forceps
Narrow Point serrated
05-55 12.7cm
05-56 14 cm
05-57 15.2cm
05-58 17.8cm
05-59 20.3cm
05-60 25.4cm
05-61 30.5cm



Semken Dressing Forceps
Straight Serrated
05-62 12.5cm
05-63 15.5cm



Semken Dressing Forceps
Curved Serrated
05-64 12.5cm
05-65 15.5cm



Potts-Smith Dressing Forceps
Straight Serrated
05-66 17.5cm
05-67 21cm
05-68 25 cm



Potts-Smith Dressing Forceps
Curved Serrated
05-69 17.5cm
05-70 21cm
05-71 25cm

TISSUE & DRESSING FORCEPS



Pinces a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección

Pinzas para tejidos, Pinces à dissection
Pinze per tessuti, Pinze per dissezione



Potts-Smith Tissue Forceps
1x2 teeth Straight
05-72 17.5cm
05-73 21.0cm
05-74 25.0cm



Potts-Smith Tissue Forceps
1x2 teeth Curved
05-75 17.5cm
05-76 21.0cm
05-77 25.0cm



Brophy Dressing Forceps
Straight Serrated
05-78 20cm



Brophy Dressing Forceps
Curved Serrated
05-79 20cm



Brophy Tissue Forceps
1x2 teeth
Straight Serrated
05-80 20cm




Brophy Tissue Forceps
1x2 teeth
Curved Serrated
05-81 20cm




Wangensteen Dressing Forceps
05-82 23cm




Stone Tissue Forceps
4X5 Teeth
05-83 15.2cm



Lerche Tissue Forceps
5X6 Teeth
05-84 15cm



Adlerkreutz Tissue Forceps



Gillies Tissue Forceps
1X2 Teeth



Treves Tissue Forceps
1x2 Teeth
05-91 12.5cm
05-92 15cm
05-93 18cm
Dissecting Forceps
teeth, on smooth jaws.



FIMA INSTRUMENTS

TISSUE & DRESSING FORCEPSPinzas para tejidos, Pincers à dissection
Pinze per tessuti, Pinze per dissezionePincers a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección

Tissue Forceps
2x3 Teeth
05-94 11.4cm
05-95 13.0cm
05-96 14.5cm
05-97 16.0cm
05-98 17.5cm
05-99 20.0cm
05-100 25.0cm
05-101 30.0cm



Tissue Forceps
3x4 Teeth
05-102 11.4cm
05-103 13.0cm
05-104 14.5cm
05-105 16.0cm
05-106 17.5cm
05-107 20.0cm
05-108 25.0cm
05-109 30.0cm



Tissue Forceps
4x5 Teeth
05-110 11.4cm
05-111 13.0cm
05-112 14.5cm
05-113 16.0cm
05-114 17.5cm
05-115 20.0cm
05-116 25.0cm
05-117 30.0cm



Iris Tissue Forceps
1x2 Teeth, Straight
05-118 10.2 cm



Iris Tissue Forceps
1x2 Teeth, Half Curved
05-119 10.2 cm



Iris Tissue Forceps
1x2 Teeth, Full Curved
05-120 10.2 cm



Iris Forceps
05-121
10cm



Iris Tissue Forceps
1x2 teeth
05-122 10cm



Waugh Dressing Forceps
05-123 15 cm
05-124 18 cm
slightly heavier
non-toothed,
serrated jaws.



Waugh Tissue Forceps
1x2 Teeth
05-125 15 cm
05-126 18 cm
as above but one two
tooth on serrated jaws.
Known as Scott.



Leedham Tissue Forceps
1x2 Teeth
05-127 14.5 cm
05-128 18 cm
Green Dissecting
Forceps two three
teeth, smooth jaws.



Charnley Tissue Forceps
1x2 Teeth
05-129 14.5cm
05-130 18.0cm
Dissecting Forceps
one into two teeth,
block on jaws.

TISSUE & DRESSING FORCEPS



Pinzas para tejidos, Pincas à dissection
Pinze per tessuti, Pinze per dissezione

Pincas a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección



Ramsay Dissecting Forceps
05-131 18cm
serrated jaws.



05-132
18cm
1x2 teeth

Ramsay Tissue Forceps
Dissecting Forceps
one into teeth on
serrated jaws.



05-133
18cm
2x3 teeth



Dissecting Dressing Forceps
05-134 12.5cm
05-135 15.0cm
05-136 18.0cm
turn over end,
serrated jaws.



Dissecting Dressing Forceps
05-137 12.5cm
05-138 15.0cm
05-139 18.0cm
05-140 20.0cm
block end, serrated
jaws.



Ferris-Smith Tissue Forceps
1x2 Teeth
05-141 18cm



Ferris-Smith Tissue Forceps
2x3 Teeth
05-142 18cm



Lane Tissue Forceps
1x2 Teeth
05-143 12.5cm
05-144 15cm
05-145 17cm
05-146 18cm
05-147 20cm



Lane Tissue Forceps
2x3 Teeth
05-148 12.5cm
05-149 15cm
05-150 17cm
05-151 18cm
05-152 20cm



Bonney Tissue Forceps
1x2 Teeth
05-153 18cm



Bonney Tissue Forceps
2x3 Teeth
05-154 18cm



Durante Tissue Grasping Forceps
05-155 14.5cm
05-156 16cm
05-157 20cm



Nelson Tissue Forceps
05-158 23cm
6x7 Teeth



TISSUE & DRESSING FORCEPS

Pinzas para tejidos, Pinces à dissection
Pinze per tessuti, Pinze per dissezione


Pinces a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección



**Russian-Modell
Tissue Forceps**
05-159 15cm
05-160 20cm
05-161 25cm



**Singley-Tuttle
Forceps**
05-162 23cm



**Collin-Duval
Forceps**



**Stille-Barraya Tissue
Forceps 3x4 Teeth**
05-166 18cm
05-167 20cm
05-168 25cm



**Duehrssen Dressing
Forceps**
Serrated
05-169 20cm



**Duehrssen Tissue
Forceps**
1x2 Teeth
05-170 20cm



**Duehrssen Tissue
Forceps**
2x3 Teeth
05-171 20cm



MC INDOE

05-172
150mm
1x2 teeth



Ligature Forceps

Ochsner
05-173
145mm



SWEDISH PATTERN

05-174
13cm

05-175
15cm

05-176
17cm

05-177
20cm



ADLERKREUTZ

05-178
15cm
2x3 teeth

05-179
20cm
2x3 teeth

05-180
15cm
4x5 teeth

05-181
20cm
4x5 teeth








05-182
25cm
4x5 teeth

TISSUE & DRESSING FORCEPS



Pinzas para tejidos, Pincas à dissection
Pinze per tessuti, Pinze per dissezione

Pincas a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección



TISSUE & DRESSING FORCEPS

Pinzas para tejidos, Pincas à dissection
Pinze per tessuti, Pinze per dissezione

Pincas a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección



LUCAE / ENT Dressing Forceps



Paper Dressing Forceps
05-213 15cm
Curved



TISSUE & DRESSING FORCEPS



Pinzas para tejidos, Pincas à dissection
Pinze per tessuti, Pinze per dissezione

Pincas a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección



Yasargil
05-231
20cm



Hunt
05-232
20cm



Hunt
05-233
20cm



Plain Forceps
05-234 18.5cm



Jansen
05-235 16 cm
05-236 18 cm
05-237 20 cm



Plain Forceps
05-238 16cm




Tissue Forceps
05-239 16cm



Gruenwald Forceps
05-240 14cm
05-241 20cm



Gruenwald Tissue Forceps
05-242 20cm



Hoen
05-243
20cm



D'Errico
05-244
22cm



D'Errico
05-245
22cm


TISSUE & DRESSING FORCEPS

Pinzas para tejidos, Pincers à dissection
Pinze per tessuti, Pinze per dissezione

Pincers a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección



Adson
05-246
23cm



Troeltsch
05-247 10cm
05-248 12cm
05-249 15cm
05-250 17cm



Troeltsch
05-251
13cm



Wilde
05-252
11cm




Troeltsch
05-253
11.5cm



Landolt
05-254
21.5cm



Yasargil



05-255 18cm
05-256 20 cm
05-257 22 cm
05-258 24 cm
05-259 18 cm
05-260 20 cm
05-261 22 cm
05-262 24 cm
05-263 18 cm
05-264 20 cm
05-265 22 cm
05-266 24 cm
1x2 teeth



YASARGIL
22 cm



05-267



05-268



YASARGIL
22 cm



05-269



05-270



YASARGIL
24 cm



05-271



05-272






05-273

TISSUE & DRESSING FORCEPS



Pinces a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección


Pinzas para tejidos, Pinces à dissection
Pinze per tessuti, Pinze per dissezione




TISSUE & DRESSING FORCEPS

*Pinces a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección*

*Pinzas para tejidos, Pinces à dissection
Pinze per tessuti, Pinze per dissezione*



Adson Brown Tissue



**Adson Brown
Tissue Forceps**
Fenestrated Handle

TISSUE & DRESSING FORCEPS



Pinces a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección

Pinzas para tejidos, Pinces à dissection
Pinze per tessuti, Pinze per dissezione



Graefe Forceps

05-324
7cm
05-325
10cm
straight
serrated



05-326
7cm
05-327
10cm
half curved
serrated



05-328
7cm
05-329
10cm
full curved
serrated



05-330
7cm
05-331
10cm
1x2 teeth
straight



05-332
7cm
05-333
10cm
1x2 teeth
half curved



05-334
7cm
05-335
10cm
1x2 teeth
full curved



Stevens Tissue Forceps
1x2 Teeth



05-336
11cm
straight



Stevens Tissue Forceps
1x2 Teeth



05-337
11cm
Curved



Graefe Dressing Forceps

05-338 0.5 mm
05-339 0.7 mm
05-340 1.0 mm
7cm
Straight Serrated



Graefe Dressing Forceps

05-341 0.5 mm
05-342 0.7 mm
05-343 1.0 mm
7cm
curved serrated



Graefe Tissue Forceps
1x2 Teeth

05-344 0.5 mm
05-345 0.7 mm
05-346 1.0 mm
7cm Straight



Graefe Tissue Forceps
1x2 Teeth

05-347 0.5 mm
05-348 0.7 mm
05-349 1.0 mm
7cm Curved



Graefe Fixation Forceps

05-350 11.2 cm
standard jaws 4.5mm wide
with fine teeth



05-351
90mm
05-352
105mm
05-353
115mm
05-354
125mm

JEWELER TYPE, DISSECTING FORCEPS

Chirurgische pinzetten, Uhrmacher-Pinzetten
Pinzas de disección con dientes, Pinzas de relojero

Pinces à dissection à groffes, Pinces d'horloger
Pinze per dissezione con denti, Pinze de orologiaio



Dumont Forceps
05-355 11.5cm
curved



05-356
10.5cm
full curved



05-357
10.5cm
half curved



05-358
11cm



05-359
11cm
fine point



05-360
11cm
extra fine
point



05-361
12cm
fine point



05-362
12cm
extra fine
point



05-363
12cm



Dumont Forceps
05-364
style 1, fine point,
12cm



Dumont Forceps
05-365
style 2, strong
fine point,
12cm



Dumont Forceps
05-366
style 3, narrow
fine point,
12cm



Dumont Forceps
05-367
style 3C, narrow
fine point,
11cm



Dumont Forceps
05-368
style 4, fine point,
11cm



Dumont Forceps
05-369
style 5, super fine
11cm point



Dumont Forceps
05-370
style 5, micro fine
11cm point



Dumont Forceps
05-371
style 7, micro fine
11.5 cm curved



Dumont Forceps
05-372
style 7, fine curved,
11.5 cm



Dumont Forceps
05-373
style 6, pick up,
11.5 cm

CILIA FORCEPS



Cilienpinzetten

Pinzas depilatorias

Pinces à épiler les cils

Pinze depilatoria



Ziegler Cilia Forceps
05-374 9cm



Beer Cilia Forceps
05-375 9cm



Douglas Cilia Forceps
05-376 9cm



Gradle Cilia Forceps
05-377 10.5cm



Kuehne Cilia Forceps
05-378 10.5cm



Kuehne Cilia Forceps
05-379 10.5cm
angled sideways



Barraquer Cilia Forceps
05-380 10.5cm



Littauer Cilia Forceps
05-381 8.5cm



Littauer Cilia Forceps
05-382 9cm



Splinter Forceps
05-383 8 cm
05-384 9 cm
05-385 10.5cm
05-386 11.5cm
05-387 12.5cm
05-388 14.5cm
Straight Plain



Splinter Forceps
05-389 8cm
05-390 9cm
05-391 10.5cm
05-392 11.5cm
05-393 12.5cm
05-394 14.5cm
Curved Plain



Splinter Forceps
05-395 8cm
05-396 9cm
05-397 10.5cm
05-398 11.5cm
05-399 12.5cm
05-400 14.5cm
Straight serrated



Splinter Forceps
05-401 8cm
05-402 9cm
05-403 10.5cm
05-404 11.5cm
05-405 12.5cm
05-406 14.5cm
Curved serrated



Feilchenfeld Forceps

05-407 7.5cm
05-408 9cm
05-409 11cm
Straight

SPLINTER FORCEPS

Splitterpinzetten

Pinzas Para espinas

Pinces à échardes

Pinze per scheggie



05-410
10.5cm
Straight

Hunter Forceps



05-411
10.5cm
Curved

Hunter Forceps



05-412
10.5cm
Straight

Carmalt Forceps



05-413
10.5cm
Curved

Carmalt Forceps



Splinter Forceps
05-414 8cm
05-415 9cm
05-416 10.5cm
05-417 11.5cm
05-418 12.5cm
05-419 14.5cm
Straight serrated
with guide pin



Splinter Forceps
05-420 8cm
05-421 9cm
05-422 10.5cm
05-423 11.5cm
05-424 12.5cm
05-425 14.5cm
Curved serrated with
guide pin



Splinter Forceps
05-426 8cm
05-427 9cm
05-428 10.5cm
05-429 11.5cm
05-430 12.5cm
05-431 14.5cm
Straight plain with
guide pin



Splinter Forceps
05-432 8cm
05-433 9cm
05-434 10.5cm
05-435 11.5cm
05-436 12.5cm
05-437 14.5cm
Curved plain with
guide pin



**Stieglitz Forceps
with Magnifying
Mirrors**
05-438 12.5cm



Ochsner Forceps
05-439 14cm



**Lange Tissue
Forceps**
05-440 10cm



**Childe Tissue
Forceps**
05-441 18cm



**Michel Tissue
Forceps**
05-442 12cm



**Wachenfeldt
Tissue Forceps**
05-443 12.5cm



**Hegenbarth
Tissue Forceps**
05-444 12.5cm

SPLINTER FORCEPS



Splitterpinzetten
Pinzas Para espinas


Pinces à échardes
Pinze per scheggie



**Wachenfeldt
Tissue Forceps**
05-445 12.5cm



Stieglitz Splinter Forceps
05-446 14cm Str.
05-447 14cm 45° Angled




Stieglitz Splinter Forceps
05-448 14cm, 90° Angled

90°



Peet Splinter Forceps
05-449 11.5cm



Ralk Splinter Forceps
05-450 15cm, Str.
05-451 15cm, Cvd.



Michel Forceps
05-452 12cm



Michel Clips
05-453
sizes; 11mm x 3mm
& 18mm x 3mm.



Mod. Michel Forceps
05-454 12cm



Collin Forceps
05-455 13cm



Richter-Heath Forceps
05-456 14cm



05-457 14cm
05-460 20cm



05-458 16cm
05-461 25cm



05-459 18cm
05-462 30cm

T.C Oehler Tissue Forceps



05-463 15cm
Semken
with dissector
end



05-464 16cm
05-465 18cm
Semken
1x2 teeth



05-466 15cm
05-467 19cm
05-468 24cm
**T.C De Bakey
Dressing Forceps**



05-469 15cm
**T.C Gillies
Tissue Forceps**






FIMA INSTRUMENTS

TISSUE & DRESSING FORCEPS T.C

Pinces a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección

Pinzas para tejidos, Pinces à dissection
Pinze per tessuti, Pinze per dissezione






TISSUE & DRESSING FORCEPS T.C



FIMA INSTRUMENTS

Pinces a tissu, Anatomische pinzetten
Gewebe, Pinzas de disección

Pinzas para tejidos, Pinces à dissection
Pinze per tessuti, Pinze per dissezione



HAEMOSTATIC FORCEPS



HAEMOSTATIC FORCEPS



Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche



Mosquito Haemostatic Forceps

Halsted-Mosquito
Haemostatic Forceps



Halsted-Mosquito
Haemostatic Forceps

Kelly Haemostatic Forceps

HAEMOSTATIC FORCEPS

Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche



**Dandy Haemostatic
Fine point Forceps**
06-18 14cm
Curved sideways



**Dandy Haemostatic Fine
point Tissue Forceps**
06-20 14cm
(rat tooth 1x2 teeth)
Curved sideways




06-21
14cm
06-22
16cm
straight

06-23
14cm
06-24
16cm
curved

Dandy Haemostatic Forceps
06-17 14cm
Curved sideways

**Dandy Haemostatic
Tissue Forceps**
06-19 14cm
(rat tooth 1x2 teeth)
Curved sideways


Crile Haemostatic Forceps



06-25
14cm
(rat tooth
1x2 teeth)
(kocher)
straight

06-26
14cm
(rat tooth
1x2 teeth)
(kocher)
curved


**Crile Haemostatic
Tissue Forceps**



06-27
14cm
06-28
16cm
straight

06-29
14cm
06-30
16cm
curved

**Crile Haemostatic
Fine Point Forceps**



06-31
14cm
straight
(rat tooth
1x2 teeth)
Crile Haemostatic
Fine Point
Tissue Forceps


06-32
14cm
curved
(rat tooth
1x2 teeth)
Crile Haemostatic
Fine Point
Tissue Forceps

HAEMOSTATIC FORCEPS



Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche



Halsted Haemostatic Forceps



06-33
18.5cm
straight




06-35
18.5cm
(rat tooth
1x2 teeth)
straight



06-34
18.5cm
curved



06-36
18.5cm
(rat tooth
1x2 teeth)
curved



Halsted Haemostatic
Fine point Forceps



06-37
18.5cm
straight




06-39
18.5cm
(rat tooth
1x2 teeth)
straight



06-38
18.5cm
curved



06-40
18.5cm
(rat tooth
1x2 teeth)
curved




Packer Mosquito Forceps
06-41 12.7cm
Curved, flat jaws




Coller Haemostatic Forceps



06-42
14cm
06-43
16cm
Straight,
serrated



06-44
14cm
06-45
16cm
Curved,
serrated




Christophe Haemostatic Forceps




06-46
13cm
06-47
15cm
Curved,
Serrated



Saroot Haemostatic
Forceps



06-48
24cm
Straight,
Serrated



06-49
24cm
Curved
Serrated

HAEMOSTATIC FORCEPS

Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche



06-50
18.5cm
straight

06-51
18.5cm
curved


Adson Haemostatic Fine point Forceps



06-52
18.5cm
straight
(rat tooth
1x2 teeth)

06-53
18.5cm
Curved
(rat tooth
1x2 teeth)

Adson Haemostatic Tissue Forceps



06-54
15cm
06-55
20cm
Curved

06-56
15cm
06-57
20cm
Curved
Fine Point

O'Shaughnessy Haemostatic Forceps



06-58 14.5cm
Straight,
box joint, as
above but
curved on flat.

Cairns Haemostatic Forceps



06-59
13cm
06-60
14cm
straight

06-61
13cm
06-62
14cm
curved

Pean Haemostatic Forceps Slender Type



06-63
14.5cm
Straight


Kocher Haemostatic Forceps Strong Pattern

HAEMOSTATIC FORCEPS



Arterienklemmen
Pincas Hemostaticas


Pinzas Hemostaticas
Pinze Emostatiche



06-64
13cm
06-65
15cm
06-66
17.5cm
06-67
20cm
straight




06-68
13cm
06-69
15cm
06-70
17.5cm
06-71
20cm
curved



06-72
13cm
06-73
14cm
straight



06-74
13cm
06-75
14cm
curved



06-76
16cm
06-77
20cm
Straight




06-78
16cm
06-79
20cm
Curved

Spencer-Wells Haemostatic Forceps
Slender Type

Kocher Haemostatic Forceps
1x2 teeth

Rochester-Carmalt Haemostatic
Forceps




06-80
16cm
06-81
20cm
Straight



06-82
14cm
1x2 teeth
Straight



06-83
14cm
1x2 teeth
Curved




06-84
23cm



06-85
24cm



06-86
19.5cm



06-87
19.5cm

Rochester-Ochsner
Haemostatic Forceps


Moynihan Haemostatic Forceps

Negus Haemostatic Forceps

HAEMOSTATIC FORCEPS

Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche




06-88 14cm
06-89 16cm
06-90 18.5cm
06-91 20cm
06-92 22.5cm
06-93 24cm
06-94 25.4cm
06-95 28cm
06-96 30.5cm
06-97 35.6cm
06-98 40.5cm
06-99 45.5cm
06-100 50.8cm
06-101 55.6cm
06-102 60cm
straight

Rochester-Pean Haemostatic Forceps
Straight



06-103 14cm
06-104 16cm
06-105 18.5cm
06-106 20cm
06-107 22.5cm
06-108 24cm
06-109 25.4cm
06-110 28cm
06-111 30.5cm
06-112 35.6cm
06-113 40.5cm
06-114 45.5cm
06-115 50.8cm
06-116 55.6cm
06-117 60cm
Curved

Rochester-Pean Haemostatic
Forceps Curved




06-118
18.5cm
1x2 teeth

Mikulicz Haemostatic Tissue Forceps



06-119 16cm
06-120 18.5cm
06-121 20cm
06-122 22.5cm
06-123 24cm
06-124 26cm
Straight
1x2 teeth

Kocher-Ochsner Haemostatic
Forceps 1x2 teeth



06-125 16cm
06-126 18.5cm
06-127 20cm
06-128 22.5cm
06-129 24cm
06-130 26cm
Curved
1x2 teeth

Rochester-Ochsner
Haemostatic Forceps
1x2 teeth



06-131 14cm
06-132 15cm
06-133 16cm
06-134 18cm
06-135 20cm
06-136 22cm
06-137 24cm
06-138 26cm
06-139 28cm
06-140 25cm

06-141 14cm
06-142 15cm
06-143 16cm
06-144 18cm
06-145 20cm
06-146 22cm
06-147 24cm
06-148 26cm
06-149 28cm
06-150 25cm

HAEMOSTATIC FORCEPS





Arterienklemmen
Pinces Hemostatiques



Pinzas Hemostaticas
Pinze Emostatiche




Halsted Haemostatic Forceps




Bengolea Haemostatic Forceps



Bengolea Haemostatic Tissue Forceps
(rat tooth 1x2 teeth)



Mixer-O'Shaugnessy
Haemostatic Forceps




Mixer-O'Shaugnessy
Haemostatic Forceps

06-162
14cm
06-163
16.5cm
06-164
19cm
06-165
21.5cm
06-166
24cm
micro point



Craford Haemostatic Forceps



HAEMOSTATIC FORCEPS

Arterienklemmen


Pinces Hemostatiques

Pinzas Hemostaticas



Pinze Emostatiche




Herrick Haemostatic Forceps
06-168 23cm




Roberts Haemostatic Forceps




**Dixon-Lovelace
Haemostatic Tissue Forceps**



**Crafoord-Sellors
Haemostatic Forceps**




**Heaney / Ballantine
Haemostatic Forceps**



HAEMOSTATIC FORCEPS

Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche




Heaney / Rezak
Haemostatic Forceps

06-178
21cm
Straight, Cross
Serrated
06-179
21cm
Slight Curved,
Cross
Serrated



Stille Haemostatic Forceps
06-180 23cm



Lower Haemostatic Forceps
06-181 18cm



Lahey (Sweet)
Haemostatic Forceps

06-182
19cm
06-183
22cm




Desjardins Haemostatic
Forceps
06-184 22cm



Mixer Haemostatic
Gall Duct Forceps

06-185
22cm




06-186
25cm


HAEMOSTATIC FORCEPS

Arterienklemmen
Pinces Hemostatiques


Pinzas Hemostaticas
Pinze Emostatiche



Mixer Forceps,
06-188 18.4cm
Fully curved,
original model



06-189
19.5cm
06-190
21.5cm
06-191
24cm




06-192
19.5cm
06-193
21.5cm
06-194
24cm


Mixer Forceps, Lee Modif
06-187 13.3cm
Curved jaws and angled
shanks, delicate

Heaney Haemostatic Forceps

**Heaney Modif. Haemostatic
Forceps**




06-195
21cm



06-196
24cm
Straight

06-197
24cm
Curved



06-198
1x2 teeth
20cm

**Heaney-Kanter
Haemostatic Forceps**

**Segond-Landau
Haemostatic Forceps**

**Schindler Haemostatic Tissue
Forceps**


HAEMOSTATIC FORCEPS




FIMA INSTRUMENTS

Arterienklemmen
Pincas Hemostaticas


Pinzas Hemostaticas
Pinze Emostatiche



Parametrium Haemostatic Forceps



06-199
21cm
06-200
24cm
Straight



06-201
21cm
06-202
24cm
25° Angled



06-203
21cm
06-204
24cm
45° Angled




06-205
21cm
06-206
24cm
90° Angled



Bozemann Haemostatic Forceps



06-207
26cm
Curved
Serrated



Bozemann Haemostatic Forceps

06-208
26cm
Straight



Westermarck Haemostatic Forceps



06-209
23cm
Straight



Phaneuf Haemostatic Tissue Forceps



06-210
1x2 teeth
21.5cm
Straight



06-211
1x2 teeth
21.5cm
Hard Curved



Long Haemostatic Tissue Forceps



06-212
1x2 teeth
19cm
Curved



06-213
1x2 teeth
19cm
Straight

HAEMOSTATIC FORCEPS

Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche



06-214
1x2 teeth
14cm

**Baby-Mikulicz Haemostatic
Tissue Forceps**




06-215
1x2 teeth
18.5cm

Mikulicz Haemostatic Tissue Forceps



06-216
1x2 teeth
20.5cm

Mikulicz Haemostatic Tissue Forceps



06-217
1x2 teeth
20.5cm

Mikulicz Haemostatic Tissue Forceps



06-218
1x2 teeth
20.5cm

Mikulicz Haemostatic Tissue Forceps



06-219
1x2 teeth
20.5cm

Faure Haemostatic Tissue Forceps

HAEMOSTATIC FORCEPS



Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche



06-220
1x2 teeth
23.5cm

Wertheim Haemostatic
Tissue Forceps




06-221
1x2 teeth
24.5cm

Wertheim Haemostatic
Tissue Forceps



06-222
25.5cm

Wertheim Hysterectomy
Forceps



06-223
26cm

Wertheim Parametrium
Clamp



Heaney Atraumata
Hysterectomy Forceps



06-224
22.5cm
Straight



06-225
22.5cm
25° Angled



06-226
21.5cm
45° Angled



06-227
20.5cm
90° Angled



Wertheim Parametrium Clamps



06-228
25cm
Straight



06-229
24.5cm
25° Angled



06-230
24cm
45° Angled



06-231
22cm
90° Angled

HAEMOSTATIC FORCEPS

Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche



Gemini Mixer Hemostatic Forceps
06-232 14cm
06-233 17.75cm
full curve, delicate



Dick Haemostatic Forceps



06-234
25.5cm

06-235
26.5cm

06-236
27.5cm



06-237
21.5cm



**Wertheim-Cullen
Haemostatic Forceps**



**Wertheim-Cullen Haemostatic
Forceps**
23cm



06-238
60mm

06-239
80mm



06-240
25cm



Wertheim Haemostatic Forceps

TOWEL FORCEPS



Tuchklemmen | Pincès à champs
Pinzas de campo | Pinze fissatele



**Jones / Perforating
Towel Forceps**
06-241 9 cm
Slight Curve

**Jones-Martin
06-242**
Towel Clamp sharp
8.5 cm

**Engl. Modell
Towel Forceps**
06-243 9 cm


**Schaedel Towel
Forceps**
06-244 9 cm

**Thomson-Walker's
Penile Clamps**
06-245 9 cm
Cross-Action model.

**Doyen
06-246**
11cm

**Doyen Towel
Forceps**
06-247
18 cm

**Fricke Artery
Forceps**
06-248
13cm



06-249 8cm
06-250 11cm
06-251 13cm
06-252 15cm
Curved

Backhaus Towel Forceps



Backhaus Clip Towel Forceps
06-253 11cm
with clip for tubing or cable



Roeder Towel Forceps
06-254 13cm
06-255 15cm
Curved, Serrated



Lorna Towel Forceps
06-256 9cm
06-257 11.4cm
06-258 13.5cm



Ball & Socket Towel Forceps
06-259 13cm
Non-perforating



Tohoku Towel Forceps
06-260 10cm
06-261 13cm
Non-Perforating



Non-Perforating Towel Forceps
06-262 11cm
06-263 13cm



Peers-Bertman Towel Forceps
06-264 14.5cm
Serrated

HAEMOSTATIC FORCEPS

Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche




Bernhard Towel Forceps
06-265 16.5cm

Presbyterian Tubing Clamp
06-266 16.5cm
06-267 18cm
Straight, Smooth


Tubing Clamp
06-268 18cm
06-269 20cm
Cross Serrated

Klinikum Berlin Tubing Clamp

06-270
14.5cm
with long
ratchet




Tubing Clamp with Protector
Cross Serrated




06-271
15cm

06-272
18cm

06-273
20cm



Fergusson Tubing Clamp



06-274
16cm

06-275
20cm

06-276
16cm

06-277
20cm

HAEMOSTATIC FORCEPS



Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche



Tuffier Haemostatic Forceps
06-278 13cm



Chaput Haemostatic Forceps
06-279 13cm




Pratt Haemostatic Forceps
06-280 15cm




Pennington / Siezing Haemostatic Forceps
06-281 15.5cm
06-282 20cm
Straight



Martel Haemostatic Forceps
06-283 14cm



Willet-Martel Haemostatic Forceps
06-284 19cm



Lister Sinus Forceps



06-285
13cm

06-286
15.5cm

06-287
17.5cm

HAEMOSTATIC FORCEPS

Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche




06-288
12.5cm



06-289
14cm

**Collin Dressing
Haemostatic Forceps**




06-290
16.5cm



06-291
18.5cm


**Colling Dressing
Haemostatic Forceps**



Gross Maier Dressing Forceps
06-292 14.5cm straight
Without Ratchet




Gross Maier Dressing Forceps
06-293 14.5cm straight
With Ratchet




Gross Maier Dressing Forceps
06-294 14.5cm, curved
With Ratchet



Gross Maier Dressing Forceps
06-295 14.5cm, curved
With Ratchet



**Gross Maier Dressing
Sponge & Dressing Forceps**
06-296 20.5cm, straight
Without Ratchet




**Gross Maier Dressing
Sponge & Dressing Forceps**
06-297 17.8cm, straight
With Ratchet

HAEMOSTATIC FORCEPS




Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche



**Gross Maier Dressing
Sponge & Dressing Forceps**
06-298 20.5cm, curved
Without Ratchet
(bullet holding forceps)



**Gross Maier Dressing
Sponge & Dressing Forceps**
06-299 17.8cm, curved
With Ratchet




Gross Maier Dressing Forceps
06-300 26.5cm, straight
Without Ratchet
(bullet holding forceps)



Gross Maier Dressing Forceps
06-301 26.5cm, straight
With Ratchet
(bullet holding forceps)



Gross Maier Dressing Forceps
06-302 26.5cm, curved
Without Ratchet
(bullet holding forceps)




Gross Maier Dressing Forceps
06-303 26.5cm, curved
With Ratchet
(bullet holding forceps)



**Collin Sponge &
Dressing Forceps**

06-304
24cm




**Duplay Sponge &
Dressing Forceps**
06-305 19cm
Slight Curve

HAEMOSTATIC FORCEPS

Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche



06-306
28cm
Straight

06-307
28cm
Curved


Sims-Maier Sponge & Dressing Forceps



06-308
22cm
Straight

06-309
25cm
06-310
27cm
Straight

Ulrich Sponge Forceps



06-311
22cm
Curved

06-312
25cm
06-313
27cm
Curved


Ulrich Sponge Forceps



06-314
22cm
Straight

06-315
22cm
Curved

Aesculap Sponge & Dressing Forceps



06-316
24cm
Straight

06-317
24cm
Slight
Curve

Fletcher Sponge Forceps



06-318
25cm
Straight

**Pelkmann Sponge
.Forceps**

HAEMOSTATIC FORCEPS

Arterienklemmen

Pinzas Hemostaticas


Pinces Hemostatiques

Pinze Emostatiche



06-319
25cm
Curved

Pelkmann Sponge Forceps



06-320
24cm

Cheron Forceps
Handles Laterally
Angled



06-321
20.5cm
Straight

Duplay Sponge Forceps




06-322
20.5cm
Curved

Duplay Sponge Forceps



**Heywood Smith
Sponge Forceps**
06-323 25cm



06-324
21cm

**Emmet Sponge &
Dressing Forceps**



06-325
15.5cm

Tendon Seizing Forceps



06-326
15.5cm
coarse
serration

Herz Tendon Seizing Forceps



06-327
15.5cm
fine
serrations

HAEMOSTATIC FORCEPS

Arterienklemmen
Pinces Hemostatiques

Pinzas Hemostaticas
Pinze Emostatiche



Noto Sponge & Dressing Forceps
06-328 28cm

Kelly Sponge & Dressing Forceps

06-329
32cm
Curved

06-330
32cm
Straight

Heywood-Smith Sponge & Dressing Forceps
06-331 20cm


Collin Sponge & Dressing Forceps

06-332
16cm
06-333
17cm

06-334
16cm



Foerster-Ballenger Sponge Forceps
& (bullet holding forceps)




06-335 18cm
06-336 20cm
06-337 24cm
Straight Smooth



06-338 18cm
06-339 20cm
06-340 24cm
Straight Serrated



06-341 18cm
06-342 20cm
06-343 24cm
Curved Smooth



06-344 18cm
06-345 20cm
06-346 24cm
Curved Serrated


NEEDLE HOLDERS



NEEDLE HOLDERS

Nadelhalter
Prote-aiguilles

Porta-agujas
Porta aghi



Mayo-Hegar
Needle Holders

07-01
12.5 cm
07-02
14 cm
07-03
16 cm
07-04
18 cm
07-05
20 cm
07-06
24 cm
07-07
26 cm
07-08
30 cm
Serrated-
Grooved


07-09
12.5 cm
07-10
14 cm
07-11
16 cm
07-12
18 cm
07-13
20 cm
07-14
24 cm
07-15
26 cm
07-16
30 cm
Serrated

07-17
16 cm
07-18
18 cm
07-19
20 cm
Serrated-
Grooved
left handed

07-20
12.5 cm
07-21
14 cm
07-22
16 cm
07-23
18 cm
07-24
20 cm
07-25
24 cm
07-26
26 cm
07-27
30 cm
Serrated-
Grooved
Curved

07-28
12.5 cm
07-29
14 cm
07-30
16 cm
07-31
18 cm
07-32
20 cm
07-33
24 cm
07-34
26 cm
07-35
30 cm
Serrated
Curved


Also available in
Ultra Smooth



Baumgartner
Needle Holder

07-36
12.5 cm
07-37
14 cm
07-38
16 cm
07-39
18 cm
07-40
20 cm
07-41
24 cm
07-42
26 cm
07-43
30 cm
Serrated-Solid


Also available in
Ultra Smooth



Webster
Needle Holder

07-44
11.5 cm
07-45
12 cm
07-46
12.5 cm
07-47
13 cm
07-48
14 cm
07-49
15 cm
07-50
18 cm
07-51
20 cm
Smooth


Also available in
Ultra Smooth



Brown
Needle Holder

07-52
13 cm
07-53
18 cm
Serrated-Solid

Also available in
Ultra Smooth




Baby Crile-Wood
Needle Holder

07-54
15 cm

NEEDLE HOLDERS


Nadelhalter	Porta-agujas
Prote-aiguilles	Porta aghi



**Crile-Wood
Needle Holder**


07-55
13 cm
07-56
15 cm
07-57
18 cm
07-58
20 cm
07-59
22 cm
07-60
25 cm
grooved
serrated

07-61
13 cm
07-62
15 cm
07-63
18 cm
07-64
20 cm
07-65
22 cm
07-66
25 cm
solid
serrated



Sarot Needle Holders

07-67
15 cm
07-68
18 cm
07-69
22 cm
07-70
26 cm




Bozemann Needle Holders

07-71
20 cm
07-72
24 cm
07-73
26 cm
07-74
30 cm




Also available in
Ultra Smooth



Jameson Needle Holders

07-75
23 cm




Finocchietto Needle Holders

07-76
27 cm



Also available in
Ultra Smooth



Stratte Needle Holders


07-77
20 cm

07-78
23 cm

NEEDLE HOLDERS

Nadelhalter
Prote-aiguilles


Porta-agujas
Porta aghi



07-79
18 cm

Also available in
Ultra Smooth

Jacobson Needle Holders



07-80
20 cm
Box Joint

Also available in
Ultra Smooth

**Thomson Walker
Needle Holders**



07-81 18 cm
07-82 25 cm
Pound Action Joint.

Naunton Morgan Needle Holders



07-83
17 cm

**Mass. General Hospital
Needle Holders**



07-84
18 cm

Adson Needle Holders



07-85
18 cm

New Orleans Needle Holders

NEEDLE HOLDERS




Nadelhalter
Prote-aiguilles

Porta-agujas
Porta aghi




Heaney Needle Holders



07-86
20 cm
Cvd.

07-87
27 cm
Cvd.




Senning Needle Holders
(Schwed. Model)




07-88
15 cm

07-89
17 cm

Also available in
Ultra Smooth



Wertheim Needle Holders
07-90 20 cm
07-91 24 cm



Gillies Needle Holders
07-92 16 cm



A

B




Masson-Luethy
Needle Holders



07-93
27 cm
Serrated



Olsen-Hegar Needle Holders



07-94
11.5 cm
07-95
14 cm
07-96
17 cm
07-97
19 cm
serrated

07-98
11.5 cm
07-99
14 cm
07-100
17 cm
07-101
19 cm
smooth

NEEDLE HOLDERS

Nadelhalter
Prote-aiguilles


Porta-agujas
Porta aghi



Halsey Needle Holders



Derf Needle Holders



Kilner Needle Holders



Neivert Needle Holders



Collier Needle Holders



Doyen Needle Holders

NEEDLE HOLDERS




Nadelhalter Porta-agujas
Prote-aiguilles Porta aghi



07-111
15 cm
Slender
pattern

**Crile-Murray
Needle Holders**



07-112
18.5 cm
One grooved and
fenestrated jaw

**Metzenbaum
Needle Holders**




07-113
27 cm
Serrated-Solid

**Wangensteen
Needle Holders**



07-114
27 cm

Masson Needle Holder



07-115
27 cm

Johnson Needle Holders



De Bakey Needle Holders
07-116 18 cm
07-117 20 cm
07-118 22 cm
07-119 25 cm




Intracardiac Needle Holders
07-120 18 cm
07-121 20 cm
07-122 22 cm
07-123 25 cm



Microvascular Needle Holders
07-124 18 cm
07-125 20 cm
07-126 22 cm
07-127 25 cm

NEEDLE HOLDERS

Nadelhalter | Porta-agujas
Prote-aiguilles | Porta aghi



Sarot Needle Holders



07-128
18 cm




07-129
26.5 cm



Micro Ryder Needle Holders



07-130
15 cm
07-131
18 cm
07-132
20 cm



Ryder Needle Holders
Narrow jaws.



07-133
13 cm
07-134
15 cm
07-135
18 cm
07-136
20 cm
07-137
23 cm

Narrow jaws
Smooth inserts.



07-138
13 cm
07-139
15 cm
07-140
18 cm
07-141
20 cm
07-142
23 cm

Narrow jaws
serrated



French-Eye Needle Holder

07-143 16 cm
07-144 18 cm
07-145 20 cm
07-146 22 cm
07-147 25 cm
Narrow Tip



French-Eye Needle Holder

07-148 16 cm
07-149 18 cm
07-150 20 cm
07-151 22 cm



**Cooley Intracardiac
Needle Holders**
Straight, narrow jaw.



07-152
18 cm
07-153
20 cm
07-154
23 cm
07-155
26 cm



**Debakey Cardiovascular
Needle Holders**
Delicate patterns.
Straight jaw, serrated.



07-156
18 cm
07-157
23 cm
07-158
27 cm
07-159
30 cm

NEEDLE HOLDERS




Nadelhalter	Porta-agujas
Prote-aiguilles	Porta aghi



Finochietto Needle Holders
 07-160 13 cm
 07-161 16 cm
 07-162 19 cm
 07-163 21 cm



Finochietto Needle Holders
 07-164 13 cm Angled
 07-165 16 cm Angled
 07-166 19 cm Angled



07-167
 20 cm
 Curved



07-168
 14 cm
07-169
 17 cm
07-170
 20 cm




07-171
 14 cm
07-172
 16.5 cm
07-173
 20 cm
 Narrow Jaws

Mathieu Needle Holders



07-174
 18 cm
07-175
 23 cm

**Mayo Rochester
 Needle Holders**
 button release, screw joint.



07-176
 15 cm
07-177
 18 cm
07-178
 20 cm
 lined with
 copper.

Macphail Needle Holders



07-179
 15 cm

Crile Needle Holders

NEEDLE HOLDERS

Nadelhalter
Prote-aiguilles

Porta-agujas
Porta aghi



07-180
18.5 cm

Toennis Needle Holders



07-181
13 cm
with lock
straight



07-182
13 cm
with lock
curved



07-183
13 cm
without lock
straight



07-184
13 cm
without lock
curved

Castroviejo Needle Holders



07-185
10 cm
07-186
14 cm
Straight, Smooth
with lock Straight

07-187
10 cm
07-188
14 cm
Curved, Smooth
with lock curved

**Troutman
Needle Holders**



**Barraquer Micro
Needle Holders**
07-189 11cm
07-190 13cm



**Barraquer
Needle Holders**



07-191
13 cm
Straight, Smooth
without lock
Straight



07-192
13 cm
Curved, Smooth
without lock
curved



**Castroviejo
Needle Holders**



07-193
14 cm
Straight
Smooth
with lock



07-194
14 cm
Straight
Smooth
without lock



07-195
14 cm
Curved
Smooth
with lock



07-196
14 cm
Curved
Smooth
without lock



07-197
18 cm
Straight
Smooth
with lock



07-198
18 cm
Straight
Smooth
without lock



07-199
18 cm
Curved
Smooth
with lock



07-200
18 cm
Curved
Smooth
without lock

NEEDLE HOLDERS

Nadelhalter | Porta-agujas
Prote-aiguilles | Porta aghi



Castroviejo Olsen Hegar /
Needle Holders / Scissors



07-201
13 cm
07-202
15 cm



Boynton Needle Holders



07-203
12 cm



Arruga Needle Holders



07-204
16 cm
Str.



07-205
16 cm
Cvd.



Micro Needle Holders



07-206
19 cm
07-207
23 cm
Smooth
Str.



07-208
19 cm
07-209
23 cm
Smooth
Cvd.



Micro Needle
Holders



07-210
16 cm
Smooth
Cvd.



Castroviejo Needle
Holders



07-211
12 cm
Smooth Str
with lock



07-212
12 cm
Smooth Cvd
with lock



Barraquer-Troutman
Needle Holders



07-213
10 cm
07-214
11.5 cm
Smooth
Str.



07-215
10 cm
07-216
11.5 cm
Smooth
Cvd.



Kalt Needle Holders



07-217
14 cm
Serrated



Arruga Needle Holders



07-218
14 cm
Serrated
Str.



07-219
14 cm
Serrated
Cvd.


NEEDLE HOLDERS T.C.

Nadelhalter	Porta-agujas
Prote-aiguilles	Porta aghi



07-220 14 cm
07-221 18 cm
07-222 20 cm
07-223 23 cm
07-224 25 cm

Microvascular
Needle Holders T.C.



07-225 15 cm
07-226 18 cm
07-227 20 cm

Micro Ryder Needle
Holders T.C.



07-228 13 cm
07-229 15 cm
07-230 18 cm
07-231 20 cm
07-232 22 cm
07-233 25 cm

Crile-Wood Needle
Holders T.C.




07-234 18 cm
07-235 20 cm
07-236 22 cm
07-237 25 cm

Intracardiac Needle Holders T.C.



07-238 18 cm

New Orleans Needle
Holders T.C.



07-239 18 cm
Serrated-Solid



07-240 15 cm
07-241 17 cm

Senning Needle Holders T.C.
(Schwed. Model)



07-242 15 cm
Serrated-Solid

Carroll Needle Holders T.C.

NEEDLE HOLDERS T.C.

Nadelhalter
Prote-aiguilles

Porta-agujas
Porta aghi




Wertheim Needle Holders T.C
07-243 20 cm
07-244 24 cm



Kilner Needle Holders T.C



Jacobson Needle Holders T.C



Bozemann Needle Holders T.C



Baumgartner Needle Holders T.C



Olsen-Hegar Needle Holders T.C

07-247
20 cm
07-248
24 cm
07-249
26 cm
07-250
30 cm

07-251
13 cm
07-252
15 cm

07-253
14 cm
07-254
17 cm
07-255
19 cm
Serrated-
Solid

07-256
14 cm
07-257
17 cm
07-258
19 cm
Serrated-
Solid

NEEDLE HOLDERS T.C.

Nadelhalter
Prote-aiguilles

Porta-agujas
Porta aghi



French-Eye Needle Holders T.C

French-Eye Needle Holders T.C

Sarot Needle Holders T.C

Wangensteen Needle Holders T.C



Finochietto Needle Holders T.C

Finochietto Needle Holders T.C

Mayo-Hegar Needle Holders T.C
Also Available Right Handle

NEEDLE HOLDERS T.C.

Nadelhalter Porta-agujas
Prote-aiguilles Porta aghi



07-298
20 cm
Curved



07-299
18 cm



07-300
20 cm

07-301
23 cm

07-302
25 cm
Curved



07-303
20 cm

07-304
23 cm

07-305
25 cm
Curved

Finochietto Needle Holders T.C

Rochester Needle Holders T.C

Stratte Needle Holders T.C



07-306
16 cm



07-307
12.5 cm
Smooth

07-309
12.5 cm
Tungsten
Carbide
Smooth



07-308
12.5 cm
Serrated-Solid

07-310
12.5 cm
Tungsten
Carbide
Serrated-Solid



07-311
14 cm
Serrated-Delicate

Gillies Needle Holders T.C


Halsey Needle Holders T.C

Eufrete Pasque T.C
Needle Holders


NEEDLE HOLDERS T.C.

Nadelhalter
Prote-aiguilles


Porta-agujas
Porta aghi



Ryder Needle Holders T.C
Narrow jaws.



Cooley Intracardiac Needle Holders T.C
Narrow Jaw.



Baby Crile-Wood Needle Holders T.C




Ryder French-Eye/Standard/Micro Needle Holders T.C



Converse Needle Holders



Derf Needle Holders T.C



NEEDLE HOLDERS T.C.

Nadelhalter
Prote-aiguilles

Porta-agujas
Porta aghi



De Bakey Needle Holders T.C
07-341 18 cm
07-342 20 cm
07-343 22 cm
07-344 25 cm



Debaquey Cardiovascular Needle Holders T.C
Delicate patterns.
serrated.



Debaquey-Vascular/Extra Light Needle Holders T.C



Gillies Slim T.c Needle Holders with Scissors



Wire Twisters Needle Holders T.C



Sternum Needle Holders T.C



Cooley-Baumgaften Wire Twister Needle Holders T.C
Extra heavy needle holder designed for twisting of wire and placement of retention suture. Heavy, blunt jaws.

07-355
15 cm
Delicate.

07-356
15 cm
07-357
18 cm
07-358
20 cm
Heavy

07-359
18 cm
Wire Twister

07-360
20 cm

NEEDLE HOLDERS T.C.


Nadelhalter	Porta-agujas
Prote-aiguilles	Porta aghi




Rubio Wire Twister
Needle Holders T.C



Neivert Needle Holders T.C



Webster Needle Holders T.C



Heaney-Gyne Needle Holders T.C



Wire Tightening
Needle Holders T.C
07-370 18 cm



Julian Needle Holders T.C

NEEDLE HOLDERS T.C.



Nadelhalter
Prote-aiguilles

Porta-agujas
Porta aghi



07-373
27 cm

07-374
30 cm

07-375
15 cm
07-376
17 cm
07-377
20 cm

07-378
20 cm
07-379
25 cm

Masson Needle Holders T.C O'Brien Needle Holders T.C

Martin
Needle Holder mini profile elastic type

Thomson-Walker
Needle Holders T.C



07-382
18 cm
07-383
22 cm

07-384
14 cm
Straight

07-385
14 cm
Curved

Hepp-Scheidel
07-380
Micro Needle Holder smooth jaw
18 cm,

Lichtenberg
07-381 20 cm
Needle Holder lightweight model
standard profile

Castroviejo Needle
Holders T.C

Castroviejo Needle
Holders T.C

NEEDLE HOLDERS T.C.



Nadelhalter | Porta-agujas
Prote-aiguilles | Porta aghi



**Jacobson
Needle Holders T.C**



**Castroviejo
Needle Holders T.C**



Micro Needle Holder
07-398 15 cm
07-399 18 cm
07-400 21 cm
07-401 23 cm




Toennis Needle Holders T.C
07-402 18.5 cm



Boynton Needle Holders T.C
07-403 12 cm



Arruga Needle Holders T.C
07-404 16 cm Str.
07-405 16 cm Cvd.



**Spencer Coronary
Needle Holders T.C**
07-406 18 cm



**Jaboma
07-407**
Micro Needle Holder
18 cm



NEEDLE HOLDERS T.C.

Nadelhalter Porta-agujas
Prote-aiguilles Porta aghi



Castroviejo Micro Needle Holder
07-408 14.5 cm
07-409 18 cm
07-410 21.5 cm



Barraquer
07-411
Micro Needle Holder
13 cm



Barraquer
07-412
Micro Needle Holder
13 cm



Mathieu Needle Holder
07-413 14 cm
07-414 17 cm
07-415 20 cm



Mathieu Needle Holder
07-416 14 cm
07-417 17 cm
07-418 20 cm



Mathieu Needle Holder
07-419 17 cm
07-420 20 cm
07-421 23 cm
standard profile with interior
ratchet



Castroviejo
07-422
Micro Needle Holder
12.5 cm



Kalt Needle Holders T.C
07-423 14 cm Str.



**Mathieu-Ryder Needle
Holders T.C**



07-424
17 cm
07-425
20 cm
07-426
24 cm



07-427
17 cm
07-428
20 cm
07-429
24 cm

**Mathieu-Ryder Vascular
Needle Holders T.C**

NEEDLE HOLDERS T.C.

Nadelhalter
Prote-aiguilles

Porta-agujas
Porta aghi



07-430
15 cm



07-431
15 cm



07-432
17 cm



07-433
20 cm



**Langenbeck
Needle
Holders T.C**
07-434
16 cm
07-435
18 cm
07-436
20 cm



**Langenbeck-
Delicate
Needle
Holders T.C**
07-437
18 cm
07-438
20 cm
07-439
23 cm



**Langenbeck-
Ryder
Needle
Holders T.C**
07-440
18 cm
07-441
20 cm
07-442
23 cm



Corwin Needle Holders T.C

Lichtenberg Needle Holders T.C



07-443
20 cm



07-444
23 cm



07-445
18.5 cm

Eiselsberg-Matheu
Needle Holders T.C



07-446
20 cm

Luethy-Beck
Needle Holders T.C

Zweifel Needle Holders T.C

REVERDIN NEEDLES



Reverdin nadeln
aguja de reverdin

Aguilles de reverdin
aghi di reverdin



Reverdin Needle
short



07-447
14.5 cm



07-448
14.5 cm



07-449
14.5cm



07-450
14.5cm



21cm



2
07-451



1
07-452



4
07-453



3
07-454



Reverdin Needle
19.5 cm



07-455
19.5cm
07-456
23cm



07-457
19.5cm
07-458
23cm



07-459
19.5cm
07-460
23cm



07-461
19.5cm
07-462
23cm



Reverdin Needle
23 cm



07-463
25cm



07-464
23cm
07-465
28cm



Konig
07-466
20 cm



Payr
07-467
22 cm



07-468
19.5cm
/ 3 mm

Konig Ligature
Conductor







Michel
07-469

REVERDIN NEEDLES

Reverdin nadeln
aguja de reverdin

Aguilles de reverdin
aghi di reverdin



07-483
20cm
sharp
for left
hand

Deschamps Needle




07-485
27cm
blunt
for left
hand

Deschamps Needle



07-487
27cm
sharp
for left
hand

Deschamps Needle



07-489
28cm
blunt
for left
hand

Deschamps Needle

LIGATUTR NEEDLES



Unterbindungsnadeln
agujas para ligaduras

Aiguilles à ligatures
aghi di legatura




Deschamps Needle
07-490 20cm sharp for right hand
07-491 20cm sharp for left hand



Deschamps Needle
07-492 20cm blunt for right hand
07-493 20cm blunt for left hand



Deschamps Needle
07-494 21cm blunt for right hand
07-495 21cm blunt for left hand



Deschamps Needle
07-496 20cm blunt for right hand
07-497 20cm blunt for left hand



Deschamps Needle
07-498 20cm sharp for right hand
07-499 20cm sharp for left hand



Deschamps Needle
07-500 20cm blunt for right hand
07-501 20cm blunt for left hand

REVERDIN NEEDLES

Reverdin nadeln
aguja de reverdin

Aguilles de reverdin
aghi di reverdin




SURGEONS NEEDLES

- 07-502 Fistula Needle, 1/2 Circle, Cutting Edge, #00
- 07-503 Fistula Needle, 1/2 Circle, Cutting Edge, #0
- 07-504 Fistula Needle, 1/2 Circle, Cutting Edge, #1
- 07-505 Fistula Needle, 1/2 Circle, Cutting Edge, #2
- 07-506 Fistula Needle, 1/2 Circle, Cutting Edge, #3
- 07-507 Fistula Needle, 1/2 Circle, Cutting Edge, #4
- 07-508 Fistula Needle, 1/2 Circle, Cutting Edge, #5
- 07-509 Fistula Needle, 1/2 Circle, Cutting Edge, #6
- 07-510 Fistula Needle, 1/2 Circle, Cutting Edge, #7
- 07-511 Fistula Needle, 1/2 Circle, Cutting Edge, #8
- 07-512 Fistula Needle, 1/2 Circle, Cutting Edge, #9



- 07-513 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #0
- 07-514 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #1
- 07-515 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #2
- 07-516 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #3
- 07-517 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #4
- 07-518 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #5
- 07-519 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #6
- 07-520 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #7
- 07-521 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #8
- 07-522 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #9
- 07-523 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #10
- 07-524 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #11
- 07-525 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #12
- 07-526 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #13
- 07-527 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #14
- 07-528 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #15
- 07-529 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #16
- 07-530 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #17
- 07-531 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #18
- 07-532 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #19
- 07-533 Regular Surgeons Needle, 1/2 Circle, Cutting Edge, #20



3/8 Circle - Cutting Edge

SURGEONS NEEDLES

- 07-534 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #1
- 07-535 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #2
- 07-536 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #3
- 07-537 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #4
- 07-538 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #5
- 07-539 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #6
- 07-540 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #7
- 07-541 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #8
- 07-542 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #9
- 07-543 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #10
- 07-544 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #11
- 07-545 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #12
- 07-546 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #13
- 07-547 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #14
- 07-548 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #15
- 07-549 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #16
- 07-550 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #17
- 07-551 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #18
- 07-552 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #19
- 07-553 Regular Surgeons Needle, 3/8 Circle, Cutting Edge, #20




Keith Abdominal Straight, Triangular Point

- 07-566 Keith Abdominal Needles, Straight, Triangular Point, 4"
- 07-567 Keith Abdominal Needles, Straight, Triangular Point, 3-1/2"
- 07-568 Keith Abdominal Needles, Straight, Triangular Point, 3"
- 07-569 Keith Abdominal Needles, Straight, Triangular Point, 2-7/8"
- 07-570 Keith Abdominal Needles, Straight, Triangular Point, 2-1/8"
- 07-571 Keith Abdominal Needles, Straight, Triangular Point, 2"
- 07-572 Keith Abdominal Needles, Straight, Triangular Point, 1-3/4"
- 07-573 Keith Abdominal Needles, Straight, Triangular Point, 1-1/2"

SURGEONS NEEDLES

- 07-554 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #0
- 07-555 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #1
- 07-556 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #2
- 07-557 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #4
- 07-558 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #6
- 07-559 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #8
- 07-560 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #10
- 07-561 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #12
- 07-562 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #14
- 07-563 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #16
- 07-564 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #20
- 07-565 Regular Surgeons Needles, 1/2 Curved, Cutting Edge, #22




half Curved - Cutting Edge

REVERDIN NEEDLES

Reverdin nadeln
aguja de reverdin

Aguilles de reverdin
aghi di reverdin




- 07-574** Regular Surgeons Needles, 1/2 Curved, Tapered Point, #2
07-575 Regular Surgeons Needles, 1/2 Curved, Tapered Point, #4
07-576 Regular Surgeons Needles, 1/2 Curved, Tapered Point, #6
07-577 Regular Surgeons Needles, 1/2 Curved, Tapered Point, #8
07-578 Regular Surgeons Needles, 1/2 Curved, Tapered Point, #10
07-579 Regular Surgeons Needles, 1/2 Curved, Tapered Point, #12
07-580 Regular Surgeons Needles, 1/2 Curved, Tapered Point, #14
07-581 Regular Surgeons Needles, 1/2 Curved, Tapered Point, #16
07-582 Regular Surgeons Needles, 1/2 Curved, Tapered Point, #18
07-583 Regular Surgeons Needles, 1/2 Curved, Tapered Point, #20
07-584 Regular Surgeons Needles, 1/2 Curved, Tapered Point, #22



SURGEONS NEEDLES

- A 07-585** Veterinary Needles, 3/8 Circle, Cutting Edge, #1
07-586 Veterinary Needles, 3/8 Circle, Cutting Edge, #2
07-587 Veterinary Needles, 3/8 Circle, Cutting Edge, #3
07-588 Veterinary Needles, 3/8 Circle, Cutting Edge, #4
07-589 Veterinary Needles, 3/8 Circle, Cutting Edge, #5
07-590 Veterinary Needles, 3/8 Circle, Cutting Edge, #6
07-591 Veterinary Needles, 3/8 Circle, Cutting Edge, #7
07-592 Veterinary Needles, 3/8 Circle, Cutting Edge, #8



- B 07-593** Post Mortem Needles, Double Curve Cutting Edge, Size 1
07-594 Post Mortem Needles, Double Curve Cutting Edge, Size 2
07-595 Post Mortem Needles, Double Curve Cutting Edge, Size 3
07-596 Post Mortem Needles, Double Curve Cutting Edge, Size 4
07-597 Post Mortem Needles, Double Curve Cutting Edge, Size 5
07-598 Post Mortem Needles, Double Curve Cutting Edge, Size 6
07-599 Post Mortem Needles, Double Curve Cutting Edge, Size 7
07-600 Post Mortem Needles, Double Curve Cutting Edge, Size 8

SURGEONS NEEDLES

- A 07-601** Martin Uterine Needles, 1/2 Curved Cutting Edge, #1
07-602 Martin Uterine Needles, 1/2 Curved Cutting Edge, #2
07-603 Martin Uterine Needles, 1/2 Curved Cutting Edge, #3
07-604 Martin Uterine Needles, 1/2 Curved Cutting Edge, #4
07-605 Martin Uterine Needles, 1/2 Curved Cutting Edge, #5
07-606 Martin Uterine Needles, 1/2 Curved Cutting Edge, #6
07-607 Martin Uterine Needles, 1/2 Curved Cutting Edge, #7
07-608 Martin Uterine Needles, 1/2 Curved Cutting Edge, #8

- B 07-609** Lane Cleft Palate, 1/2 Circle, Cutting Edge, #1
07-610 Lane Cleft Palate, 1/2 Circle, Cutting Edge, #2
07-611 Lane Cleft Palate, 1/2 Circle, Cutting Edge, #3
07-612 Lane Cleft Palate, 1/2 Circle, Cutting Edge, #4
07-613 Lane Cleft Palate, 1/2 Circle, Cutting Edge, #5

- C 07-614** Post Mortem Needles, Half Curved Cutting Edge, Size 1
07-615 Post Mortem Needles, Half Curved Cutting Edge, Size 2
07-616 Post Mortem Needles, Half Curved Cutting Edge, Size 3
07-617 Post Mortem Needles, Half Curved Cutting Edge, Size 4
07-618 Post Mortem Needles, Half Curved Cutting Edge, Size 5
07-619 Post Mortem Needles, Half Curved Cutting Edge, Size 6



Martin, 1/2 Circle Cutting Edge

A



Left Cleft Palate, 1/2 Circle Cutting Edge


B



Post Mortem - Half Curved Cutting Edge

C

RETRACTORS INSTRUMENTS



HOOKLET & SKIN HOOKS



Wundhäkchen, Mikrohäkchen
Ganchito, Micro Ganchitos

Crochet, Micro crochetes
Divaricatore, Micro Divaricatori



**GUTHRIE
HOOKLET**



**GILLIES
SKIN HOOK**



KILNER HOOK



**FISCH HOOKLET
& Skin Hook**



**FRAZIER HOOKLET
& Skin Hook**



**DOUBLE PRONG MICRO
HOOK & Skin Hook**
08-012
Sharp, 2mm, wide jaws,
14.5 cm, 6 3/4"



**WIENER HOOK
& Skin Hook**
08-013
Sharp,
14 cm, 5 1/2"



**TOENNIS HOOK
& Skin Hook**
08-014
15 cm, 6"
Sharp



ROLLET HOOK



**FREER SKIN
HOOK**

HOOKLET & SKIN HOOKS

Wundhäkchen, Mikrohäkchen
Ganchito, Micro Ganchitos

Crochet, Micro crochets
Divaricatore, Micro Divaricatori



**CRILE HOOKLET
& SKIN HOOKS**
08-020 20.4 cm, 8"



**VON GRAEFE Strabismus
& SKIN HOOKS**
08-021 13.0 cm, 5" - A
08-022 14.0 cm, 5 1/2" - B
08-023 15.0 cm, 6" - C



**WAGENER HOOKLET
& SKIN HOOKS**
08-024 14 cm, 5 1/2" - A
08-025 14 cm, 5 1/2" - B
08-026 14 cm, 5 1/2" - C
08-027 14 cm, 5 1/2" - D




**DAY HOOKLET
& SKIN HOOKS**
08-028 16 cm, 6 1/4" - A
08-029 16 cm, 6 1/4" - B



**LUCAE HOOKLET
& SKIN HOOKS**
08-030 14 cm, 5 1/2" - A
08-031 14 cm, 5 1/2" - B
08-032 14 cm, 5 1/2" - C



Converse Hook
08-033
12 cm, 5"



MANNERFELT-LANGENBECK HOOK
08-034 15 x 6 mm 15.5 cm, 6" Fig.1
08-035 20 x 6 mm 15.5 cm, 6" Fig.2
08-036 25 x 6 mm 15.5 cm, 6" Fig.3
08-037 30 x 8 mm 15.5 cm, 6" Fig.4



MANNERFELT HOOK
15.5 cm, 6"
large curve, single
hook



08-038
Sharp



08-039
2 prong
Blunt



08-040
3 prong
Blunt



**Brom Hooklet
Retractors**
08-041
17 cm, 6 3/4"




**Hooklet
Retractors**
08-042
13 cm, 5"

HOOKLET, RETRACTORS



Wundhäkchen, Wundhaken
Ganchito, Separadores


Crochet, Ecarteurs
Divaricatore, Divaricatori




Tracheal Retractor
18 cm



Tracheal Retractor
18 cm



Trachea Retractor
16.5 cm Sharp



Tracheal Retractor
16.5 cm Blunt



Rollet Hooklet
Retractors
08-057
12.5cm



Rollet (Blair)
Hooklet Retractors
08-058
13.5cm




Paul Hooklet
Retractors
08-059
12.5cm



Paul Hooklet
Retractors
08-060
16cm



Miller Retractor
08-061
13.5cm
Sharp 6 prong




Miller Retractor
13.5 cm

HOOKLET, RETRACTORS

Wundhäkchen, Wundhaken
Ganchito, Separadores

Crochet, Ecarteurs
Divaricatore, Divaricatori




Hooklet Retractors
17 cm

Finger Rake, Retractor
16.5 cm

Finger Rake, Retractor
16.5 cm

Hooklets, Retractors
flexible,
16 cm



Peck Rake Retractor
16.5 cm

Kocher Hooklet,
Retractors

Kocher Hooklet,
Retractors

Hooklet,
Retractors
08-090
17 cm
Sharp

HOOKLET, RETRACTORS



Wundhäkchen, Wundhaken
Ganchito, Separadores

Crochet, Ecarteurs
Divaricatore, Divaricatori



Volkman Retractors
08-091 23 cm
Sharp



Volkman Retractors
08-092 23 cm
Blunt



Lambotte Retractors
20 cm

08-093
13 x 235 mm
08-094
20 x 255 mm
08-095
27 x 270 mm



Volkman Retractors
11.5 cm



08-096
2 prong

08-100
2 prong



08-097
3 prong

08-101
3 prong



08-098
4 prong

08-102
4 prong



08-099
6 prong

08-103
6 prong



Murphy Retractors
19.1 cm 7 1/2"



Senn Retractor
08-108
16 cm 6 1/4"
Blunt



Senn Retractor
08-109
16 cm 6 1/4"
Sharp



Senn-Mueller Retractor
08-110
18 x 6 mm
16 cm, 6 1/4"
sharp



Volkman Retractor
08-111
22 cm 8 3/4"
4 Prongs Sharp



Volkman Retractor
08-112
22 cm 8 3/4"
4 Prongs Blunt

HOOKLET, RETRACTORS

Wundhäkchen, Wundhaken
Ganchito, Separadores

Crochet, Ecarteurs
Divaricatore, Divaricatori



Sternberg Retractor
08-113 12 cm, 4 ¾"
08-114 14 cm, 5 ½"
08-115 16 cm, 6 ¼"



Collin-Hartmann Retractor
08-116 38 x 28 / 32 x 21 mm 14 cm, 5 ½"
08-117 52 x 34 / 36 x 27 mm 16 cm, 6 ¼"



Herzfeld Retractor
08-118 20 cm, 8"



Strandell-Stille Retractor
08-119 19 cm, 7 ½"



Strandell-Stille Retractor
08-120 19 cm, 7 ½"



Schoenbore Retractor
08-121 21 cm, 8 ¼"



Konig Retractor
17 cm, 6 ¾"



08-122
13 x 13 mm



08-123
12 x 16 mm



Senn-Green Retractor
15 cm, 6"



08-124
10 x 6 mm



08-125
20 x 6 mm



Cottle Retractor
08-126 15 cm, 6"
10 x 12 mm





Ragnell-Davis Retractor
08-127 14 cm, 5 ½"
08-128 15 x 5.5 mm

RETRACTORS




Wundhaken | Ecarteurs
Separadores | Divaricatori





RETRACTORS

Wundhaken | Ecarteurs
Separadores | Divaricatori



- 08-174** Blunt, 2 prongs.
- 08-175** Blunt, 3 prongs.
- 08-176** Blunt, 4 prongs.
- 08-177** Blunt, 5 prongs.
- 08-178** Blunt, 6 prongs.
- 08-179** Sharp, 2 prongs.
- 08-180** Sharp, 3 prongs.
- 08-181** Sharp, 4 prongs.
- 08-182** Sharp, 5 prongs.
- 08-183** Sharp, 6 prongs.

Volkmann Retractors
Long Pattern
22.2 cm, 8 1/4"



Koerte Retractor
08-195 28cm
25 x 40 mm
semi-sharp

- 08-184** 22 cm, 8 3/4" 15 x 10 mm 2 prong Fig. 1
- 08-185** 22 cm, 8 3/4" 15 x 15 mm 3 prong Fig. 2
- 08-186** 22 cm, 8 3/4" 15 x 20 mm 4 prong Fig. 3
- 08-187** 22 cm, 8 3/4" 15 x 30 mm 6 prong Fig. 4
- 08-188** 22 cm, 8 3/4" 15 x 40 mm 8 prong Fig. 5

- 08-189** 22 cm, 8 3/4" 15 x 10 mm 2 prong Fig. 1
- 08-190** 22 cm, 8 3/4" 15 x 15 mm 3 prong Fig. 2
- 08-191** 22 cm, 8 3/4" 15 x 20 mm 4 prong Fig. 3
- 08-192** 22 cm, 8 3/4" 15 x 30 mm 6 prong Fig. 4
- 08-193** 22 cm, 8 3/4" 15 x 40 mm 8 prong Fig. 5

Korte-Wagner Retractor
08-194 28cm
Sharp

RETRACTORS



Wundhaken | Ecarteurs
Separadores | Divaricatori





FIMA INSTRUMENTS

RETRACTORS

Wundhaken
Separadores

Ecarteurs
Divaricatori



Kocher Retractor
28 cm 11"



08-211



08-213



08-212



08-214



Kocher Retractor
08-215
22.5 cm 8 3/4"
70 x 25 mm



**Kocher-Wagner
Retractor**
28 cm, 11"



08-216
22 x 15 mm



08-217
38 x 20 mm



08-218
67 x 25 mm



Kocher Retractor
23 cm, 9"

08-219
40 x 18 mm

08-220
60 x 20 mm

08-221
60 x 25 mm

08-222
75 x 30 mm

08-223
75 x 40 mm

RETRACTORS



Wundhaken | Ecarteurs
Separadores | Divaricatori



Kocher-Langenbeck Retractor
21.5 cm, 8 1/2"



08-224
25 x 6 mm



08-225
35 x 8 mm



08-226
35 x 11 mm



08-227
35 x 15 mm



08-228
40 x 11 mm



08-229
55 x 11 mm



08-230
70 x 14 mm



08-231
80 x 12 mm



08-232
80 x 16 mm



Langenbeck Retractor
24 cm, 8 1/4"



08-233
35 x 7mm



08-234
35 x 15mm



08-235
35 x 11mm



08-236
35 x 20mm



Langenbeck Retractor 08-237
21.5 cm 8 1/2"
Handle vertical to the blade



Langenbeck Retractor 08-238
21 cm 8 1/4"
38 x 9mm



Caspar Retractor
21 cm, 8 1/4"

RETRACTORS

Wundhaken
Separadores

Ecarteurs
Divaricatori



Langenbeck Retractor
20.4 cm, 8"



08-241
20 x 7 mm



08-242
30 x 7 mm



08-243
30 x 11 mm



08-244
30 x 15 mm




08-245
30 x 20 mm



08-246
44 x 10 mm




08-247
44 x 15 mm



08-248
45 x 20 mm



08-249
50 x 20 mm



08-250
60 x 20 mm



Meyerding Finger Retractor
18 cm, 7"



08-251
long blunt blade,
6 x 17 mm



08-254
long toothed blade,
5 x 16mm



08-252
short toothed blade,
3 x 6 mm



08-255
short blunt blade,
3 x 10mm



08-253
curved blunt blade,
3 x 8mm
Saddle Shaped



08-256
curved three prong
sharp blade,
10 x 4mm



Langenbeck Retractor
22 cm, 8 1/2"



08-257
25 x 6 mm



08-260
35 x 15 mm



08-264
25 x 6 mm



08-267
35 x 15 mm



08-258
35 x 8 mm



08-261
41 x 11 mm



08-265
35 x 8 mm



08-268
41 x 11 mm



08-259
35 x 11 mm



08-262
55 x 11 mm



08-266
35 x 11 mm



08-269
55 x 11 mm



08-263
70 x 14 mm



Langenbeck Retractor
22 cm, 8 1/2"



08-270
70 x 14 mm

RETRACTORS



Wundhaken | Ecarteurs
Separadores | Divaricatori




Blount Knee Retractor
08-271
17.5 cm, 6 3/4"
42 x 4 mm



Taylor Spinal Retractors
08-272 16.5 cm, 6 1/2"
75 x 30 mm / 100 x 30 mm



Luer Retractor
08-273 11 cm, 4 1/4"
10 x 8 mm / 12 x 11 mm



Rose Retractor
08-274 13 cm, 5"
12 x 12 mm /
12 x 15 mm



Crile Retractor
08-275 11 cm, 4 1/4"
12 x 4 mm /
20 x 6 mm



Mathieu Retractor
Set 2 pcs.
08-276 20 cm, 8"
47 x 13 / 46 x 26 mm
08-277 20 cm, 8"
42 x 13 / 42 x 26 mm



Us-Army Retractor
Set 2 pcs.
08-278 21.5 cm, 8 1/2"
28 x 15 / 42 x 15 mm
08-279 21.5 cm, 8 1/2"
25 x 15 / 38 x 15 mm



Richardson Eastman Retractor
08-280
26 cm, 10 1/4"
small



Richardson Eastman Retractor
08-281
28 cm, 11"



Sofield Retractor
08-282
set of 6 pec.

RETRACTORS

Wundhaken
Separadores

Ecarteurs
Divaricatori



Roux Retractor

Set 3 pcs. consisting of

08-283 14 cm, 5 1/2" 24 x 20 / 30 x 28 mm

08-284 15 cm, 6" 29 x 26 / 36 x 36 mm

08-285 17 cm, 6 3/4" 33 x 29 / 38 x 43 mm



Baby-Roux Retractor

08-286 13 cm, 5"

15 x 15 mm/

18 x 21 mm



Roux Retractor

08-287 16.5 cm, 6 1/2"

28 x 23 mm/

43 x 23 mm



Roux Retractor

08-288 18.5 cm, 7 1/4"

28 x 27 mm/

45 x 27 mm



Mini-Roux Retractor

08-289

14 cm, 5 1/2"




Roux-Magnus Retractor

08-290 15 cm, 6"

20 x 22 mm/

25 x 38 mm




Roux-Magnus Retractor

08-291 16 cm, 6 1/4"

23 x 26 mm/

29 x 38 mm



Roux-Magnus Retractor

08-292 17 cm, 6 3/4"

29 x 29 mm/

34 x 44 mm



08-293

Baby

(7x100mm)

08-294

infantil

(10x120mm)

08-295

Adult

(13x125mm)

08-296

Adult

(15x150mm)

08-297

Extra Large

(20x180mm)

Farabeuf Retractor

RETRACTORS



Wundhaken | Ecarteurs
Separadores | Divaricatori



Farabeuf Retractor
Set 2 Pcs.
08-298 13 cm, 5"
30 x 10 / 34 x 13 mm
08-299 13 cm, 5"
26 x 10 / 30 x 13 mm




Farabeuf Retractor
Set 2 Pcs.
08-300 15 cm, 6"
27 x 16 / 30 x 16 mm
08-301 15 cm, 6"
23 x 16 / 26 x 16 mm



Parker Retractor
Set 2 pcs.
08-302 13.5 cm, 5 1/4"
22 x 15 mm / 22 x 15 mm
08-303
18.4 cm 7 1/4"



Parker Retractor
08-304 13.5 cm, 5 1/4"
18 x 15 mm / 19 x 15 mm



Collin Retractor
Set 2 pcs.
08-305 16.5 cm, 6 1/2"
25 x 12 / 35 x 16 mm
08-306 16.5 cm, 6 1/2"
20 x 12 / 30 x 16 mm



Mayo-Collin Retractor
Set 2 pcs.
08-307 15 cm, 6"
28 x 18 / 25 x 18 mm



Parker-Langenbeck Retractor
Set of 2 Pcs.
08-308 25 x 15 mm fig. 2
08-309 21 x 15 mm fig. 1
08-310 40 x 15 mm fig. 1
08-311 46 x 15 mm fig. 2



Parker-Mott Retractor
08-312
16 cm 6 1/4"
68 x 24 / 74 x 24 mm

RETRACTORS

Wundhaken
Separadores

Ecarteurs
Divaricatori

08-313
115 x 65 mm




08-314
115 x 38 mm



08-315
82 x 22 mm



08-316
130 x 38 mm



Coryllos Retractor
22.5 cm, 8 3/4"

08-317
72 x 20 mm

08-318
55 x 23 mm

08-319
45 x 15 mm

Sauerbruch Retractor
25 cm, 10"



Davidson Retractor
08-320 16 cm, 6 1/4" 75 x 50 mm
08-321 16 cm, 6 1/4" 75 x 90 mm

RETRACTORS




Wundhaken | Ecarteurs
Separadores | Divaricatori



Bruner Retractor
25 cm, 10"

- 08-322**
180 x 30 mm
- 08-323**
160 x 30 mm
- 08-324**
140 x 30 mm
- 08-325**
130 x 30 mm
- 08-326**
120 x 25 mm
- 08-327**
100 x 25 mm
- 08-328**
80 x 20 mm



Bruner Retractor
24 cm, 9 1/2"

- 08-329** 80 x 20 mm
- 08-330** 100 x 20 mm
- 08-331** 120 x 25 mm
- 08-332** 140 x 25 mm
- 08-333** 160 x 30 mm
- 08-334** 180 x 30 mm



Bruner Retractor
24 cm, 9 1/2"

- 08-335** 80 x 20 mm
- 08-336** 100 x 20 mm
- 08-337** 120 x 25 mm
- 08-338** 140 x 25 mm
- 08-339** 160 x 30 mm
- 08-340** 180 x 30 mm

RETRACTORS

Wundhaken
Separadores

Ecarteurs
Divaricatori



RETRACTORS

Wundhaken | Ecarteurs
Separadores | Divaricatori



Mikulicz Retractor
26 cm, 10 1/4"
08-351
85 x 50 mm



Kelly Retractor
27 cm, 10 1/2"
08-352
190 x 38 mm
08-353
190 x 57 mm



Mikulicz Retractor
26 cm, 10 1/4"
08-354
95 x 35 mm

Mikulicz Retractor
26 cm, 10 1/4"
08-355
120 x 50 mm
08-356
150 x 50 mm



Cheek-And Lip Retractor
19 cm, 7 1/2"
08-357
Movable



Tuffier Retractor
23 cm, 9"
08-358 55 x 63 mm



Bennet Retractor
08-359
24cm



Doyen Retractor
25 cm, 10"
08-360 50 x 85 mm



Fritch Retractor
08-361
2.75 x 2 mm



Fritsch Retractor
24 cm, 9 1/2"
08-362 35 x 40 mm

RETRACTORS

Wundhaken
Separadores

Ecarteurs
Divaricatori



Fritsch Retractor

08-363 24 cm, 9 1/2" 30 x 30 mm
08-364 24 cm, 9 1/2" 32 x 40 mm
08-365 24 cm, 9 1/2" 40 x 50 mm
08-366 24 cm, 9 1/2" 40 x 60 mm
08-367 24 cm, 9 1/2" 45 x 75 mm
08-368 24 cm, 9 1/2" 65 x 85 mm



Mikulicz Retractor
08-369

29 cm 11.5"



Tuffier Retractor

08-370 23 cm 9"
64 x 60mm



Harrington Retractor

32 cm, 12 3/4"
08-371
127 x 40 mm
08-372
127 x 62 mm



Harrington Extra Deep Retractor

08-373 30cm 11 3/4"
Hollow handle, 5.7cm
width x 20.5cm 8" depth blade



Harrington Splanchnic Retractor

08-374 31.5cm 12 1/2"
Hollow handle, 3.8cm 1 1/2"
width x 12.0cm 4 3/4"
depth blade



Debakey-Cooley Retractor

08-375 31.0cm 12 1/4"
Hollow handle, 5.7cm 2 1/4"
width x 18.4cm 7 1/4"
depth blade




Meyerdig Retractor

08-376
22.5 cm, 9"
75 x 25mm



08-377
50 x 16mm



08-378
75 x 15mm



Goelet Retractor

08-379
19 cm 7 1/2"
Double Ended

RETRACTORS, VEIN HOOKS



Wundhaken, Venenhäkchen

Separadores, Ganchitos Para venas

Ecarteurs, Ecarteurs sous

Divaricatori, Divaricatori per veni



08-380



Kaye-Millard Retractor
08-381
8 cm, 3 1/4"



Richardson/Kelly
Retractor
24 cm, 9 1/2"



Ross Retractor
25.5 cm, 10"




Weinberg Vagotomy
Retractor (joe's Hoe)
08-392 20 cm, 8"
Blade 102 x 165 mm wide

RETRACTORS

Wundhaken
Separadores


Ecarteurs
Divaricatori



Turner-Doyen Retractors

23.5 cm
08-393
Blade 140 x 60 mm wide
08-394
Blade 76 x 140 mm wide
08-395
Blade 165 x 76 mm wide

08-396 55 x 35 mm
08-397 85 x 35 mm
08-398 115 x 35 mm
22.5 cm, 9"
08-399 55 x 45 mm
08-400 85 x 45 mm
08-401 115 x 45 mm
08-402 160 x 45 mm
25 cm, 10"
08-403 55 x 60 mm
08-404 85 x 60 mm
08-405 115 x 60 mm
08-406 160 x 60 mm
23.4 cm, 9 1/4"



Doyen Vaginal Retractors



08-407 23cm
Small Blade 70x19mm wide
08-408 25cm
Medium Blade 82x29mm wide
08-409 28cm
Large Blade 22X64mm
08-410 28cm
Large Blade 38X89mm

Sawyer Retractors



08-411 23cm
Blade 60x24mm wide
08-412 25cm
Blade 75x35mm wide

Ferguson-Moon Rectal Retractors



08-413 23cm
Blade 64x24mm wide
08-414 23cm
Blade 76x32mm wide
08-415 25cm
Blade 76x35mm wide

Hill Ferguson Rectal Retractors



08-416 24cm
50x40mm wide
08-417 24cm
65x50mm wide
08-418 25cm
75x65mm wide
08-419 25cm
90x75mm wide

Kelly Retractors



Richardson Retractors



08-420 24cm
Blade 25x20mm wide
08-421 24cm
Blade 36x28mm wide
08-422 24cm
Blade 44x28mm wide
08-423 25cm
Blade 52x22mm wide
08-424 25cm
Blade 70x52mm wide

RETRACTORS




Wundhaken | Ecarteurs
Separadores | Divaricatori




Kelly Retractor
08-425 25cm
55x50mm



Kelly Retractor
26cm



08-426
70 X 65 mm



08-427
64 X 54 mm



08-428 70x40mm
08-429 70x50mm
08-430 70x65mm

Morris Retractor
24.5 cm 9 3/4"



Richardson Retractor
08-431 24cm



Richardson Retractor
24.0 cm 9.5"



08-432
23 x 20 mm




08-433
30 x 29 mm



08-434
37 x 37 mm




08-435
64 x 43 mm



38 x 37 mm



22 x 21 mm



08-436
39 x 30 mm

Richardson-Eastman Retractor
25.0 cm 10", Set of 2 piece

RETRACTORS

Wundhaken
Separadores

Ecarteurs
Divaricatori



Malleable Retractor

Item	Description	Tip
08-437	20 cm	12 mm
08-438	20 cm	17 mm
08-439	20 cm	25 mm
08-440	33 cm	25 mm
08-441	33 cm	32 mm
08-442	33 cm	40 mm
08-443	33 cm	50 mm



**Ribbon Retractor
08-444**
20 cm, 8"
Width 4 cm



**Ribbon Retractor
08-445**
20 cm, 8"
Width 2.5 cm



**Tuffier / Haberer
Retractor
08-446** 20 cm, 8"
17 x 25 mm
08-447 28 cm,
11" 25 x 30 mm
08-448 30 cm,
11 3/4" 37 x 45 mm
malleable



**Bruening
Spatula
08-449**
19 cm



**Flexival Abdominal
Spatula
08-450** 300x20mm
08-451 300x25mm
08-452 300x30mm
08-453 300x35mm
08-454 300x40mm
08-455 300x45mm
08-456 300x50mm



**Reverdin
Retractor
08-457**
28cm
46x60mm



**Reverdin
Retractor
08-458**
29cm



**Blood Scoop
08-459**
15cm



**Heiss Retractor
08-460** 10cm
Sharp 4:4
08-461 10cm
Blunt 4:4




**Finsen Retractor
08-462** 7 cm, 2 3/4"
3 x 4 teeth, sharp
08-463 7 cm, 2 3/4"
3 x 4 teeth, Blunt



**Finsen Retractor
08-464** 5 cm, 2"
2 x 3 teeth sharp/blunt
08-465 7 cm, 2 3/4"
2 x 3 teeth sharp/blunt

RETRACTORS

Wundhaken | Ecarteurs
Separadores | Divaricatori



Buckley Retractor
08-466 6cm



Lancaster Eye Speculum
08-467 8.3cm
blades 15 X 5 mm



Castroviejo Eye Speculum
08-468 9.5cm
medium blades
15 X 5 mm O.D.



ALM Retractor
08-469 7 cm, 2 3/4"
4 x 4 teeth sharp/blunt
08-470 10 cm, 4"
4 x 4 teeth sharp/blunt



1
1x1
Sharp Prong

Wire Retaining Retractor
08-471 8.5 cm, 3 1/2" 3 x 3 Prong
08-472 8.5 cm, 3 1/2" 1 x 1 Prong Sharp



Farr Wire Retractor
08-473 14 cm 5 1/2"
prongs 1.9 cm 3/4"

RETRACTORS

Wundhaken
Separadores

Ecarteurs
Divaricatori



Jansen Retractor
08-474 10 cm, 4"
Sharp 3x3
08-475 10 cm, 4"
Blunt 3x3



Houzel Retractor
08-476 29 cm, 11 1/2"



Mollison Retractor
08-477 15cm, 6"
sharp



Milligan Retractor
08-478 13 cm, 5"
18.5 mm x 20 mm
Blunt, 3x3



Markham/Meyerding Retractor
08-479 18cm, Blade 55 x 25 mm Left
08-480 18cm, Blade 55 x 25 mm Right



Darling Retractor
08-481 18 cm, 7"
40 x 76/25 x 54 mm
08-482 18 cm, 7"
25 x 64/40 x 76 mm
Laminectomy Retractors
1 blade toothed



Debaquey Gelpi Retractor
08-483 13.5cm 5 3/8"
Spread 3 1/8" 8cm, Sharp,

RETRACTORS



Wundhaken | Ecarteurs
Separadores | Divaricatori



Gelpi-Loktite Retractor

08-484 9cm
6 mm Sharp
08-485 11cm
10 mm Sharp
08-486 14cm
10 mm Sharp
08-487 16cm
10 mm Sharp
08-488 19cm
10 mm Sharp



Gelpi-Perineal Retractor
08-489 17.8cm
ball stops on sharp tenaculum
prongs, arm length 12cm



Gelpi-Loktite Retractor
08-490 9 cm, 6mm Sharp Cross Over
08-491 11 cm, 10mm Sharp Cross Over
08-492 14 cm, 10mm Sharp Cross Over
08-493 16 cm, 10mm Sharp Cross Over
08-494 19 cm, 10mm Sharp Cross Over



Seletz-Gelpi Retractor
08-495 17.8cm, 7"
modified with long
curve in shanks



Schuhknecht Retractor
08-496 Sharp
08-497 Blunt
14 cm, 5 1/2"



Wullstein Retractor
08-498 08-501 10.5 cm, 4"
08-499 08-502 13.0 cm, 5"
08-500 08-503 16.0 cm, 6 1/4"
Sharp Blunt



Travers Retractor
08-504
21 cm, 8 1/4" Blunt

RETRACTORS


Wundhaken | Ecarteurs
Separadores | Divaricatori



Righby Retractor
08-505 35 x 20 mm
08-506 60 x 20 mm
08-507 65 x 15 mm
18 cm, 7"



Weitlaner Retractor
08-508
13.0 cm 5.25"
Sharp



Weitlaner Retractor
08-509
16.5 cm 6.5"
Sharp



Plester Retractor
08-510
11.0 cm 4.25"
Blunt



Perkins Retractor
08-511 13cm



Weitlaner Retractor
08-512 11 cm
Sharp




Weitlaner Retractor
08-513 11 cm
Semi-Sharp

RETRACTORS



Wundhaken	Ecarteurs
Separadores	Divaricatori



Weitlaner Retractor
08-514 11cm
Blunt

Weitlaner Retractor
08-515 13cm
Sharp

Weitlaner Retractor
08-516 13cm
Semi-Sharp

Weitlaner Retractor
08-517 13cm
Blunt



Kolbel Soft Tissue Retractor
08-518 18.0cm 7"
Angled, 2 x 2 Pronged


Meyerding Retractor
08-519 16.5cm 6 1/2"
08-520 24.1cm 9 1/2" Large Blade
8.9 cm 3 1/2" X 5.1 cm 2"
08-521 24.1 cm 9 1/2" Small Blade
5.1 cm 2" X 1.6 cm 5/8"

D'Errico-Adson Retractor
08-522 19.0cm 7 3/8"
4 X 4-pronged,
2.0cm 3/4" depth blade
Sharp


RETRACTORS

Wundhaken
Separadores

Ecarteurs
Divaricatori




Anderson-Adson Retractor
08-523 19.0 cm 7 1/2"
Sharp



Norfolk & Norwich Retractor
08-524 22.0 cm 8 3/4"
Blunt



Adson Retractor
08-525
19.1 cm 7-1/2" sharp,
4 X 4 teeth straight Arms



Adson Retractor
08-526 26.5 cm 10 1/4"
Sharp



Beckman-Weitlaner Retractor
08-527 16.5cm 6 1/2"
3 x 4 Sharp prongs



Beckman-Weitlaner Retractor
08-528 16.5cm 6 1/2"
3 x 4 Blunt prongs



Cone Retractor

08-529 31 cm, 12"
3x4 Blunt prongs,
13 mm x 14.5 mm
08-530 16.5 cm, 6 1/2"
17 mm x 19 mm
3x4 Sharp prongs,
08-531 16.5 cm, 6 1/2"
17 mm x 19 mm
3x4 Blunt prongs,
08-532 20 cm, 8"
21 mm x 210 mm
3x4 Blunt prongs,

RETRACTORS



Wundhaken	Ecarteurs
Separadores	Divaricatori



Beckman Retractor
08-533 31.7cm 12 1/2"
Spread 3.5cm 5 3/8" 4 x 4
Sharp prongs, Arm length
8.5cm 3 3/8"



Beckman Retractor
08-534 31.1cm 12 1/4"
4 x 4 Sharp prongs 27mm,
Side Blade 1.7cm 5/8" deep,



Beckman Retractor
08-535 31.1cm 12 1/4"
4 x 4 Blunt prongs 27mm,
Side Blade 1.7cm 5/8" deep,



Joll (Friedman)
08-536 15.5 cm, 6 1/2"
Thyroid- and Vaginal Retractor



Jackson Burrows Retractor
08-537 18.5cm
Blunt



Gelpi Vaginal Retractor
08-538 13.5cm



Gelpi Skin & Vaginal Retractor
08-539 17.5cm



De Bakey Femur Retractor
08-540 14.5cm


RETRACTORS

Wundhaken
Separadores


Ecarteurs
Divaricatori




**Gosset Retractor
08-541**
Center Blades
40 x 30 mm




**Gosset Retractor
08-542**
Lateral Blades
58 mm deep
Center blade
59 x 57 mm




**Gosset Retractor
08-543**
Blades:
55mm deep




**Gosset Retractor
08-544**
Blades:38 mm deep




**Hertzier Retractor
08-545**
Blades 20 x 25 mm



**Balfour-Baby
Retractor
08-546**
Lateral Blades
27 mm deep
Center blade
21 x 24 mm




**Balfour Retractor
08-547**
Lateral blades
60 mm deep
Center blade
47 x 80 mm



**Balfour Retractor
08-548**
Lateral blades
105 mm deep
Center blade
107 x 59 mm



08-549
Center blade only
54 x 80 mm



**Balfour
Retractor
08-550**
Spreader only

- 08-551**
50 mm
- 08-552**
60 mm
- 08-553**
75 mm
- 08-554**
90 mm
- 08-555**
80 x 80 mm
- 08-556**
80 x 80 mm
- 08-557**
80 x 80 mm
- 08-558**
80 x 80 mm



Finochietto Retractors
08-555 18x22 8cm
08-556 25x40 13cm
08-557 38x45 16cm
08-558 40x45 18cm



Gosset Abdominal Retractors
08-559 12cm
08-560 16cm
08-561 24cm
08-562 30cm

RETRACTORS



Wundhaken | Ecarteurs
Separadores | Divaricatori



O'Sullivan/O'Connor Retractor
08-563 (Set includes frame & 3 blades)
08-564 Frame Only (not pictured)



08-565
Single Retractor Blade
110 x 70 mm



08-566
Set of 2 Retractor Blade
80 x 45 mm



**ABRE BOCA MC IVOR C/3
LAMINAS S/ SUCCAO
ABRE BOCA
08-567**



08-568
Center Blade
45 x 80 mm
70 x 85 mm


Sizes of These Blades to
Same as about Blades.



08-569
Center Blade
100 x 60 mm



Collin Retractor
08-570 Complete
08-571 Spreader only
08-572 Blades 55 x 25 mm
08-573 Blades 80 x 25 mm
08-574 Center Blade
115 x 45 mm



Scoville-Richter Retractor
08-575 Complete with 1 pair blades.
08-576 Frame only.
The Retractor bar and
shanks may be easily
rotates 180 when
blades are in situ.



08-577
62 x 52 mm



08-578
62 x 37 mm



08-579
47 x 51 mm



08-580
30 x 40 mm





08-581
37 x 62 mm

RETRACTORS

Wundhaken
Separadores

Ecarteurs
Divaricatori



Denis-Browne Retractor
08-594 Complete with 4 blades

08-595 Frame 190 x 160 mm
08-596 2 blades only 40 x 30 mm
08-597 2 blades only 40 x 40 mm

08-598
Blades only 20 x 30 mm

08-599
Blades only 20 x 40 mm



RETRACTORS





Wundhaken | Ecarteurs
Separadores | Divaricatori

Abdominal Retractors
well suitable for gynaecologists

Semm Abdominal Retractors
consisting of:

- 08-600** 1 frame 225 x 220 mm
- 08-601** 2 blade 46 x 64 mm
- 08-602** 1 blade 53 x 54 mm
(for bladder)
- 08-603** 1 blade 82 x 97 mm
(for intestine)
- 08-604** 1 blade 98 x 51 mm
(for intestine)
- 08-605** 1 blade, narrow
107 x 36 mm



Kirschner Retractor
08-606
with folding frame &
2 blades of each **F8-584** &
F8-585

Kirschner Retractor
08-607
with vaulted fixed frame &
2 blades of each **08-608** &
08-609

08-610
300 x 380 mm
Frame only, over-sized

08-611
320 x 300 mm
Frame only

Additional valves for corpulent persons



Kirschner Retractor
08-612
64 x 80 mm



Kirschner-Balfour
Retractor
08-613
80 x 82 mm



Kirschner Retractor
08-614 90 x 69 mm
08-615 92 x 95 mm
08-616 78 x 58 mm



Kirschner Retractor
08-617 108 x 42 mm
08-618 91 x 55 mm
08-619 121 x 55 mm

PLASTIC SURGERY INSTRUMENTS





09-001
Grossmann Areola Marker,
40mm, 44mm, 48mm, 52mm



Freeman Areola Marker,
2 Sizes,
09-002 34mm / 36mm
09-003 38mm / 42mm
09-004 42mm / 45mm
09-005 45mm / 50mm



Freeman Areola Marker,
09-006 34mm
09-007 36mm
09-008 38mm
09-009 42mm
09-010 45mm
09-011 50mm



Freeman Areola Marker,
09-012 38 mm
09-013 42 mm
09-014 50 mm
09-015 45 mm
09-016 34 mm
09-017 36 mm



Solz Areola Marker, Blunt,
09-018 36mm Diameter
09-019 40mm Diameter
09-020 42mm Diameter
09-021 45mm Diameter



Solz Atraumatic
Breast Dissector,
09-022 34cm
09-023 37cm



Solz Breast Hook
Dissector,
09-024 36cm



Agris Dingman Breast
Dissector, Set of 2,
Left And Right,
09-025 36cm



Maccollum Dingman
Submammary
Dissector,
09-026 30cm



Maccollum Dingman
Breast Dissector,
09-027 31.5cm
09-028 41cm



Reynolds Transauxillary
Breast
Dissector, Fenestrated,
09-029 33cm
09-030 40cm



Breast Dissector,
Angulated Blade,
09-031 33cm
09-032 42cm



Breast Dissector,
Spatulated Blade,
09-033 33cm
09-034 42cm



Breast Dissector,
Oval Spatulated
Blade,
09-035 33cm
09-036 42cm



Spatula Breast
Retractor
Dissector,
Serrated,
09-037 1/2",
Semi Malleable



Trocarn Breast Tube,
09-038 356mm
Total Length,
Large, 17.3mm Id,
19mm Od
09-038 343mm
Small, 14mm Id, 16mm Od



Penn Breast
Dissector,
Straight,
09-039
25mm Width
324mm Length



Maccollum
Dingman
Breast
Dissector
09-040
38mm Width,
432mm Length



Mcinnis
Polycurette
Breast
Instrument
09-041
30mm Width
09-042
318mm Length



Iconoclast,
09-043 Large
292mm Length
146mm Blade
09-044 Small
222mm Length,
83mm) Blade



Silverstein Breast
Dissector,
09-045
356mm Length



Breast Elevator For
Mammary Plasty
09-046



Temporal T
Dissector
09-047
23.5cm
12mm Tip,
Straight



Flap Dissector
09-048 23.5cm
5mm Tip, Straight



Midface Facial
Dissector
09-049 24cm
Multiple Curves,
12mm Wide Tip



Transoral
Dissector
09-050
23.5cm
Curved Tip



Frontotemporal Dissector
09-051 22.5cm
5mm Tip, Slightly Curved
09-052 23cm
7mm Tip, Slightly Curved
09-053 23cm
10mm Tip, Slightly Curved



Scalp Elevator 24cm
09-054 12.5mm Tip, Slightly Curved
09-055 15mm Tip, Slightly Curved
09-056 7mm Tip, Straight
09-057 12.5mm Tip, S Shaped
09-058 12.5mm Tip, Slightly Curved
09-059 7mm Tip, Curved Tip
09-060 7mm Tip, Strongly Curved
09-061 7mm Tip, Slightly Curved



Orbital Rim Dissector
09-062 22.5cm
Slightly Curved



Frontoglabellar Dissector
09-063 26.5cm
7mm Tip, S Shaped



Arcus Marginalis Dissector
09-064 24cm
10mm Tip, S Shaped



Periosteal Elevator
09-065 25cm
S Shaped



Periosteal Spreader
09-066 24cm
8mm Tip, S Shaped



Periosteal Spreader
09-067 24cm
10mm Wide Blade, S Shaped




Zygomatic Arch Elevator
09-068 20.5cm
Straight Blade, 5mm



Scalp Elevator With Suction
09-069 24cm
7mm Wide Tip, Slightly Curved



Scalp Elevator With Suction
09-070 26cm
8.5mm Wide Tip, Curved Tip



Langenbeck Bone Elevator
09-071 20cm, 7mm, Blunt
09-072 20cm, 10mm, Blunt



Suction Elevator With Guide
09-073 21cm, 4mm



Cottle Periosteal Elevator,
09-074 20cm,
2/4mm, Double Ended, Fine
09-075 23cm,
3/4mm, Double Ended



FIMA INSTRUMENTS

DISSECTORS AND ELEVATORS

Masing Septum Elevator,
09-076 22cm
5mm, Double
Ended



Nagata Auricular
Septum Elevator,
09-077 16cm, Blunt
09-078 17cm, Blunt



Halle Septum
Elevator,
Narrow, 3 mm
Width
09-079 178mm
Length



Freer Septum Elevator,
Double Ended, Sharp /
Blunt, Wide, 6.5mm Width,
09-080 197mm Length



Hajek Ballenger Septum
Elevator, Double Ended,
8.5mm Width
09-081 191mm Length



Freer Septum Elevator,
Curved, Double Ended,
5mm Width
09-082 191mm
Length



Killian Septum Elevator,
Double Ended, 4 mm
With 3.4 Mm Taper, And
4.3 Mm With 3mm Taper,
09-083 184mm Length



Killian Septum Elevator,
Double Ended
09-084 229mm Length,
6mm With 4.5 Mm Taper
And 4.3mm With 3mm Taper



Howarth Septum Elevator,
Double Ended, 5mm Width
09-085 216mm Length



Gruber Combination
Knife Handle
/ Perichondrium Elevator #7
09-086 165mm) Length,
10/15c Blade Size



Obwegeser Rib Raspatory,
09-087 30cm, Left
09-088 30cm, Right



Joseph Masing
Septum Raspatory,
09-089 22cm, 4mm,
Sharp



Cottle Periosteal
Elevator/raspatory,
09-090 18cm,
8mm, Sharp

BREAST RETRACTORS



Marx Breast Retractor,
09-091 27cm, 125 X 30mm
09-092 27cm, 140 X 30mm

Ferriera Style Retractor, With Fiber Optic Light Guide

With Smoke Tube

09-093	70 X 16 mm
09-094	80 X 16 mm
09-095	90 X 27 mm
09-096	110 X 25 mm
09-097	120 X 25 mm
09-098	130 X 16 mm
09-099	130 X 25 mm
09-100	150 X 25 mm
09-101	150 X 30 mm
09-102	150 X 36 mm
09-103	180 X 16 mm
09-104	180 X 25 mm
09-105	190 X 25 mm
09-106	190 X 36 mm
09-107	220 X 27 mm

Without Smoke Tube

09-108	70 X 16 mm
09-109	80 X 16 mm
09-110	90 X 27 mm
09-111	110 X 25 mm
09-112	120 X 25 mm
09-113	130 X 16 mm
09-114	130 X 25 mm
09-115	150 X 25 mm
09-116	150 X 30 mm
09-117	150 X 36 mm
09-118	180 X 16 mm
09-119	180 X 25 mm
09-120	190 X 25 mm
09-121	190 X 36 mm
09-122	220 X 27 mm



Breast Spatula With Handle
09-123 31.5cm, 25mmX18cm
Blade, Serrated Tip, Semi Malleable



Wieder Breast And Facial Retractor /
Tongue Depressor, Serrated,
09-124 14cm, 28mm, With Fiber Optic Light Guide
09-125 14cm, 38mm, With Fiber Optic Light Guide
09-126 15cm, 45mm, With Fiber Optic Light Guide



Wieder Breast And Facial Retractor
/ Tongue Depressor, Serrated,
09-127 14cm, 28mm
09-128 14cm, 38mm
09-129 15cm, 45mm



Breast Retractor, Double Ended,
09-130 31cm,
Blades 22mm X 55mm and 38mm x 70mm
09-131 24cm,
Blades 22mm x 55mm and 38mm x 70mm



Biggs Mammaplasty Retractors,
09-132 30mm
09-133 50mm
With Fiber Optic Light Guide
09-134 30mm,
09-135 50mm,
With Fiber Optic Light Guide And Suction



Biggs Mammaplasty Retractor,
09-136
50mm Width, Standard,
20cm Length



09-137 Biggs Mammaplasty Retractor,
With Fiber Optic Light Guide 5.0cm
Blade Width, 20cm Length



09-138 Biggs Mammaplasty Retractor,
3cm Blade Width, Narrow,
171mm Length



09-139 Biggs Mammaplasty Retractor,
With Fiber Optic Light Guide,
Narrow, 3.0cm Blade Width,
15cm Length



09-140 Ferreira Breast Augmentation Retractor,
25x171cm Blade, 130mm Handle Length



09-141 Ferreira Breast Augmentation Retractor, With Fiber Optic Light Guide, 25 X 171 Mm Blade Width, 130 Mm Handle Length Width



Tebbetts Retractor With Fiber Optic Light Guide, 18.5cm

With Irrigation Tube

09-142 70 X 16 mm
09-143 80 X 16 mm
09-144 90 X 27 mm
09-145 110 X 25 mm
09-146 120 X 25 mm
09-147 130 X 16 mm
09-148 130 X 25 mm
09-149 150 X 25 mm
09-150 150 X 30 mm
09-151 150 X 36 mm
09-152 180 X 16 mm
09-153 180 X 25 mm
09-154 190 X 25 mm
09-155 190 X 36 mm
09-156 220 X 27 mm

Without Irrigation Tube

09-157 70 x 16 mm
09-158 80 x 16 mm
09-159 90 x 27 mm
09-160 110 x 25 mm
09-161 120 x 25 mm
09-162 130 x 16 mm
09-163 130 x 25 mm
09-164 150 x 25 mm
09-165 150 x 30 mm
09-166 150 x 36 mm
09-167 180 x 16 mm
09-168 180 x 25 mm
09-169 190 x 25 mm
09-170 190 x 36 mm
09-171 220 x 27 mm



Tebbetts Fiber Optic retractor, With Double Handle and Teeth at Tip

With Smoke Tube

09-172 70 x 16 mm
09-173 80 x 16 mm
09-174 90 x 27 mm
09-175 110 x 25 mm
09-176 120 x 25 mm
09-177 130 x 16 mm
09-178 130 x 25 mm
09-179 150 x 25 mm
09-180 150 x 30 mm
09-181 150 x 36 mm
09-182 180 x 16 mm
09-183 180 x 25 mm
09-184 190 x 25 mm
09-185 190 x 36 mm
09-186 220 x 27 mm

Without Smoke Tube

09-187 70 X 16 mm
09-188 80 X 16 mm
09-189 90 X 27 mm
09-190 110 X 25 mm
09-191 120 X 25 mm
09-192 130 X 16 mm
09-193 130 X 25 mm
09-194 150 X 25 mm
09-195 150 X 30 mm
09-196 150 X 36 mm
09-197 180 X 16 mm
09-198 180 X 25 mm
09-199 190 X 25 mm
09-200 190 X 36 mm
09-201 220 X 27 mm


Tebbetts Breast Augmentation Fiber Optic Retractor, With Double Handle And Smooth End

With Smoke Tube

09-202 70x16mm
09-203 80x16mm
09-204 90x27mm
09-205 110x25mm
09-206 120x25mm
09-207 130x16mm
09-208 130x25mm
09-209 150x25mm
09-210 150x30mm
09-211 150x36mm
09-212 180x16mm
09-213 180x25mm
09-214 190x25mm
09-215 190x36mm
09-216 220x27mm

Without Smoke Tube

09-217 70x16mm
09-218 80x16mm
09-219 90x27mm
09-220 110x25mm
09-221 120x25mm
09-222 130x16mm
09-223 130x25mm
09-224 150x25mm
09-225 150x30mm
09-226 150x36mm
09-227 180x16mm
09-228 180x25mm
09-229 190x25mm
09-230 190x36mm
09-231 220x27mm



Tebbetts Breast Augmentation Fiber Optic Retractor, Insulated Blade With Smooth End, 18.5 Cm

With Smoke Tube

09-232 70x16mm
09-233 80x16mm
09-234 90x27mm
09-235 110x25mm
09-236 120x25mm
09-237 130x16mm
09-238 130x25mm
09-239 150x25mm
09-240 150x30mm
09-241 150x36mm
09-242 180x16mm
09-243 180x25mm
09-244 190x25mm
09-245 190x36mm
09-246 220x27mm

Without Smoke Tube

09-247 70x16mm
09-248 80x16mm
09-249 90x27mm
09-250 110x25mm
09-251 120x25mm
09-252 130x16mm
09-253 130x25mm
09-254 150x25mm
09-255 150x30mm
09-256 150x36mm
09-257 180x16mm
09-258 180x25mm
09-259 190x25mm
09-260 190x36mm
09-261 220x27mm

Tebbetts Fiber Optic retractor, Insulated Blade With Teeth 18.5cm

With Smoke Tube

09-262 70x16mm
09-263 80x16mm
09-264 90x27mm
09-265 110x25mm
09-266 120x25mm
09-267 130x16mm
09-268 130x25mm
09-269 150x25mm
09-270 150x30mm
09-271 150x36mm
09-272 180x16mm
09-273 180x25mm
09-274 190x25mm
09-275 190x36mm
09-276 220x27mm

Without ISmoke Tube

09-277 70x16mm
09-278 80x16mm
09-279 90x27mm
09-280 110x25mm
09-281 120x25mm
09-282 130x16mm
09-283 130x25mm
09-284 150x25mm
09-285 150x30mm
09-286 150x36mm
09-287 180x16mm
09-288 180x25mm
09-289 190x25mm
09-290 190x36mm
09-291 220x27mm



Sculpo Breast Retractors, 180 X 25mm
Blade Size, With
09-292 5mm Endoscopic Scope Sheath
09-293 10mm Endoscopic Scope Sheath



09-294 Ring Breast Retractor
With Fiber Optic Illumination
And Irrigation Tube



Alach Forceps,
09-295 16cm,
5mm Tip

Dingman Cartilage
Abraders,
Set of Two
09-296

Rubin Septal Morselizer
Forceps,
09-297 21cm

Rubin Septal Morselizer
Forceps,
T.c. (tungsten Carbide)
09-298 20cm

Mayo Scissors,
Supercut Plus T.c.
(tungsten Carbide)
09-299 17cm, Curved



Cottle Cartilage Crusher, With Clamp, For
Taking Out The Prepared Material,
09-300 65x30x30mm



Cartilage Cutting Block, Teflon,
09-302 85x55x10mm




Gubisch Cartilage Crusher
09-301 65x30mm



Bone Mill, With Assembling Wrench
09-303




FIMA INSTRUMENTS

DISSECTORS AND ELEVATORS


Aufricht Nasal Retractor,
09-304 18cm
45x7-10mm




Aufricht Nasal Retractor
Double Ended, Blunt,
09-305 13.5cm
282-38x7mm



Aufricht Walter Nasal Retractor
Double Ended, Blunt,
09-306 13.5cm
283-38x7mm



Aufricht Walter Nasal Retractor
Double Ended, Sharp,
09-307 13.5cm
284-42x7mm



Cottle Alar Retractor,
09-308 15cm,
22x13mm



Walter Alar And Nasal Retractor
With Fiber Optic Illumination,
09-309 18.5cm



Gubisch Aufricht
Nasal Retractor,
09-310 17cm,
50x6-12mm

RHINOPLASTY RETRACTORS



Aufrecht Nasal Retractor
09-311 16cm
7mm Tip, 45mm
Blade, With Fiber Optic
Light Guide
09-312 18cm 10mm



Converse Blade Nasal Retractor,
09-313 20mm X 10mm Wide,
With Fiber Optic Light Guide
09-314 30mm X 12mm Wide,
09-315 40mm X 12mm Wide,
09-316 50mm X 12mm Wide,



Fomon Nasal Retractor,
09-317 1 Ball End Prong
16cm
09-318 2 Ball 16cm




Killian Nasal Speculum
With Fiber
Optics Illumination
09-319




Joseph Nasal Saw,
09-320 17cm,
Straight



Joseph Nasal Saw,
09-321 17cm, Left
09-322 17cm, Right



Joseph Nasal Saw,
09-323 Angled To Left
09-324 Angled To Right



Tuppers Hand Retractor Set
09-325 Medium
09-326 For Infants
09-327 For Adults



Obwegeser Alveolar
(wisdom Tooth) Retractor,
09-328 22cm, 13/23 X 70mm




Obwegeser Nasal Spine Retractor,
09-329 22cm, 10/23 X 45mm



Obwegeser Iliac
Crest Retractor
09-330 22cm
30/60 X 60mm



Obwegeser Iliac
Crest Retractor
09-331 22cm
30/60 X 80mm



Obwegeser Mandibular Channel
Retractor,
09-332 17.5cm, 8mm
09-333 17.5cm, 10mm



Obwegeser Chin
Retractor,
09-334 16cm



Saldanha Solz
Abdominal
Retractor,
09-335 23.5cm



Killian Nasal
Speculum With
Fiber Optics
Illumination
09-336



Deaver Retractor
With Fiber Optics
09-337
Illumination, 16mm
09-338
Illumination, 25mm



Ferreira Face Lift Retractor
With Fiber Optic
Illumination 22 X 108mm
09-339



Freeman Flap Retractor
With Fiber Optics Illumination,
09-340 20cm, 13mm
09-341 20cm, 25mm



09-342
Mandibular
Channel
Retractor
8mm,
17.5cm

09-343
Bremerich
Mandibular
Channel
Retractor
10mm, 18cm

09-344
Ramus Retractor 20x80mm
16cm

09-345
Ramus Retractor 25x80mm
16cm

09-346
Ramus Retractor 30x80mm
16cm

09-347
Ramus Retractor 35x80mm
16cm



09-348
Medial Ramus
Retractor
left 19.5cm

09-349
Medial Ramus
Retractor
right 19.5cm



09-350
Bremerich Channel
Retractor 10mm,
22.5cm w/suction



09-351
Extra Oral
Retractor
24cm



09-352
Intra Oral
Retractor
24cm



09-353
Intra Oral
Retractor
L and M
w/o F/O
fitting.
25cm



09-354
Intra Oral
Retractor
right w/o
F/O fitting.
21cm



09-355
Intra Oral
Retractor
left w/o
F/O fitting.
21cm



09-356
Intra Oral Retractor
w/o F/O
Fitting 23cm



09-357
Intra Oral Retractor
With
F/O fitting 23cm



09-358
Intra Oral Retractor
Left
with F/O fitting
21cm



09-359
Intra Oral Retractor
right
w/F/O fitting. 21cm



**Levasseur Merrill Fiber
Optic Retractor for
Vertical Osteotomies**
09-360 Length 25cm
Working Tip 27mm



**Shea with Fiber Optic
(Multipurpose Retractor)**
09-361 Length 25cm
Working Tip 14.5mm



**Brunner Retractors
with Fiber Optic**
09-362 Length 24cm
Blades 145 x 25mm



09-363
Cheek and Lip
Retractor D/E
21x30mm,
26x35mm, 18cm

RETRACTORS



09-364
Schwenzer Cheek
Retractor 65x50mm,
17cm



09-365
Reuther Lliac
Crest
Retractor 24cm



09-366
Tessier Rib Graft
Retractor adult
120x28mm, 20.5cm



09-367
Lester
Retractor



09-368
Richardson Retractor
28x20mm



09-369
Richardson Retractor
38x28mm



09-370
Richardson Retractor
38x37mm



09-371
Richardson Retractor
52x22mm



**Aufrich Nasal Retractors
with Fiber Optic**
09-372 16cm 7mm
09-373 18cm 10mm



**Harrington Retractors
with Fiber Optic**
09-374 40 x 127mm
09-375 60 x 127mm

for Right cheem



for Left cheem



**Bauer Retractor with
Fiber Optic**
09-376 26.5cm / 16mm
09-377 26.5cm / 16mm



**Petri Pterygoid Retractor
with Fiber Optic**
09-378 23cm, 25.5mm



**Sigmoid Nptch Retractor
with Fiber Optic**
09-379 22cm, 16mm



**Maliniac Nasal Retractor
with Fiber Optic**
09-380 9mm x 45mm



09-381
Yancoskie Retractor
18cm, 6 sharp prongs



09-382
Coronal Brow Lift Retractor
18cm, 7 Sharp prongs



Freeman Retractor

09-383 18cm,
4 Sharp prongs
22mm, wide
09-384 18cm,
4 Blunt prongs
22mm, wide



**Hook Retractor in-line
prongs / four prongs 18cm**
09-385 straight
09-386 narrow head
09-387 angled shaft
09-388 angled shaft 4 prong
with shorter inner prongs 18cm



**Anderson Bear
Claw Retractor**
09-389 12cm



**Miniature model /
with thumb ring /
in-line prongs /
four prongs**
09-391 11.5cm


**Miniature model /
with thumb ring /
in-line prongs /
four prongs**
09-390 6cm



**Freeman Hook
Retractor
four Prongs /
with shorter
inner prongs /
seven prongs**
09-392 18 cm



**Freeman Hook
Retractor
Seven Prongs /
with shorter
inner prongs /
seven prongs**
09-393 18 cm
concave



**Freeman Hook
Retractor
Seven Prongs /
with shorter
inner prongs /
seven prongs**
09-394 18 cm
convex



Rhinoplasty set

09-395

01. Scalpel Handle, #7k,
02. Surgical Operating Scissors, 14.5cm, Curved - Sharp/Blunt
03. Ragnell Kilner Scissors, 15cm, Curved - Flattened At Distal End
04. Fox T.C. Gum And Tissue Scissors, 13cm, Curved, 1 Serrated Blade
05. Adson Brown Delicate Tissue Forceps, 12cm, 7x7 Teeth, Straight
06. Adson Brown Delicate Tissue Forceps, 12cm, 9x9 Teeth, Straight
07. Semken Delicate Tissue Forceps, 12.5cm, 2x3 Teeth, Straight
08. Debakey At Vascular Forceps, 15cm - 1.5mm Tips
09. Bayonet Grasping Forceps, 16cm
10. Bayonet Grasping Forceps, 16cm, 1x1 Teeth
11. Troeltsch Wilde Angular Forceps, 13cm, Adult Size - Serrated
12. Grooved Director And Tongue Tie, 11.5cm
13. Frazier Fergusson Suction Tube, 18cm, Angled, 10 French e
14. Senning T.C. Needle Holder, 17cm, Serrated
15. Halsted Micro Mosquito Forceps, 12.5cm, Curved
16. Backhaus Towel Clamp, 8.5cm
17. Gross Dressing Forceps, 20cm, Straight - Without Catch, Box Lock
18. Walter Cope Retractor, 17.5cm
19. Joseph Skin Hook, 16cm, Single Prong, Sharp
20. Guthrie Double Retractor Hook, 13cm, 2mm Wide, Sharp
21. Lempert Raspatory, 17cm, 4mm, Delicate, Sharp
22. Supercut Walter Dorsal Scissors, 14cm, Angular
23. Supercut Walter Cottle Dorsal Scissors, 16cm, Curved
24. Medicon Mallet, 18.5cm, Plastic Facings, 200 Grams, 25mm Diameter
25. Walter Williger Bone Curette, 13.5cm, 2/2mm, Double Ended, Sharp, Angled
26. Defourmental Bone Rongeur, 19cm, Straight, Slender Beaks - Jaws 3mm Wide
27. Weil Blakesley Wilde Ethmoid Bone Forceps #3, 120mm Shaft, 4.2mm
28. Weil Blakesley Wilde Ethmoid Bone Forceps #2, 45 Degree Angle, 120mm Shaft, 3.6mm
29. Myles Nasal Through Cutting Forceps #4, 120mm Shaft
30. Cottle Nasal Septum Speculum, 14cm, 50mm
31. Cottle Elevator, 16cm, 2.5mm
32. Masing Septum Elevator, 22cm, 5/5mm, Double Ended
33. Aufricht Walter Nasal Retractor, 13.5cm, Double Ended, Sharp, 42x7mm
34. Aufricht Walter Nasal Retractor, 13.5cm, Double Ended, Blunt, 38x7mm
35. Walter Chisel/Osteotome, 19cm, 3mm, Modified
36. Walter Chisel/Osteotome, 19cm, 5mm, Modified
37. Walter Osteotome, 19cm, 2mm, Straight
38. Walter Osteotome, 19cm, 12mm, Straight
39. Walter Chisel/Osteotome, 19cm, 12mm, Modified
40. Walter Osteotome, 19cm, 6mm, Curved
41. Walter Cantanares Lifting Scissors, 16cm, Curved
42. Fomon Rhinoplasty Scissors, 12.5cm
43. Cottle Columella Forceps, 11cm
44. Fomon Bone And Nasal Rasp, 20cm
45. Asch Septum Straightening Forceps, 22cm, Angled
46. Brunings Tongue Depressor, 19cm, 14mm



Rhinoplasty set

09-396

Steel Ruler 6"	1
Metz Lahey Crv Del 5 3/4" Tc	1
Metz Sciss Cvd Del 7 Tc	1
Crile Wood Nh Serr 7" Tc	1
Mallet W/Repl Plst Discs	1
Maltz Rasp 7" Bkd Cut	1
Joseph Rasp Fine Crss Serr 6 1/4"	1
Aufricht Glabellar Rasps 8" Bkd	1
Cottle Dorsal Sciss Ang 6 1/2"	1
Black Cotted Cottle Septal Spec 2"	2
Cottle Elev Cm-Grad D/E 9	1
Joseph Perio Elev 7" X6.5Mm	1
Joseph Perio Elev 6 1/4" X4.5Mm	1
Cottle Nasal Knife Sh 5 1/2" 4Mm	1
Cottle Nasal Osteo 6 1/2" Cvd	1
Cinelli Osteo Grd 2Sd Str 10Mm	1
Cot Knife Grd/Retr 8" 5Mm Hook	1
Neiv Knife Gd/Retr 6 1/2" D-End	1
Cottle Ret 5 1/2" 4Pmg 2Reces	1
Fomon Retr Ball End Dbl 6 1/4"	1
Converse Alar Retr D/E 4 1/2"	1
Iris Sciss Str S/S 3 1/2"	1
Freer Elevator De Sh/Bl 7 1/4"	1
Joseph Hook 1/Pr Sh 6 1/4"	1
Joseph Hook 2/Pr Sh 6 1/4" 5Mm	2
Gillies Skin Hook 7" 3Mm	1
Adson Dressing Fcp Serr 4 3/4"	1
Strabismus Sciss Str 4 1/2"	1
Strabismus Sciss Crv 4 1/2"	1
Desmarres Lid Retr 12Mm	1
Castroviejo Caliper 3 1/4"	1
Cottle Columella Fcp 4 1/2"	1
Hartman Ear Fcp Light Serr Jaw	1
Jansen Nasal Fcp Bayonet 6 1/2"	1
Sheehan Osteo Str 6 1/4" 4Mm	1
Sheehan Osteo Str 6 1/4" 8Mm	1
Sheehan Osteo Str 6 1/4" 10Mm	1
Sheehan Osteo Str 6 1/4" 12Mm	1
Killian Med W/Set Screw 3"	1
Foman Nas Low Lat Scis Cv 5 3/4"	1
Ash Sept Straightening Fcp 9"	1
Ballenger Swivel Knife 4Mm	1
Joseph Bone Saw Str 7 1/2"	1
Schmidt Half Crv Tons Fcp 7 1/2"	2
Maltz Rasp Se Fine Tc	1
Maltz Rasp Se Medium Tc	1



Rhinoplasty set

09-397

01. Scalpel Blades #15 - Sterile,
02. Scalpel Handle, #7
03. Surgical Operating Scissors, 13cm, Straight - Sharp/Blunt
04. Mayo Stille T.C. Operating Scissors, 17cm, Curved
05. Reynolds Jameson Scissors, 14cm, Curved - With Fine, Rounded Points
06. Adson T.C. Tissue Forceps, 12cm
07. Micro Adson Forceps, 12cm, 1x2 Teeth - 0.8mm Tips
08. Adson Brown Delicate Tissue Forceps, 12cm, 7x7 Teeth, Straight - Atraumatic
09. Dressing Forceps, 20cm - Modified Swedish Pattern
10. Nelson Tissue Forceps, 15cm
11. Gruenwald Jansen Nasal Dressing Bayonet Forceps, 16cm, Serrated
12. Gubisch T.C. Needle Holder, 12.5cm, Extra Fine
13. Crile Wood T.C. Needle Holder, 15cm
14. Halsted Micro Mosquito Forceps, 12.5cm, Curved
15. Leriche Artery Forceps, 15cm, 1x2 Teeth, Straight
16. Pean Hemostatic Forceps, 14cm, Straight
17. Backhaus Towel Clamp, 13cm
18. Dressing Forceps, 17.5cm, Curved, Delicate
19. Gross Dressing Forceps, 20cm, Curved - With Catch, Box Lock
20. Barsky Skin Hook, 15cm, Single Prong, Sharp
21. Joseph Delicate Retractor, 17cm, 2 Prong, 5mm, Sharp
22. Guthrie Double Retractor Hook, 13cm, Sharp
23. Supercut Plus T.C. Mayo Scissors, 14cm, Straight
24. Supercut Plus T.C. Mayo Scissors, 17cm, Curved
25. Supercut Plus Metzenbaum T.C. Scissors, 18cm, Curved
26. Supercut Plus Kelly T.C. Scissors, 16cm, Curved
27. Langenbeck Bone Elevator, 20cm, 7mm, Blunt
28. Freer Raspatory, 19cm, 5/5mm, Double Ended, Sharp / Sharp
29. Williger Raspatory, 16cm, 5mm, Sharp
30. McKenty Raspatory, 15cm, 4mm, Sharp, Light Curve
31. Supercut Heymann Nasal Scissors, 18cm
32. Supercut Plus Fomon T.C. Dorsal Scissors, 13cm, Angular
33. Mallet, 26cm, Novotex Plastic, 240 Grams, 45mm Diameter
34. Round Bowl For Solutions, 40x20mm, 20cc
35. Round Bowl For Solutions/Marking Color, 60x30mm, 60cc
36. Round Bowl For Solutions, 80x40mm, 160cc
37. Round Bowl For Solutions, 110x55mm, 400cc
38. Castroviejo Caliper, 7cm, 20mm Measuring Range
39. Desmarres Lid Retractor #2, 14cm, Medium, 10x14mm
40. Killian Asept Nasal Speculum, #4, 90mm (Modified)
41. Weil Blakesley Wilde Ethmoid Bone Forceps #3, 120mm Shaft, 4.2mm
42. Freer Septum Knife (Mucosa Scalpel), 16cm
43. Cottle Nasal Speculum, Delicate, 15cm, 55mm
44. Gubisch Suction Raspatory With Guide And Bypass Hole, 20cm, 4mm
45. Gubisch Suction Raspatory With Guide, 17cm, 4mm
46. Cartilage Cutting Block, Teflon, 85x55x10mm
47. Gubisch Cartilage Crusher, 65x30mm
48. Gubisch Kilner Double Nasal Hook (Alar Retractor), 8cm, Sharp, 5x10mm
49. Gubisch Aufricht Nasal Retractor, 17cm, 50x6-12mm
50. Cottle Gubisch Chisel, 18cm, 5mm, Graduated, Hollow Cut
51. Cottle Gubisch Chisel, 18cm, 12mm, Graduated, Hollow Cut
52. Gubisch Osteotome, 18.5cm, 3mm, Aluminium Handle
53. Fanous Osteotome, 19.5cm, 4mm, Straight, With Guide Thorn - 1)
54. Caplan Nasal Septum Scissors, 20cm, Serrated
55. Gubisch Fomon Nasal Rasp, 20cm
56. Rubin Septal Morselizer Forceps, 20cm, With Protective Guard
57. Buchwald Tongue Spatula, 18cm, 18x23mm




Rhinoplasty set

09-398

- Adson Forceps TC, Serrated, 12cm 2
- Aston SuperCut Scissors, 5 1/2" (14cm), Angled, Serrated 1
- Aufrecht Nasal Retractor with Fiber Optic Light Guide, 7" (18cm), 1.75in Tapered 1
- Backhaus Towel Forceps, 3 1/2" (9cm) 2
- Beasley-Babcock Forceps, 4 3/4" (12cm), 3.5mm tips 1
- Bishop-Harmon Forceps, 3 1/2" (9cm), 1x2 teeth, 0.6mm tips 2
- Boies Nasal Fracture Elevator, 7 3/4" (19.5cm) 1
- Brown-Adson Forceps, 4 3/4" (12cm), Brown teeth 1
- Caplan Nasal Shears Scissors, 7 3/4" (19.5 cm), Angled, Serrated 1
- Castroviejo Caliper, 3 1/2" (9cm) 1
- Cottle Cartilage Crusher Morselizer 1
- Cottle Double-Ended Septum Elevator, 8 1/2" (21cm), 4mm/3mm, Sharp/Blunt 1
- Cottle Modified Nasal Speculum, 5 1/2" (14cm), 50mm 1
- Cottle Retractor, 5 1/2" (14cm), 2 Sharp Prongs 1
- Cushing Bayonet Insulated Bipolar Forceps, 7 1/4" (18.5cm), 1mm tips, American Pattern 1
- Foerster Sponge Holding Forceps, 7" (18cm) 1
- Fomon TC Rasp, 8 1/2" (21.5cm), #3 & #4, Double 1
- Fomon TC Rasp, 8 1/2" (21.5cm), #5 & #6, Double 1
- Fomon TC Rasp, 8 1/2" (21.5cm), #7 & #8, Double 1
- Frazier Suction Tubes 10 French 2
- Freer Chisel, 6 1/2" (17cm), 4mm 1
- Gerzog Stainless Steel Mallet (Lead filled) 1
- Goldman Septum Elevator, 7 1/4" (18.4cm), 6mm 1
- Gunter Lower Lateral Longitudinal Jaw Morselizer, 7" (18cm) 1
- Guthrie Double Hook, 4 3/4" (12cm), 2.5mm 2
- Halsey Needle Holder, 5" (13cm), Delicate, Smooth Jaw 2
- Halstead Mosquito Forceps, Curved, 5" (12.5cm) 1
- Iris SuperCut Scissors, 4 1/2" (11.5 cm), Curved, Serrated 1
- Jansen Middleton Forceps, 7 3/4" (19.5cm), Large Fenestrated jaw 1
- Joseph Nasal Saw Left Bayonet, 7 1/4" (18.5 cm) 1
- Joseph Nasal Saw Right Bayonet, 7 1/4" (18.5 cm) 1
- Joseph Periosteal Elevator, 6 1/2" (16.5cm) 1
- Joseph Skin Hook, 6 1/4" (16cm), 10mm Double Hook 1
- Joseph Skin Hooks, 6 1/4" (16cm), 5mm Double Hook 1
- Joseph Skin Hooks, 6 1/4" (16cm), Single Hook 2
- Knife Handle #3, 5" (12.5 cm) for B-P (10-15) blades 4
- Lempert Ronqueur, 6 1/4" (16cm), 2.5mm tips 1
- Lucae Nasal Dressing Forceps, 5 1/2" (14cm), Serrated jaw 1
- Mitchells Trimmer Double-Ended Elevator, 6 5/8" (17cm) 1
- Operating Sharp/Blunt TC Scissors, 5 1/2" (14 cm), Straight 1
- Pean Hemostatic Forceps, Straight, 6 1/4" (16cm) 2
- Pennington Double-Ended Elevator, 8 1/2" (22cm), 4mm/4.5mm 1
- Ribbon Suture Scissors, 3 3/4" (9.5 cm), Curved 1
- Senn Double-Ended Blunt Retractor, 6 1/4" (16cm), blade 1/4in x 3 sharp prongs 2
- Sheen Cartilage Grid 1
- Sheen Fragment Forceps, 7" (18cm) 1
- Silver Osteotome, 7" (18cm) Left 1
- Silver Osteotome, 7" (18cm) Right 1
- Silver Osteotome, 7" (18cm) Straight 1
- Spencer Stitch Removal Scissors, 3 1/2" (9 cm), Straight 1
- Strabismus SuperCut Scissors, 4 1/2" (11.5 cm), Curved, Serrated 1
- SuperCut Scissors, 7 1/2" (19cm), Angled, Serrated 1
- Walsham Septal Forceps, 9" (23cm), 7mm x 33mm Jaws 1
- Webster Glabella TC Rasp, 8 1/2" (21.5cm), Coarse, Push/Pull 1
- Webster TC Needle Holder, 4 1/2" (11.5cm), Delicate, Smooth Jaw 1
- Woodson Double-Ended Elevator, 7" (18cm), straight/angled 1

BONE SURGERY INSTRUMENTS



AMPUTATING & RESECTION SAWS

amputations-und resektiosagen
Scies a amputation et resection

sierra para amputacion y reseccion
seghe per amputazione



Satterlee Saw
10-001 29cm
Stainless Steel Ring Handle
Straight tip, Straight shaft



Charriere Bone Saw
10-002 27cm
10-003 30cm
10-004 34cm
Straight tip, Straight shaft



Charriere Saw
10-005 35cm
Spare saw blades
10-006 12mm
10-007 15mm



Charriere Bone Saw Complete
10-008 35.5cm
Straight tip,
Straight shaft



Helferich-Bier Saw
10-009 39cm
Spare saw blades
10-010 3mm
10-011 5mm
10-012 5mm



Charriere Saw
10-013 25 cm
10-014 27 cm
Double saw



Engel Plaster Saw
10-015 15.5 cm



Serra Gigli
10-016 30 cm
10-017 40 cm
10-018 50 cm



Percy Amputation Shield
10-019 30.5cm
Handles Detachable,
Simple Pattern



Putti Bone Rasp
10-020 27cm
Double ended,
One end curved up
one end curved sideways



Putti Rasp
10-021 30cm
Double Ended



Putti Bone File
10-022 28cm
semicircular



Putti Bone File
Bayonet Shape
10-023 24.5cm
X=13 mm
semicircular



Putti Bone File
10-024 22cm
X=20 mm
flat



Joseph Bone Saw / Rasp
10-025 19cm
Single Sided



Linderman Amputation Saw
10-026 19cm
20cm Blade



Wellenschliff Nasal Saw
10-027 19cm
Straight



Langenbeck Metacarpal Saw
10-028 23cm
11.5cm Blade

BONE FILES



Knochenfeilen
Limas para huesos

Limes à os
Lime per ossa



Patella Bone Saw
10-029 Includes 1 blade
10-030 Replacement blade only
Delrin Handle



Foman Rasp
10-031 21.5cm
Double ended, Fine & coarse serrations 4-sided convex and flat blades



Foman Rasp
10-032 21cm



Miller-Colburn Bone File
18 cm



Fig. 1
10-033



Fig. 2
10-034




Fig. 3
10-035



Fig. 4
10-036

cross cut



Kleinert-Kutz Bone Rasp
10-037 22.5cm
Double ended,
Fine & coarse serrations



Miller Alveotomy File
10-038 19cm



Miller Alveotomy File
10-039 19cm



Gallaher Bone Rasp
10-040 18cm
10-041 25cm
Curved



Gallaher Bone Rasp
10-042 18cm
10-043 25cm
Convex Curve



Gallaher Bone Rasp
10-044 18cm
10-045 25cm
Straight



Good Bone File / Rasp
10-046 16.5cm
Straight



Kistler Ear Cartilage Rasp
10-047 left
10-048 right
16.5cm Curved



Micro Rasp
10-049 14cm
Single ended
double sided



Nail Rasp
10-050 17cm
Single sided



Maltz Rasp
10-051 18cm
Single-Ended
Single sided



Improved Maltz Rasp
10-052 18cm
Single-Ended,
up & down rasping



Putti Rasp
10-053 23cm
falt / flat
10-054 30.5cm
Round / flat

RASPS, FILES

Raspeln, Feilen
Rapes, Limes

Escofinas, limas
Raspe, lime



Miller Rasp
10-055 21.5cm
Double Ended straight
single Sided



Fomon Rasp
10-056 21.5cm
Double Ended straight
single Sided
one side concave



Fomon Rasp
10-057 21.5cm
Double Ended straight
single Sided
one side convex



Fomon Rasp
10-058 21.5cm
Double Ended straight
single Sided



Fomon Rasp
10-059 21.5cm
Double Ended straight
single Sided
one side concave



Fomon Rasp
10-060 21.5cm
Double Ended straight
single Sided
one side convex



Chisel Blade Handle
10-061 18.5cm



10-062 8.5cm
5mm
10-063 8.5cm
10mm
10-064 8.5cm
16mm
10-065 8.5cm
25mm



10-066 8.5cm
5mm
10-067 8.5cm
10mm
10-068 8.5cm
16mm
10-069 8.5cm
25mm



Chisel Blade key
10-070 9cm



Osteotome
13.5cm
10-071 8mm
10-072 10mm
10-073 12mm
10-074 14mm
10-075 16mm



Chisel
13.5cm
10-076 8mm
10-077 10mm
10-078 12mm
10-079 14mm
10-080 16mm



Gouge
13.5cm
10-081 8mm
10-082 10mm
10-083 12mm
10-084 14mm
10-085 16mm



Lucas Chisel
16cm
10-086 3mm
10-087 4mm
10-088 5mm
10-089 6mm



Lucas Gouge
16cm
10-090 3mm
10-091 4mm
10-092 5mm
10-093 6mm



**U.S. Army Pattern
Chisel**
16.5 cm
10-094 6 mm wide
10-095 12 mm wide
10-096 18 mm wide



Parts of Chisel	Parts of Chisel	Parts of Chisel
11.5cm	13.5cm	17.5cm
10-097 2mm	10-103 2mm	10-109 2mm
10-098 3mm	10-104 3mm	10-110 3mm
10-099 4mm	10-105 4mm	10-111 4mm
10-100 6mm	10-106 6mm	10-112 6mm
10-101 7mm	10-107 7mm	10-113 7mm
10-102 8mm	10-108 8mm	10-114 8mm

BONE CHISELS & GOUGES



Knochenmeißel

Ciseaux-burins, gouges

Cineles, gubias

Scalpelli, sgorbi



Partsch Gouge	Partsch Gouge	Partsch Gouge
11.5cm	13.5cm	17.5cm
10-115 2mm	10-121 2mm	10-127 2mm
10-116 3mm	10-122 3mm	10-128 3mm
10-117 4mm	10-123 4mm	10-129 4mm
10-118 6mm	10-124 6mm	10-130 6mm
10-119 7mm	10-125 7mm	10-131 7mm
10-120 8mm	10-126 8mm	10-132 8mm



Alexander Chisel
18cm
10-133 4mm
10-134 6mm
10-135 8mm
10-136 10mm
10-137 12mm
10-138 14mm



Alexander Gouge
18cm
10-139 4mm
10-140 6mm
10-141 8mm
10-142 10mm
10-143 12mm
10-144 14mm



Smith Peterson Gouge
20cm Straight
10-145 09 mm
10-146 14 mm
10-147 19 mm
10-148 25 mm
Curved
10-149 09 mm
10-150 14 mm
10-151 19 mm
10-152 25 mm



Straight
10-153 06 mm
10-154 09 mm
10-155 13 mm
10-156 16 mm
10-157 19 mm
10-158 25 mm
10-159 32 mm
Curved
10-160 06 mm
10-161 09 mm
10-162 13 mm
10-163 16 mm
10-164 19 mm
10-165 25 mm
10-166 32 mm

Smith Peterson Osteotome
20 cm



Straight
10-167 10 mm
10-168 15 mm
10-169 20 mm
10-170 25 mm
Curved
10-171 10 mm
10-172 15 mm
10-173 20 mm
10-174 25 mm

Stille Gouge
20 cm



Stille
20 cm




Stille Chisel
20 cm



10 mm
10-185



15 mm
10-186



20 mm
10-187



25 mm
10-188

Chisels
10-175 10 mm
10-176 12 mm
10-177 15 mm
10-178 20 mm
10-179 25 mm
Osteotome
10-180 10 mm
10-181 12 mm
10-182 15 mm
10-183 20 mm
10-184 25 mm



Stille Osteotome
20 cm




10 mm
10-189



15 mm
10-190



20 mm
10-191



25 mm
10-192



Lexel Osteotome
22 cm
10-193 7 mm
10-194 10 mm
10-195 12 mm
10-196 15 mm
10-197 20 mm
10-198 25 mm
10-199 30 mm
Lexel Modif. Osteotome
29 cm
10-200 7 mm
10-201 10 mm
10-202 12 mm
10-203 15 mm
10-204 20 mm
10-205 25 mm
10-206 30 mm



Lexel Modif. Chisel
29 cm
10-207 7 mm
10-208 10 mm
10-209 12 mm
10-210 15 mm
10-211 20 mm
10-212 25 mm
10-213 30 mm
Lexel Chisel
22 cm
10-214 7 mm
10-215 10 mm
10-216 12 mm
10-217 15 mm
10-218 20 mm
10-219 25 mm
10-220 30 mm



Lexel Gouge
25 cm
10-221 7 mm
10-222 10 mm
10-223 15 mm
Straight

BONE CHISELS & GOUGES

Knochenmeißel
Ciseaux-burins, gouges

Cineles, gubias
Scalpelli, sgorbi



Lexer Gouge
25 cm
10-224 7 mm
10-225 10 mm
10-226 15 mm
Curved



Lexer Osteotome Flat
22 cm
10-227 10 mm
10-228 15 mm
10-229 20 mm
10-230 25 mm



Lexer Gouge Tuffnel Handle
22 cm
10-231 10 mm
10-232 15 mm
10-233 20 mm
10-234 25 mm



Mini Lexer Osteotome
10-235
18 cm



Mini Lexer Chisel
10-236
18 cm



Capener Osteotome With Tuffnel Handle
25cm
10-237 10mm
10-238 12mm
10-239 15mm
10-240 20mm
10-241 25mm
10-242 30mm
Straight



Capener Osteotome With Tuffnel Handle
25cm
10-243 10mm
10-244 12mm
10-245 15mm
10-246 20mm
10-247 25mm
10-248 30mm
Curved



Capener Chisels With Tuffnel Handle
25cm
10-249 6mm
10-250 8mm
10-251 10mm
10-252 12mm
10-253 15mm
10-254 20mm
10-255 25mm
10-256 30mm
Straight



Capener Osteotome With Tuffnel Handle
25cm
10-257 6mm
10-258 8mm
10-259 10mm
10-260 12mm
10-261 15mm
10-262 20mm
10-263 25mm
10-264 30mm
Curved



Raspatories
10-265 20cm
Medium straight blade
round edge, width 14mm



Raspatories
10-266 20cm
Large straight blade
round edge, width 20mm



Raspatories
10-267 20cm
Medium curved blade
round edge, width 14mm



Raspatories
10-268 20cm
Large curved blade
round edge, width 20mm



Mini Lexer Osteotome
10-269 20cm
18mm



Raspatories
10-270
28cm



Raspatories
10-271
28cm



Chisel
10-272
115mm/280mm
for extraction of cement

BONE CHISELS & GOUGES



Knochenmeißel

Ciseaux-burins, gouges

Cineles, gubias

Scalpelli, sgorbi



Chisel
10-273
for extraction
of cement 9mm
Positive 400mm



Chisel
10-274
with 2 prong



**Lambotte
Osteotome**
25cm
10-275 6mm
10-276 8mm
10-277 10mm
10-278 13mm
10-279 15mm
10-280 20mm
10-281 25mm
10-282 30mm
10-283 38mm
10-284 44mm
10-285 50mm
Straight



**Lambotte
Osteotome**
25cm
10-286 6mm
10-287 8mm
10-288 10mm
10-289 13mm
10-290 15mm
10-291 20mm
10-292 25mm
10-293 30mm
10-294 38mm
10-295 44mm
10-296 50mm
Curved



**Mini-Lambotte
Osteotome**
12.5cm
10-297 3mm
10-298 4mm
10-299 6mm
10-300 8mm
10-301 10mm
10-302 12mm
10-303 15mm
Straight



**Mini-Lambotte
Osteotome**
17cm
10-304 3mm
10-305 4mm
10-306 6mm
10-307 8mm
10-308 10mm
10-309 12mm
10-310 15mm
Straight



Lambotte Osteotome
12.5cm, Straight
10-311 2mm tip width
10-312 4mm tip width
10-313 6mm tip width
10-314 8mm tip width
10-315 10mm tip width
10-316 12mm tip width
10-317 15mm tip width



Hex Handle Osteotome
15cm, Straight
10-318 2mm tip width
10-319 3mm tip width
10-320 4mm tip width
10-321 5mm tip width
10-322 6mm tip width
10-323 7mm tip width
10-324 9mm tip width
10-325 10mm tip width



Cottle Osteotome
18.5cm, Straight
10-326 4mm tip width
10-327 6mm tip width
10-328 9mm tip width
10-329 12mm tip width



**Anderson-Nievert
Osteotome**
20.4cm
10-330 Curved right
10-331 Curved left
10-332 Straight



Freer Gouge
10-333 15cm
06mm



Freer Chisel
10-334 16cm
04mm



**Shoulder
Penetrating Gouge**
10-335 21.6cm
Straight, 3mm distal tip



**Ballenger
Chisel** 16cm
10-336
4mm
10-337
6mm
10-338
8mm



**Mannerfelt
Modif.**



Osteotomes
10-339 4mm
10-340 6mm
10-341 8mm
15.5cm



Gouge
10-342
4mm
15.5cm



**Kilian,
Septum Gouge**
10-343 16cm
5mm



**Kilian-Claus
Septum Gouge**
10-344 16cm
5mm



FIMA INSTRUMENTS

BONE CHISELS & GOUGES

Knochenmeißel

Ciseaux-burins, gouges

Cineles, gubias

Scalpelli, sgorbi



Stille
20cm



Stille Chisel
10-345 10mm
10-346 15mm
10-347 20mm
10-348 25mm
Straight



Stille Osteotome
10-349 10mm
10-350 15mm
10-351 20mm
10-352 25mm
Straight



Stille Osteotome
10-353 10mm
10-354 15mm
10-355 20mm
10-356 25mm
Curved



Stille Gouge
10-357 10mm
10-358 15mm
10-359 20mm
10-360 25mm
Straight



Stille Gouge
10-361 10mm
10-362 15mm
10-363 20mm
10-364 25mm
Curved



Smith-Petersen
20cm 8"



Smith-Petersen osteotome
10-365 6mm
10-366 9mm
10-367 13mm
10-368 16mm
10-369 19mm
10-370 25mm
10-371 32mm
Straight



Smith-Petersen Osteotome
10-372 6mm
10-373 9mm
10-374 13mm
10-375 16mm
10-376 19mm
10-377 25mm
10-378 32mm
Curved



Smith-Petersen Gouge
10-379 9mm
10-380 16mm
10-381 19mm
10-382 25mm
Straight



Smith-Petersen Gouge
10-383 9mm
10-384 16mm
10-385 19mm
10-386 25mm
Curved



Hoke Osteotome
14cm
10-387 3mm
10-388 5mm
10-389 7mm
10-390 8.5mm
10-391 10mm
10-392 12mm



Hibbs
24cm



Hibbs Osteotome
10-393 6mm
10-394 13mm
10-395 19mm
10-396 25mm
10-397 32mm
Straight



Hibbs Chisel
10-398 6mm
10-399 13mm
10-400 19mm
10-401 25mm
10-402 32mm
Straight



Hibbs Osteotome
10-403 6mm
10-404 13mm
10-405 19mm
10-406 25mm
10-407 32mm
Curved



Hibbs Gouge
10-408 6mm
10-409 13mm
10-410 19mm
10-411 25mm
10-412 32mm
Straight



Hibbs Gouge
10-413 6mm
10-414 13mm
10-415 19mm
10-416 25mm
10-417 32mm
Curved

BONE CHISELS & GOUGES



Knochenmeißel

Ciseaux-burins, gouges

Cineles, gubias

Scalpelli, sgorbi



Cottle Osteotome
10-418 7mm
10-419 12mm
15cm



Cottle Osteotome
10-420 19cm
04mm



Cottle Osteotome with depth markings
10-421 18cm
06mm



Fomon Chisel
17cm



10-422
04mm



10-423
07mm



10-424
05mm




10-425
08mm



10-426
06mm



Cinelli Chisel
18cm



10-427 10mm
10-428 12mm
10-429 14mm
10-430 16mm



10-431
04mm



10-432
09mm



10-433
07mm



10-434
12mm



Cottle Osteotomes
10-435 18cm
16mm Protected
and depth markings



10-436
04mm



10-437
09mm



10-438
07mm



10-439
12mm

Cottle Chisel
18cm



Cottle Osteotomes
10-440 18cm
06mm Cross bar
and depth markings



Cottle Osteotomes
10-441 18cm
16mm Cross bar and
depth markings



10-442
3mm tip width
10-443
4mm tip width
10-444
6mm tip width
10-445
8mm tip width
10-446
10mm tip width
10-447
13mm tip width
10-448
16mm tip width
10-449
19mm tip width

Hoke Osteotome
13.5cm
Curved



10-450
3mm



10-451
4.5mm



10-452
4mm



10-453
9mm



10-454
7mm



10-455
12mm

Cottle Osteotome
20cm



10-456
10mm



10-457
12mm



10-458
14mm

Rubin osteotome
with stabilizer to reduce
chance of twisting for
more control
18cm



10-459
18cm
Straight



10-460
18cm
right



10-461
18cm
Left

Silver Rhinoplastic Chisel



10-462
20cm
Straight



10-463
20cm
right





10-464
20cm
Left

Neivert-Anderson osteotome



FIMA INSTRUMENTS

CHISELS, OSTEOTOMES FOR RHINOPLASTIC*Meißel zur nasenplastik**Osteotomes, ciseaux-burins pour rhinoplastique**Cinceles para rinoplastica**Scalpelli per plastica nasale*

Hibbs Chisel
22.5 cm straight
10-485 6 mm
10-486 9 mm
10-487 12 mm
10-488 15 mm
10-489 18 mm
10-490 20 mm
10-491 25 mm
10-492 30 mm
10-493 35 mm
10-494 40 mm

22.5 cm curved
10-495 6 mm
10-496 9 mm
10-497 12 mm
10-498 15 mm
10-499 18 mm
10-500 20 mm
10-501 25 mm
10-502 30 mm
10-503 35 mm
10-504 40 mm

Osteotomes
22.5 cm straight
10-505 6 mm
10-506 9 mm
10-507 12 mm
10-508 15 mm
10-509 18 mm
10-510 20 mm
10-511 25 mm
10-512 30 mm
10-513 35 mm
10-514 40 mm

22.5 cm curved
10-515 6 mm
10-516 9 mm
10-517 12 mm
10-518 15 mm
10-519 18 mm
10-520 20 mm
10-521 25 mm
10-522 30 mm
10-523 35 mm
10-524 40 mm

Hibbs Chisel / Osteotomes



Cobb Chisel
27.5 cm straight
10-525 6 mm
10-526 9 mm
10-527 12 mm
10-528 15 mm
10-529 18 mm
10-530 25 mm
10-531 35 mm

27.5 cm Curved
10-532 6 mm
10-533 9 mm
10-534 12 mm
10-535 15 mm
10-536 18 mm
10-537 25 mm
10-538 35 mm

Osteotomes
27.5 cm Straight
10-539 6 mm
10-540 9 mm
10-541 12 mm
10-542 15 mm
10-543 18 mm
10-544 25 mm
10-545 35 mm

27.5 cm curved
10-546 6 mm
10-547 9 mm
10-548 12 mm
10-549 15 mm
10-550 18 mm
10-551 25 mm
10-552 35 mm

Cobb Chisel / Osteotomes



Hibbs Chisel Endoscopic
10-553 46cm
10 mm tip width straight
Hibbs Osteotomes Endoscopic
10-554 46cm
10 mm tip width straight

MALLETS



Hammer | Martillos
Maillets | Martelli



Doyen Mallets
25cm
10-555
25mmØ - 120g
10-556
30mmØ - 180g
10-557
35mmØ - 260g



Collin Mallets
10-558
20cm
30mmØ - 220g



Hajek Mallets
10-559
21.5cm
27mmØ - 140g



Kirk Mallets
10-560 19cm
38mmØ - 750g
10-561 20.4cm
Head diameter
38.0mm,
Weight: 1 lb. 4 oz.
570 grams



Partsch Mallets
10-562 18cm
22mmØ - 200g



Lucae Mallets
10-563
19mmØ - 220g
10-564
25mmØ - 260g
19cm



Mallet Mallets
10-565 19cm
250mmØ - 205g
Solid steel, with
plastic faces



Cottle Mallets
10-566
18.5cm
30mmØ - 240g
1 side flat /
1 side convex



Williger Hammer
10-567 250g
10-568 350g



Neufield Hammer
10-569



**Solid Copper
Head Mallet**
10-570 16.5cm
Solid, Head
diameter
20mm



Nylon Head Mallet
10-571 19cm
Head diameter 25.4mm,
7.76 oz.
10-572 19.5cm
Head diameter 3cm
Weight 11 oz.



**Extra 3cm Nylon Face
FOR 10-572**
10-573 Head diameter
3cm



Combination Mallet
10-574 19.5cm
Solid, Head
diameter 32.0mm,
Weight 12 oz 340
grams, Head weight:
8 oz 227 grams



**Bergmann
Mallets**
10-575
24.5cm
45mmØ - 300g

MALLETS

Hammer	Martillos
Maillets	Martelli



Ombredanne Mallets
10-576
24cm
40mmØ - 550g



Heath Mallets
10-577 18cm
40mmØ - 700g
10-578 20cm
45mmØ - 1000g



Williger Mallets
10-579
28mmØ - 340g
10-580
30mmØ - 400g
24cm
Lead filled



Hajek Hammer
10-581 250gm
10-582 500gm



Universal Hammer
10-583 500gm
10-584 750gm



Mallet, Lightweight
10-585 19cm
Stainless Steel and Aluminum
Head diameter 25.0mm
Weight 7 oz 200 grams
Head Weight: 4 oz 113 grams



Nylon Caps Left & Right
FOR 10-585 (PAIR)
10-586
2.5cm



Mallet #49
10-587 19.5cm
Head diameter
22.0mm
Weight 7 oz 200
grams
Head Weight:
5oz/150 grams



Mead Mallet
10-588 19cm
Head diameter
26.0mm
Weight 11 oz 315
grams
Head Weight:
6 oz 170 grams



Heath Mallet, Heavyweight
10-589 19cm
Head diameter
40.0mm
Weight 2 lbs.
(907 grams)
Head Weight:
20 oz (567 grams)



Orthopedic Mallet
10-590 27.5cm
Lightweight
Solid, Head
diameter 35.0mm
Weight 2 lbs.
(907 grams)
Head Weight:
1 lb. (454 grams)



Orthopedic Mallet
10-591 28cm
Heavyweight
Solid, Head
diameter 50.0mm
Head Weight:
2 lb. (900 grams)



Ortho Mallets
10-592
26.5cm
(450g)



Nylon Cap
10-593
For Phenolic Handle Mallet

Phenolic Handle Mallet
10-593 28 cm
Head diameter 32mm
Head weight 16 oz (454 grams)
10-594 Head diameter 45mm
Weight 11 oz (312 grams);
Head weight: 2 lbs (907 grams),



Nylon Cap
10-593
For Phenolic Handle Mallet

Phenolic Handle Mallet
10-593 28 cm
Head diameter 32mm
Head weight 16 oz (454 grams)
10-594 Head diameter 45mm
Weight 11 oz (312 grams);
Head weight: 2 lbs (907 grams),

MALLETS



Hammer | Martillos
Maillets | Martelli



**Reprocussion
Free Mallet**
Cap For 10-596

Reprocussion Free Mallet
10-595 28cm
Head diameter 30mm
Weight 1lb 5 oz (595 grams)
Head weight: 10 oz (283 grams)



Mallet
10-597 26cm
Plane (Nylon)
convex (steel)
Head 30mm Ø
weight 675 grams.



Mallet
10-598 26cm
Head 42/30mm Ø
weight 530 grams



Mallet
10-599 26cm
Plastic Handle
head 30mm Ø
weight 510 grams



Vickers Mallets
10-600 19 cm
25mmØ - 300g
nylon disks plastic



Herbert Mallets
10-601 19 cm
25mmØ - 300g
for ulnar head
prosthesis plastic
handle



Mallet
10-602 18cm
Made of Plastic
Head 30mm Ø
weight 125 grams



Mallet
10-603 26 cm
Made of Plastic
with lead filling
Head 42mm Ø
weight 250 grams



Mallet
10-604 22.5 cm
Disks & Plastic Handle
Head 37mm Ø
weight 370 grams



Spare Facing (pair)
10-605 Plastic



Mallet
26 cm
Made of Plastic
with lead filling
10-606
Head 42 mm Ø
weight 400 grams
10-607
Head 42 mm Ø
weight 250 grams



Mallet
10-608 24cm
Plastic Handle
Head 30mm Ø
weight 480
grams



Mallet
10-609



Mallet
10-610







Mallet
10-611



Mallet
10-612



FIMA INSTRUMENTS

PERIOSTEAL ELEVATORS, TAMPERS*Raspatorien, Stobel
Periostótomos, Pilons**Rugines-raspatoires, Punzones
Periostotomi, Pestoni***Cobb Spinal Elevator**
23.8 cm sharp
with hexagonal handle**10-613**
12 mm**10-614**
19 mm**10-615**
25 mm**10-616**
32 mm**Harrington Spinal Elevator**
10-617 30cm
19 mm tip width
Sharp**Cobb Modif.**
28 cm
Round Handle**10-618**
7 mm**10-619**
9 mm**10-620**
11 mm**Cobb Modif.**
28 cm
Round Handle**10-621**
4 mm**10-622**
6 mm**10-623**
8 mm**10-624**
17.5 cm**10-625**
15.5cm
Head diameter
2mm**10-626**
15.5cm
Head diameter
3mm**10-627**
15.5cm
Head diameter
4mm**10-628**
15.5cm
Head diameter
5mm**10-629**
15.5cm
Head diameter
6mm**10-630**
15.5cm
Head diameter
8mm**10-631**
15.5cm
Head diameter
10mm**10-632**
2.5mm x 2.5mm**10-633**
3.5mm x 3.5mm**Passow
Bone Tamper**
12cm**10-634**
3mm**10-635**
5mm**10-636**
8mm**Caspar
Bone Tamper**
20cm**10-637**
12mm**10-638**
14mm**10-639**
16mm**Caspar
Bone Tamper**
20cmTip Size
2 mm
3 mm
4 mm
5 mm
6 mm
8 mm
10 mm**Bone Tamps**
10-640
15.4 cm

BONE CURETTES



Scharfe Löffel

Curettes tranchantes

cucharillas cortantes

Cucchiai taglienti



**Volkmann
Bone Curette**

10-641
17cm
0000

10-644
17cm
0

10-647
17cm
3

10-642
17cm
000

10-645
17cm
1

10-648
17cm
4

10-643
17cm
00

10-646
17cm
2

10-649
17cm
5

10-650
17cm
6



**Schede
Bone Curette**
17cm

10-651
17cm
000

10-653
17cm
0

10-655
17cm
2

10-652
17cm
00

10-654
17cm
1

10-656
17cm
3

10-657
17cm
4



**Bruns
Bone Curette**

10-658
17cm
000

10-660
17cm
0

10-662
17cm
2

10-664
17cm
4

10-659
17cm
00

10-661
17cm
1

10-663
17cm
3

10-665
17cm
5

10-666
17cm
6



**Spratt Bone Curette
/ Hollow Handle**

10-667
16.5cm
00000

10-670
16.5cm
00

10-673
16.5cm
2

10-668
16.5cm
0000

10-671
16.5cm
0

10-674
16.5cm
3

10-669
16.5cm
000

10-672
16.5cm
1

10-675
16.5cm
4

10-676
16.5cm
5



**Halle Bone Curette
Flexible**

10-677
21cm
1

10-678
21cm
2

10-679
21cm
3



**Uffenorde Bone
Curette Left**
10-680 17.5cm



**Uffenorde Bone
Curette Right**
10-681 17.5cm



FIMA INSTRUMENTS

BONE CURETTES

Scharfe Löffel


Curettes tranchantes

cucharillas cortantes


Cucchiai taglienti




Bruns Bone Curette
23 cm




Bone Curette
Straight




Bone Curette
Curved



Schroeder Bone Curette
Sharp / rigid



Simon Bone Curette
Sharp / rigid



BONE CURETTES



Scharfe Löffel

cucharillas cortantes

Curettes tranchantes

Cucchiari taglienti



Bone Curette
25 cm
Straight




Daubenspeck Bone Curette
20 cm



Scoville Bone Curette
25 cm




Lempert Bone Curette
21.5 cm



Cushing Bone Curette
23.5 cm
Pituitary Spoon,
malleable, silver-plated



Scoville Bone Curette
25 cm
Ruptured Disc Curette



Tay-Bronson Bone Curette
22 cm
Pituitary Curette, double-edged



BONE CURETTES

Scharfe Löffel
Curettes tranchantes


cucharillas cortantes
Cucchiai taglienti



Ray Bone Curette, 19 cm
Pituitary Curette

Volkmann
Curette
Oval / Oval

Volkmann
Curette
Round / Oval



Curette
10-774
21cm
round / oval

Martini
Curette
13.5cm

Martini
Curette
15.5cm

Williger
Curette
13.5cm

Williger
Curette
15.5cm

BONE CURETTES



Scharfe Löffel

Curettes tranchantes

cucharillas cortantes

Cucchiari taglienti



Jansen Curette
10-793 16cm



Jansen Curette
10-794 15cm



Williger Curette
10-795 13.25cm
Double Ended 00/0



10-796 13cm
Double Ended
6.2 mm & 6.7 mm

10-797 16cm
Double Ended
7 mm & 9 mm

10-798 20cm
Double Ended
9.1 mm & 11.1 mm

Curette



Lucas Curette
double-ended
18cm

Fig. 1 R/L 10-799
Fig. 2 R/L 10-800
Fig. 3 R/L 10-801
Fig. 4 R/L 10-802

PERIOSTEAL ELEVATORS



Obwegeser Elevator



10-803
17cm
6 mm



10-804
17cm
7 mm



10-805
17cm
9 mm



10-806
17cm
11 mm



10-807
20cm
5 mm



10-808
20cm
8 mm



10-809
20cm
9 mm




Obwegeser-Freer Elevator
10-810 21cm



Obwegeser-Freer Elevator
10-811 21cm



Sachs Nerve elevator
10-812
18.5cm
3 mm blunt







Kleinert Kutz Elevator and Dissector
10-813
19.5cm
2 mm x 3 mm
Sharp ends



Freer Elevator
10-814 18.5cm
double ended
5 mm sharp end
5 mm blunt end



FIMA INSTRUMENTS

PERIOSTEAL ELEVATORS*Raspatorien*
*Periostótomos**Rugines-raspatoires*
Periostotomi**Woodson Elevator/
Spatula
10-815**
22.5cm
double ended
3 mm semi-
blunt ends**Woodson Dural
Separator/Packer
10-816** 17.5cm
double ended with
grooved
blade 3 mm blunt
end**Penfield
Elevators
10-817**
style #1
17.5 cm
double ended**Penfield
Elevators
10-818**
style #2
18.5 cm
double ended**Penfield
Elevators
10-819**
style #3
28 cm
double ended**Penfield
Elevators
10-820**
style #4
28 cm
single
ended**Penfield
Elevators
10-821**
style #5
18.5 cm
double
ended**Alderdyce
Elevators
10-822**
21.5 cm
double
ended**Mc Donald
Elevators
10-823**
19 cm
double
ended**Watson-
Cheyne
Elevators
10-824**
19 cm
double
ended**Obwegeser
Needle
10-825**
17.5cm**Obwegeser
Needle
10-826**
23cm**Obwegeser
Needle
10-827**
18cm**Obwegeser
Needle
10-828**
14cm**10-829**
Right Crved
For Adults

10-830
Right Crved
For Children**Doyen Rib Raspatory**
17cm**10-831**
Left Crved
For Adults

10-832
Left Crved
For Children**Doyen Rib Raspatory**
17cm**Sedillot
Elevator
10-833**
18.5cm
10-834
21cm

PERIOSTEAL ELEVATORS



Raspatorien
Periostótomos

Rugines-raspatoires
Periostotomi



**Crego
Periosteal
Elevator
10-835**
15.5cm



**Williger
Elevator
10-836**
13cm
5mm




**Williger
Elevator
10-837** 4mm
10-838 6mm
16cm



**Williger
Elevator
10-839**
16cm
5mm




**Williger
Elevator
10-840**
16cm
6.5mm



**Farabeuf-
Collin
Elevator
10-841**
16cm
11mm
Straight



**Farabeuf-
Collin
Elevator
10-842**
16cm
11mm
Curved




**Farabeuf
Elevator
10-843**
15cm
Straight



**Farabeuf
Elevator
10-844**
15cm
Curved




Lambotte Elevator
10-845 5mm
10-846 10mm
10-847 15mm
10-848 20mm
10-849 25mm
21cm Straight



Lambotte Elevator
10-850 5mm
10-851 10mm
10-852 15mm
10-853 20mm
10-854 25mm
21cm Curved



Adson Elevator



10-855
17cm
1



10-856
17cm
2



10-857
17cm
3



10-858
17cm
4



Quervain Elevator
10-859 06mm
19.5cm




Trelat Elevator
10-860 8mm
15cm






FIMA INSTRUMENTS

PERIOSTEAL ELEVATORSRaspatorien
PeriostótomosRugines-raspatoires
Periostotomi

10-861 4.5mm
15.5cm
Half Curved
König Elevator



10-862 4.5mm
15.5cm
Full Curved



Obwegeser Elevator



10-863
15cm
2.5mm



10-864
15cm
2.5mm



Obwegeser Elevator
10-865
18.5cm




Obwegeser - Freer Osteotome Elevator
10-866
18.5cm



Howarth Periosteal Elevator
10-867
21.5cm
Sharp / Blunt



Langen-Beck Periosteal Elevator
10-868
18cm



Langen-Beck Periosteal Elevator
10-869 19cm



Alexander Periosteal Elevator
10-870 20cm



Periosteal Elevator
17 cm



10-872



Semb Periosteal Elevator



10-873
17.5cm
12mm



10-874
21cm
12mm



10-875
17.5cm
12mm



10-876
21cm
13mm



10-877
22.5cm
13mm




10-878
22.5cm
15mm

PERIOSTEAL ELEVATORS




Raspatorien
Periostótomos


Rugines-raspatoires
Periostotomi




10-879
14mm
Straight
1



10-880
20mm
Straight
3



10-881
14mm
Curved
2




10-882
20mm
Curved
3

Ferrozell Periosteal Elevator
20cm



10-883
6mm
1



10-884
3mm
2



10-885
6mm
3

Ferrozell Periosteal Elevator
18cm



Periosteal Raspatories
10-886
18.5cm
cvd



Periosteal Elevator
10-887 18cm
Fiber/Phenolic Handle
blade 6mm




Hoen Periosteal Elevator
10-888
16.5cm



Chandler Bone Elevator
10-889 19cm
13 mm
10-890 20cm
20 mm
10-891 23cm
25 mm
10-892 25cm
31 mm



Key Periosteal Elevator 19cm
10-893
6 mm
10-894
12 mm
10-895
20 mm
10-896
25 mm
sharp,
slightly
curved



Pennybacker Periosteal Elevator
10-897
6 mm
17 cm sharp,
curved



Langenbeck Periosteal Elevator
10-898
15 mm
19 cm
sharp, slightly
curved



Cushing Periosteal Elevator 15 mm
10-899
Round Edge
10-900
Square Edge
18 cm



FIMA INSTRUMENTS

PERIOSTEAL ELEVATORSRaspatorien
PeriostótomosRugines-raspatoires
Periostotomi

Hoen Periosteal Elevator
10-901 7mm, 16cm
10-902 15mm, 19cm
10-903 20mm, 19cm
 Semisharp



Sayre Periosteal Elevator
10-904 17cm
 one straight blunt
 and one curved
 sharp




Langenbeck Periosteal Elevator
10-905 7mm
10-906 10mm
 20 cm Blunt

7 mm

10 mm




Cushing Periosteal Elevator
 19cm
10-907
 Left
10-908
 Right



Cobb Periosteal Elevator
10-909 10mm
10-910 13mm
10-911 19mm
 29cm
 Spinal Raspatory



10-912
 18cm
 8mm



10-915
 20cm
 19mm

Key Periosteal Elevator



10-913
 19cm
 10mm



10-914
 19cm
 12mm



10-916
 21cm
 25mm



Cobb Periosteal Elevator
10-917 23.4cm
 Hex Handle



Darrach Elevator
10-918 26cm
 19mm
10-919 35cm
 25.4mm
 Blunt, Slightly
 Curved



Carroll Periosteal Elevator
10-920 14cm
 sharp, double ended
 3 mm round end,
 10 mm flat end



Langenbeck Periosteal Elevator
 19.5cm
10-921
 blunt edge
10-922
 sharp edge



Key Periosteal Elevator sharp edge
10-923 17.7 cm 6.5 mm width
10-924 18.5 cm 9.5 mm width
10-925 18.5 cm 12.7 mm width
10-926 20.0 cm 19.0 mm width
10-927 20.9 cm 25.0 mm width



Lane Bone Elevator
10-928
 16cm
 Serrated



Lane Bone Elevator
10-929
 27cm



Trethowan's Bone Elevator
10-930
 22cm



Bristow's Bone Lever
10-931
 24cm

PERIOSTEAL ELEVATORS, MENISCOTOMY KNIVES



FIM INSTRUMENTS

Raspatorien, Meniskotome

Rugines-raspatoires, Bisturies para meniscos

Periostótomos, Bistouries a menisques

Periostotomi, Bisturi per menisco



**St. Thomas
10-932**
16cm
Pointed



**St. Thomas
10-933**
16cm
Rounded



**Robinson's
10-934**
17cm



**Freer D
Septum
Knife
10-935**
15cm




**Bircher
Meniscotomy
Knife
10-936**
19cm



**Matson
Meniscotomy
Knife
10-937**
22cm



**Downing
Cartilage
Knife
10-938**
25.5cm
Curved tip,
Straight shaft



**Smillie
Meniscotomy
Knife
10-939** 16cm
Straight Shaft
Blade tip width
is 7.8mm
Straight tip




**Smillie
Meniscotomy
Knife
10-940** 16cm
Curved Shaft
Blade tip width
is 7.8mm
Curved left tip



**Kanadisches
Modell
Meniscotomy
knife
10-941**
20cm



**Kanadisches
Modell
Meniscotomy
knife
10-942**
20cm



**Meniscotomy
Knife
10-943** 18cm



**Neff Meniscotomy
Knife**



10-944
22cm
left



10-945
22cm
right



**Bunnell Tendon
Strippers**
23cm
Malleable shaft



10-946
4mm




10-947
4.5mm



10-948
5.5mm



**Bunnell
Tendon
Strippers**
15.2cm




10-949
Fig. 1
3mm



10-951
Fig. 3
5mm



10-950
Fig. 2
4mm



10-952
Fig. 4
6mm



10-953
Fig. 5
7mm



**Fischer
Tendon
Strippers
10-954**
31cm

TENDON STRIPPERS, BONE LEVERS

Sehnenstripper, Knochenhebel

Instruments a ecorcher les tendons, Leviers a os

enucleadores de tendones, Elevadores para huesos

Scollatendini, Elevatori per ossa



Bunnell Tendon Strippers
10-955 23cm
Tendon passer



Kleinert-Kutz Tendon Retriever
10-956 20cm
Serrated



Kleinert-Kutz Tendon Retriever
10-957 20cm
Serrated, flexible shaft



Brand Tendon Tunelling Forceps
10-958 15cm
10-959 19cm
10-960 24cm
Straight, serrated
2 x 3 teeth



Hohmann Bone Elevator
10-961 43mm
24.5 cm



Hohmann Bone Elevator
10-962 22mm
26cm



Hohmann Bone Elevator
10-963 33mm
28cm



Hohmann Bone Elevator
10-964 24mm
26.5cm



Hohmann Bone Elevator
10-965 22mm
26cm



Capsule Retractor
10-966 25.5cm
3 sharp-pronged,
22mm width blade



Bristow-Bankart Humeral Retractor
10-967 19cm
Double-pronged,
25mm width x 102mm
depth blade




Bristow-Bankart Soft Tissue Retractor
10-968 20cm
3-pronged
21mm width blade,



Humeral Head Retractor
10-969 22cm
18mm width blade



Hohmann Bone Elevator
10-970 22cm
08mm width blade
10-971 23cm
10mm width blade



Hohmann Bone Elevator
10-972 24cm
18mm



Hohmann Bone Elevator
10-973 24cm 70mm
10-974 28.5cm 35mm width blade
10-975 33cm 40mm width blade

BONE LEVERS



Knochenheber
Leviers a os

Elevadores para huesos
Elevatori per ossa



**Hohmann Long
Shank Retractor**
10-976 23cm
43mm width blade



**Mini-Hohmann
Bone Elevator**

10-977 16cm
6mm width
blade



10-978 16cm
8mm width
blade



**Hohmann
Bone Lever**
10-979
15cm



**Hohmann
Bone Lever**
10-980
15.4cm
Short tip



**Hohmann
Retractor**
10-981
15cm
15mm width
blade



**Hayes Hand
Retractor**
10-982
16cm
3mm width
blade



**Hayes Hand
Retractor**
10-983
15cm
6mm width
blade



10-984
14cm
2.8mm



10-985
14cm
4.5mm



Bone Lever
10-986
14cm
6mm




**Bristow Bone
Lever**
10-987
23cm



**Lane Bone
Lever**
10-988
21cm



Periosteal Elevator 17.5cm
For Metacarpals Small
10-989 6 mm width
For Metacarpals Medium
10-990 10 mm width
Light Curve, Rounded End,
10-991 5 mm width
Contra Angled Curve
10-992 4.5 mm width



Hohmann Retractor
10-993 18cm
Small, width 6 mm
short, narrow tip



Hohmann Retractor
10-994 18cm
Small, width 8 mm
short, narrow tip



Hohmann Retractor
10-995 12.5cm
Finger Ring handle,
16mm width blade



FIMA INSTRUMENTS

BONE LEVERS / FORCEPS

Knochenhebel / Pinzette

Palancas de hueso / Pinzas

Leviers / pincers à os

Leve / pinze ossee




**Hohmann-Williger
Retractor
10-996**
16cm



**Blount Retractor
10-997**
17.5cm



**Taylor Retractor
16.5cm
10-998**
a=75mm - b=30mm
10-999
a=100mm - b=30mm



**Cushing Vein
Retractor
10-1000** 22.5cm
with fenestrated
handle
Blade 12 mm
width x 10 mm



**Cushing
Decompression
Retractor
10-1001** 21.5cm
with fenestrated
blade and handle
Blade 13 mm
width x 22 mm



**Scoville Nerve
Retractor
10-1002**
24cm
7 mm
straight



**Scoville Nerve
Retractor
10-1003**
24 cm, 7 mm
angled



**Taylor Retractor
17cm
10-1004** 30 x 75 mm
10-1005 30 x 100 mm



**Hibbs Retractor
10-1006**
25 x 75 mm, 23cm
10-1007
25 x 80 mm, 24cm



**Meyerding Retractor
10-1008**
16 X 50 mm, 20cm
10-1009
25 X 75 mm, 23cm
10-1010
50 X 90 mm, 24cm



**Nail
Extracting
Forceps
10-1011**
14cm



**Nail Extracting
Forceps
10-1012**
14.5cm

ORTHOPAEDICS, TENDON STRIPPERS




Orthopädie, Sehnenstripper

Orthopédie, enucleadores de tendones

Ortopedia, Instruments a ecorcher les tendons

Ortopedia, Scollatendini





Flat Nose Pliers
10-1033
13.5cm
Broad, 4mm tip,



Flat Nose Pliers
10-1034 16.5cm
Combination plier and
cutter



Needle Nose Pliers
10-1035 13.5cm



Flat Nose Pliers
10-1036
15cm Delicate,
1.5mm tip



K-Wire Extraction Pliers
W/Screw Lock
10-1037
13cm, 3mm tip



Flat Nose
K-Wire Pliers
10-1038
14cm
10-1039
17cm



Nose Pliers
10-1040
17cm



Utility Technical
Forceps Pliers
10-1041 19cm
adjustable



Wire Bending
Forceps
10-1042 18.5cm



Flat Nose Plier
10-1043 18.5cm
parallel grip, incorporated
cutter, for hard wires up to
1.6mm

TECHNICAL FORCEPS

Technische Zangen
Pinces techniques

Pinzas técnicas
Pinze tecniche



Locking Plier
10-1044 18cm
10-1045 20cm
10-1046 24cm



Plate Bending Plier
10-1047 20cm
for plates up to 2.7mm
10-1048 24cm
for plates up to 4mm



Locking Plier
10-1049 22cm
10-1050 25cm
10-1051 30cm
Locking plier, only
10-1052 30cm
with holding pins
10-1053 46cm
with extraction device



Needle-Nose Locking Pliers
10-1054



Needle-Nose Locking Pliers
10-1055



Needle-Nose Locking Pliers
10-1056




Flat Nosed Plier
10-1057 18cm
parallel, with lateral
wire cutter



**Cast Spreader
Model Usa**
10-1058 23cm
3 Prong




Nail Extractor
10-1059 16cm
with Longitudinal Groove
(for Intramedullary Gliding Nails)




10-1060
15cm
10-1061
18cm
Fig.1



10-1062
15cm
10-1063
18cm
Fig.2



10-1064
17cm
Fig. 1



10-1065
17cm
Fig. 2




Friedmann Bone Rongeur
10-1066 14.5cm
Single action, Curved,
4 x 11mm


BONE RONGEURS

Hohlmeißelzangen
Pinzas gubias


Pinces-gouges
Pinze-sgorbie




Blumenthal Bone Rongeur




Jasen Bone Rongeur



Jansen Bone Rongeur
10-1072 19cm



Stellbrink Bone Rongeur



Lempert Bone Rongeur



Rongeur #4
10-1079 14cm
Rongeur #4A
10-1080 16.5cm
Single action,
Curved, 4 x 9mm



Chichirelli (Blumenthal) Rongeur
10-1081 15cm
Single action, Curved 45°,
4 x 13mm



Dean (Blumenthal) Rongeur
10-1082 15cm
Single action,
Curved 90°, 4 x 11mm



Dean Rongeur
10-1083 14.5cm
Single action
Curved, 3mm

BONE RONGEURS



Hohlmeißelzangen
Pinzas gubias

Pinces-gouges
Pinze-sgorbie



Hartmann Bone Rongeur
10-1084 18.5cm
Single action
Curved, 4 x 14mm



Mead Bone Rongeur
10-1085 16.5cm
Single action,
Curved, 5 x 13mm



Cleveland Bone Rongeur
10-1086 14cm
10-1087 17cm



Universal Bone Rongeur

10-1088
13.25cm
Angled
17 mm x 3.5 mm



Bane Bone Rongeur



10-1089
18cm
Single action
Curved,
3 x 16mm



Bohler Bone Rongeur



10-1090
15c
Straight



10-1091
15cm
Curved



Zaufal-Jansen Bone Rongeur
10-1092 18cm
Double action
Curved, 5x15mm



Marquardt Bone Rongeur
10-1093 18cm



Beyer Bone Rongeur
18cm
10-1094 3 x 14mm
10-1095 4 x 15mm
Double action, Curved



Smith-Peterson Laminectomy Rongeur
10-1096 24cm
Double action,
Straight 3 x 16mm



Smith-Peterson Laminectomy Rongeur
10-1097 24cm
Double action
Curved 3 x 16mm

BONE RONGEURS

Hohlmeißelzangen
Pinzas gubias

Pinces-gouges
Pinze-sgorbie



Ruskin Rongeur
10-1098 18cm 4x16mm
10-1099 18.5cm 5x10mm
Double action, Curved on side



Ruskin Rongeur
10-1100 18.5cm
Double action, 6 x 16mm,
Straight, Curved



Ruskin Rongeur
10-1101 15.5cm Straight, 3.2mm
10-1102 15cm Curved, 2 x 15mm
10-1103 15cm Curved, 3 x 15mm
Double action,



Riskin Bone Rongeur



10-1104 15.5cm
10-1105 18cm
Double action
Straight, 2 x 14mm



10-1106 15.5cm
10-1107 18cm
Double action
Curved 2 x 14mm



Mayfield Bone Rongeur
10-1108 18cm



Stille-Ruskin Bone Rongeur



10-1109
23cm
Curved



10-1110
23cm
Straight



Leksell Bone Rongeur
10-1111 24.5cm 5 x 13mm
10-1112 24cm 8 x 14.5mm
Double action, Curved 50°



Leksell Rongeur
10-1113 24cm
6 x 15mm
10-1114 24cm
7 x 13mm
Double action, Curved



Leksell Rongeur
10-1115 24cm
3 x 18mm
Double action, Curved



Leksell Rongeur
10-1116 23cm
4 x 13mm
Double action
Fully curved



Stille Bone Rongeur



10-1117
23cm
4mm



10-1118
23cm
6mm



10-1119
23cm
8mm



Stille-Luer Rongeur Forceps



10-1120
23cm
10-1121
27cm
straight



10-1122
23cm
10-1123
27cm
curved



Echlin Duckbill Bone Rongeur 23cm



10-1124
2 x 12mm
Double action
Curved on side



10-1125
3 x 12mm
Double action
Curved on side



10-1126
4.5 x 11mm
Double action
Curved on side

BONE RONGEURS



Hohlmeißelzangen
Pinzas gubias

Pinces-gouges
Pinze-sgorbie



Stille-Leur Rongeur
10-1127 21.5cm
Curved, Straight
Double action, 10 x 18mm



Stille-Leur Rongeur
10-1128 26.5cm
Straight, Curved
Double action, 10 x 18mm



**Stille-Luer Echlin
Duckbill Rongeur**
23cm
10-1129 6 x 13mm
10-1130 6 x 15mm
Double action Curved on side



Stille-Leur Rongeur
10-1131 21.5cm
Double action
Curved 10 x 16mm




**Rottgen-Ruskin
Bone Rongeur**
10-1132 23cm



Semb Bone Rongeur
10-1133 24cm



Frykholm Bone Rongeur
10-1134 24cm



Sauerbruch Bone Rongeur
30cm

10-1135
Round
Fig. 1

10-1136
Straight
Fig. 2



Luer-Stille Bone Rongeur
27cm

10-1137
Fig. 1



10-1138
Fig. 2



Leksell-Stille Bone Rongeur

10-1139
24cm

10-1140
24cm

BONE RONGEURS

Hohlmeißelzangen
Pinzas gubias

Pinces-gouges
Pinze-sgorbie



Leksell Bonerongeur



10-1141 25cm 3 x 16mm
10-1142 23cm 5 x 16mm
10-1143 25cm 4.5 x 16mm
10-1144 25cm 6 x 16mm
10-1145 25cm 7 x 16mm
10-1146 25cm 8 x 16mm
Double action, Curved 15°



Hakansson Bone Rongeur
10-1147 20cm



Kleinert-Kutz Bone Rongeur



10-1148
15.25cm
2 x 13 mm



Ruskin-Mini Bone Rongeur



10-1149
15.25cm
Curved
18 x 5 mm



Sybert Rongeur
10-1150 35.5cm 4 x 16mm
10-1151 35.5cm 8 x 16mm
10-1152 35.5cm 12 x 16mm
10-1153 37cm 5 x 16mm
Double action, Angled




Cottle-Kazanjian Bone
Cutting Rongeur Forceps
10-1154 19cm



Littauer / Liston Bone
Cutting Rongeur Forceps

10-1155
15cm
Straight,



Cleveland Bone Cutting
Rongeur Forceps



10-1156
15cm
curved



10-1157
17cm
curved



Liston Bone Cutting
Rongeur Forceps



10-1158
14cm
Angled



10-1159
17cm
Angled



10-1160
19cm
Angled



10-1161
22cm
Angled

BONE CUTTING FORCEPS



Knochensplitterzangen
Pinces coupantes a os

Pinzas cortantes para huesos
Cesoie per ossa



Liston Bone Cutting
Rongeur Forceps



10-1162 14 cm
10-1163 17 cm
10-1164 19 cm
10-1165 20 cm
10-1166 22 cm
Straight



Stille-Liston Bone Cutting
Rongeur Forceps



10-1167
27 cm
Straight



10-1168
27cm
Angled 45°



10-1169
27 cm
Angled 90°



Semb Bone Cutting Rongeur
Forceps
10-1170 24cm



Liston Key Horsley Bone
Cutting Rongeur Forceps
10-1171 27cm



Horsley Bone Cutting Rongeur
Forceps 10-1172 25cm



Ruskin Bone Cutting
Rongeur Forceps
10-1173 20cm
Double-Action




Cleveland-Gardner Bone
Cutting Rongeur Forceps

10-1174
16cm



Markwalder Bone Cutting
Rongeur Forceps

10-1175
20cm



Bone Cutting Rongeur Forceps
10-1176 15cm

10-1177
15cm
Straight

10-1178
15cm
Angled




Bohler Bone Cutting
Rongeur Forceps

BONE CUTTING FORCEPS

Knochensplitterzangen
Pinces coupantes a os

Pinzas cortantes para huesos
Cesoie per ossa



10-1179
18.5cm
Straight



10-1180
18.5cm
Angled



10-1181
15cm
Straight



10-1182
15cm
Angled



10-1183 Straight
10-1184 Curved
10-1185 Angled
19cm delicate



10-1186
17.5cm
Straight



10-1187
17.5cm
Curved



Mc Indoe Bone
Cutting Forceps
10-1188 18cm
Angled, delicate,



Kleinert-Kutz Bone
Cutting Forceps
10-1189 15cm
Straight, delicate,



Kleinert-Kutz
Synovectomy Rongeur
10-1190 15cm 2 x 10mm
Double action,
Curved on side, delicate,



10-1191
27cm
Straight



Stille-Horsley Bone Cutting
Forceps
10-1192 24cm



10-1193
23cm
Angled



10-1194
27cm
Angled




Liston-Stille Bone Cutting
Rongeur Forceps
10-1195 26cm

WIRE CUTTING INSTRUMENTS



Drahtschneidezangen
Pinces coupantes

Alicates para cortar alambre
Pinze taglia filo



Vital Wire Cutter
10-1196 22cm
for hard wires up to 2.2mm,
for soft wires up to 3mm



Vital Wire Cutter
10-1197 22cm
for hard wires up to 2.2mm,
for soft wires up to 2.8mm



Vital Wire Cutter
10-1198 23cm
for hard wires up to 2.2mm,
for soft wires up to 3mm



Vital Wire Cutter
10-1199 18cm
for hard wires up to 1.5mm,
for soft wires up to 2mm



**Pin Puller W/
Groove TC**
10-1200 14cm



Vital Wire Cutter
10-1201 18cm
for hard wires up to 1.7mm,
for soft wires up to 2mm



Vital Wire Cutter
10-1202 15cm
for hard wires up to 1mm,
for soft wires up to 1.5mm



Vital Wire Cutter
10-1203 14.5cm
for hard wires up to 1mm,
for soft wires up to 1mm



**Cannulated Wire
/ Pin Cutter**
10-1204 19cm
Cuts wire upto 1.6 mm



Hard Wire Cutter Cuts
10-1205 26cm
for hard wires up to 3mm



**Kirschner Wires &
Steinmann
Pins Cutter**
10-1206 57cm
up to 6mm



Side Pin Cutter
10-1207
57.5cm
Side cutting w/
carbide inserts
max. cap: 5/32"



**Vital Heavy
Duty Pliers TC**
10-1208
22cm
Grip Forceps,
Combined
with cutter



**Vital Wire
Extraction
Pliers T.C**
10-1209 5mm
10-1210 7mm
18cm

WIRE CUTTING INSTRUMENTS

Drahtschneidezangen
Pinces coupantes

Alicates para cortar alambre
Pinze taglia Filo



Side Cutter
10-1211 18.5cm
Welded carbide jaws
max capacity. 1.6mm



Delicate Wire Cutter
10-1212 12.5cm
Side cutting w/ carbide inserts
max capacity. 0.7mm



Front & Side Wire Cutter
10-1213 17cm
Front and side cutting w/ carbide inserts,
double action. Max capacity. 1.6mm



Double Action Wire Cutter
Angled, Side cutting w/
carbide inserts
10-1214 17cm
max cap 1.6mm
10-1215 21.5cm
max cap. 2.4mm



Flush Wire Cutter
10-1216 17cm
Front and side cutting w/ carbide inserts,
double action. Max capacity. 1.5mm



Wire Forceps
10-1217 17cm
Straight, serrated 5.4 mm tip



Pin Cutter Small
10-1218 26cm
Max cap: 3.2mm



Pin Cutter Medium
10-1219 36.5cm
Max cap: 4.7mm



Large Pin Cutter
10-1220 57.5cm
Large, Straight, Max.
capacity 6.3mm



Stille Plaster Shear
10-1221 37cm



**Schumacher Sternum
Shears**
10-1222 21cm

BONE HOLDING CLAMPS - BONE HOLDING FORCEPS

Knochenhalteklammern - Knochenhaltezangen
Daviers a os - Daviers porte-so

Pinzas porta huesos - Pinzas para huesos
Tenaglie per ossa




Gerster-Lowmann Bone Clamp
10-1223 12cm
10-1224 14.5cm
10-1225 18cm, 28mm
10-1226 20cm, 35mm
10-1227 22cm, 40mm
1x2 Jaws



Lowmann Bone Clamp
10-1228 12cm
10-1229 15cm
10-1230 17cm
10-1231 18.5cm
10-1232 20cm




Lambert-Lowman Bone Clamp
10-1233 17.5cm
10-1234 20.3cm
2x2 Jaws



Gerster Bone Holding Clamp
10-1235 15cm



Gerster Bone Holding Clamp
10-1236



Gerster Bone Holding Clamp
10-1237



Parham-Martin Bone Clamp
10-1238 21cm
Straight



Parham Martin Bands
10-1239 19cm

FIMA INSTRUMENTS

BONE HOLDING FORCEPS - REPOSITION FORCEPS

Knochenhaltezaugen - Repositionszangen
Daviers porte-so-pinces de reposition

Pinzas para huesos - pinzas para reposicion
Tenaglie per ossa - Pinze per riposizioni



Bone Reduction Forceps
10-1240 13.5cm



Locke Phalangeal Forceps
10-1241 14cm



10-1242
 17.5cm
 serrated



Stagbeetle Bone Reduction Forceps
10-1243 Pointed tips, 13cm



10-1244
 13.5cm
 Angled
 Serrated



Shoulder Percussion Awl
 Pointed, 21.6cm
10-1245 Small
10-1246 Medium
10-1247 Large
10-1248 Pointed, 23.5cm



Meniscus Clamp
10-1249 13.5cm
 Serrated



Meniscus Clamp
10-1250 13.5cm
 Angled 20°, serrated



Preston Logamentum Flavum Forceps
10-1251 19cm
 Curved, 1 x 2 teeth



Meniscus Clamp
10-1252 18cm
 Straight, Curved shaft



Dingman Bone Clamp
10-1253 19cm
 Straight, Serrated,
 1 1/2" tip length, Curved shaft



Fergusson Bone Holding Forceps
10-1254 21cm 2 x 2 teeth

BONE HOLDING FORCEPS



Knochenhaltezaugen
Daviers porte-so

Pinzas para huesos
Tenaglie per ossa



Van Buren Sequestrum Forceps
10-1255 23cm
Curved on side, serrated



Langenbeck Bone Holding Forceps
10-1256 21cm



Semb Bone Holding Forceps
10-1257 19cm without ratchet
curved sideways deep teeth



Semb Bone Holding Forceps
10-1258 20.3cm, with ratchet
curved sideways deep teeth



Sequestrum (Nircosis) Bone Holding Forceps
10-1259 17.8cm
10-1260 20cm
10-1261 23cm straight




Farabeuf Bone Holding Forceps
10-1262 17.8cm
10-1263 20cm
10-1264 23cm angled on flat



Farabeuf Bone Holding Forceps
with adjustable jaw
10-1265 23cm
10-1266 26cm



Farabeuf-Lambotte Bone Holding Forceps
with adjustable jaw
10-1267 26cm



Farabeuf-Lambotte Bone Holding Forceps
with adjustable jaw
10-1268 26cm



Bishop Bone Clamp
10-1269 22.5cm
Straight



Kern Bone Holding Forceps
10-1270 17cm
10-1271 24cm without ratchet



Kern Bone Holding Forceps
10-1272 17cm
10-1273 24cm with ratchet



Hey Groves Bone Holding Forceps
10-1274 21cm
10-1275 27cm
10-1276 30cm



FIMA INSTRUMENTS

BONE HOLDING FORCEPSKnochenhaltezaugen
Daviers porte-soPinzas para huesos
Tenaglie per ossa

10-1277
33cm 11 mm tip W/ Ratchet
10-1278
33cm 13 mm tip W/ Ratchet
10-1279
43cm 13.5 mm tip W/ Ratchet
10-1280
13cm 10 mm tip W/o Ratchet


Lane Bone Holding Forceps
Light Model




Lane Bone Holding Forceps
10-1281 33cm 13 mm tip
10-1282 43cm 12 mm tip
without ratchet



Ulrich Bone Holding Self Retaining Forceps



10-1283 18 cm mini
10-1284 22 cm
10-1285 26 cm
10-1286 28 cm
Straight with speed lock



10-1287 18cm
10-1288 22cm
10-1289 26cm
10-1290 28cm
Angular, with
speed lock

Ulrich Bone Holding Forceps



10-1291
28cm

Ulrich Bone Holding Forceps
Angular With Speed Lock
& swivel jaw



10-1292
22cm
10-1293
27cm
10-1294
30cm

Lambotte Bone Holding Forceps
One Jaw Mobile
adjustable jaw



Lane Bone Holding Forceps
10-1295 31.5cm
10-1296 39cm



10-1297
15 cm, 2.5mm
10-1298
19 cm, 6.5mm
10-1299
24 cm, 9.5mm
10-1300
26 cm, 10mm
10-1301
28 cm, 11mm

Lambotte Bone Holding Forceps
Self-Centering, with speed lock



Subacromion Spreader
10-1302 26.75cm Right
10-1303 26.75cm Left



Bone Reduction Forceps
10-1304 14cm
Serrated with pointed tips

BONE HOLDING FORCEPS



Knochenhaltezaugen
Daviers porte-so

Pinzas para huesos
Tenaglie per ossa



Bone Reduction Forceps
10-1305 14cm Curved
10-1306 15cm Serrated



10-1307
14cm Serrated
10-1308
14cm Curved
10-1309
17cm Curved

Bone Reduction Forceps,
Pointed Tips




Semb Bone Holding Forceps
Curved on side
10-1310 19cm 4 x 5 teeth
10-1311 21cm 5 x 5 teeth



Martin Cartilage Bone Holding
Clamp Forceps
10-1312 19 cm Str.



10-1313
20cm
Str.




10-1314
20cm
Cvd.

Bircher-Ganske
Bone Holding Forceps




10-1315
17.5cm
5.0mm


Verbrugge Bone
Holding Forceps




10-1316
25cm
9.5mm



10-1317
26cm
10.5mm




10-1318
27cm
12.0mm



10-1319
18cm
10-1320
26cm
10-1321
27cm
10-1322
28cm

Verbrugge Bone
Holding Forceps
W/ Long Ratchet



10-1323
17.2cm

Verbrugge Bone Holding
Forceps
with short ratchet



10-1324 15 cm
10-1325 19 cm
10-1326 24 cm
10-1327 26 cm
10-1328 29 cm

Verbrugge Self-Centering
Bone Holding Forceps
with long ratchet



Lane Bone Plate Holding
Forceps
10-1329 22cm

BONE HOLDING FORCEPS

Knochenhaltezaugen
Daviers porte-so

Pinzas para huesos
Tenaglie per ossa



Verbrugge Bone Holding
Forceps
10-1330 18cm




Verbrugge Bone Holding
Forceps
10-1331 18cm



10-1332
25cm
10-1333
26cm
10-1334
27cm

Verbrugge Bone Holding
Forceps



10-1335
25cm
10-1336
26cm
10-1337
27cm

Verbrugge Bone Holding
Forceps



Verbrugge Bone Holding
Forceps
10-1338 18cm



10-1339
25cm
10-1340
26cm
10-1341
27cm

Verbrugge Bone Holding
Forceps



Bone Reduction Forceps
10-1342 13.5cm



Reill Bone Reduction Forceps
10-1343 17cm
Curved, pointed tips



Bone Reduction Forceps
10-1344 20.4cm
Pointed tips



Bone Reduction Forceps
10-1345 13.5cm
for maxillary reposition



Haase
10-1346 16.5cm
for mandibular fractures



Bone Holding forceps



10-1347
9.5cm



10-1348
9.5cm

BONE HOLDING FORCEPS



Knochenhaltezaugen
Daviers porte-os

Pinzas para huesos
Tenaglie per ossa



Fibula Bone Holding Forceps
10-1349 22cm



Glenoid Perforating Forceps
10-1350 17cm
Lightly curved, pointed tips
Strongly curved, pointed tips



Reil Bone Holding Forceps
10-1351 17cm



Bone Reduction Forceps
10-1352 20cm
with Points



Pelvic Reduction Forceps Asymmetric
10-1353 40cm



Malleolar Bone Holding Forceps
10-1354 21cm
Straight, 1 x 1 teeth



Patella Malleolar Forceps
10-1355 21cm



Reduction Forceps with Thread Lock

10-1356 16cm
10-1357 24cm



Reposition Forceps
10-1358 23cm



Pelvic Reposition Forceps
10-1359 19cm
10-1360 23cm



Pelvic Reposition Forceps
10-1361 25cm



Pelvic Reposition Forceps
10-1362 25cm
3 Balls Spikes



Pelvic Reduction Forceps Asymmetric with 3 ball Spikes
10-1363 25cm



Pelvic Reduction Forceps Asymmetric
10-1364 25cm



Pelvic Reposition Forceps
10-1365 40cm



Pelvic Reposition Forceps
10-1366 40cm

HAND DRILLS

Handbohrapparate
Perforateurs a main

Perforadores de mano
Perforatori a mano



Bunnell Drill
10-1367 16cm
Opening up to 4mm
10-1368 spare key
light metal housing
longitudinally cannulated



Ralk Drill
10-1369 16cm
Opening up to 6.35mm
light metal housing
longitudinally cannulated



10-1372



Plaster Saw Electronic Power
Volt #110
10-1370



10-1373
Spare Blade

Volt #220
10-1371



Moore Drill And Bit Set Complete
10-1374



Fig. 1 10-1376 2.5mm
Fig. 2 10-1377 3.0mm
Fig. 3 10-1378 3.5mm
Fig. 4 10-1379 4.0mm
Fig. 5 10-1380 5.0mm
Fig. 6 10-1381 8.0mm
Fig. 7 10-1382 12.0mm
Fig. 8 10-1383 16.0mm
Fig. 09 10-1384
Fig. 10 10-1385 1.0mm
Fig. 11 10-1386 1.5mm
Fig. 12 10-1387 2.0mm

Stille Pattern Drill Set
10-1375
Complete set of stainless
drill bits and burrs

HAND DRILLS



Handbohrapparate
Perforateurs a main

Perforadores de mano
Perforatori a mano



10-1388
Hudson Brace 26.7cm
with standard thumb screw chuck

- C 10-1390 Hudson Bur, 9mm Diameter
- D 10-1391 Hudson Bur, 14mm Diameter
- E 10-1392 Hudson Bur, 16mm Diameter
- F 10-1393 Hudson Bur, 22mm Diameter
- G 10-1394 Cushing Perforator, for Hudson Brace
- H 10-1395 Hudson Cerabellar Extension, 4" for Hudson Brace
- I 10-1396 McKenzie Perforator for Hudson Brace

10-1389
Hudson Brace 27.9cm
with Improved Snap-Lock chuck



Jacobs
10-1397 28cm
opening up to 6.35mm



Hand Drill
10-1398
opening up
to 6.35mm



T-Handle
10-1399
AO Coupling



Jacob's Chuck
10-1400 T-Handle and key
10-1401 only
10-1402 Key




Key Handle
10-1403



T-Handle
10-1404
for AO Coupling



Memoral Head Extractor
10-1405



HEXAGONAL HEAD

	Size of screws	Head Size	Fiber-Handle	AO Coupling	Dental Coupling
A	1.5mm	1.5mm	10-1406	10-1409	10-1412
B	2.0mm				
C	2.7mm	2.5mm	10-1407	10-1410	
D	3.5mm				
E	4.0mm				
F	4.5mm	3.5mm	10-1408	10-1411	
G	6.5mm				



Parallel
Target
Device
10-1413
for guide
wires 2mm



Mortising
Chisel
10-1414



Wire Guide
10-1415
45mm




Wire Guide
10-1416
70mm



Countersink
10-1417



Countersink
10-1418



Countersink
10-1419



Perthes
10-1420 18cm
10-1421 21cm



10-1422
14cm



10-1423
14cm



Mini Drill Sleeve
10-1424
1.1 mm and 1.5 mm drill bits
10-1425
1.5 mm and 2.0 mm drill bits



Drill Guide
10-1426
drill sleeve
for 2.0 mm
drill bit

DRILL SLEEVES



10-1427 Tap Sleeve, 2.0 mm $\emptyset$
10-1428 Drill Sleeve, 2.7 mm $\emptyset$
10-1429 Tap Sleeve, 2.5 mm $\emptyset$
10-1430 Drill Sleeve, 3.5 mm $\emptyset$



Neutral & Load Drill Guide
10-1431 1.5 mm for 2.0 mm Screws
10-1432 2.0 mm for 2.7 mm Screws
10-1433 2.5 mm for 3.5 mm Screws
10-1434 3.2 mm for 4.5 mm Screws



Tap Sleeve
10-1435
 3.5 mm



Tap Sleeve
10-1436
 3.5 mm also for used
 as 3.2 mm drill Guide



Tap Sleeve
10-1437
 4.5 mm



Drill Sleeve
10-1438 3.2 mm
 for Round Hole Plates



Drill Guide
10-1439 4.5 mm Screws
 Drill Bit 4.5 mm & 3.2 mm



Drill Guide
10-1440 6.5 mm Screws
 Drill Bit 3.2 mm



Drill Guide
10-1441 3.5 mm for neutral and
 compression Positioning Drill Bit 2.5 mm



C - Clamp
10-1442
 Drill Guide with thread
 for Drill Bit 4.5 mm



Insert Drill Sleeve
10-1443 2.7 / 2.0 mm Drill Bit 2.0 mm
10-1444 3.5 / 2.5 mm Drill Bit 2.5 mm
10-1445 4.5 / 3.2 mm Drill Bit 3.2 mm



Drill Small Pointed
10-1446 Drill Bit 3.5 mm



Drill Pointed Guide
10-1447 6.5 mm



Drill Pointed Guide
10-1448 6.5 mm

PLATE BENDING DEVICES

Plattenbiegevorrichtungen

Dispositivos de doblado de placas

Dispositifs de cintrage de plaques

Dispositivi di piegatura della piastra



**Bending Device
10-1449**
21.5cm
to bend straight
and angled plates



**Bending Device
10-1450**
21.5cm
to bend straight
and angled plates



**Rod Bender
10-1451**
to bend straight
and angled



**Universal Rod Bender
10-1452** 28.5cm
to bend straight
and angled



**Rod Bender
10-1453**
to bend straight
and angled



**Slotted Hammer
10-1454** 14mm



**Plate Bender
10-1455**



**Plate Bender
10-1456**



**Plate Bender
10-1457**



**Plate Bender
10-1458**



**Plate Bender
10-1459**



**Plate Bender
10-1460**



**Plate Bender
10-1461**

EXTENSION BOWS



Extensionsbugel
Etriers pour extension

Estribos para extension
Archi per estensione



10-1462 18cm
Wire strainer for
Kirschner wires




Kirschner Extension Bows
10-1463 x=95mm y=70mm
10-1464 x=100mm y=100mm
10-1465 x=120mm y=120mm
10-1466 x=155mm y=155mm
10-1467 x=200mm y=155mm
with 3 traction hooks




Kirschner Extension Bows
10-1468 x=130mm
10-1469 x=185mm
10-1470 x=220mm
with 3 traction hooks



Wrench
Kirschner
Extension Bows
10-1471
14mm/17mm



Kirschner Extension Bows
10-1472 x=180mm y=100mm
10-1473 x=270mm y=160mm
with 3 traction hooks



Bohler Kirschner Extension Bows
10-1474 x=75mm y=90mm
10-1475 x=90mm y=160mm
10-1476 x=110mm y=210mm
10-1477 x=150mm y=210mm
with 3 traction hooks



Charnley-mueller
10-1478
Compression Clamp



Screw Depth Gauges

- 10-1479** Mini Gauge, Longmeasures up to 30mm
- 10-1480** Small Gauge, Longmeasures up to 50mm
- 10-1481** Large Gauge, Longmeasures up to 100mm
- 10-1482** X-Large Gauge, Longmeasures up to 120mm




Screw Depth Gauges

- 10-1483** Mini Gauge, Longmeasures up to 30mm
- 10-1484** Small Gauge, Longmeasures up to 50mm
- 10-1485** Large Gauge, Longmeasures up to 100mm
- 10-1486** X-Large Gauge, Longmeasures up to 120mm




**Measuring Gauge
10-1487
for nail**

10-1488 Tri Angular Positioning Plate




NEUROSURGERY INSTRUMENTS



HAND DRILLS

Handbohrapparate
Perforateurs a main

Perforadores a mano
Perforatori a mano



Hubson, Brace Hand Drill
11-004 27cm
complete with 4 twist drills
(Fig. 1 - Fig. 4) and extension piece

Hudson Brace
11-005
28cm with Improved
Snap-Lock chuck

FIMA INSTRUMENTS



Huson Hand
drill bit
Fig.1
11-006
9mm

Huson Hand
drill bit
Fig.2
11-007
14mm

Huson Hand
drill bit
Fig.3
11-008
16mm

Huson Hand
drill bit
Fig.4
11-009
22mm

Cushing
Hand drill bit
11-010
14mm

Hudson
Cerabellar
Extension
10cm
11-011
Hudson Brace

McKenzie
Hand drill bit
11-012
13mm

CERVICAL TRACTION TONGS




Schadelhalter


Pinzas de traccion cranial

Pinces de traction Cranienne

Pinze per trazione craniale



**Crutchfield
Traction Tongs
11-013 13cm
small**



Crutchfield drill bit 3.5mm
Crutchfield drill bit 5mm
Set of 2 Pieces for
11-013 and 11-014



**Crutchfield
Traction Tongs
11-014 18.5cm**



11-015
Wrench double-ended for



**Crutchfield Skull
11-016**
large type, cervical traction tongs -
pins adjustable no drilling to be done
11cm



**Cloward Vertebral Cervical
Spreader
11-017**



**Spinal Distractor
11-018 24cm**
45° angle / small style
6mm rods



**Spinal Distractor
11-019 25.5cm**
small style



**Inge Retractor
11-020 20.5cm**
toothed w/spindle

RETRACTORS

Wundhaken | Ecarteurs
Separadores | Divaricatori



Supra Cut Distraction Screw

- 11-021 12mm Non Sterile
- 11-022 14mm Non Sterile
- 11-023 16mm Non Sterile
- 11-024 18mm Non Sterile
- 11-025 12mm Sterile
- 11-026 14mm Sterile
- 11-027 16mm Sterile
- 11-028 18mm Sterile



Supra Cut
Distraction
Screw Driver
11-029



Cloward Retractor
11-030 15.5 cm




Cloward Retractor
11-031 13.5 cm



Inge
11-032
Lamina Spreader
16 cm



Inge Retractor
11-033 17 cm



Inge Retractor
11-034 27.5 cm

INSTRUMENTS FOR CRANIAL OPENING



Instrumente zur Schädelöffnung

Instrumentos para apertura craneana

Instruments pour l'ouverture de crâne

Strumenti per l'apertura del cranio



De Martel
Saw Conductor
11-035 33cm



Adson
Drill guide
11-036
15cm




De Vilbiss
Skull Forceps
11-037
20cm



Dahlgren
11-038
Cranial Rongeur
20 cm



Gigle-Strully
Handle for
wire saw
11-039



Olivecrona
Wire Saw
11-040 30cm
11-041 40cm
11-042 50cm
11-043 60cm
11-044 70cm



Gigle, Wire Saw
11-045 30cm
11-046 40cm
11-047 50cm
6 wires twisted



Bailey-Gigli
11-048
Saw Guide 30.5cm
flexible



Davis
11-049
Guide for Gigli Saws,
20cm



Michele Trephine,
11-050 16cm shaft, 0.47cm
O.D graduated 0.64cm To 3.2cm
11-051 16cm shaft, 0.79cm
O.D graduated 0.64cm To 3.2cm
11-052 16cm shaft, 0.56cm
O.D graduated 0.64cm To 3.2cm

11-050

11-051

11-052


SELF RETAINING BRAIN RETRACTOR

selbsthaltender Hirnretractor


Retracteur cerebrale autostatique

Separador cerebral autostatico


Retratore celebrale autostatico



11-053
wire cable only



11-054
Support only,
flat brain spatulas,



11-055
Support only,
for brain spatulas
with round shaft
up to 5.5 mm diam.
Fixing for flexible




Riechert
11-056
Fixation base, for skull
mounting, holding 1
flexible arms




Riechert
11-057
Fixation base, for skull
mounting, holding 2
flexible arms



11-058
Fixation base, for skull
mounting, holding 1
flexible arms



Leyla
Self-Retaining Brain Retractor
with 1 flexible arm, complete consisting of
11-059
11-060
11-061



Leyla
Self-Retaining Brain
Retractor
with 2 flexible arms, complete
consisting of
11-062
11-063
11-064



11-065
Fixation base only
holding 2 flexible arms
11-066
Support only for flat
brain spatulas

11-067
Fixation base only
holding 1 flexible arms

11-068
Support only for brain spatulas with round shaft of
4mm and 5mm to fix on the flexible arm

11-069
Flexible arm only for attaching
coupling heads and for fixation bases



11-070
Brain retractor self retaining, with 1 flexible arm,
complete consisting of fixation base for 1 flexible
arm flexible arm only support for flat brain spatulas


11-071
Brain retractor self retaining, with 2 flexible arm,
complete consisting of fixation base for 2 flexible
arm flexible arm only support for flat brain spatulas

McKENZIE INSTRUMENTS



McKenzie Instrumentarium
McKenzie instruments

Instrumentos McKenzie
Strumenti McKenzie



McKenzie Clip Applying Forceps

11-072
15cm
11-073
19cm
Straight

11-074
15cm
11-075
19cm
Curved



McKenzie 11-076
148mm x 40mm x 6mm
Brain clip rack for 20 clips



McKenzie Clip Applying Forceps 11-077
Straight Handle, Curved Jaws, 15cm



McKenzie Clip Applying Forceps 11-078
Curved Handle, Curved Jaws, 15cm



Smithwick Clip Applying Forceps 11-079 23 cm
11-080 26.5 cm
Curved Jaws
Straight Handle



Schwartz Clip Applying Forceps 11-081 20cm
Schwartz Vessel Clip



11-082
2.5cm Straight
Jaws 1.7x 8mm




11-083
2.5cm Slightly
Angled Jaws
1.7x 8mm



11-084
2.5cm Strong
Angled Jaws
1.7x 8mm




Olivecrona Toennis (Scoville), Clip Forceps 11-085 14cm



11-087
Raney Hemostasis Clip



11-088
Kolner Scalp Clip



Raney Applyig & Removing Forceps 11-086 15cm


CLIPS & APPLYING FORCEPS

clips und anlegezangen

Clips y aplicadores

Clips et pinces a appliquer

Clips e pinze per mettere clips



Example
For Rotation

Mini

11-089 7cm
11-090 9cm
11-091 11cm
11-092 13cm

non rotating

Standard

11-093 7cm
11-094 9cm
11-095 11cm
11-096 13cm

non rotating

Mini

11-097 7cm
11-098 9cm
11-099 11cm
11-100 13cm

rotating

Standard

11-101 7cm
11-102 9cm
11-103 11cm
11-104 13cm

rotating



Mini
11-105 18cm
Straight with ratchet
Standard
11-106 18cm
Straight with ratchet



Mini
11-107 17cm
11-108 21cm
11-109 23cm
bayonett shaped, with ratchet
Standard
11-110 17cm
11-111 21cm
11-112 23cm
bayonett shaped, with ratchet

Applying Forceps



Cushing, brain spatula
11-113 7mm x 9mm
11-114 11mm x 13mm
11-115 15mm x 18mm
11-116 18mm x 22mm
18cm



Olivecrona
11-117
brain spatula
18cm



Brain Spatula
11-118 7mm x 9mm
11-119 11mm x 13mm
11-120 15mm x 18mm
11-121 18mm x 22mm
18cm Convex



Brain Spatula
11-122 7mm x 9mm
11-123 11mm x 13mm
11-124 15mm x 18mm
11-125 18mm x 22mm
18cm Concave



Heifetz, brain spatula
11-126 8mm
11-127 11mm
11-128 14mm
11-129 17mm
11-130 20mm
20cm Malleable

BRAIN SPATULAS



Hirnspatel

Espatulas cerebrales

Spatules cerebrales

Spatole cerebrali

NERVE HOOKS



**Scoville, Nerve
Root Retractor**
11-131 23cm
8mm



**Holsher Nerve
Root Retractor**
11-132 24cm
7mm



**Martin Nerve
Root Retractor**
11-133 21cm
7mm



**Toennis Nerve
Dura Hook**




**Halle Nerve
Dura Hook**
11-136 20cm



**Smithwick Nerve
Dura Hook**
11-137 17cm
11-138 21cm



**Hoen Nerve
Dura Hook**
11-139 22cm



**Cushing Nerve
Dura Hook**



Olivecrona
11-144
Dura Dissector
18 cm



Olivecrona
11-145
Dura Dissector
20 cm



Milligan
11-146
Dura Dissector
21.5 cm



Davis
11-147
Dura Dissector sharp
21.5 cm


NERVE HOOKS - DURA HOOKS

Nervenhooken - Durahaken

Crochets a nerfs - Crochets pour duremers

Ganchos para nervios - Ganchos para duramadre

Uncini per nervi - Uncini per duramadre



**Graham Nerve
Dura Hook**
11-148
18.5cm



**Frazier Nerve
Dura Hook**
11-149
17cm



**Dandy Nerve
Dura Hook**
11-150
20cm



**Dandy Nerve
Dura Hook**
11-151
20cm



**Dandy Nerve
Dura Hook**
11-152
20cm



**Adson Nerve
Dura Hook**



**Yasargil Nerve Dura
Hook Retractor**
11-155 21cm
11-156 41cm
Galea retractor




**Love Nerve Dura
Hook Retractor**
11-157 19cm
Straight



**Love Nerve Dura
Hook Retractor**
11-158 19cm
45° Angled



**Love Nerve Dura
Hook Retractor**
11-159
19cm
90° Angled



**Cushing Nerve Dura
Hook Retractor**
11-160 14cm



**Frazier Nerve Dura
Hook Retractor**



11-161
13cm
11-162
18cm
Blunt



**Lahey Nerve Dura
Hook Retractor**
11-163
13cm

DURA DISSECTORS



Dura-Dissektoren

Disectores de duramadre

Decolleurs duremere

Dissetori per duramadre



Toennis Dura Dissector
11-164
25cm



Davis Dura Dissector
11-165
24.5cm



Dura Dissector
11-166 22cm




Sachs Dura Dissector
11-167 21cm




Sachs Dura Dissector
11-168 21cm



Sachs Dura Dissector
11-169
21cm



Horsley Dura Dissector
11-170
20cm



Smithwich Dura Dissector
11-171
30.5cm



Groningen
11-172
Nerve Hook
26.5 cm



Groningen
11-173
Nerve Hook
26.5 cm



11-174
A




11-175
B

Cushing
Nerve Hook blunt
19 cm


DURA DISSECTORS

Dura-Dissektoren
Decolleurs duremere

Disectores de duramadre
Dissetori per duramadre



NEURO-DISSECTORS AND CURETTES



NEURO-MICRO DISSECTORS



Neuro-Mikro-Dissektoren
Neuro-micro dissecteurs

Disectores para micro-neurocirugia
Strumento per micro-neuro chirurgia



11-191
19cm
Straight
11-192
19cm
bayonet



11-193
19cm
Straight
11-194
19cm
bayonet



11-195
19cm
Straight
11-196
19cm
bayonet



Kraus
11-211
19cm
bayonet



Yasargil
11-212
19cm
bayonet



Koos
11-213
19cm
bayonet



11-197
19cm
Straight
11-198
19cm
bayonet



11-199
19cm
Straight
11-200
19cm
bayonet



Koos
11-214
19cm
bayonet



Koos
11-215
19cm
bayonet



11-201
19cm
Straight
11-202
19cm
bayonet



11-203
19cm
Straight
11-204
19cm
bayonet



11-205
19cm
Straight
11-206
19cm
bayonet



11-207
19cm
Straight
11-208
19cm
bayonet



11-209
19cm
Straight, with
Twist-fixationhook



11-210
19cm
Straight, with
Twist-fixationhook



**Neuro-Micro
Dissector**

**Yasargil Neuro-Micro
Dissector**

LAMINECTOMY PUNCHES

Laminektomiestanzen

Pinces emporte-pieces a laminectomie

Pinzas cortantes para laminectomia

Pinze taglianti per laminectomia

Cushing Intervertebral Disc Rongeurs

2 X 10mm, 15cm
11-216 straight
11-217 UP
11-218 Down

1.5 X 10mm, 18cm
11-219 straight
11-220 Up
11-221 Down

2 X 10mm, 18cm
11-222 straight
11-223 UP
11-224 Down

1.5 X 10mm, 23cm
11-225 straight
11-226 Up
11-227 Down

2 X 10mm, 23cm
11-228 straight
11-229 UP
11-230 Down

2 X 10mm, 25.5cm
11-231 straight
11-232 UP
11-233 Down

2 X 10mm, 28cm
11-234 straight
11-235 UP
11-236 Down

2 X 10mm, 30.5cm
11-237 straight
11-238 UP
11-239 Down

Love-Gruenwald Intervertebral Disc Rongeurs

3 X 10mm, 15cm
11-240 straight
11-241 Up
11-242 Down

3 X 10mm, 18cm
11-243 straight
11-244 Up
11-245 Down

3 X 10mm, 23cm
11-246 straight
11-247 Up
11-248 Down

3 X 10mm, 25.5cm
11-249 straight
11-250 Up
11-251 Down

3 X 10mm, 28cm
11-252 straight
11-253 Up
11-254 Down

3 X 10mm, 30.5cm
11-255 straight
11-256 Up
11-257 Down

Spurling Intervertebral Disc Rongeurs

4 X 10mm, 15cm
11-258 straight
11-259 Up
11-260 Down

4 X 10mm, 18cm
11-261 straight
11-262 Up
11-263 Down

4 X 10mm, 23cm
11-264 straight
11-265 Up
11-266 Down

4 X 10mm, 25.5cm
11-267 straight
11-268 Up
11-269 Down

4 X 10mm, 28cm
11-270 straight
11-271 Up
11-272 Down


4 X 10mm, 30.5cm
11-273 straight
11-274 Up
11-275 Down

FIMA INSTRUMENTS



Cushing

Laminectomy Rongeur jaw 2 x 10 to 4 x 10
mm straight shaft length 180 mm



Cushing

Laminectomy Rongeur jaw 2 x 10 to 4 x 10
mm upwards 45 degree

LAMINECTOMY PUNCHES



Laminektomiestanzen

Pinzas cortantes para laminectomia

Pinces emporte-pieces a laminectomie

Pinze taglianti per laminectomia

Wilde Intervertebral Disc Rongeurs

18cm Straight
11-286 2.5mm
11-287 3mm
11-288 4mm
11-289 5mm

18cm Angled up
11-290 2.5mm
11-291 3mm
11-292 4mm
11-293 5mm

23.5cm Straight
11-294 3mm
11-295 4mm
11-296 5mm

Caspar Type Intervertebral Disc Rongeurs

16cm Straight
11-297 2 x 12mm
11-298 3 x 12mm
11-299 4 x 14mm
11-300 5 x 14mm
11-301 6 x 16mm

16cm Up
11-302 3 x 12mm
11-303 4 x 14mm
11-304 5 x 14mm

16cm Down
11-305 3 x 12mm
11-306 4 x 14mm

Caspar Type Intervertebral Disc Rongeurs

18.5cm Straight
11-307 2 x 12mm
11-308 3 x 12mm
11-309 4 x 14mm
11-310 5 x 14mm
11-311 6 x 16mm

18.5cm Up
11-312 2 x 12mm
11-313 3 x 12mm
11-314 4 x 14mm
11-315 5 x 14mm

18.5cm Down
11-316 2 x 12mm
11-317 3 x 12mm
11-318 4 x 14mm

Cushing Intervertebral Disc Rongeurs

w/teeth 2 x 10mm, 18cm
11-319 Straight
11-320 Up
11-321 Down

L-G Intervertebral Disc Rongeurs

w/teeth 3 x 10mm, 18cm
11-322 Straight
11-323 Up
11-324 Down

3 x 10mm, 25.5cm
11-325 Straight

Schlesinger Intervertebral Disc Rongeurs

11-326 2 x 10mm, 13cm Straight
11-327 3 x 10mm, 13cm Straight
11-328 4 x 10mm, 13cm Straight
11-329 3 x 10mm, 13cm Angled Up

Cloward Intervertebral Disc Rongeurs

serrated 6 x 10mm, 13cm
11-330 Straight
11-331 Up
11-332 Down

Cushing Intervertebral Disc Rongeurs

2 X 10mm, 13cm
11-336 Straight
11-337 UP
11-338 Down

Love-Gruenwald Intervertebral Disc Rongeurs

3 X 10mm, 13cm
11-339 Straight
11-340 Up
11-341 Down

Spurling Intervertebral Disc Rongeurs

4 X 10mm, 13cm
11-342 Straight
11-343 Up
11-344 Down

Oldberg Intervertebral Disc Rongeurs

11-333 6mm Diameter
18.5cm

Selverstone Intervertebral Disc Rongeurs

11-334 2 x 3mm, 15cm
11-335 2 x 5mm, 15cm

LAMINECTOMY PUNCHES

Laminektomiestanzen

Pinces emporte-pieces a laminectomie

Pinzas cortantes para laminectomia

Pinze taglianti per laminectomia



Decker Intervertebral Disc Rongeurs
2 x 5mm, 15cm
11-345 Right
11-346 Left
11-347 Down
11-348 Up
11-349 Straight



Jackson Intervertebral Disc Rongeurs
70° Up 18cm
11-350 2 x 8mm
11-351 3 x 8mm
11-352 4 x 8mm




Williams Style Micro Disc Intervertebral Disc Rongeurs
1 x 2 teeth, Straight
11-353 14cm
11-354 18cm
11-355 25.5cm



Micro Cup Forceps
14cm, 0.5mm
11-356 Straight
11-357 Right
11-358 Left



Micro Cup Forceps
14cm 1mm
11-359 Straight
11-360 Right
11-361 Left
14cm 2mm
11-362 Straight



Micro Cup Forceps
20.5cm, 2mm
11-363 Straight
11-364 Up
11-365 Down



Well Blakesley Laminectomy rongeurs



12cm
Straight

Fig. 1
3.0mm
11-366

Fig. 2
3.5mm
11-367

Fig. 3
4.5mm
11-368

Fig. 4
5.0mm
11-369

Fig. 5
5.5mm
11-370



12cm
45° up

11-371

11-372

11-373

11-374

11-375



12cm
90° up

11-376

11-377

11-378



18cm
45° up

Fig. 1
3.0mm
11-379

Fig. 2
3.5mm
11-380

Fig. 3
4.5mm
11-381

Fig. 4
5.0mm
11-382

Fig. 5
5.5mm
11-383



18cm 7"
90° up

11-384

11-385

11-386

11-387

LAMINECTOMY PUNCHES



FMA INSTRUMENTS

Laminektomiestanzen

Pinzas cortantes para laminectomia

Pinces emporte-pieces a laminectomie

Pinze taglianti per laminectomia


Ferris Smith Micro Kerrison Laminectomy Rongeurs

45° UP 9mm bite opening 18cm

- 11-388 1mm
- 11-389 2mm
- 11-390 3mm
- 11-391 4mm
- 11-392 5mm
- 11-393 6mm

90° UP 9mm bite opening 20.5cm

- 11-394 1mm
- 11-395 2mm
- 11-396 3mm
- 11-397 4mm
- 11-398 5mm
- 11-399 6mm



Ferris Smith Micro Kerrison Laminectomy Rongeurs

45° Down 9mm bite opening 18cm

- 11-400 1mm
- 11-401 2mm
- 11-402 3mm
- 11-403 4mm
- 11-404 5mm
- 11-405 6mm

90° Down 9mm bite opening, 20.5cm

- 11-406 1mm
- 11-407 2mm
- 11-408 3mm
- 11-409 4mm
- 11-410 5mm
- 11-411 6mm



Ferris Smith Micro Kerrison Laminectomy Rongeurs UPWARD

40° UP

9mm bite opening, 18cm

- 11-412 1mm
- 11-413 2mm
- 11-414 3mm
- 11-415 4mm
- 11-416 5mm
- 11-417 6mm

40° UP

9mm bite opening, 23cm

- 11-435 1mm
- 11-436 2mm
- 11-437 3mm
- 11-438 4mm
- 11-439 5mm
- 11-440 6mm

40° Down

9mm bite opening, 18cm

- 11-418 2mm
- 11-419 3mm
- 11-420 4mm
- 11-421 5mm
- 11-422 6mm

40° UP

9mm bite opening, 25.5cm

- 11-441 2mm
- 11-442 3mm
- 11-443 4mm
- 11-444 5mm
- 11-445 6mm

40° UP

9mm bite opening, 20.5cm

- 11-423 1mm
- 11-424 2mm
- 11-425 3mm
- 11-426 4mm
- 11-427 5mm
- 11-428 6mm

40° UP

9mm bite opening,

- 11-446 2mm
- 11-447 3mm
- 11-448 4mm
- 11-449 5mm
- 11-450 6mm

40° Down

9mm bite opening, 20.5cm

- 11-429 1mm
- 11-430 2mm
- 11-431 3mm
- 11-432 4mm
- 11-433 5mm
- 11-434 6mm

40° UP

9mm bite opening, 30.5cm

- 11-451 2mm
- 11-452 3mm
- 11-453 4mm
- 11-454 5mm
- 11-455 6mm



Ferris Smith Micro Kerrison Laminectomy Rongeurs DOWNWARD

130° Down

9mm bite opening, 18cm

- 11-456 1mm
- 11-457 2mm
- 11-458 3mm
- 11-459 4mm
- 11-460 5mm
- 11-461 6mm

130° Down

9mm bite opening, 28cm

- 11-480 1mm
- 11-481 2mm
- 11-482 3mm
- 11-483 4mm
- 11-484 5mm
- 11-485 6mm

130° Down

9mm bite opening, 20.5cm

- 11-462 1mm
- 11-463 2mm
- 11-464 3mm
- 11-465 4mm
- 11-466 5mm
- 11-467 6mm

130° Down

9mm bite opening, 30.5cm

- 11-486 1mm
- 11-487 2mm
- 11-488 3mm
- 11-489 4mm
- 11-490 5mm
- 11-491 6mm

130° Down

9mm bite opening, 23cm

- 11-468 1mm
- 11-469 2mm
- 11-470 3mm
- 11-471 4mm
- 11-472 5mm
- 11-473 6mm

130° Down

9mm bite opening, 25.5cm

- 11-474 1mm
- 11-475 2mm
- 11-476 3mm
- 11-477 4mm
- 11-478 5mm
- 11-479 6mm


LAMINECTOMY PUNCHES

Laminektomiestanzen

Pinces emporte-pieces a laminectomie


Pinzas cortantes para laminectomia

Pinze taglianti per laminectomia



**Micro Ejector Kerrison
Laminectomy Rongeurs**

40° Up 20.5cm
11-492 2mm
11-493 3mm
11-494 4mm
11-495 5mm
11-496 6mm



**Wide Mouth Ejector Kerrison
Laminectomy Rongeurs**

40° Up 18cm
11-497 2mm
11-498 3mm
11-499 4mm
11-500 5mm
11-501 6mm

40° Up 20.5cm
11-502 2mm
11-503 3mm
11-504 4mm
11-505 5mm
11-506 6mm



**Spurling Kerrison
Laminectomy Rongeurs
Large Handle**

40° UP,
15mm Bite Opening, 18cm
11-507 2mm
11-508 3mm
11-509 4mm
11-510 5mm
11-511 6mm



**Spurling Kerrison
Laminectomy Rongeurs
Large Handle**

40° UP,
15mm Bite Opening, 25.5cm
11-522 2mm
11-523 3mm
11-524 4mm
11-525 5mm
11-526 6mm



**Spurling Kerrison
Laminectomy Rongeurs
Large Handle**

90° Down,
15mm Bite Opening, 18cm
11-542 2mm
11-543 3mm
11-544 4mm
11-545 5mm
11-546 6mm

40° UP,
15mm Bite Opening, 20.5cm
11-512 2mm
11-513 3mm
11-514 4mm
11-515 5mm
11-516 6mm

40° UP,
15mm Bite Opening, 28cm
11-527 2mm
11-528 3mm
11-529 4mm
11-530 5mm
11-531 6mm

90° UP,
15mm Bite Opening, 20.5cm
11-532 2mm
11-533 3mm
11-534 4mm
11-535 5mm
11-536 6mm

90° Down,
15mm Bite Opening, 20.5cm
11-547 2mm
11-548 3mm
11-549 4mm
11-550 5mm
11-551 6mm

40° UP,
15mm Bite Opening, 23cm
11-517 2mm
11-518 3mm
11-519 4mm
11-520 5mm
11-521 6mm

LAMINECTOMY PUNCHES



Laminektomiestanzen

Pinzas cortantes para laminectomia

Pinces emporte-pieces a laminectomie


Pinze taglianti per laminectomia



Foramenotomy Kerrison

Laminectomy Rongeurs, Curved

- 11-552 UP, 9mm Bite opening, 20.5cm 2mm
- 11-553 UP, 9mm Bite opening, 20.5cm 3mm
- 11-554 UP, 9mm Bite opening, 20.5cm 4mm
- 11-555 UP, 9mm Bite opening, 20.5cm 5mm
- 11-556 UP, 9mm Bite opening, 20.5cm 6mm
- 11-557 Left, 9mm Bite opening, 20.5cm 6mm
- 11-558 Right, 9mm Bite opening, 20.5cm 6mm



Bayonet Micro Kerrison

Laminectomy Rongeurs

- 40° Up, 150 mm Working Length 20.5cm
- 11-559 1mm
- 11-560 2mm
- 11-561 3mm
- 11-562 4mm
- 11-563 5mm
- 11-564 6mm



Spurling Laminectomy Rongeurs

40° UP, 15mm Bite Opening, 18cm

- 11-583 2mm
- 11-584 3mm
- 11-585 4mm
- 11-586 5mm
- 11-587 6mm

40° UP, 15mm Bite Opening, 23cm

- 11-593 2mm
- 11-594 3mm
- 11-595 4mm
- 11-596 5mm
- 11-597 6mm

40° UP, 15mm Bite Opening, 20.5cm

- 11-588 2mm
- 11-589 3mm
- 11-590 4mm
- 11-591 5mm
- 11-592 6mm

40° UP, 15mm Bite Opening, 25.5cm

- 11-598 2mm
- 11-599 3mm
- 11-600 4mm
- 11-601 5mm
- 11-602 6mm

40° UP, 15mm Bite Opening, 28cm

- 11-603 2mm
- 11-604 3mm
- 11-605 4mm
- 11-606 5mm
- 11-607 6mm

90° Up, 150 mm Working Length 20.5cm

- 11-565 1mm
- 11-566 2mm
- 11-567 3mm
- 11-568 4mm
- 11-569 5mm
- 11-570 6mm

Bayonet Spurling Kerrison

Laminectomy Rongeurs

- 40° Up, 150 mm Working Length 20.5cm
- 11-571 1mm
- 11-572 2mm
- 11-573 3mm
- 11-574 4mm
- 11-575 5mm
- 11-576 6mm

90° Up, 150 mm Working Length 20.5cm

- 11-577 1mm
- 11-578 2mm
- 11-579 3mm
- 11-580 4mm
- 11-581 5mm
- 11-582 6mm

LAMINECTOMY PUNCHES

Laminektomiestanzen

Pinces emporte-pieces a laminectomie

Pinzas cortantes para laminectomia

Pinze taglianti per laminectomia

90° UP, 15mm Bite Opening, 18cm

11-608 2mm
11-609 3mm
11-610 4mm
11-611 5mm
11-612 6mm



Spurling Laminectomy Rongeurs

90° UP, 15mm Bite Opening, 20.5cm

11-613 2mm
11-614 3mm
11-615 4mm
11-616 5mm
11-617 6mm

90° Down, 15mm Bite Opening, 18cm

11-618 2mm
11-619 3mm
11-620 4mm
11-621 5mm
11-622 6mm



Spurling Laminectomy Rongeurs

90° Down, 15mm Bite Opening, 20.5cm

11-623 2mm
11-624 3mm
11-625 4mm
11-626 5mm
11-627 6mm

Micro Ejector Kerrison
Laminectomy Rongeurs

40° Up 20.5cm

11-628 2mm
11-629 3mm
11-630 4mm
11-631 5mm
11-632 6mm




Wide Mouth Ejector Kerrison
Laminectomy Rongeurs

40° Up 18cm

11-633 2mm
11-634 3mm
11-635 4mm
11-636 5mm
11-637 6mm

40° Up 20.5cm

11-638 2mm
11-639 3mm
11-640 4mm
11-641 5mm
11-642 6mm



Ferris Smith Micro Kerrison
Laminectomy Rongeurs

130° Down Bite, 28cm

11-643 4mm
11-644 5mm
11-645 6mm

130° Down Bite, 30.5cm

11-646 1mm
11-647 2mm
11-648 3mm
11-649 4mm
11-650 5mm
11-651 6mm



Kerrison Cervical Rongeur

11-657

17.8 cm shaft, 3 mm
40 degree forward angle

Spurling Intervertebral
Disc Rongeurs

w/teeth 4 x 10mm, 18cm

11-652 Straight
11-653 Up
11-654 Down

4 x 10mm, 25.5cm

11-655 Straight

5 x 10mm, 25.5cm

11-656 Straight



Schlesinger Cervical Rongeur

11-658

thin tips, 15 cm shaft, up-bite, 3 mm



LAMINECTOMY PUNCHES



Laminektomiestanzen

Pinzas cortantes para laminectomia

Pinces emporte-pieces a laminectomie

Pinze taglianti per laminectomia



Ferris-Smith Intervertebral Disc Rongeurs

18cm straight

11-659 2 x 10mm


11-660 3 x 10mm

11-661 4 x 8mm

11-662 4 x 10mm

11-663 6 x 10mm

11-664 6 x 12mm



Ferris-Smith Intervertebral Disc Rongeurs

18cm Up

11-665 2 x 10mm

11-666 3 x 10mm

11-667 4 x 10mm

11-668 6 x 10mm

18cm Down

11-669 2 x 10mm

11-670 3 x 10mm

11-671 4 x 10mm

11-672 6 x 10mm



Peapod Intervertebral

Disc Rongeurs, 14cm

11-673 Angle up 1.5mm, x-small

11-674 Angle up 2mm, small

11-675 Angle up 2.5mm, medium

11-676 Angle up 3mm, large



Peapod Intervertebral

Disc Rongeurs,

11-677 Angle up 2.5mm, 18cm

11-678 Angle up 3mm, 18cm

11-679 Angle up 4mm, 18cm

11-680 Straight 2mm, 14cm

11-681 Straight 2.5mm, 14cm

11-682 Straight 3mm, 14cm


LAMINECTOMY PUNCHES

Laminektomiestanzen

Pinzas cortantes para laminectomia

Pinces emporte-pieces a laminectomie

Pinze taglianti per laminectomia



Kerrison Rongeur

Up-cutting, 90° angled jaws,
with ejector

Ferris-Smith Kerrison handle	Kerrison Small handle	Kerrison Large handle	Jaw size (width x height) Length*	
				
		11-707	2,0 x 2,0 mm	18 cm
11-683	11-695	11-708	2,5 x 2,5 mm	18 cm
11-684	11-696	11-709	2,0 x 3,0 mm	18 cm
11-685	11-697	11-710	3,0 x 4,0 mm	18 cm
11-686	11-698	11-711	4,0 x 5,0 mm	18 cm
11-687	11-699	11-712	5,0 x 5,0 mm	18 cm
11-688	11-700	11-713	6,0 x 5,0 mm	18 cm
		11-714	7,0 x 6,0 mm	18 cm
		11-715	2,0 x 2,0 mm	20 cm
11-689	11-701	11-716	2,5 x 2,5 mm	20 cm
11-690	11-702		2,0 x 3,0 mm	20 cm
11-691	11-703	11-717	3,0 x 4,0 mm	20 cm
11-692	11-704	11-718	4,0 x 5,0 mm	20 cm
11-693	11-705	11-719	5,0 x 5,0 mm	20 cm
11-694	11-706	11-720	6,0 x 5,0 mm	20 cm
		11-721	7,0 x 6,0 mm	20 cm
Maulgrößen / Jaw sizes:  2,0 x 2,0 mm 2,5 x 2,5 mm 2,0 x 3,0 mm 3,0 x 4,0 mm 4,0 x 5,0 mm 5,0 x 5,0 mm 6,0 x 5,0 mm 7,0 x 6,0 mm		11-722	2,0 x 2,0 mm	22 cm
		11-723	2,5 x 2,5 mm	22 cm
		11-724	2,0 x 3,0 mm	22 cm
		11-725	3,0 x 4,0 mm	22 cm
		11-726	4,0 x 5,0 mm	22 cm
		11-727	5,0 x 5,0 mm	22 cm
		11-728	6,0 x 5,0 mm	22 cm
		11-729	7,0 x 6,0 mm	22 cm
		11-730	2,0 x 2,0 mm	25 cm
		11-731	2,0 x 3,0 mm	25 cm
		11-732	3,0 x 4,0 mm	25 cm
		11-733	4,0 x 5,0 mm	25 cm
		11-734	5,0 x 5,0 mm	25 cm
		11-735	2,0 x 2,0 mm	30 cm
* Lengths are measured with closed jaws		11-736	3,0 x 4,0 mm	30 cm
		11-737	4,0 x 5,0 mm	30 cm
		11-738	5,0 x 5,0 mm	30 cm

LAMINECTOMY PUNCHES




Laminektomiestanzen

Pinzas cortantes para laminectomia

Pinces emporte-pieces a laminectomie

Pinze taglienti per laminectomia



Kerrison Rongeur

Up-cutting, 50° angled jaws,
with ejector

Ferris-Smith Kerrison handle	Kerrison Small handle	Kerrison Large handle	Jaw size (width x height)	Length*
		11-763	2,0 x 2,0 mm	18 cm
11-739	11-751	11-764	2,5 x 2,5 mm	18 cm
11-740	11-752		2,0 x 3,0 mm	18 cm
11-741	11-753	11-765	3,0 x 4,0 mm	18 cm
11-742	11-754	11-766	4,0 x 5,0 mm	18 cm
11-743	11-755	11-767	5,0 x 5,0 mm	18 cm
11-744	11-756	11-768	6,0 x 5,0 mm	18 cm
		11-769	7,0 x 6,0 mm	18 cm
		11-770	2,0 x 2,0 mm	20 cm
11-745	11-757	11-771	2,5 x 2,5 mm	20 cm
11-746	11-758	11-772	2,0 x 3,0 mm	20 cm
11-747	11-759	11-773	3,0 x 4,0 mm	20 cm
11-748	11-760	11-774	4,0 x 5,0 mm	20 cm
11-749	11-761	11-775	5,0 x 5,0 mm	20 cm
11-750	11-762	11-776	6,0 x 5,0 mm	20 cm
		11-777	7,0 x 6,0 mm	20 cm
		11-778	2,0 x 2,0 mm	22 cm
		11-779	2,5 x 2,5 mm	22 cm
		11-780	3,0 x 4,0 mm	22 cm
		11-781	4,0 x 5,0 mm	22 cm
		11-782	5,0 x 5,0 mm	22 cm
		11-783	6,0 x 5,0 mm	22 cm
		11-784	7,0 x 6,0 mm	22 cm
		11-785	2,0 x 2,0 mm	25 cm
		11-786	3,0 x 4,0 mm	25 cm
		11-787	4,0 x 5,0 mm	25 cm
		11-788	5,0 x 5,0 mm	25 cm
		11-789	2,0 x 2,0 mm	30 cm
		11-790	3,0 x 4,0 mm	30 cm
		11-791	4,0 x 5,0 mm	30 cm
		11-792	5,0 x 5,0 mm	30 cm

Maulgrößen / Jaw sizes:

- 2,0 x 2,0 mm
- 2,5 x 2,5 mm
- 2,0 x 3,0 mm
- 3,0 x 4,0 mm
- 4,0 x 5,0 mm
- 5,0 x 5,0 mm
- 6,0 x 5,0 mm
- 7,0 x 6,0 mm

* Lengths are measured with closed jaws


LAMINECTOMY PUNCHES

Laminektomiestanzen

Pinzas cortantes para laminectomia

Pinces emporte-pieces a laminectomie

Pinze taglianti per laminectomia



Kerrison Rongeur

Down-cutting, 90° angled jaws,
with ejector

Ferris-Smith Kerrison handle	Kerrison Small handle	Kerrison Large handle	Jaw size (width x height)	Length*
				
		11-817	2,0 x 2,0 mm	18 cm
11-793	11-805	11-818	2,5 x 2,5 mm	18 cm
11-794	11-806		2,0 x 3,0 mm	18 cm
11-795	11-807	11-819	3,0 x 4,0 mm	18 cm
11-796	11-808	11-820	4,0 x 5,0 mm	18 cm
11-797	11-809	11-821	5,0 x 5,0 mm	18 cm
11-798	11-810	11-822	6,0 x 5,0 mm	18 cm
		11-823	7,0 x 6,0 mm	18 cm
		11-824	2,0 x 2,0 mm	20 cm
11-799	11-811	11-825	2,5 x 2,5 mm	20 cm
11-800	11-812	11-826	2,0 x 3,0 mm	20 cm
11-801	11-813	11-827	3,0 x 4,0 mm	20 cm
11-802	11-814	11-828	4,0 x 5,0 mm	20 cm
11-803	11-815	11-829	5,0 x 5,0 mm	20 cm
11-804	11-816	11-830	6,0 x 5,0 mm	20 cm
		11-831	7,0 x 6,0 mm	20 cm
		11-832	2,0 x 2,0 mm	22 cm
		11-833	2,5 x 2,5 mm	22 cm
		11-834	2,0 x 3,0 mm	22 cm
		11-835	3,0 x 4,0 mm	22 cm
		11-836	4,0 x 5,0 mm	22 cm
		11-837	5,0 x 5,0 mm	22 cm
		11-838	6,0 x 5,0 mm	22 cm
		11-839	7,0 x 6,0 mm	22 cm
		11-840	2,0 x 2,0 mm	25 cm
		11-841	2,0 x 3,0 mm	25 cm
		11-842	3,0 x 4,0 mm	25 cm
		11-843	4,0 x 5,0 mm	25 cm
		11-844	5,0 x 5,0 mm	25 cm
		11-845	2,0 x 2,0 mm	30 cm
		11-846	2,0 x 3,0 mm	30 cm
		11-847	3,0 x 4,0 mm	30 cm
		11-848	4,0 x 5,0 mm	30 cm
		11-849	5,0 x 5,0 mm	30 cm

Maulgrößen/ Jaw sizes:

-  2,0 x 2,0 mm
-  2,5 x 2,5 mm
-  2,0 x 3,0 mm
-  3,0 x 4,0 mm
-  4,0 x 5,0 mm
-  5,0 x 5,0 mm
-  6,0 x 5,0 mm
-  7,0 x 6,0 mm

* Lengths are measured with closed jaws

LAMINECTOMY PUNCHES




Laminektomiestanzen

Pinzas cortantes para laminectomia

Pinces emporte-pieces a laminectomie

Pinze taglienti per laminectomia



Kerrison Rongeur

Down-cutting, 50° angled jaws,
with ejector

Ferris-Smith Kerrison handle	Kerrison Small handle	Kerrison Large handle	Jaw size (width x height)	Length*
		11-874	2,0 x 2,0 mm	18 cm
11-850	11-862	11-875	2,5 x 2,5 mm	18 cm
11-851	11-863	11-876	2,0 x 3,0 mm	18 cm
11-852	11-864	11-877	3,0 x 4,0 mm	18 cm
11-853	11-865	11-878	4,0 x 5,0 mm	18 cm
11-854	11-866	11-879	5,0 x 5,0 mm	18 cm
11-855	11-867	11-880	6,0 x 5,0 mm	18 cm
		11-881	7,0 x 6,0 mm	18 cm
		11-882	2,0 x 2,0 mm	20 cm
11-856	11-868	11-883	2,5 x 2,5 mm	20 cm
11-857	11-869		2,0 x 3,0 mm	20 cm
11-858	11-870	11-884	3,0 x 4,0 mm	20 cm
11-859	11-871	11-885	4,0 x 5,0 mm	20 cm
11-860	11-872	11-886	5,0 x 5,0 mm	20 cm
11-861	11-873	11-887	6,0 x 5,0 mm	20 cm
		11-888	7,0 x 6,0 mm	20 cm
		11-889	2,0 x 2,0 mm	22 cm
		11-890	2,5 x 2,5 mm	22 cm
		11-891	2,0 x 3,0 mm	22 cm
		11-892	3,0 x 4,0 mm	22 cm
		11-893	4,0 x 5,0 mm	22 cm
		11-894	5,0 x 5,0 mm	22 cm
		11-895	6,0 x 5,0 mm	22 cm
		11-896	7,0 x 6,0 mm	22 cm
		11-897	2,0 x 2,0 mm	25 cm
		11-898	2,0 x 3,0 mm	25 cm
		11-899	3,0 x 4,0 mm	25 cm
		11-900	4,0 x 5,0 mm	25 cm
		11-901	5,0 x 5,0 mm	25 cm
		11-902	2,0 x 2,0 mm	30 cm
		11-903	2,0 x 3,0 mm	30 cm
		11-904	3,0 x 4,0 mm	30 cm
		11-905	4,0 x 5,0 mm	30 cm
		11-906	5,0 x 5,0 mm	30 cm

Maulgrößen Jaw sizes:

- 2,0 x 2,0 mm
- 2,5 x 2,5 mm
- 2,0 x 3,0 mm
- 3,0 x 4,0 mm
- 4,0 x 5,0 mm
- 5,0 x 5,0 mm
- 6,0 x 5,0 mm
- 7,0 x 6,0 mm

* Lengths are measured with closed jaws


LAMINECTOMY PUNCHES

Laminektomiestanzen

Pinces emporte-pieces a laminectomie

Pinzas cortantes para laminectomia


Pinze taglianti per laminectomia



Link Rongeur

11-907

Small handle, gold-plated, up-cutting, 50° angled jaws, bayonet shaped upward, jaw width 3 mm, with ejector, 20 cm



Link Rongeur

Up-cutting, straight jaws, 20 cm

Curved to left

Curved to right



11-908

11-911

11-909

11-912

11-910

11-913

Jaw size
(width x height)



3,0 x 4,0 mm




4,0 x 4,0 mm



5,0 x 4,0 mm

Lengths are measured with closed jaws



LINK-Beatty Rongeur

Shank curved 70° / 90° upward, up-cutting, 20 cm



Jaw size
(width x height)

11-914

11-918

11-915

11-919

11-916

11-920

11-917

11-921



2,0 x 3,0 mm



3,0 x 4,0 mm




4,0 x 4,0 mm



5,0 x 4,0 mm

Lengths are measured with closed jaws




The angled LINK-Beatty Rongeur allows to perform foraminotomies while preserving dorsal processes and facet joints.

LAMINECTOMY PUNCHES



Laminektomiestanzen	Pinzas cortantes para laminectomia
Pinces emporte-pieces a laminectomie	Pinze taglianti per laminectomia



50° angled jaws thin footplate			50° angled jaws		
	Jaw width	Length*		Jaw width	Length*
11-922	2 mm	18 cm	11-934	2 mm	18 cm
11-923	3 mm	18 cm	11-935	3 mm	18 cm
11-924	4 mm	18 cm	11-936	4 mm	18 cm
11-925	5 mm	18 cm	11-937	5 mm	18 cm
			11-938	6 mm	18 cm
11-926	2 mm	20 cm	11-939	2 mm	20 cm
11-927	3 mm	20 cm	11-940	3 mm	20 cm
11-928	4 mm	20 cm	11-941	4 mm	20 cm
11-929	5 mm	20 cm	11-942	5 mm	20 cm
			11-943	6 mm	20 cm
11-930	2 mm	23 cm	11-944	2 mm	23 cm
11-931	3 mm	23 cm	11-945	3 mm	23 cm
11-932	4 mm	23 cm	11-946	4 mm	23 cm
11-933	5 mm	23 cm	11-947	5 mm	23 cm
			11-948	6 mm	23 cm

Lengths are measured with closed jaws



The SimpleClean™ mechanism makes cleaning of the Rongeurs an easy matter

To open the SimpleClean™ mechanism pull the handle, turn the release and lift up the pusher.

The coloured release help to recognize more easily the individual widths of the jaw.

LAMINECTOMY PUNCHES

Laminektomiestanzen

Pinces emporte-pieces a laminectomie

Pinzas cortantes para laminectomia

Pinze taglienti per laminectomia



40° UPWARD

Black Line' Kerrison Laminectomy Punch
40° UPWARD Jaw angle
with special titanium coating on surface
Bone ejector



40° UPWARD

Thin Footplate Kerrison Laminectomy Punch
40° UPWARD Jaw angle
Thin Footplate For cervical surgery

40° UPWARD JAW

18cm/7" Shaft	20cm/8" Shaft	23cm/9" Shaft	Bite Size
11-949	11-955	11-961	1mm
11-950	11-956	11-962	2mm
11-951	11-957	11-963	3mm
11-952	11-958	11-964	4mm
11-953	11-959	11-965	5mm
11-954	11-960	11-966	6mm

40° UPWARD JAW

18cm/7" Shaft	20cm/8" Shaft	Bite Size
11-967	11-972	2mm
11-968	11-973	3mm
11-969	11-974	4mm
11-970	11-975	5mm
11-971	11-976	6mm

A reduced footplate thickness useful for a restricted operating area

40° UPWARD



40° UPWARD



Ferris-Smith Kerrison Laminectomy Punch
40° UPWARD Jaw angle
Modified ring handle



Long Anterior Kerrison Laminectomy Punch
40° UPWARD Jaw angle
with special Titanium coating on surface
Extra long shaft for Thoracic approach
Bone ejector

40° UPWARD JAW

18cm/7" Shaft	20cm/8" Shaft	Bite Size
11-977	11-982	2mm
11-978	11-983	3mm
11-979	11-984	4mm
11-980	11-985	5mm
11-981	11-986	6mm

40° UPWARD JAW

28cm/11" Shaft	Bite Size
11-987	2mm
11-988	3mm
11-989	4mm
11-990	5mm
11-991	6mm


BRAIN RETRACTORS SET



<i>Brain Retractors Set</i>	<i>Conjunto de retractores cerebrales</i>
<i>Ensemble d'écarteurs cérébraux</i>	<i>Set di divaricatori per cervello</i>



11-992
Book Walter Retractor
Complete Set with Box



11-993
Book Walter Retractor
Complete Set with Wooden Box



11-994
Leyla
Self-Retaining Brain Retractor

CARDIOVASCULAR & THORACIC INSTRUMENTS



VARICOSE VEIN PROBE SET

Krampfaderbesteck
Trousse a varices

Juego para varices
Astuccio per varici



Nabatoff Varicose Vein Stripper

12-001

135mmx100mmx25mm
complete

12-002


Nabatoff
Metal Case

12-003

Nabatoff
Handle

12-004

Nabatoff
Cable 900mm



3mm



6mm



9mm



12mm

12-005

Nabatoff

Probe tip in plastic



Nabatoff
Metal Probe
12-006
3mm



Nabatoff
Metal Olive
12-007
6mm



Nabatoff
Metal Olive
12-008
9mm



Nabatoff
Metal Olive
12-009
12mm



Nabatoff
Metal Olive
12-010
15mm



FIMA INSTRUMENTS

ATRAUMA MICRO BULLDOG CLAMPS

Atrauma Mikro-Bulldogklemmen

Micro-Clamps Bulldog con bocas Atrauma

Micro-Clamps Bulldog à mors Atrauma

Attrae i morsetti Micro Bulldog



De Bakey

12-011
7.5 cm



12-012
8.5 cm



12-013
10 cm



12-014
12 cm



De Bakey



12-015
7.5 cm



12-016
8.5 cm



12-017
9.5 cm



12-018
11 cm



Glover



12-020
7 cm



12-021
8 cm



12-022
9.5 cm



12-023
11 cm



Glover



12-024
6.5 cm



12-025
7.5 cm



12-026
9 cm



12-027
10.5 cm



Applying forceps
for bulldog clamps
12-019
24 cm

VASCULAR FORCEPS



Gefäßklemmen | Pinzas vasculares
Pincers vasculaires | Pinze vascolari



De Bakey
2.7 mm



12-028
16 cm



12-029
20 cm



12-030
25 cm



12-031
30 cm



De Bakey
1.5 mm



12-032
16 cm



12-033
20 cm



12-034
25 cm



De Bakey
2 mm



12-035
16 cm



12-036
20 cm



12-037
24 cm



12-038
30 cm



De Bakey
2 mm



12-039
16 cm



12-040
19.5 cm



12-041
25 cm



12-042
30 cm



De Bakey
3.5 mm



12-043
16 cm



12-044
20 cm



12-045
25 cm




12-046
30 cm



Cooley
2.0 mm



12-047
16 cm



12-048
20 cm

VASCULAR FORCEPS

Gefäßklemmen | Pinzas vasculares
Pincers vasculaires | Pinze vascolari

Hemostatic forceps




De Bakey-Rankin
16 cm



12-049
Straight



12-050
Curved



De Bakey-Mosquito
15 cm



12-051
Straight



12-052
Curved



Ochsner
12-053
22 cm



Diethrich
12-054
21 cm

COARCTATION CLAMPS



FIM INSTRUMENTS

Coarctationsklemmen | Pincers-Clamps pour coarctation
Pinzas para coarctación | Morsetti di Coartazione




Potts



Potts



Potts



Potts
12-061
21.5 cm
55 mm

ANASTOMOSIS CLAMPS, LIGATURE FORCEPS

Anastomosenklemmen, Ligaturklemmen

Pinzas para anastomosis, Pinzas para ligaduras

Pinces-clamps à anastomose, Pinces à ligatures

Morsetti per anastomosi, pinze per legature



Potts
12-062
22 cm
55 mm



De Bakey



12-063
19 cm




12-064
24 cm



De Bakey
12-065
16 cm



De Bakey
12-066
19 cm



PERIPHERAL VASCULAR CLAMPS

Peripheral-Gefäßklemmen
Pinzas para vasos periféricos

Pinces-clamps vasculaires périphériques
Morsetti vascolari periferici



Multi purpose Vascular Clamps



De Bakey
12-067
20 cm



De Bakey
12-068
angled, 20 cm



Leland-Jones

FIMA INSTRUMENTS

12-069
18 cm





12-070
19 cm



12-071
20 cm



MULTI PURPOSE VASCULAR CLAMPS



TANGENTIAL OCCLUSION FORCEPS



Tangential-Verschlußklemmen
Pinzas tangentes de oclusión

Pinces vasculaires pour occlusion tangentielle
Morsetti di occlusione tangenziale




Aneurysm Clamps



Satinsky



De-Bakey-Bahnson



MULTI PURPOSE VASCULAR CLAMPS

Mehrzweck-Gefäßklemmen

Pinzas vasculares de usos múltiples

Pinces vasculaires universelles

Morsetti vascolari multiuso



**Potts
straight**



**12-088
21.5 cm**



**12-089
23.5 cm**



**12-090
27 cm**



**12-091
31.5 cm**



**Potts
curved**



**12-092
21 cm**



**12-093
23 cm**



**12-094
26 cm**




**12-095
31 cm**



De Bakey



**12-096
19 cm**



**12-097
20 cm**



**12-098
23 cm**



**12-099
27.5 cm**

MULTI PURPOSE VASCULAR CLAMPS



Mehrzweck-Gefäßklemmen

Pinces vasculaires universelles

Pinzas vasculares de usos múltiples

Morsetti vascolari multiuso



De Bakery



12-100
21 cm

12-101
22 cm

12-102
26 cm

12-103
31.5 cm



De Bakery
12-104
straight
17.5 cm




De Bakery



12-105
strongly
angled
15 cm



De Bakery Beck
12-106
17 cm



MULTI PURPOSE VASCULAR CLAMPS

Mehrzweck-Gefäßklemmen

Pinzas vasculares de usos múltiples

Pinces vasculaires universelles

Morsetti vascolari multiuso



PEDIATRIC VASCULAR CLAMPS



Kinder-Gefäßklemmen

Pinces vasculaires pédiatriques

Pinzas vasculares pediátricas

Morsetti vascolari pediatrici



De Bakey

12-114
90° angled
10 cm




12-115
curved
13 cm



12-116
straight
13 cm




De Bakey



12-117
s-shape
12.5 cm

12-118
angled
12.5 cm

Anastomosis Clamps



Cooley

12-119
straight
14.5 cm



12-120
angled
14.5 cm



12-121
angled
14.5 cm



12-122
double curved
14 cm



12-123
spoon shape
14.5 cm



Cooley
12-124
angled
22 cm

MULTI PURPOSE VASCULAR CLAMPS

Mehrzweck-Gefäßklemmen

Pinzas vasculares de usos múltiples



Pinces vasculaires universelles

Morsetti vascolari multiuso

Neonatal Vascular Clamps



Castaneda Multipurpose Clamp & Atraumatic Cardiovascular Forceps



Cooley 12-133
curved sideways
17 cm



Cooley 12-134
double curved
16 cm



12-135
straight angled
14.5 cm



12-136
medium angled
16 cm



12-137
slight angled
16.5 cm

RIB CONTRACTORS



Rippenkontraktoren

Aproximadores para costillas

Rapprocheurs des côtes

Approssimatori costali



**Bailey-Gibbon
Rib Contractors
12-138**
20 cm
For Adult



**Bailey-Gibbon
Rib Contractors
12-139**
20 cm
For children




**Bailey-Baby
Rib Contractors
12-140**
15 cm



**Sellers Rib Contractor
12-141** 20 cm



**Blalock
Rib Contractors**



12-142
15 cm



12-143
20 cm



12-144
20 cm



**Johns-Hopkins
(Potts Modif.)
Rib Contractors**



12-145
20 cm



12-146
20 cm



12-147
20 cm


LUNG SPATULAS, AORTIC DILATORS

Lungenspatel, Aorta-Dilatatoren
Espátulas pulmonares, Dilatadores aórticos


Spatules pulmonaires, Dilatateurs pour aorte
Spatole per pulmone, Dilatori per aorta



**Allison Lung
Spatula Retractor**
12-148 25 cm
12-149 30 cm
12-150 32 cm
wide 3.8cm
Small Size



**Allison Lung Spatula
Retractor**
12-151 32.4 cm
wide 5.2cm
large size



**Tubbs Mitral Valve Aortic
Dilator**
12-152
Blade Opening
From 8 To 42 mm



**Cooley Blade Aortic
Dilator**
12-153
Blade Opening
From 8 To 15 mm

RIB RETRACTORS



Rippensperrer	Ecarteur de côtes
Separador costal	Divaricatore costale




Morse 12-154
Sternal Retractor 4 movable
blades 22 x 27mm spreading 200 mm
Size 1,



DeBakey 12-155
Rib Retractor spreading 150 mm
Size 1,



DeBakey 12-156
Rib Retractor spreading 200 mm
Size 1,



Cherry 12-157
mm mm mm
44 70 135



Cooley Rib Retractors 12-158
lateral blades 12 x 16 mm spreading 95 mm



Cooley Rib Retractors 12-159 lateral blades 15 x 30 mm, spreading 140 mm
12-160 lateral blades 21 x 44 mm, spreading 150 mm

RIB RETRACTORS

Rippensperrer

Ecarteur de côtes

Separador costal

Divaricatore costale



Cooley Rib Retractors

12-161

lateral blades 35 x 50 mm, spreading 200 mm



Finocchietto, Rib Retractors

12-162 lateral blades 15 x 15 mm, spreading 75 mm, infant

12-163 lateral blades 31 x 31 spreading 100 mm, baby - height

12-164 lateral blades 40 x 45 spreading 135 mm, baby - height

12-165 lateral blades 52 x 62 spreading 175 mm, baby - height

12-166 lateral blades 84 x 65 spreading 175 mm, baby - height



Haight Rib Retractor

12-167

lateral blades 30 x 30 mm, spreading 120 mm



Morse-Favaloro Sternal Retractor

12-168

4 movable blades 17 x 22 mm, spreading 160mm for children



Finocchietto (Modif.) Rib Retractors

12-169

lateral blades 14 x 15 mm, spreading 120 mm



Tuffier Rib Retractor

12-170

lateral blades 45 x 45 mm, spreading 160 mm

12-171


lateral blades 22 x 27 mm, spreading 200 mm

RIB RETRACTORS




Rippensperrer
Separador costal

Ecarteur de côtes
Divaricatore costale



Hertzer Rib Retractors
12-172
blades 20 x 25 mm, spreading 105 mm



Tuffier Rib Retractors
12-173
blades 43 x 53 mm, spreading 130 mm



Thoracic Rib Retractors
12-174 Retractor only
12-175 pair of blades, 48 x 68 mm
12-176 1 pair of blades, 70 x 75 mm



Thoracic & I.M.A Rib Retractors
12-177 Complete set
12-178 lateral blades for I.M.A Retractor
left 40 x 50 x 80 mm
12-179 lateral blades for I.M.A Retractor
Right 40 x 50 x 80 mm



Burford Rib Retractors
12-180
18 cm For children
spread with 2 sets of Blades
12-181 48 mm x 42 mm spreading 135 mm
12-182 38 mm x 42 mm spreading 135 mm



Burford Rib Retractors
12-183
18 cm For Adults
spread with 2 sets of Blades
12-184 with two pair of blade 47 x 62 mm + 65 x 62 mm spreading 200 mm
12-185 with 2 pair of blade 47 x 62 mm + 65 x 62 mm spreading 260 mm
12-186 with 2 pair of blade 47 x 62 mm + 65 x 62 mm spreading 300 mm

RIB RETRACTORS

Rippensperrer

Ecarteur de côtes

Separador costal

Divaricatore costale



Finocchietto/Baby Rib Retractors

12-187 7.5 cm

14 mm x 15 mm spreading 75 mm

Finocchietto, Rib Retractors

20 cm

12-188 30 mm x 45 mm spreading 150 mm

12-189 46 mm x 62 mm spreading 190 mm

12-190 80 mm x 65 mm spreading 190 mm



Castaneda Rib Retractors

12-191 10.5 mm x 32 mm

spreading 65 mm



Castaneda Rib Retractors

12-192 12 mm x 45 mm

spreading 70 mm



Castaneda Rib Retractors

12-193 15 mm x 60 mm

spreading 120 mm

RIB SHEARS, STERNUM CHISEL



Rippenscheren, Sternummeißel
Costotomos, Cinceles

Costotomes, ciseau-Burin à sternum
Costotomi, Scalpello per sterno



Sauerbruch-Frey Rib Shears
12-194
for the first rib
37 cm



Schumacher Rib Shears
12-195
sternum shears
21 cm



Bone Rib Shear
12-196
23 cm



Rib Shear
12-197
23 cm



Gluck Rib Shear
12-198
22 cm



Lebsche Sternum Chisel
12-199 24.5 cm



Semb Rib Shears
12-200 20 cm



Semb Rib Shears
12-201 20 cm


RIB SHEARS, STERNUM CHISEL

Rippenscheren, Sternummeißel
Costotomos, Cinceles

Costotomes, ciseau-Burin à sternum
Costotomi, Scalpello per sterno



Giertz-Stille Rib Shears
12-202 27 cm



Brunner Rib Shears
12-203 33.5 cm



Maurer-Collin Rib Shears



12-204
Hooked Left
23 cm



12-205
Hooked Right
23 cm



Sauerbruch Rib Shears
12-206 24 cm, For Children
12-207 27 cm, For Adults



Schuardt Rib Shears
12-208 17 cm



Collin Rib Shears
12-209 19 cm



Doyen Rib Shears
12-210 19 cm

OPHTHALMIC INSTRUMENTS



Micro Forceps



Chalazion Clamp
13-001
11cm



Chalazion Forceps
13-002
Oval 12 X 14mm, 10 cm long



13-003
Round 8mm,
10 cm long



13-004
Round 15mm,
10 cm long



13-005
Oval 17 X 9mm,
10 cm long



Chalazion Forceps
13-006



Cauer
13-007
10 cm
Chalazion Forceps



Suture Tying Forceps
13-008 9cm
13-009 11cm



McPherson Forceps
13-010 9cm
13-011 11cm



Corneal Forceps
1x2 (0.12)
13-012 9cm
13-013 11cm



Bonn Forceps
13-014
11cm



Utrata
Capsulorhexis
Forceps
13-015 11cm
Curved



Utrata
Capsulorhexis
Forceps
13-016 11cm
Straight



Colibri Forceps
1x2 (0.12)
13-017



Holding & Folding
Forceps
for Acrylic Lenses
13-018
11cm



Holding & Folding
Forceps
for Acrylic Lenses
13-019
11cm



Morefeldt Forceps
13-020
11cm

Round Handle Available



Suture Tying Forceps
13-021 9cm
13-022 11cm

Round Handle
Available



McPherson Forceps
13-023 9cm
13-024 11cm



Corneal Forceps
1x2 (0.12)
13-025 9cm
13-026 11cm



Bonn Forceps
13-027 9cm
13-028 11cm



Utrata Capsulorhexis Forceps
13-029 9cm
13-030 11cm
Curved



Utrata Capsulorhexis Forceps
13-031 9cm
13-032 11cm
Straight



Holding & Folding Forceps
for Acrylic Lenses
13-033 9cm
13-034 11cm



Holding & Folding Forceps
for Acrylic Lenses
13-035 9cm
13-036 11cm



Colibri Forceps
with tying platform
13-037



Gerl
Tying Forceps
curved shanks
6mm long platform
For 8-0 to 11-0 sutures.
titanium
13-038



Tennant
Tying Forceps
extra delicate smooth jaws
with 6mm long platform
round handle with guide pin
For 9-0 to 11-0 sutures.
13-039 straight
13-040 straight, titanium
13-041 curved
13-042 curved, titanium



13-043
10 cm
sharp



13-044
10 cm
blunt



Berke
13-045
11 cm
Ptosis Forceps



Förster
13-046
10 cm
Conjunctiva
Forceps



Fechtner Conjunctiva
Forceps
13-047 11 cm long
Very delicate ring tips,
with tying
platforms



Scleral Plug Forceps
13-048
Designed to securely
hold scleral plugs during
insertion and removal.
The cross-action handle
reduces hand fatigue
while grasping plug,
11.5 cm long

Round Handle Available



McPherson
Forceps
13-049



13-050



Suture Tying
Forceps
13-051



13-052



Utrata Capsulorhexis
Forceps
13-053



13-054



Corneal Forceps
1x2 (0.12)
13-055



13-056



Colibri
Forceps
extra delicate
1x2 teeth, 0.12mm
with tying platform
13-057



13-058
Cooley
18 cm



13-059
Castroviejo TC
18cm



13-060
Noyes Micro Scissors
CVD 18cm



Needle Holder
delicate, curved, without lock,
overall length: 115mm, titanium
Ideal for 9-0 to 10-0 sutures
13-061 with lock
13-062 without lock
13-063 titanium, with lock
13-064 titanium, without lock



Barraquer
Needle Holder
delicate jaws, straight
13-065 with lock
13-066 without lock
13-067 titanium, with lock
13-068 titanium, without lock

Micro Scissors



IOL Holding Forceps
13-069 11.5cm



Micro Needle Holders
13-070 11.5cm



Corneal Scissors
13-071 11.5cm



Cephalotomy Scissors
13-072 11.5cm



Vannas Scissors
13-073 11.5cm
Straight



Vannas Scissors
13-074 11.5cm
Curved



Westcott Scissors
13-075 11.5cm



Iris Scissors
13-076 11.5cm



Vannas Scissors
13-077

Sharp pointed tips. Curved 7.5mm blades, 9.5mm from pivot to tip, 8 cm long



Vannas Scissors
13-078

Sharp tips. Angled to side 7mm blades, 9mm from pivot to tip, 8 cm long



13-079



13-080



13-081



13-082



13-083



13-084



13-085

Capsulotomy Scissors Vannas

5mm long blades, sharp tips

13-079 straight

13-080 curved

13-081 angled forward

13-082 angled to side

Gills-Vannas

7mm long blades, sharp tips

13-083 straight

13-084 curved

13-085 angled forward

Mcpherson-Vannas

8mm long blades, sharp tips

13-086 curved



Barraquer
Needle Holder
delicate jaws, curved
13-087 with lock
13-088 without lock



McPherson
Needle Holder
very delicate jaws
13-089 with lock straight
13-090 without lock straight
13-091 with lock curved
13-092 without lock curved



Micro Needle Holder
Needle Holder
delicate jaws
strongly curved
13-093 with lock
13-094 without lock



Castroviejo
Needle Holder
delicate jaws, curved
13-095 with lock
13-096 without lock



Castroviejo
Needle Holder
standard jaws
straight, with lock
13-097
13-098 TC Tungsten Tip



Kalt
Needle Holder
delicate smooth jaws
straight, with lock
13-099



Micro Incision Forceps
13-100



Micro Incision Scissors
13-101



Endothelial Grasping Forceps
curved, 23 gauge
delicate end-gripping
platforms jaws 90°
off-set from "squeeze handle"
mechanism
13-102



Vitreoretinal Forceps
Lens Cartridge
Loading Forceps
straight, 20 gauge
13-103



Vitreotomy Scissors
T20 gauge shank
Stainless steel
construction with
flat handle.
Autoclavable.



Scissors
horizontal opening
straight
13-104



Scissors
horizontal opening
angled 60°
13-106



Scissors
horizontal opening
angled 30°
13-105



Scissors
horizontal opening
angled 80°
13-107



Troutman-Castroviejo
delicate patterns, with stop



13-108
10,5 cm
right




13-109
10,5 cm
left



Wecker
13-110
11 cm
pointed-pointed

Cannulas



Sterilizing Boxes (with silicon Pads)





**Curette
13-124**



**Wire Vectus
13-125**



**Strabismus Hook
13-126**



**Lense Expressor
13-127**



**Bimanual Hand Irrigation
13-128**



**Bimanual Hand Aspiration
13-129**



**Rosen
Phco Splitter
60° angled wedge
shaped inferior edge
with blunt tip
13-130
Titanium**



**Nagahara
Chopper
13-131**



**Sinskey
Hook
13-132**



**Sinskey
Hook
13-133**



**Nagahara
Chopper
13-134**



**Lachrymal Dilator
Double Ended
13-135**



**Lachrymal Dilator
13-136 Fine Taper
13-137 Medium Taper
13-138 Heavy Taper**



**Lachrymal Probes
13-139**





Sinskey II
Lens Manipulating Hook
0.2mm diameter blunt tip
13-140 straight
13-141 angled



Lens Manipulating Hook
0.15mm diameter tip
13-142
Fenzl, round handle
straight with guard
13-143
Fenzl, round handle, angled
13-144
Maltzman-Fenzl, round handle,
13-145
Maltzman-Fenzl, round handle,
angled
13-146
Jaffe-Maltzman, flat handle, angled



Fixation Hook
double, sharp
13-147
15 cm
13-148
18 cm



VON GRAEFE Knife
Graefe
1,2,3,4,5,



13-149



13-152

KNIVES



13-150



13-153



13-151



13-154



Jaeger Lid Plate
13-155
Jaeger
4 1/2", 20mm, 23mm,



13-156
Fig. 1



13-157
Fig. 2



13-158
Fig. 3

Eye Knives

Jaeger



Barraquer Wire Speculums
13-159
heavy wire, 15mm blades



Barraquer Solid Blade
Wire Speculum
13-160
Barraquer, 10mm blades




Kratz - Barraquer
Wire Speculum
15mm blades
one solid and one open wire
13-161
for right eye



Kratz - Barraquer
Wire Speculum
15mm blades
one solid and one open wire
13-162
left eye



Wire
Speculums
13-163



Wire
Speculums
13-164



Laster-Burch Speculum
eye Speculum with flanges on
blades for retaining eye lashes
adjustable mechanism
with locking screw
13-165
11mm blades



Lieberman
Speculum
13-166



Clark
Speculums
13-167



Jameson
Caliper chrome
plated
13-168



Braunstein Caliper
3.5/4.0mm
Stainless Steel
13-169



13-170
Calliper
Fixed



Fixed
Caliper
13-171



Castroviejo Caliper
20mm spread
13-172 straight
13-173 angled
13-174 curved



Maloney
Intra-Operative
Keratometer titanium
For qualitative
measurement
of astigmatism.
13-175

ENT SURGERY



EAR SPECULUM



Ohrspekulum

Espéculo de oreja

Spéculum auriculaire

Speculum auricolare



14-001
4.5 mm Ø



14-002
5.5 mm Ø



14-003
6.5 mm Ø



14-004
7.5 mm Ø



14-007
2 mm Ø



14-008
3 mm Ø



14-009
4 mm Ø



14-010
5 mm Ø



Toynbee Ear Speculum
3.6 cm 1 1/2"
Set for adults

14-005
Set of 3
5.5 mm
6.5 mm
7.5 mm

14-006
Set of 4
4.5 mm
5.5 mm
6.5 mm
7.5 mm



Poltzer Ear Speculum
Set for adults

14-011
Set of 3
2 mm
3 mm
4 mm

14-012
Set of 4
2 mm
3 mm
4 mm
5 mm



Gruber Ear Speculum
Set for adults



14-013
a - 4.5 mm Ø



14-014
a - 5.5 mm Ø



14-015
a - 6.5 mm Ø



14-016
a - 7.5 mm Ø



Gruber Ear Speculum
Set for adults



14-017
a - 4.5 mm Ø



14-018
a - 5.5 mm Ø



14-019
a - 6.5 mm Ø



14-020
a - 7.5 mm Ø



Boucheron Ear Speculum
Set for adults



14-021
4.5 mm Ø



14-022
5.5 mm Ø



14-023
6.5 mm Ø



14-024
7.5 mm Ø



Boucheron Ear Speculum
Set for children



14-025
3 mm Ø



14-026
4.5 mm Ø



14-027
5.5 mm Ø

EAR SPECULUM

Ohrspekulum

Espéculo de oreja

Spéculum auriculaire

Speculum auricolare



Hartmann Ear Speculum
for children

- 14-028 2 mm Ø
- 14-029 3.5 mm Ø
- 14-030 5 mm Ø



Hartmann Ear Speculum
Set for adults

- 14-031 4.5 mm Ø
- 14-032 5.5 mm Ø
- 14-033 6.5 mm Ø
- 14-034 7.5 mm Ø



Zollner Ear Speculum
Set for adults, Black

- 14-035 4.5 mm X 6.5 mm
- 14-036 6.5 mm X 7.5 mm
- 14-037 7.5 mm X 8.5 mm
- 14-038 8.5 mm X 9.5 mm
- 14-039 9.5 mm X 10.5 mm
- 14-040 10.5 mm X 11.5 mm



Hartmann Ear Speculum
Set for adults, Black


- 14-041 2 mm
- 14-042 3.5 mm
- 14-043 5 mm
- 14-044 6.5 mm
- 14-045 7.5 mm




Gruber Aural Specula Ear Speculum
Black finish, 38 mm Ø
14-046 Set of 4
14-047 3 mm
14-048 4 mm
14-049 5 mm
14-050 7 mm



Rosen Ear Speculum
Black finish, Slotted round ends 35 mm Ø
14-051 Set of 4
14-052 3 mm
14-053 4 mm
14-054 5 mm
14-055 7 mm



Rosen Ear Speculum
Black finish, 35 mm Ø
14-056 Set of 4
14-057 4 mm
14-058 5 mm
14-059 6 mm
14-060 7 mm



Tumarkin Ear Speculum
Tapered slot set of 8
14-061 4.5 mm
14-062 6.5 mm
14-063 7.0 mm
14-064 7.5 mm
14-065 8.0 mm
14-066 8.5 mm
14-067 9.0 mm
14-068 11 mm

EAR SPECULUM



Ohrspekulum

Espéculo de oreja

Spéculum auriculaire

Speculum auricolare



14-069
6.5 mm Ø



14-070
7 mm Ø



14-071
8 mm Ø

Siegle Magnifier Ear Speculum
14-072
Complete Set
with 3 Ear Specula & Otoscope



Siegle Magnifier
Otoscope



Bruenings complete
14-073
Bruenings Ear Speculum
Complete Set
with 4 Ear Specula & Otoscope



14-074
6.5 mm Ø



14-075
7.5 mm Ø



14-076
8.5 mm Ø



14-077
9.5 mm Ø



Bruening Magnifier
Otoscope

Ohrenspekulum, Zerstäuber und Ohrenspritzen
Spéculum auriculaire, atomiseur et seringues auriculaires

Espéculo de oreja, atomizador y jeringas de oreja
Speculum auricolari, atomizzatori e siringhe auricolari



**Richnau Holmgren
Ear Speculum**
3.5 cm

14-078
3 mm Ø


14-079
4 mm Ø

14-080
5 mm Ø

14-081
6 mm Ø



Vacher Nasal Speculum
14-082



Voltolini Nasal Speculum
14-083 Fig. 1
14-084 Fig. 2
14-085 Fig. 3



Yankauer
14-086
19cm

14-087
2 mm Ø

14-088
2.5 mm Ø

14-089
3 mm Ø

14-090
4 mm Ø

14-091
5 mm Ø



Hartmann
Eustachian Catheters



Reiner Ear syringe
14-092
100CC Complete
14-093
50CC Complete
14-093
75CC Complete



**Reiner Ear
Syringe
Tip Only**
14-096
Tip Only



**Reiner Ear
Syringe
Tip Only**
14-097
Tip Only



**Reiner Ear
Syringe Shield**
14-098

14-094
100CC Complete
14-095
150CC Complete
Reiner Ear syringe
With Tip & Shield


TUNING FORKS



Stimmgabeln	Diapasones
Diapasons	Diapasoni



Hartmann
Tuning Fork
14-099
C 64



Hartmann
Tuning Fork
14-100
c 128



Hartmann
Tuning Fork
14-101
c1 256



Hartmann tuning
fork with or without
plastic base a1 440
14-102
a1 440



Tuning Forks neurol
with hammer
14-103



Barany Noise Drum
14-104

NOISE DRUM

Lärmtrommel	Aparato
Assourdisseur	Assordatore



Lucae
Tuning Fork
14-105
c 128
Variable up to h (ehn, b)



Rydel-Seiffer
Tuning Fork
14-106
C 64 / c 128
With clamps C 64,
Without clamps c 128



French
Tuning Fork
14-107
C2 16



French
Tuning Fork
14-108
C1 32

TUNING FORKS

Stimmgabeln | Diapasones
Diapasons | Diapasoni



14-109
c 128



14-110
c1 256



14-111
c2 512



14-112
c3 1024



14-113
c4 2048




14-114
c5 4096

Hartmann Tuning Fork



14-115
Hartmann tuning Forks Set case/5



14-116
c 128



14-117
c1 256



14-118
c2 512



14-119
c3 1024



14-120
c4 2048



14-121
TUNING FORKS Set in case of 5

EAR POLYPUS SNARES




Ohrpolypenschnürer

Serre noeuds auriculaires

Polipótomos auriculares

Serranodi auricolari



**Nasal Polypus
Snare
14-122**

0.3 mm Ø X 100 mm
In bundles of 12 pieces

**Krause Nasal
Polypous Snare
14-123**

25.4 cm

**Kraus -Vosse
Nasal Polypous Snare
14-124**

25.4 cm
With 2 inserts

**Kraus -Wilde
Nasal Polypous Snare
14-125**

24 cm

**Krause
14-126**

26 cm


EVE Tonsil Snare

28 cm

14-127 without ratchet
14-128 with ratchet



Eves Tonsil Snare



14-129
28 cm

14-130
28 cm




Eves Tonsil Snare



14-131
28 cm



**Bruenings
Tonsil Snare**



14-132
28 cm

TONSIL INSTRUMENTS

Tonsilleninstrumente

Instruments pour amygdales

Instrumentos para amígdalas


Strumenti per tonsille



Bruenings
14-133
10 m 33 ft
0.3 mm Ø
14-134
10 m 33 ft
0.4 mm Ø



Bruenings
14-135
40 mm
Package of 100 pieces
14-136
50 mm
Package of 100 pieces




Tyding Tonsil Snare
14-137
With one straight tip

NASAL SPECULUM



Thudichum (Goldsmith) Nasal Speculum
6.5cm
14-138 Fig. 1
14-139 Fig. 2
14-140 Fig. 3
14-141 Fig. 4
14-142 Fig. 5
14-143 Fig. 6
14-144 Fig. 7



ST-CLAIR THOMSON
14-145 Fig 1
14-146 Fig 2
14-147 Fig 3

NASAL SPECULUM



Nasenspekulum
Spéculum nasal

Espéculo nasal
Speculum nasale




**Thudichum
(Goldsmith)
Nasal Speculum**



14-148
6.5 cm
21 mm X 9.5 mm



14-149
6.5 cm
21 mm X 10.5 mm




14-150
6.5 cm
21 mm X 11.5 mm




**Hartmann
Nasal Speculum**



14-151
16 cm
Tip Size
27mm



14-152
16 cm
Tip Size
31mm



14-153
16 cm
Tip Size
35mm



**Hartmann-Halle
Nasal Speculum**



14-154
15 cm



14-155
15 cm



14-156
15 cm

NASAL SPECULUM

Nasenspekulum
Spéculum nasal

Espéculo nasal
Speculum nasale




**Tieck-Halle
Nasal Speculum**



14-157
13.5 cm
For babies



Schwarz Nasal Speculum



14-158
14 cm



Vienna Speculum
14-159 13cm
14-160 14.6cm
medium standard patten



14-161
Nose Speculum,
thigh 30mm




Killian Nasal Speculum




14-162
13 cm
X = 35 mm



14-163
13 cm
X = 50 mm



14-164
13 cm
X = 75 mm



14-165
13 cm
X = 90 mm

NASAL SPECULUM



Nasenspekulum
Spéculum nasal

Espéculo nasal
Speculum nasale



**Killian Septum Speculum
14-166**
15 cm
2.51 cm blades
with set screw



**Killian Septum Speculum
14-167**
15 cm
8.9cm blades
set screw




**Killian Septum Speculum
14-168**
Fiber Optic 15 cm
length 2.51mm deep

FIMA INSTRUMENTS



**Beckmann
Nasal Speculum**



14-169
15 cm
25mm



14-170
15 cm
28mm



14-171
15 cm
32mm




**Mod. Wien
Nasal Speculum**



14-172
14.5 cm
25cm



14-173
14.5 cm
27cm



14-174
14.5 cm
30cm


NASAL SPECULUM

Nasenspekulum
Spéculum nasal


Espéculo nasal
Speculum nasale




Cottle Nasal Speculum



14-175
15 cm
X = 35 mm



14-176
15 cm
X = 50 mm




14-177
15 cm
X = 75 mm

14-178
15 cm
X = 90 mm


FIMA INSTRUMENTS



Cottle Nasal Speculum



14-179
14 cm
X = 45 mm




14-180
14 cm
X = 75 mm

EAR CURETTES



Ohrküretten	Curetas auriculares
Curettes auriculaires	Curette auricolari



Pictures show the outer diameter

EAR CURETTES AND SHARP SPOONS

Ohrküretten und scharfe Löffel


Curettes auriculaires et cuillères tranchantes

Curetas y cucharillas auriculares cortantes

Curette auricolari e cucchiari taglienti



EAR LOOPS - HOOKS



TONSIL KNIVES



NASAL KNIVES



Nasenmesser

Cuchillos nasales

Bistouris pour nez

Bisturi per naso



**Bruenings
Tonsil Knife**
14-238
23 cm 9"



**Abraham
Tonsil Knife**
14-239
22 cm 8 3/4"



**Fischer Tonsil
Knife and Dissector**
14-240
21.6 cm 8 1/2"
semi-sharp tip, with teeth



**Fomon
Nasal Knife**
14-241
16.5 cm



**Joseph
Nasal Knife**

**14-242
Fig.1**
16.5 cm
S/S Str.



**14-244
Fig.1**
16.5 cm
S/S Cvd.



**14-243
Fig.2**
16.5 cm
S/S Str.



**14-245
Fig.2**
16.5 cm
S/S Cvd.



**Joseph
Nasal Knife**
14-246
14 cm



**Cottle
Nasal Knife**
14-247
15 cm



**Joseph
Nasal Knife**
14-248
17.5 cm



**Converse
Nasal Knife**
14-249
16 cm



**Freer
Nasal Knife**
14-250
15 cm



**Freer
Nasal Knife**
14-251
15 cm



**Freer
Nasal Knife**
14-252
15 cm



**Freer
Nasal Knife**
14-253
19 cm



**Freer
Nasal Knife**
14-254
19 cm



**Bruenings
Septatome**
14-255
16.5 cm
14-256
Bruenings
Spare blade

ANTRUM CURETTES

Antrumküretten | Curetas para el antrum
Curettes d'antrum | Curette per antrum



**Coakley
Antrum Curette**
14-257
21 cm



14-258
20 cm



14-259
20 cm



Antrum Curette
14-260
20 cm
7 mm X 11 mm



**Pratt Antrum
Curette**
14-261
21 cm



**Faulkner
Antrum Curette**
14-262
21 cm



**Coakley
Antrum Curette**



14-263
17 cm



14-264
17 cm



14-265
17 cm



14-266
17 cm



14-267
17 cm



14-268
17 cm

Antrum Curette
14-269
17 cm
Complete set



Hypophyseal
14-270
Curette 19cm



Hypophyseal
14-271
Curette 19cm

SEPTUM ELEVATORS - RASPATORIES



Septum elevatorien - Raspatorien

Elevadores - leguas para el tabique nasal

Elevateurs - rugines pour la cloison nasale

Elevatori - periostotomi per setto



**Lempert
Septum
Elevators
Raspatory
14-272**
20.5 cm 8"
X = 6 mm



**Lempert
Septum
Elevators
Raspatory
14-273**
18 cm 7"
X = 6 mm



**Cottle
Septum
Elevators
Raspatory
14-274**
19 cm 7 1/2"
X = 8 mm
Sharp



**Joseph
Septum
Elevators
Raspatory
14-275**
16 cm 6 1/4"
X = 5 mm
Periosteal elevators



**Joseph-Masing
Septum Elevators
Raspatory
14-276**
21.5 cm 8 1/2"
Sharp



**Pennington
Septum
Elevators
Raspatory
14-277**
21.5 cm 8 1/2"
Blunt



**Freer
Septum
Elevators
Raspatory
14-278**
19 cm 7 1/2"
Blunt



**Freer
Septum
Elevators
Raspatory
14-279**
18 cm 7"
Sharp/Blunt



**Freer
Septum
Elevators
Raspatory
14-280**
20.5 cm 8"
Sharp/Blunt



**Cottle
Septum
Elevators
Raspatory
14-281**
23.5 cm 9 1/4"
X = 5 mm
Sharp

Periostelevatorien, Septumelevatorien

Legras para periostio, Elevadores para el tabique nasal

Elevateurs pour la périoste, Elevateurs pour la cloison nasale

Elevatori per il periosto, Elevatori per setto




Masing Periosteal Septum Elevators
14-282 20 cm 8"



Ballenger-Hajek Periosteal Septum Elevators
14-283 19 cm 7 1/4"
14-284 22 cm 8 3/4"



Cottle Periosteal Septum Elevators
14-285 22 cm 8 3/4"



Cottle Periosteal Septum Elevators
14-286 21 cm 8 1/4"



2.5 mm 2.0 mm
14-287
2.7 mm 2.0 mm
14-288

Watson-Cheyne
13 cm 5" Raspatory



Boies
14-289 19 cm 7 1/2"
Nasal fracture elevator, Blunt



Goldman
14-290 17 cm 6 3/4"
Septum displacer



Henke Tonsil Dissector
14-291 23 cm



Henke Tonsil Dissector
14-292 23 cm



Hurd Tonsil Dissector
14-293 22 cm



Carpenter Enucleator
14-294 25 cm

NASAL SAWS / FILE



Nasensägen / Feile
Scies nasales / Lime

Sierras nasales / archivo
Seghe / file nasali



**Luniatschek
Tonsil
Tampon
Plugger
14-295**
19 cm



**Lewis
Tonsil
Screw
14-296**
19.5 cm



**Joseph
Nasal
Saw
14-297**
19 cm 7 1/2"



**Joseph
Nasal
Saw
14-298**
19 cm 7 1/2"



**Joseph
Nasal
Saw
14-299**
19 cm 7 1/2"



**Maltz Rasp
14-300**
21 cm 8-1/4"
improved double end
pattern, forward and
backward cutting



**Joseph
Rasp
14-301**
15.9 cm
fine cross
serrations



**Joseph Nasal Saw
14-302**
17 cm 6 3/4"

**Joseph Nasal Saw
14-303**
17 cm 6 3/4"

RASPS WITH TUNGSTEN CARBIDE




**Fomon Rasp
with Tungsten
Carbide TC
20cm 8"**



Teeth 12
Teeth 14

14-304




Teeth 16
Teeth 20

14-305




Teeth 24
Teeth 30

14-306



Teeth 36
Teeth 44

14-307



Teeth 52
Teeth 62

14-308


NASAL SAWS / FILE

Nasensägen / Feile

Sierras nasales / archivo

Scies nasales / Lime


Seghe / file nasali




Cottle Nasal Saw
14-309
21 cm 8 1/4"



Maltz Nasal Saw



14-310
17 cm 6 3/4"
34 mm x 7 mm
Pushing cut



14-311
17 cm 6 3/4"
34 mm x 7 mm
Drawing cut




Aufricht Nasal Saw
Glabella Rasp,
backward cutting



14-312
21 cm 8 1/4"
Pushing cut



14-313
21 cm 8 1/4"
Drawing cut



14-314
17.5 cm 7"
40 mm X 9 mm

Maltz Nasal Saw



Wiener Nasal Saw



14-315
18 cm 7"
3.5 mm Ø
Antrum rasps



14-316
18 cm 7"
3.5 mm Ø
Antrum rasps



Gallaher Nasal Saw



14-317
16 cm 6 1/4"



14-318
16 cm 6 1/4"



14-319
16 cm 6 1/4"



Fomon Nasal Saw
14-320
21 cm 8 1/4"



RINGMESSER



adenoid currettes

Couteaux a vegetations adenoides

Cuchillos para vegetaciones adenoides

Adenotomi



St. Clair-Thomson
Adenoid Curette



14-321
21 cm 8 1/4"



14-322
21 cm 8 1/4"



Barnhill
Adenoid Curette



14-325
22 cm 8 3/4"



14-326
22 cm 8 3/4"



14-323
21 cm 8 1/4"



14-324
21 cm 8 1/4"



14-327
22 cm 8 3/4"



14-328
22 cm 8 3/4"



Beckmann
Adenoid Curette



14-329
22 cm 8 3/4"
Straight



14-330
22 cm 8 3/4"
Straight



14-331
22 cm 8 3/4"
Straight



Beckman
Adenoid Curette



14-335
22 cm 8 3/4"
Curved



14-336
22 cm 8 3/4"
Curved



14-337
22 cm 8 3/4"
Curved



14-332
22 cm 8 3/4"
Straight



14-333
22 cm 8 3/4"
Straight



14-334
22 cm 8 3/4"
Straight



14-338
22 cm 8 3/4"
Curved




14-339
22 cm 8 3/4"
Curved



14-340
22 cm 8 3/4"
Curved




La Force




La Force
adenotome
unperforated
Fig 1
14-341
21.6 cm 8 1/2"
25 x 16.0mm
Infant




La Force
adenotome
perforated
Fig 1
14-342
21.6 cm 8 1/2"
25 x 16.0mm
Infant




La Force
adenotome
unperforated
Fig 2
14-343
24.1 cm 9 1/2"
26 x 19.3mm
Medium



La Force
adenotome
perforated
Fig 2
14-344
24.1 cm 9 1/2"
26 x 19.3mm
Medium



La Force
adenotome
unperforated
Fig 3
14-345
25 cm 10"
30 x 22.2mm
Large



La Force
adenotome
perforated
Fig 3
14-346
25 cm 10"
30 x 22.2mm
Large



14-347
20 cm 8"
X = 3 mm
Straight
14-348
20 cm 8"
Straight
X = 4 mm
14-349
20 cm 8"
X = 5 mm
Straight

Ballenger Knife



Ballenger Knife




14-350
20 cm 8"
X = 3 mm
Bayonet
14-351
20 cm 8"
Bayonet
X = 4 mm
14-352
20 cm 8"
X = 5 mm
Bayonet



Wieder Tongue Depressor
14-353 large blade 36mm
wide serrated



Wieder Tongue Depressor
14-354 43.5mm
width blade



Blade
Fig 1
14-355

Blade
Fig 2
14-356

Blade
Fig 3
14-357

Sluder-Ballenger Tonsil Guillotines & Punch Forceps
14-358
complete set with 3 blades


Sluder-Ballenger Tonsil Guillotines & Punch Forceps
Handle
14-359 Handle



14-360
X = 12 cm
13 mm Ø



14-361
X = 15 cm
13 mm Ø



14-362
X = 17 cm
13 mm Ø

Hartmann Tonsil Guillotines & Punch Forceps



NASAL POLYPUS AND SEPTUM FORCEPS



Nasenpolypenzangen-Septumzangen
Pinces à polypes et la cloison nasale

Pinzas para pólipos y el tabique nasal
Pinze per polipi e setto nasale



Westmacott Nasal Polypus and Septum Forceps
14-363 20 cm 8"



Lubet-Barbon Nasal Polypus and Septum Forceps
14-364 20 cm 8"



Hartmann Nasal Polypus and Septum Forceps
14-365 16 cm 6 1/4"
14-366 20 cm 8"



Tilley Nasal Polypus and Septum Forceps
14-367 14 cm 5 1/2"



Heymann-Knight Nasal Polypus and Septum Forceps
14-368 18 cm 7"



Craig Nasal Polypus and Septum Forceps
14-369 18 cm 7"



Knight Nasal Polypus and Septum Forceps
14-370 18 cm 7"



Littauer Nasal Polypus and Septum Forceps
14-371 18 cm 7"

EAR POLYPUS FORCEPS

Ohrpolypenzangen

Pinces pour polypes auriculaires

Pinzas para pólipos auriculares

Pinze per polipi auricolari



**Littauer Ear
Polypus Forceps
14-372**
12 cm 4 3/4"



**Hartmann Ear
Polypus Forceps
14-373**
12 cm 4 3/4"
Delicate jaw



**Hartmann Ear
Polypus Forceps
14-374**
14 cm 5 1/2"
Normal jaw



**Quire Ear
Polypus Forceps
14-375**
10.5 cm 4 1/4"
Foreign body lever



**Buck Ear
Polypus Forceps
14-376**
11.5 cm 4 1/2"
Foreign body forceps



**Wilde Nasal
Dressing Forceps
14-377**
16.5cm 6 1/2"
angular

FIMA INSTRUMENTS

TONSIL INSTRUMENTS



**Troeltsch
Dressing Forceps
14-378**
11cm, 4 1/2"



**Troeltsch
Dressing Forceps
14-379**
13 cm, 5"



**Lucae Dressing
Forceps
14-380**
14 cm, 5 1/2"



**Jansen Dressing
Forceps
14-381**
16 cm, 6 1/4"

Nasal Polypus & Nasal Dressing Forceps



Nasenpolypenzangen, Nasentamponzange


Pinzas para pólipos nasales, Pinza nasale para tamponiamento

Pinces à polypes nasales, Pince à pansement pour le nez


Pinze per polipi nasali, Pinza per tamponamento nasale



Heymann Nasal Dressing Forceps




Bruenings Nasal Cutting Forceps



Takahashi Nasal Cutting Forceps



Noyes Nasal Dressing Forceps
14-391
X = 11 cm 3 1/2"
1 X 2



Luc Nasal Cutting Forceps



Watson-Williams Nasal Cutting Forceps
14-398
X = 12 cm 4 3/4"


NASAL CUTTING FORCEPS

Scharfe Nasenzangen

Pinces coupantes pour le nez


Pinzas nasales cortantes

Pinze taglienti per naso



Weil-Blakesley
Nasal Cutting Forceps

FIMA INSTRUMENTS



Myles Nasal
Cutting Forceps

Hartmann Nasal
Cutting Forceps

PUNCHES SPHENOID BONE

Keilbeinstanzen

Pinzas para esfenoides

Pinces pour le sinus sphénoïdal

Pinze sfenoidali taglienti




**Citelli Bone
Ear Rongeurs**
14-427 17 cm 6 3/4"
X = 60 mm



**Kerrison Bone
Ear Rongeurs**
14-428 20 cm 8"
X = 80 mm



Beyer Bone Ear Rongeurs
14-429 14 cm 1 mm
14-430 19 cm 2 mm



**Hajek Rongeur
punch**
14-431
X - 14 cm
4 mm X 4 mm
14-432
X - 14 cm
4 mm X 4 mm

**Hajek-Kofler Rongeur
punch**
14-433
X - 14 cm
3.5 mm X 3.5 mm
14-434
X - 14 cm
3.5 mm X 3.5 mm



**Hajek-Claus
Rongeur punch**
14-435 20 cm



**Hajek-Kofler
Rongeur punch**
14-436
X - 18 cm
3.5 mm
Shaft rotates 360°



Bruenings-Brunk
14-437
19 cm 7 1/2"
Antrum Retractor
Serrated



Ritter-Halle
14-438 16 cm 6 1/4"
2.5 mm Ø
14-439 16 cm 6 1/4"
3.0 mm Ø
14-440 16 cm 6 1/4"
4.0 mm Ø
Forntal sinus bougies



Cottle
14-441
1.2 cm 4 3/4"
Handle



Cottle
14-442
7 mm
Guide support
for osteotomes



Cottle
14-443
65 mm X 30 mm X 30 mm
Cartilage crusher



Cottle-Walsham
Septum Forceps



14-444
23 cm 9"
Straight



14-445
23 cm 9"
Left



14-446
23 cm 9"
Right



Ash Septum
Forceps
14-447
23 cm 9"

SEPTUM FORCEPS



Septumzangen

Pinzas para el tabique nasal

Pinces pour la cloison nasale


Pinze per il setto nasale



**Middleton-Jansen
Septum Forceps
14-448 19 cm**




**Middleton-Jansen
Septum Forceps
14-449 19 cm**



**Caplan Septum
Forceps
14-450 20 cm
Toothed**



**Rubin Septum
Forceps
14-451
21 cm
Protection cap**



**Killian
14-452
20 cm
Laryngotomy scissors**



**Reichert Epiglottis
Retractor
14-453
23 cm**



**Aufrecht Retractor
Speculum
14-454
16.5cm
solid blade**

RETRACTORS

Wundhaken	Ecarteurs
Separadores	Divaricatori



Tobold
14-455
Tongue Depressor
14.5 cm



Hartmann
14-456
15 cm



Cheek Retractor
14-457
14.5 cm



Cawood-Minnesota
14-458
Cheek Retractor
15 cm



Sternberg
14-459
Lip and Cheek Retractor
14 cm



Obwegeser
14-460
Cheek Retractor
12 cm



Vestibulum Retractor
14-461
14 cm



Oringer
Lip and Cheek Retractor
14-462 9.5cm
14-463 10 cm
14-464 10.5cm
14-465 11 cm

SPATULA, RETRACTORS



<i>Spachtel, Retraktoren</i>	<i>Espátula, retractores</i>
<i>Spatule, écarteurs</i>	<i>Spatola, divaricatori</i>



**Obwegeser
14-466**
Channel Retractor
17 cm

**Obwegeser
14-467**
Chin Retractor
16.5 cm



**Retractor
14-468**
18.5 cm



**Tongue Spatula
14-469**
16 cm



**Martin
14-470**
Tongue Spatula
18.5 cm



**Haberer flexible
Abdominal and Intestinal
14-471**
spatulas 33cm (40mm)



**Bruenings Tongue
Depressors
14-472**
19 cm



**Obwegeser
14-473**
Soft Tissue Retractor
21.5 cm



**Obwegeser
14-474**
Nasal Spine Retractor
22 cm



**Obwegeser
14-475**
Ramus Retractor
22 cm



**Doyen
14-476**
Tongue Spatula
13.5 cm

MOUTH GAGS

Mundsperrer	Abrebocas
Ouvre-bouches	Apribocca




**Whitehead
Mouth Gag**
14-477 11 cm 4 1/4"
14-478 13 cm 5"
14-479 15 cm 6"



**Heister
Mouth Gag**
14-484
13 cm 5"



**Maunder
Mouth Gag**
14-485
10 cm 4"
4 mm Ø
Wood



**Jennings
Mouth Gag**
14-480 9 cm 3 1/2"
14-481 11 cm 4 1/4"
14-482 13 cm 5"
14-483 15 cm 6"



Molt Mouth Gag
14-486
11 cm 4 1/4"
For children



Molt Mouth Gag
14-487
14 cm 5 1/2"
For adults



**Denhart
Mouth Gag**
14-488
12 cm 4 3/4"



**Fergusson-Ackland
Mouth Gag**
14-489
14.5 cm 5 3/4"



**Roser-Konig
Mouth Gag**
14-490
16 cm 6 1/4"
14-491
19 cm 7 1/2"

MOUTH GAGS

Mundsperrer | Abrebocas
Ouvre-bouches | Apribocca



Complete with frame, 3 tongue depressors and 2 adjustable lateral blades



MOUTH GAGS

Mundsperrer

Ouvre-bouches

Abrebocas

Apribocca

**Davis-boyle Mouth
Gag Complete Set
with five blades
14-509**



14-510
Mouth Gag
Blade
20 mm X 26 mm



14-511
Mouth Gag
Blade
23 mm X 65 mm



14-512
Mouth Gag
Blade
24 mm X 80 mm



14-513
Mouth Gag
Blade
25 mm X 93 mm



14-514
Mouth Gag
Blade
26 mm X 108 mm



**Davis-Boyle
Mouth Gag
Handle
14-515**



14-517
Mouth Gag Blade
23 mm X 70 mm



14-518
Mouth Gag Blade
25 mm X 86 mm



14-519
Mouth Gag Blade
26 mm X 100 mm

**McLvor Mouth Gag
Complete Set
14-516**



**McLvor Mouth Gag
Handle
14-520**



**Roser-Knig
14-521
Mouth Gag
13.5 cm**

TONSIL INSTRUMENTS



Tonsilleninstrumente
Instruments pour amygdales

Instrumentos para amígdalas
Strumenti per tonsille



Jackson
14-522
22 cm 8 3/4"
Laryngeal forceps



Marschik
14-523
22 cm 8 3/4"
Tonsil seizing forceps




Blohmke
Tonsil
Forceps
14-524
20 cm 8"



Colver Tonsil
Forceps
14-525
19 cm 7 1/2"
Straight
14-526
19 cm 7 1/2"
Curved



White Tonsil
Seizing Forceps
14-527
19.1 cm 7 1/2"
14-528
23 cm 9"




Museux
Tonsil Holding
Forceps
14-529 20 cm 8"

LARYNGEAL MIRRORS

Kehlkopfspiegel
Miroirs laryngiens

Espejos laringeos
Specchi laringei



14-530 8 mm Ø
14-531 10 mm Ø
14-532 12 mm Ø
14-533 14 mm Ø
14-534 16 mm Ø
14-535 18 mm Ø
14-536 20 mm Ø
14-537 22 mm Ø
14-538 24 mm Ø
14-539 26 mm Ø
14-540 28 mm Ø
14-541 30 mm Ø

Laryngeal mirrors



Laryngeal Mirrors
Handle
14-541
10 cm



Laryngeal Mirrors
with Handle

14-542 8 mm Ø
14-543 10 mm Ø
14-544 12 mm Ø
14-545 14 mm Ø
14-546 16 mm Ø
14-547 18 mm Ø
14-548 20 mm Ø
14-549 22 mm Ø
14-550 24 mm Ø
14-551 26 mm Ø
14-552 28 mm Ø
14-553 30 mm Ø



Michel Moveable
Mirror

14-554 10 mm Ø
14-555 12 mm Ø
14-556 14 mm Ø
14-557 16 mm Ø
14-558 18 mm Ø

GOITRE GRASPING FORCEPS - GOITRE ENUCLEATOR



Mastin Goitre
Grasping Forceps
14-559
17 cm



Lahey Goitre
Grasping Forceps
14-560 15 cm
14-561 20 cm



Kocher Goitre
Enucleator
14-562
16 cm



Kocher Goitre
Enucleator
14-563
16 cm
With Hole

TRACHEA DILATORS, TRACHEA TUBES



Trachealdilatatoren, Trachealkanülen

Dilatadores para la traqueotomia, Cabykas oara ka traqyea

Dilatateurs pour la tracheotomie, Canules pour la tracheotomie

Dilatatori tracheali, Cabyke tracheali



**Laborde Trachea
Dilator 3 Prong**
14-564 14 cm



**Trousseau Trachea
Dilator 2 Prong**
14-565 14 cm



Bowlby Trousseau Tracheal Dilator
Bowlby Paediatric
14-566 15.5 cm

FIMA INSTRUMENTS



**Uckermann
Trachea Tubes**
14-567
13 cm
10 mm Ø



**Luer-Standard
Trachea Tubes**
German silver

14-568 4 mm Ø
14-569 5 mm Ø
14-570 6 mm Ø
14-571 7 mm Ø
14-572 8 mm Ø
14-573 9 mm Ø
14-574 10 mm Ø
14-575 11 mm Ø
14-576 12 mm Ø
14-577 13 mm Ø



**Chevalier-Jackson
Trachea Tubes**
German silver

14-578 4 mm Ø
14-579 5 mm Ø
14-580 6 mm Ø
14-581 7 mm Ø
14-582 8 mm Ø
14-583 9 mm Ø
14-584 10 mm Ø
14-585 11 mm Ø
14-586 12 mm Ø
14-587 13 mm Ø

RETRACTORS

Wundhäkchen, Mikrohäkchen
Ganchito, Micro Ganchitos

Crochet, Micro crochets
Divaricatore, Micro Divaricatori



Bruenings Retractor

14-588	20 cm
14-589	30 cm
14-590	40 cm
14-591	50 cm
14-592	60 cm





Bruenings Retractor 14-593 20 cm



Huber Retractor 14-594

TRACHEOTOMY SET



14-595

Tracheotomy Broncho-Oesophagoscopy Forceps set, In metal case,
Consisting of:

- 200mm x 105mm x 30mm
- 1 Handle, 1 Blade
- 1 Handle, 1 Blade
- 1 Dressing Forceps, 130mm
- 2 Tissue Forceps, 1 x 2 teeth, 130mm
- 2 Rochester Pean Hemostatic forceps, 130mm
- 1 Grooved director, 130mm
- 1 Bose trachea dilator
- 2 iterson trachea retractors, blunt
- 1 Bose trachea retractor, 160mm, sharp, right
- 1 Bose trachea retractor, 160mm, left, right
- 1 Luer trachea tube, standard, 7mm, german silver
- 1 Luer trachea tube, standard, 8mm, german silver
- 1 Cleaning brush

EAR FORCEPS


Ohrzängchen | Pinzas auriculares
Pinces auriculaires | Pinze auricolari



**Hartman
Ear Forceps
14-596**



**Hartmann
Ear Forceps
14-597**
X = 8.5 cm
14-598
X = 14 cm




**Hartmann
Ear Forceps
14-599**
X = 8.5 cm
14-600
X = 14 cm
14-601
X = 16 cm



**Hartmann
Ear Forceps
14-602** 2 mm Ø
X = 8.5 cm
14-603 2.5 mm Ø
X = 8.5 cm
14-604 3 mm Ø
X = 8.5 cm



**Hartmann
Ear Forceps
14-606**



**Hartmann-Wullstein
Ear Forceps
14-606**
X = 8.5 cm



**Hartmann-Wullstein
Ear Forceps
14-607**
1 mm Ø
X = 8.5 cm



**Hartmann-Weingartner
Ear Forceps
14-608**
X = 8.5 cm



**Duckbill
Ear Forceps
14-609**




**Weingartner
Ear Forceps
14-610**



**Strümpel
Ear Forceps
14-611**



**Hartmann-hoffmann
Ear Forceps
14-612**



**Hartmann-Noyes
Ear Forceps
14-613**
X = 9 cm




**Struempel
Ear Forceps
14-614** 0°
X = 8.5 cm



**Struempel
Ear Forceps
14-615** 45°
X = 8.5 cm



**Struempel
Ear Forceps
14-616** 90°
X = 8.5 cm



BRONCHO OESOPHAGOSCOPY FORCEPS




MICRO INSTRUMENTS FOR LARYNGO BRONCHO OESOPHAGOSCOPY



Mikro Instrumente für Laryngo broncho Oesophagoskopie
Micro instruments pour laryngo-bronco-oesophagoscopie

micro instrumental para laringo bronco esofagoscopia
Micro strumenti per laringo bronco esofagoscopia



Grasping Broncho Oesophagocopy Forceps



14-662
Straight,
Serrated

14-663
Straight,
Fenestrated

14-664
Straight

Cup shaped Broncho Oesophagocopy Forceps O 3 mm



14-665
Straight

14-666
Angled up

14-667
Angled up,
Horizontal

14-668
Curved right

14-669
Curved left

Cup shaped Broncho Oesophagocopy Forceps




14-670
Straight

14-671
Angled up

Universal Broncho Oesophagocopy Forceps

14-672
Curved left



14-673
25.4cm

14-674
35.6cm

FIMA INSTRUMENTS

Yeoman Biopsy Forceps
with Rotating Shafts, straight
4 X 8 mm bite, with basket



14-675
X = 40 cm
14-676
X = 50 cm
14-677
X = 60 cm
Straight

14-678
X = 40 cm
14-679
X = 50 cm
14-680
X = 60 cm
Straight,
Serrated

14-681
X = 40 cm
14-682
X = 50 cm
14-683
X = 60 cm
Straight,
Fenestrated

14-684
X = 40 cm
14-685
X = 50 cm
14-686
X = 60 cm
Straight

Punch Broncho Oesophagocopy Forceps



**Cup shaped Broncho
Oesophagoscope Forceps**

14-687
X = 40 cm
4 mm
14-688
X = 50 cm
4 mm
14-689
X = 60 cm
4 mm
Straight



14-690
X = 40 cm
4 mm
14-691
X = 50 cm
4 mm
14-692
X = 60 cm
4 mm
Angled up



14-693
X = 40 cm
4 mm
14-694
X = 50 cm
4 mm
14-695
X = 60 cm
4 mm
Angled up,
horizontal



14-696
X = 40 cm
4 mm
14-697
X = 50 cm
4 mm
14-698
X = 60 cm
4 mm
Curved right



**Cup shaped Broncho
Oesophagoscope Forceps**

14-699
X = 40 cm
4 mm
14-700
X = 50 cm
4 mm
14-701
X = 60 cm
4 mm
Curved left



14-702
X = 40 cm
3 mm X 4 mm
14-703
X = 50 cm
3 mm X 4 mm
14-704
X = 60 cm
3 mm X 4 mm
Straight



14-705
X = 40 cm
4 mm X 6 mm
14-706
X = 50 cm
4 mm X 6 mm
14-707
X = 60 cm
4 mm X 6 mm
Straight



14-708
X = 40 cm
4 mm X 8 mm
14-709
X = 50 cm
4 mm X 8 mm
14-710
X = 60 cm
4 mm X 8 mm
Straight




14-711
X = 40 cm
14-712
X = 50 cm
14-713
X = 60 cm
Straight

LARYNGEAL POLYPUD FORCEPS



Kehlkopfpolypenzangen
Pincas a polypes laryngiens


Pinzas para polipos laringeos
Pinze per polipi laringei




**Fraenkel Laryngeal
Polypud Forceps**
14-714 18 cm



**Fraenkel Laryngeal
Polypud Forceps**
14-715
X = 17 cm
Both jaws articulated



**Jurasz Laryngeal
Polypud Forceps**



**Jurasz Laryngeal
Polypud Forceps**



14-716
23 cm
Right



14-717
23 cm
Left



14-720
23 cm
Right



14-721
23 cm
Left



14-718
23 cm
Downwards



14-719
23 cm
Upwards



14-722
23 cm
Downwards



14-723
23 cm
Upwards


INTESTINAL AND RECTAL INSTRUMENTS




INTESTINAL & TISSUE GRASPING FORCEPS

Darm-und Gewebefasszangen
Pinzas Para Intestinos y Tejidos

Pinces Pour Saisir Intestins et Tissus
Pinze Per Afferrare Intestino e Tessuti



Allis Tissue Forceps
15-001 15 cm
3 x 4 teeth




Allis Tissue Forceps
15-002 15.5 cm
4 x 5 teeth



Allis Tissue Forceps
15-003 16 cm
15-004 19 cm
15-005 22 cm
15-006 25.5 cm
5 x 6 teeth



Baby Allis Tissue Forceps
15-007 12 cm 4 X 5 teeth
15-008 14 cm 4 X 5 teeth
delicate jaws 4 mm wide



Judd-Allis Tissue Forceps
15-009 19.5 cm
3 x 4 teeth



Williams Grasp Forceps
15-010 16 cm



Boys-Allis Tissue Forceps
15-011 15.5 cm
5 x 6 teeth, sharp

INTESTINAL & TISSUE GRASPING FORCEPS

Darm-und Gewebefasszangen
Pinzas Para Intestinos y Tejidos

Pinces Pour Saisir Intestins et Tissus
Pinze Per Afferrare Intestino e Tessuti



Allis Tissue Grasping Forceps




**Allis-Adair
15-017**
Intestinal and Tissue Grasping Forceps
16 cm




**Thoms-Allis
15-018**
Intestinal and Tissue Grasping Forceps
20.5 cm




**Chaput Tissue Forceps
15-019** 12.5 cm
2 x 3 teeth, sharp



**Tuffier Tissue Forceps
15-020** 12.5 cm
5 x 6 teeth, sharp



**Duval-Collin Grasp Forceps
15-021** 18 cm
15-022 20 cm



**Duval Grasp Forceps
15-023** 23 cm

INTESTINAL & TISSUE GRASPING FORCEPS

Darm-und Gewebefasszangen
Pinzas Para Intestinos y Tejidos


Pinces Pour Saisir Intestins et Tissus
Pinze Per Afferrare Intestino e Tessuti



Babcock Seizing Forceps



15-024
16 cm
15-025
18 cm
15-026
20 cm
15-027
22 cm



Babcock Tissue Grasping Forceps



15-028 16 cm



15-029 21.5 cm



Babcock Tissue Grasping Forceps
15-030 18 cm
15-031 24 cm



Collin Grasp Forceps
15-032 14 cm



Collin Tissue Grasping Forceps
15-033 19.5 cm



Collin Tissue Grasping Forceps
15-034 23 cm

INTESTINAL & TISSUE GRASPING FORCEPS


Darm-und Gewebefasszangen | Pinces Pour Saisir Intestins et Tissus
Pinzas Para Intestinos y Tejidos | Pinze Per Afferrare Intestino e Tessuti



Duval Lung Tissue Forceps
15-035 1" WIDE
DeBakey type jaw, 23 cm



Lovelace Lung Tissue Forceps
15-036 1" WIDE
Angled on Side,
DeBakey type jaw, 19 cm



A
15-037
18 cm


B
15-038
20 cm

C
15-039
23 cm

Scheibe
Organ Grasping Forceps




Heywood-Smith
15-040
Polypus Grasping Forceps
21 cm



Duval Tissue Grasping Forceps

15-041
18 cm
15-042
20 cm

15-043
20 cm
15-044
23 cm



Lockwood Grasp Forceps
15-045 20 cm



Allen
Intestinal Clamp
15-046 16 cm

INTESTINAL & TISSUE GRASPING FORCEPS

Darm-und Gewebefasszangen
Pinzas Para Intestinos y Tejidos

Pinces Pour Saisir Intestins et Tissus
Pinze Per Afferrare Intestino e Tessuti



Little Wood Tissue Grasping Forceps
15-047 19 cm
Box Joint



Lane Tissue Grasping Forceps
15-048 12.5 cm
15-049 15 cm
15-050 19 cm



Duval Intestinal and Tissue Grasping Clamp Forceps



15-051
25 cm
15-052
27 cm



15-053
13 cm

15-054
18 cm
15-055
20 cm



Babcock Intestinal & Tissue Grasping Clamp Forceps T.C.




15-056
16.5 cm

15-057
20.3 cm



15-058
17.5 cm




15-059
19.5 cm



Stone-Watt
15-060
15cm
Closing Forceps



15-061
13 cm
Str.



Baby-Kocher Intestinal Forceps Clamps



15-062
13 cm
Cvd.


INTESTINAL CLAMPS

Darmklemmen

Pinces Intestinales

Pinzas Intestinales


Pinze Intestinali



15-063
16.5 cm
15-064
18 cm
Str.

15-065
16.5 cm
15-066
18 cm
Cvd.


Baby-Doyen Intestinal Forceps Clamps



15-067
21 cm
15-068
23 cm
15-069
25 cm
Str.


15-070
21 cm
15-071
23 cm
15-072
25 cm
Cvd.

Mayo-Robson Intestinal Forceps Clamps



15-073
22 cm
Straight


Kocher Intestinal Forceps Clamps



15-074
21 cm
15-075
23 cm
straight,
grooved

15-076
22 cm
Curved,
grooved

Doyen Intestinal Forceps Clamps



15-077 13 cm
15-078 15 cm
15-079 17 cm
15-080 21 cm
15-081 23 cm

Doyen-Baby Intestinal Clamp



Mayo Vessel Clamp
15-082 25.5cm
Curved, 3.8cm
tip

INTESTINAL CLAMPS

Darmklemmen | Pincas Intestinales
Pinzas Intestinales | Pinze Intestinali



15-083
25 cm
Str.

Mayo-Robson
Intestinal Clamp



15-084
25 cm
Cvd.




15-085
13 cm
15-086
15 cm
15-087
17 cm
15-088
21 cm
15-089
23 cm

Doyen-Baby
Intestinal Clamp



Hayes
Sigmoid Anastomosis Forceps

15-090
27 cm
15-091
28 cm




15-092
15 cm
15-093
18 cm
Str.

Bainbridge Intestinal
Forceps Clamps




15-094
15 cm
15-095
18 cm
Cvd.



15-096
22 cm
15-097
25 cm
15-098
28 cm
Str.

Kocher Intestinal Forceps
Clamps

15-099
22 cm
15-100
25 cm
15-101
28 cm
Cvd.



Doyen Intestinal
Forceps Clamps



15-102 13 cm
15-103 16 cm
15-104 18 cm
15-105 21 cm
15-106 23 cm
Str.

POLYOUS GRASPING FORCEPS, ANASTOMOSIS

Polypenfasszange, Darm-anastomosenklemmen

Pince Pour Saisir Polypes, Pincettes-clamps Pour Anastomose Intestinale

Pinza Para Coger Pólipos, Clamps Para Anastomosis Intestinales

Pinze Per Afferrare Polipi, Pinze Per Anastomosi Intestinali



**Doyen Intestinal
Forceps Clamps**


15-107
13 cm
15-108
16 cm
15-109
18 cm
15-110
21 cm
15-111
23 cm
Cvd.



**Scudder Intestinal Forceps
Clamps**

15-112
25 cm
15-113
33 cm
Cvd.

15-114
25 cm
15-115
33 cm
Str.



**Brunner Intestinal Polyous
Grasping Forceps**
15-116 24 cm



**Resano Intestinal Polyous
Grasping Forceps**
15-117 30 cm



**Fehland Intestinal Polyous
Grasping Forceps**
15-118 24 cm



**Nussbaum Intestinal
Forceps Clamps**

15-119
25 cm
Cvd.

15-120
25 cm
Str.



15-121
30cm
Str.



15-122
30cm
Cvd.



**Nussbaum Stomach
Intestinal
Twin Clamp Forceps**


INTESTINAL CLAMPS

Darmklemmen
Pinzas Intestinales


Pinces Intestinales
Pinze Intestinali




Intestinal & Stomach Twin Clamp Forceps



**Haberer Stomach Intestinal
Twin Clamp Forceps**
15-123 36cm



**Roosevelt Stomach Intestinal
Twin Clamp Forceps**
15-124 34cm



**Abadie Stomach Intestinal
Twin Clamp Forceps**
15-125 22cm
15-126 27cm



**Lane Stomach Intestinal
Clamp Forceps**
15-127
30 cm Str.



**Lane Stomach Intestinal
Clamp Forceps**
15-128
30 cm Cvd.




POLYOUS GRASPING FORCEPS, ANASTOMOSIS

Polypenfasszange, Darm-anastomosenklemmen

Pince Pour Saisir Polypes, Pincettes-clamps Pour Anastomose Intestinale

Pinza Para Coger Pólipos, Clamps Para Anastomosis Intestinales

Pinze Per Afferrare Polipi, Pinze Per Anastomosi Intestinali



POLYOUS GRASPING FORCEPS, ANASTOMOSIS



Polypenfasszange, Darm-anastomosenklemmen

Pince Pour Saisir Polypes, Pincés-clamps Pour Anastomose Intestinale

Pinza Para Coger Pólipos, Clamps Para Anastomosis Intestinales

Pinze Per Afferrare Polipi, Pinze Per Anastomosi Intestinali



**De Martel-Wolfson
Clamp Holder**
15-146 8.5 cm



**De Martel-Wolfson
Clamp Holder**
15-147 23 cm



**Stone-Watt
Intestinal Anastomosis Clamp**
15-148 8.5cm
15-149 10cm



De Martel-Wolfson
15-150 set of 3
anastomosis clamps




Jaw length
15-151
5 cm



Jaw length
15-152
7 cm



Jaw length
15-153
10 cm



**De Martel-Wolfson
Closing Forceps**
15-154
26 cm

INTESTINAL CLAMPS

Darmklemmen	Pinces Intestinales
Pinzas Intestinales	Pinze Intestinali



Payr-Baby Crushing Intestinal Clamp
15-155
 15 cm
 without pin
15-156
 15 cm
 with pin

Payr Crushing Intestinal Clamp
15-157
 21 cm
 without pin
15-158
 21 cm
 with pin

Payr Crushing Intestinal Clamp
15-159
 29 cm
 without pin
15-160
 29 cm
 with pin

Payr Crushing Intestinal Clamp
15-161
 35 cm
 without pin
15-162
 35 cm
 with pin

INTESTINAL RECTAL SPECULA

Mastdarmspekula
Especulos Rectales

Spécuments Rectaux
Specoli Rettali



90 x 15 mm

Mathieu Rectal Speculum
15-163 19 cm
90 x 15 mm



Sims Rectal Speculum
15-164 15 cm




Sims Rectal Speculum
15-165 15 cm



Ives-Fansler Rectal Speculum
15-166 14.5 cm
22 mm Ø



Pratt Rectal Speculum
15-167 21.5 cm
85 x 25 mm blade



Czerny Rectal Speculum
15-168
22 cm




Barr Rectal Speculum
15-169 80mm
Blade size



Smith-Blue Intestinal Rectal Speculum
15-172 15 cm
75 x 23 mm



15-170
for extra blades



15-173 22 cm
20 x 65 mm
15-174 22 cm
25 x 85 mm
15-175 22 cm
25 x 120 mm
15-176 22 cm
23 x 145 mm



Alan-Parks
15-171
Alan-Parks
14 cm



Schulze-Bergmann Intestinal Rectal Speculum

INTESTINAL RECTAL SPECULA

Mastdarmspekula
Espéculos Rectales

Spéculums Rectaux
Specoli Rettali



**Hirschmann Intestinal
Rectal Speculum**
15-177
23 mm Ø
proctoscope
working length: 145 mm



15-178
22 Ø x 50 mm
15-179
22 Ø x 100 mm
15-180
22 Ø x 140 mm
15-181
22 Ø x 200 mm

**Kelly Intestinal
Rectal Speculum**



**Barr-Shilford Intestinal
Rectal Speculum**
15-182



**Vernon-David Intestinal
Rectal Speculum**
with thumb grip handle
15-183 short handle
15-184 long (standard) handle



**Bensaude Intestinal
Rectal Speculum**
15-185 19 Ø x 70 mm
15-186 22 Ø x 70 mm



Brinkerhoff Intestinal Rectal Speculum
15-187 small 89 x 25mm tapers to 13mm
15-188 medium 101 x 25mm tapers to 15mm
15-189 large 120 x 31mm tapers to 19mm

INTESTINAL RECTAL SPECULA

Mastdarmspekula
Especulos Rectales

Spécuments Rectaux
Specoli Rettali



Fansler Operating Intestinal Rectal Speculum

with slotted / extra long handle
standard

15-190 25 mm x 65 mm

15-191 29 mm x 65 mm

15-192 38 mm x 70 mm

with slotted tube / short handle

15-193 29 mm x 60 mm

15-194 operating speculum



Chelsea-Eaton Intestinal Rectal Speculum

solid handle

15-195 small

15-196 large

with light connection available



Hirschman Anoscope Intestinal Rectal Speculum

15-197 small

57 x 14mm inside diameter

15-198 medium

52 x 17mm inside diameter

15-199 large

89 x 22mm inside diameter



Buie Hirschman Anoscope Intestinal Rectal Speculum

15-200 small

57 x 14mm inside diameter

15-201 medium

95 x 19mm inside diameter

15-202 large

98 x 25mm inside diameter



Ives-fansler Anoscope Intestinal Rectal Speculum

with hard rubber obturator /
fenestrated handle

15-203 22 mm x 85 mm



Long Chrome Handle Intestinal Rectal Speculum

15-204 81 x 29mm



Morgen Proctoscope Intestinal Rectal Speculum

with light connection

15-205 17mm

15-206 20mm

15-207 24mm



Kelly Proctoscope Intestinal Rectal Speculum

chrome plated

15-208 22mm x 50 mm

15-209 22mm x 100 mm

15-210 22mm x 145 mm


INTESTINAL & TISSUE GRASPING FORCEPS

Darm-und Gewebefasszangen
Pinzas Para Intestinos y Tejidos


Pinces Pour Saisir Intestins et Tissus
Pinze Per Afferrare Intestino e Tessuti



Hemorrhoidal Ligator-Suction
15-211




McGivney Hemorrhoidal Ligator
15-212
19cm



15-213
Rudd Hemorrhoidal Ligator with loading cone 12mmØ
31.5 cm Str.

15-214
Rudd Hemorrhoidal Ligator with loading cone 12mmØ
31.5 cm Cvd.



15-215
Lurz Hemorrhoidal Suction Ligator 8mm Ø
15-216
Lurz Hemorrhoidal Suction Ligator 10mm Ø

15-217
Ligatur rubber rings
in package of 100 pieces



15-218
McGivney Hemorrhoid Grasping Forceps
19 cm angled shank.


McGivney Hemorrhoidal Ligator
Complete Set
15-220

15-219
McGivney Hemorrhoidal Ligator
With Offset Handle 18 cm
Working Length.

INTESTINAL RECTAL SPECULA


Mastdarmspekula
Espéculos Rectales

Spéculum Rectaux
Specoli Rettali




Grasping forceps Flexible shaft
for grasping intra uterine pessary 3 mm
22 cm
15-221

Grasping forceps Flexible shaft
for grasping intra uterine pessary 3 mm
28.5 cm
15-222




Barren McGivney Hemorrhoidal Ligator forceps with stainless steel con
15-223

Barren McGivney Hemorrhoidal Ligator forceps with telon cone
15-224



McGivney Hemorrhoidal Ligator forceps with stainless steel con
15-225

Barren McGivney Hemorrhoidal Ligator forceps with telon cone
15-226



McGivney Hemorrhoidal Ligator forceps with stainless steel con
15-227 17.5 cm

Barren McGivney Hemorrhoidal Ligator forceps with telon cone
15-228 17.5 cm

UROLOGY INSTRUMENTS



KIDNEY STONE FORCEPS



FMA INSTRUMENTS

Nierensteinzangen


Pinces pour les calculs du bassinet

Pinzas para cálculoa renales


Pinze per calcoli renali




Randall Kidney Stone Clamps Forceps




16-001
23 cm
quarter curved



16-002
23 cm
half curved



16-003
23 cm
full curved



16-004
23 cm
three quarter curved

KIDNEY STONE FORCEPS

Nierensteinzangen

Pinces pour les calculs du bassin

Pinzas para cálculos renales

Pinze per calcoli renali



16-005
23 cm



16-006
24.5 cm



Randall
Kidney Stone Forceps



Heiming
Kidney Stone Forceps



16-007
16 cm



16-008
17 cm


B



Elsasser
Nephrostomy Forceps



16-009
21 cm
Fig. 1



16-010
21 cm
Fig. 2

KIDNEY STONE FORCEPS



Nierensteinzangen

Pinces pour les calculs du bassin

Pinzas para cálculos renales

Pinze per calcoli renali



**Dartigues Uterine
Elevating Forceps**
16-011 27 cm



**Collin Kidney Elevating
Stone Clamps Forceps**
16-012 28.5 cm



**Hulka-Kenwick Uterine
Elevating & Manipulating
Forceps**
16-013 27 cm




**Somer Kidney Elevating
Stone Clamps Forceps**
16-014 23.5 cm



**Strassmann Kidney Elevating
Stone Clamps Forceps**
16-015 26 cm



**Collin Kidney Elevating
Stone Clamps Forceps**
16-016 26.5 cm



**Collin Kidney Elevating
Stone Clamps Forceps**
16-017 26 cm

GALL STONE FORCEPS

Gallensteinzangen
Gallensteinlöffel

Pinces pour calculs biliaires
Curettes pour calculs biliaires



Blake Gall Stone Forceps
16-018 20 cm
Straight




Blake Gall Stone Forceps
16-019 21 cm
Curved



Mayo-Blake Gall Stone Forceps
16-020 20 cm
Straight



Mayo-Blake Gall Stone Forceps
16-021 20 cm
Curved



Uro-Tangential Gall Stone Forceps



Mixer Gall Stone Forceps
16-024 22 cm



Desjardins Gall Stone Forceps
16-025 23 cm



Desjardins Gall Stone Forceps
16-026 23 cm

CAPSULE FORCEPS, LIGATURE CARRYING FORCEPS



Kapselfasszangen, Fadenführungs-
Pinzas capsulares, Pinza porta ligaduras

Pinces pour saisir les capsules, Pince porte-fil
Pinze per afferrare capsule, Pinza porta-filo




Harris Capsule Ligature
Angled Suture Forceps



16-027
24 cm



Johnson Prostatic
Needle Holder



16-028
29 cm



Stratte Needle Holder



16-029
23 cm



Millin Prostatectomy
Forceps
16-030 24 cm



Millin Prostatectomy
Forceps
16-031 23 cm



Millin Vulsellum Forceps
16-032 23 cm



Millin Ligature Carrying
Forceps Angular
16-033 23 cm


PROSTATIC RETRACTOR

Prostatahaken


Separadores para próstata

Ecarteurs prostatiques


Divaricatori prostatici




Young Prostatic Retractor
16-034 21 cm
11 x 50 mm



Young Prostatic Retractor
16-035 22 cm
38 x 64 mm



Young Prostatic Retractor
16-036 21 cm
25 x 25 mm



Young Prostatic Retractor
16-037 22 cm
39 x 35 mm

GALL DUCT DILATOR



Gallenkanaldilatoren
Gallensteinzangen

Dilateurs pour les voies biliaires
Pincés pour calculs biliaires



**Desjardins Probe
& Catch with flexible
shaft 28 cm**



**16-038
17 F.G.**



**16-039
21 F.G.**



**16-040
24 F.G.**



**16-041
27 F.G.**



**Moynihan
16-042
34 cm**



**Bakes Common Bile
Duct Dilators, Malleable
16-043
Set of 13
32 cm**



**Fig. 1
16-044**



**Fig. 2
16-045**



**Fig. 3
16-046**



**Fig. 4
16-047**



**Fig. 5
16-048**



**Fig. 6
16-049**



**Fig. 7
16-050**



**Fig. 8
16-051**



**Fig. 9
16-052**



**Fig. 10
16-053**



**Fig. 11
16-054**



**Fig. 12
16-055**



**Fig. 13
16-056**



**Luer-Koerte
Modif With
Malleable Shaft
16-057
31 cm**



**16-058
2.8mm
Fig. 000**



**16-059
3.4mm
Fig. 00**



**16-060
4.3mm
Fig. 0**



**16-061
5.5mm
Fig. 1**



**16-062
6.7mm
Fig. 0**



**16-063
8.0mm
Fig. 3**



**16-064
11mm
Fig. 4**



**16-065
15mm
Fig. 5**

URETHRAL CATHETERS

Urethral catheter | Cateteri Uretrali
Sondes Urétrales | Catéteres Uretrales



Catheter
16-066
15 cm
No. 6-16



Catheters 29 cm
Male Catheters charr., Metal

French
16-067 8
16-068 10
16-069 12
16-070 14
16-071 16
16-072 18
16-073 20



Lister Bougies
16-074
27 cm
00 / 1 to 15 / 18
set of 17



Clutton Bougies
16-075 27 cm
06 / 10 FG to 28 / 32
set of 12



Van Buren
Urethral Sound
16-076
A to O
Set of 15






Fig.	Item	Description	Inches	Dia. Charr
A	16-077	28 cm	11"	8 French
B	16-078	28 cm	11"	10 French
C	16-079	28 cm	11"	12 French
D	16-080	28 cm	11"	14 French
E	16-081	28 cm	11"	16 French
F	16-082	28 cm	11"	18 French
G	16-083	28 cm	11"	20 French
H	16-084	28 cm	11"	22 French
I	16-085	28 cm	11"	24 French
J	16-086	28 cm	11"	26 French
K	16-087	28 cm	11"	28 French
L	16-088	28 cm	11"	30 French
M	16-089	28 cm	11"	32 French
N	16-090	28 cm	11"	34 French
O	16-091	28 cm	11"	36 French



Guyon
16-092
420 x 175 x 40 mm
Set in metal case containing:
— 21 Bougies Guyon
— 6 filiform Bougies

URETHRAL CATHETERS

Urethraalkatheter | Cateteri Uretrali
Sondes Urétrales | Catéteres Uretrales



Dittel
22.5 cm
urethral bougies and
dilating sounds

Charr. (French)	
16-093	8
16-094	9
16-095	10
16-096	11
16-097	12
16-098	13
16-099	14
16-100	15
16-101	16
16-102	17
16-103	18
16-104	19
16-105	20
16-106	21
16-107	22
16-108	23
16-109	24
16-110	25
16-111	26
16-112	27
16-113	28
16-114	29
16-115	30



Dittel
35 cm
urethral bougies and
dilating sounds

Charr. (French)	
16-116	8
16-117	9
16-118	10
16-119	11
16-120	12
16-121	13
16-122	14
16-123	15
16-124	16
16-125	17
16-126	18
16-127	19
16-128	20
16-129	21
16-130	22
16-131	23
16-132	24
16-133	25
16-134	26
16-135	27
16-136	28
16-137	29
16-138	30



33 cm
filiform guides,
olive tip

Charr. (French)	
16-139	3
16-140	4
16-141	5
16-142	6



Guyon
16-143
25 cm
Set
urethral bougies and
dilating sounds
without filiform guides,
with thread in tip

Charr. (French)	
16-144	Ben. 24 / Charr. 12 Ben. 25
16-145	Ben. 26 / Charr. 13 Ben. 27
16-146	Ben. 28 / Charr. 14 Ben. 29
16-147	Ben. 30 / Charr. 15 Ben. 31
16-148	Ben. 32 / Charr. 16 Ben. 33
16-149	Ben. 34 / Charr. 17 Ben. 35
16-150	Ben. 36 / Charr. 18 Ben. 37
16-151	Ben. 38 / Charr. 19 Ben. 39
16-152	Ben. 40 / Charr. 20 Ben. 41
16-153	Ben. 42 / Charr. 21 Ben. 43
16-154	Ben. 44 / Charr. 22 Ben. 45
16-155	Ben. 46 / Charr. 23 Ben. 47
16-156	Ben. 48 / Charr. 24 Ben. 49
16-157	Ben. 50 / Charr. 25 Ben. 51
16-158	Ben. 52 / Charr. 26 Ben. 53
16-159	Ben. 54 / Charr. 27 Ben. 55
16-160	Ben. 56 / Charr. 28 Ben. 57
16-161	Ben. 58 / Charr. 29 Ben. 59
16-162	Ben. 60 / Charr. 30

Blasenshalssperrer, Blasenwundsperrer

Ecarteur pour le col de la vessie, Ecarteurs vésicaux

Dilatador para el cuello vesical, Separadores vesicales

Divaricatore del collo vescicale, Divaricatore vescicale




Mogen Circumcision Clamp
16-165



Circumcision Shield
16-166 10cm



Winkelmann Circumcision Clamps
16-167 12 mm Ø
16-168 16 mm Ø
16-169 20 mm Ø



Gumco Circumcision Clamp
mm Ø
16-170 8
16-171 11
16-172 13
16-173 16
16-174 21
16-175 26
16-176 32



Stockmann Penile Clamp
16-177 8cm

BLADDER NECK SPREADER, BLADDER RETRACTORS

Blasenhalssperrer, Blasenwundsperrer

Ecarteur pour le col de la vessie, Ecarteurs vésicaux

Dilatador para el cuello vesical, Separadores vesicales

Divaricatore del collo vescicale, Divaricatore vescicale



Millin Bladder Neck Spreader
16-178
28 cm



Zieser Meaeus Clamp
16-179 10cm



Strauss Penile Clamp
16-180 11cm



Thomson Walker Penile Clamp
16-181



Lewkowitz Lithotomy Forceps
16-182 26cm Str.



Lewkowitz Lithotomy Forceps
16-183 26cm Cvd



Young Lobe Holding Forceps
16-184 20cm

BLADDER NECK SPREADER, BLADDER RETRACTORS

Blasenhalsperrerr, Blasenwundsperrerr

Dilatador para el cuello vesical, Separadores vesicales

Ecarteur pour le col de la vessie, Ecarteurs vésicaux

Divaricatore del collo vescicale, Divaricatore vescicale




Legueu Bladder Spatula
16-185



Gosset Bladder Retractor
16-186 160mm



Thomson Walker Bladder Retractor
16-187



Pantings Bladder Retractor
16-188



Judd Masson Bladder Retractor
16-189
w/center blade



Millin Bladder Retractor
16-190

HERNIA DIRECTOR, VASECTOMY FORCEPS, CLAMP & BLADDER RETRACTOR

Leistenbruch Director, Vasektomie Zange, Schraubstock & Blase Retractor
Directora Hernia, Pinzas vasectomía, abrazadera y la vejiga retractor

Directeur de hernie, Pincers de vasectomie, Clamp et la vessie enrouleur
Direttore Ernia, Pinze vasectomia, Clamp & vescica Divaricatore



**Child Hernia
Director**
16-191




**Key Hernia
Director**
16-192



**Teale Hernia
Director**
16-193



**Millin Suction
Tube**
16-194



Cord Hook
16-195
blunt 16cm



**Gutmann Secretions
Scoop**
16-196 26cm



**No Scalpel Vasectomy
Ringed Forceps Clamps**
16-197 14cm



**No Scalpel Vasectomy
Hemostat Forceps**
16-198 14cm
smooth sharp point curved



**Legueu Bladder
Retractor**
16-199
27cm



**Kocher's Bladder
Retractor**
16-200 23.5 cm

**RECTUM VASECTOMY SET OF 4 PCS
16-201**



**RECTUM VASECTOMY SET OF 4 PCS
16-202**



**RECTUM VASECTOMY SET OF 2 PCS
16-203**



**RECTUM VASECTOMY SET OF 3 PCS
16-204**




**RECTUM VASECTOMY SET OF 4 PCS
16-205**

GYNECOLOGY/OBSTETRICS INSTRUMENTS



Dilatation and Curettage



VAGINAL SPECULUM

Scheidenspekulum
Spéculum vaginal

Espéculo vaginal
Speculum vaginale



Cusco vaginal speculum
17-001 75 mm x 32 mm
17-002 85 mm x 35 mm
17-003 100 mm x 35 mm
17-004 75 mm x 17 mm
for vierges



Cusco vaginal speculum
17-005 75 mm x 32 mm
17-006 80 mm x 35 mm
17-007 100 mm x 37 mm



Cusco vaginal speculum
17-008 80 mm x 22 mm
17-009 100 mm x 25 mm
17-010 105 mm x 27 mm
Swiss pattern



Grave vaginal speculum
17-011 75 mm x 20 mm
17-012 95 mm x 35 mm
17-013 115 mm x 35 mm



Pederson vaginal speculum
17-014 75 mm x 13 mm
17-015 100 mm x 22 mm
17-016 120 mm x 25 mm



Winterton
17-017



Trelat vaginal speculum
17-018 85 mm x 30 mm
17-019 95 mm x 35 mm
17-020 115 mm x 43 mm



Vaginal Speculum Cusco
17-021 Extra Small
17-022 Small
17-023 Medium
17-024 Medium Extra long
17-025 Large



Semm vaginal speculum
17-026 90 mm x 14 / 18 mm
17-027 100 mm x 17 / 20 mm
17-028 110 mm x 25 / 30 mm

VAGINAL SPECULUM

Scheidenspekulum
Spéculum vaginal

Espéculo vaginal
Speculum vaginale



Collin vaginal speculum
17-029 85 mm x 30 mm
17-030 100 mm x 35 mm
17-031 110 mm x 40 mm



COLLIN
17-032
55 x 10 mm



NOTT
17-033 90 mm 30 mm
17-034 100 mm 30 mm



DE VILBISS
17-035
4 x 27 mm
102 mm x 32 mm
Medium



SEYFFERT
17-036 100 x 24 mm
17-037 110 x 32 mm
23 cm, 9"



CUSCO
17-038 75 x 32 mm
17-039 85 x 35 mm
17-040 100 x 35 mm



CUSCO
17-041 75 x 17 mm



CUSCO
17-042 90 x 22 mm
17-043 100 x 25 mm
17-044 105 x 27 mm



GRAVES
Virgin Speculaums
Color Coated, Straight
17-045 20mm x 75mm, Small
17-046 35mm x 95mm, Medium
17-047 35mm x 115mm, Large



Graves Vaginal Speculum Leep / Lietz
17-048 Small
17-049 Medium
17-050 Large
17-051 X-large



Pederson Vaginal Speculum Leep / Lietz
17-052 Small
17-053 Medium
17-054 Large
17-055 X-Large



PEDERSON
Virgin Speculaums
Color Coated
17-056 17mm x 100mm, Narrow
17-057 13mm x 75mm, Small
17-058 22mm x 100mm, Medium
17-059 25mm x 120mm, Large

VAGINAL RETRACTORS

Vaginal Retraktoren
Écarteurs vaginaux

Retractores vaginales
Divaricatori vaginali



NEW PORT
17-060
25 cm



**Lateral Vaginal
Retractor**
17-061
21 cm



**Lateral Vaginal Retractor
Blue Coated Leep / Lietz**
17-062 21 cm



Side Wall Retractors
17-063 4-Sided Pederson Speculum, Medium
17-064 4-Sided Pederson Speculum, Large
Redundant Wall Retractor, 3.5 in. x .75 in. blades
Lateral Wall Retractor, Closed Shanks, 2.25 in. blades
Lateral Wall Retractor, Open Shanks, 2.5 in. blades



**Lateral Vaginal Retractor
Blue Coated**
17-065 26cm



**McGee Lateral
Vaginal Retractor
With Ratchet Handle**
17-066 20cm

VAGINAL SPECULUM

Scheidenspekulum
Spéculum vaginal

Espéculo vaginal
Speculum vaginale



**Kogan Speculum
17-067**
24 cm
with scale-ratchet
and fixation screw



**Kogan Speculum
17-068**
24 cm



**Kogan - Modif.
Speculum
28 cm
17-069**
x=4 mm y=8mm
17-070
x=2 mm y=4mm



**Kogan Endocervical Spatula
17-071**
Ratchet, Standard Tips, 28cm
Ratchet, Narrow Tips, 28cm
Set Screw, Standard Tips, 28cm
Set Screw, Narrow Tips, 28cm



**RAMATHIBODI
17-072**
27.5 cm
Uterine Manipulator



**RAMATHIBODI
17-073**
25 cm
Tubal hook

SPECULUM RETRACTORS

Spekulumspreizer

Retractores de espéculo

Écarteurs de spéculum

Divaricatori per speculum



Kristeller Speculum
17-074 70 mm x 27 mm
17-075 80 mm x 30 mm
17-076 90 mm x 36 mm
17-077 110 mm x 27 mm
17-078 110 mm x 30 mm
17-079 110 mm x 36 mm



Kristeller Retractor
17-080 75 mm x 23 mm
17-081 85 mm x 26 mm
17-082 100 mm x 32 mm
17-083 115 mm x 23 mm
17-084 115 mm x 26 mm
17-085 115 mm x 32 mm



Kristeller Speculum & Retractor
17-086 75 mm x 23 mm
17-087 85 mm x 26 mm
17-088 100 mm x 32 mm
17-089 115 mm x 23 mm
17-090 115 mm x 26 mm
17-091 115 mm x 32 mm




Kristeller Speculum
17-092 105 mm x 18 mm
17-093 105 mm x 21 mm
 Speculum tapered, for virgins



Kristeller Retractor
17-094 110 mm x 14 mm
17-095 115 mm x 17 mm




Kristeller Speculum & Retractor
17-096 110 mm x 14 mm
17-097 115 mm x 17 mm



Jackson Retractor
17-098 75 mm x 35 mm
17-099 75 mm x 40 mm
17-100 90 mm x 35 mm
17-101 90 mm x 40 mm
17-102 100 mm x 35 mm
17-103 100 mm x 40 mm



Sims Speculum & Retractor
17-104 70 mm x 25 mm / 30 mm
17-105 75 mm x 30 mm / 35 mm
17-106 80 mm x 35 mm / 40 mm



BRIESKY
17-107
 72x25 mm / 77x30 mm
17-108
 75x30 mm / 92x38 mm
 17 cm

VAGINAL SPECULUM AND RETRACTORS



Vaginalspekulum und Retraktoren
Spéculum et écarteurs vaginaux

Espéculo vaginal y retractores
Speculum vaginale e divaricatori



Kallmorgen Speculum
17-109 70 mm x 40 mm
17-110 90 mm x 40 mm



Kallmorgen & Retractor
17-111 70 mm x 33 mm
17-112 90 mm x 33 mm



Kallmorgen Speculum & Retractor
17-113 70 mm x 33 mm
17-114 90 mm x 33 mm



MATHIEU Vaginal Speculum
17-115 Set 100x13mm-23mm
17-116 95x16mm-26mm
17-117 95x26mm-30 mm
21 cm



MATHIEU Vaginal Retractor
17-118 Set 100x17mm-30mm
17-119 100x11mm-20mm
17-120 100x14mm-26mm
22 cm



Doyen Vaginal Speculum
22 cm



Slightly Concave

- 17-121 60 x 35 mm
- 17-122 90 x 35 mm
- 17-123 120 x 35 mm
- 17-124 60 x 45 mm
- 17-125 90 x 45 mm
- 17-126 120 x 45 mm
- 17-127 160 x 45 mm
- 17-128 60 x 60 mm
- 17-129 90 x 60 mm
- 17-130 120 x 60 mm
- 17-131 160 x 60 mm



Landon Vaginal Retractor
17-132 Blade 9 cm
Deep 2.5 cm wide.



BRAUN
17-133
18 cm



BRAUN
17-134
18 cm



17-135
Handle with cold light
carrier for upper blades
28 cm



17-136
115 x 22 mm
17-137
115 x 26 mm

SIMON
115 x 22 mm
27.5 cm

VAGINAL SPECULUM AND RETRACTORS

Vaginalspekulum und Retraktoren
Spéculum et écarteurs vaginaux

Espéculo vaginal y retractores
Speculum vaginale e divaricatori



95 x 17-26 mm

**MATHIEU
17-138**
For children and
virgins fitting
21 cm



95 x 17 - 30 mm

**MATHIEU
17-139**
For children
and virgins
21 cm



LANDON
Vaginal Retractor
17-140
(2.5cm x 9cm x 20cm)
Colour coated



**Wertheim
Retractor
17-141**
(3cm x 9.5cm
x 27.5cm)



Landau
Speculum
29 cm
17-142
100x 28 mm
17-143
100x 38 mm



Breisky Speculum
30 cm
17-144
130x 20 mm
17-145
130 x 25 mm
17-146
130 x 30 mm
17-147
130 x 35 mm
17-148
130 x 40 mm
Vienna pattern



Breisky Speculum
30 cm
17-149
100 x 20 mm
17-150
100 x 25 mm
17-151
100 x 30 mm
17-152
100 x 35 mm
17-153
100 x 40 mm



Breisky Speculum
35 cm
17-154
160 x 40 mm
17-155
180 x 40 mm



Veronikis Vaginal
Retractors
17-156
(2.5 cm) blade width,
30cm length



Veronikis Vaginal
Retractors
17-157
(2.5 cm) blade width,
34 cm length



Veronikis Vaginal
Retractors
17-158
(3.8 cm) blade width,
30cm. length



Veronikis Vaginal
Retractors
17-159
(2.5 cm) Tapered blade,
30cm. length

VAGINAL SPECULUM



Scheidenspekulum
Spéculum vaginal

Espéculo vaginal
Speculum vaginale



**Scherback Speculum
17-160**
Set consisting of:
handle weight 5 blades



17-161
80 x 30 mm



17-162
85 x 35 mm



17-163
85 x 40 mm



17-164
55 x 40 - 60 mm



17-165
70 x 45 - 70 mm




**Scherback Speculum
17-166**
Set consisting of:
2 handle, 8 blades



17-167
80 x 30 mm




17-171
65 x 30 mm



17-168
85 x 33 mm



17-172
75 x 35 mm




17-169
100 x 38 mm



17-173
85 x 35 mm



17-170
110 x 42 mm



17-174
95 x 35 mm


SPECULUM RETRACTORS

Spekulumspreizer

Retractores de espéculo

Écarteurs de spéculum

Divaricatori per speculum



Doyen

speculum and retractor

24 cm slightly concave

17-175 55 mm x 35 mm

17-176 75 mm x 35 mm

17-177 115 mm x 35 mm

17-178 55 mm x 45 mm

17-179 75 mm x 45 mm

17-180 115 mm x 45 mm

17-181 160 mm x 45 mm

17-182 55 mm x 60 mm

17-183 75 mm x 60 mm

17-184 115 mm x 60 mm

17-185 160 mm x 60 mm



Doyen

speculum and retractor

24 cm medium concave

17-186 55 mm x 35 mm

17-187 75 mm x 35 mm

17-188 115 mm x 35 mm

17-189 55 mm x 45 mm

17-190 75 mm x 45 mm

17-191 115 mm x 45 mm

17-192 55 mm x 60 mm

17-193 75 mm x 60 mm

17-194 115 mm x 60 mm



DOYEN

Speculum & Retractor
23.5 cm

17-195
60 x 45 mm

17-196
90 x 45 mm

17-197
112 x 45 mm



**Garrigue Weighted
Vaginal Speculum**
17-198 4.0 cm x 8.0 cm
Weight 2 lbs



**Henrotin Weighted
Vaginal Speculum**
17-199 3.8 cm x 9.0 cm
Weight 3 lbs.



AUVARD,
speculum & retractor
17-200
80 x 38 mm - 28 mm
1,530 kg
with fixed weight
23 cm



AUVARD

speculum & retractor

17-201 75 x 42 mm

17-202 80- x 43 mm

(40 mm)

17-203 105 x 43 mm

(40 mm)

with detachable weight

23 cm



WEISSBARTH

speculum & retractor

17-204 80 x 28 mm

17-205 85 x 42 mm

21 cm

VAGINAL SPECULUM



Scheidenspekulum
Spéculum vaginal

Espéculo vaginal
Speculum vaginale



AUVARD
17-206
65 x 38 mm - 28 mm
1,530 kg
with fixed weight
23 cm



MARTIN
17-207
with hook shaped
handle for hanging
in weights
27 cm



STEINER SPECULUM
Ideal for deep vaginal surgery.
Blades may be ordered from
4" to 7", 5.5" is normal length.

17-208
1.1 Kg / 2.5 lb:
14 cm x 3 cm



HARDY-DUDDY SPECULUM
It will not fall out. Design places the 3.5 lbs
weight at a right angle to the shaft of the
speculum, preventing expulsion. Many
blade options are available.

17-209
Hardy-Duddy Speculum,
3.5 lbs., 120 mm x 35 mm
cupped blade



AUVARD SPECULUM
An assortment of weighted
Auvard specula is available.

17-210
1.1 Kg / 2.5 lb:
8.5 cm x 4 cm
17-211
1.1 Kg / 2.5 lb:
10 cm x 4.5 cm



**NICHOLS WEIGHTED
SPECULUM**
Permits retraction of posterior
vaginal wall, while providing
suction and measurement of
fluids from the surgical site.

17-212
Nichols Vaginal Weighted Speculum
with Suction, 2.5 lbs. 7cm x 3.8cm

VAGINAL RETRACTORS

Vaginal Retraktoren
Écarteurs vaginaux

Retractores vaginales
Divaricatori vaginali



O'Sullivan - O'Connor
Vaginal Retractor
17-213 Set of 3 Complete



17-214
45 x 17 mm



17-215
90 x 45 mm



17-216
90 x 37 mm



Abdominal Retractor
17-217 9 cm wide
17-218 10 cm wide
17-219 11.5 cm wide



Friedmann Retractor
17-220
10 cm



COMYN BERKELEY
Vulva Self Retaining

17-221
15.5cm



DOYEN
Supra Hollow Handle
17-222 7.5mm x 28cm
17-223 9mm x 28cm
17-224 10mm x 28cm
17-225 11.5mm x 28cm



DOYEN POZZI
Retractor
17-226 4.5mm x 5mm x 25.5cm
17-227 4.5mm x 9mm x 25.5cm
17-228 4.5mm x 11.5mm x 25.5cm

CERVICAL BIOPSY AND SPECIMEN FORCEPS

Probeexzisionszangen
Pincas emporte-pièces à biopsie

Pinzas saca-bocados para biopsia
Pinze taglienti per biopsia



Berger
Biopsy and Specimen Forceps
17-229
x = 200 mm



Eppendorfer
Biopsy and Specimen Forceps
17-230
x = 200 mm




Fig. 1
17-231
24 cm

Fig. 2
17-232
24 cm


Schumacher
Biopsy Forceps



Kevorkian
Biopsy and Specimen Forceps
17-233
x = 245 mm



Alexander
Biopsy and Specimen Forceps
17-234 28 cm
x = 200 mm




Wittner Biopsy and Specimen Forceps
17-235
x = 230 mm
straight

Wittner Biopsy and Specimen Forceps
17-236
x = 230 mm
curved



Schumacher
Biopsy and Specimen Forceps
17-237
x = 23 cm
complete

17-238
x = 23 cm
tip




Tischler - Morgan
Biopsy and Specimen Forceps
USA Modell
17-239
x = 200 mm
17-240
x = 250 mm

Tischler - Morgan
Biopsy and Specimen Forceps
17-241
x = 200 mm
17-242
x = 250 mm

Instruments as illustrated,
but without locking screw.


USA Modell



Mini - Tischler
Biopsy and Specimen Forceps
17-243
x = 200 mm
17-244
x = 250 mm

Mini - Tischler
Biopsy and Specimen Forceps
17-245
x = 200 mm
17-246
x = 250 mm

Instruments as illustrated,
but without locking screw.



Wittner Biopsy Forceps
17-247
23cm
Shaft Straight

Wittner Biopsy Forceps
17-248
23cm
Shaft Angled



Schumacher Biopsy Forceps
17-249
Small Jaw 24cm Shaft

Schumacher Biopsy Forceps
17-250
Large Jaw 24cm Shaft



Tischler Baby Biopsy Forceps
17-251
21.5cm Shaft


CERVICAL BIOPSY AND SPECIMEN FORCEPS

Probeexzisionszangen

Pinces emporte-pièces à biopsie

Pinzas saca-bocados para biopsia


Pinze taglienti per biopsia



Tischler Biopsy Forceps
17-252
21.5cm Shaft




Burke Biopsy Forceps
17-253
21.5cm Shaft



Eppendorf Biopsy Forceps
17-254
20cm Shaft



Kevorkian Young Biopsy Forceps
17-255
25.5cm Shaft



Leechwilkinson Biopsy Forceps
17-256
20cm Shaft Round



CERVICAL BIOPSY AND SPECIMEN FORCEPS

Probeexzisionszangen
Pincas emporte-pièces à biopsie

Pinzas saca-bocados para biopsia
Pinze taglienti per biopsia



Uni - Tischler
Biopsy and
Specimen Forceps
17-258
3 x 7 mm
Straight, Tip only



Uni - Tischler
Biopsy and
Specimen Forceps
17-259
3 x 7 mm
up, Tip only



Uni - Tischler
Biopsy and
Specimen Forceps
17-260
3 x 7 mm
down, Tip only



**Biopsy and Specimen
Forceps Handle only**
17-257



Uni - Townsend
Biopsy and
Specimen Forceps
17-261
2.3 x 4 mm
straight, Tip only



Uni - Townsend
Biopsy and
Specimen Forceps
17-262
2.3 x 4 mm
up, Tip only



Uni - Townsend
Biopsy and
Specimen Forceps
17-263
2.3 x 4 mm
down, Tip only



Uni - Kevorkian
Biopsy and
Specimen Forceps
17-264
3 x 9 mm
Tip only



Uni - Oval
Biopsy and
Specimen Forceps
17-265
3 x 5 mm
Tip only



Uni - Burke
Biopsy and
Specimen Forceps
17-266
3 x 5 mm
Tip only



Uni - Tischler
Biopsy and
Specimen Forceps
17-268
3 x 7 mm
Straight, Tip only



Uni - Tischler
Biopsy and
Specimen Forceps
17-269
3 x 7 mm
up, Tip only



Uni - Tischler
Biopsy and
Specimen Forceps
17-270
3 x 7 mm
down, Tip only



**Biopsy and Specimen
Forceps Handle only**
17-267



Uni - Townsend
Biopsy and
Specimen Forceps
17-271
2.3 x 4 mm
straight, Tip only



Uni - Townsend
Biopsy and
Specimen Forceps
17-272
2.3 x 4 mm
up, Tip only



Uni - Townsend
Biopsy and
Specimen Forceps
17-273
2.3 x 4 mm
down, Tip only

Uni - Kevorkian
Biopsy and
Specimen Forceps
17-274
3 x 9 mm
Tip only



Uni - Oval
Biopsy and
Specimen Forceps
17-275
3 x 5 mm
Tip only



Uni - Burke
Biopsy and
Specimen Forceps
17-276
3 x 5 mm
Tip only


POLYOUS GRASPING FORCEPS, ANASTOMOSIS

Polypenfasszange, Darm-anastomosenklemmen

Pince Pour Saisir Polypes, Pincés-clamps Pour Anastomose Intestinale

Pinza Para Coger Pólipos, Clamps Para Anastomosis Intestinales


Pinze Per Afferrare Polipi, Pinze Per Anastomosi Intestinali



YEOMANN

PILE FORCEPS

COTTON CARRIER



CERVICAL BIOPSY AND SPECIMEN FORCEPS

Probeexzisionszangen

Pinces emporte-pièces à biopsie

Pinzas saca-bocados para biopsia

Pinze taglienti per biopsia



**Schubert
Biopsy Forceps**
17-293 21 cm
17-294 28 cm



**Faure Uterine
Biopsy Forceps**
17-295 21 cm



**Schubert
Biopsy Forceps**
17-296 26 cm



**Gellhorn
Biopsy Forceps**
17-297 23 cm



**Thoms-Gaylor
Biopsy Forceps**
17-298 22 cm
17-299 24 cm



**Van Doren
Biopsy Forceps**
17-300 24 cm



**Schumacher
Biopsy Forceps**

Fig. 1
17-301
24 cm

Fig. 2
17-302
24 cm

TENACULUM FORCEPS

Hakenzangen

Pinces pour l'abaissement du col de l'utérus

Pinzas para coger el cuello del útero

Pinze per abbassare il collo dell'utero



Museux Tenaculum Forceps



24 cm
17-303
8 mm
17-304
6 mm
17-305
10 mm
Straight

24 cm
17-306
8 mm
17-307
6 mm
17-308
10 mm
Curved



Jacobs Uterine Vulsellum Forceps

17-309
21 cm



17-310
21 cm



WINTER
17-311 Fig. 1
17-312 Fig. 2
17-313 Fig. 3
28 cm



Teale Tenaculum Forceps
17-314 23 cm



Pratt Tenaculum Forceps
17-315 26 cm



Czerny Tenaculum Forceps
17-316 20.5 cm



Schroeder Vulsellum Tenaculum Forceps



17-317 24 cm
17-318 24 cm

TENACULUM FORCEPS



Hakenzangen

Pinzas para coger el cuello del útero

Pinces pour l'abaissement du col de l'utérus

Pinze per abbassare il collo dell'utero



Pozzi Tenaculum Forceps
17-319 25.5 cm




Schroeder Tenaculum Forceps
17-320 25 cm



Duplay Tenaculum Forceps
17-321 28 cm



Iowa Tenaculum Forceps
17-322 26 cm



Skene Tenaculum Forceps
17-323 24 cm



Beacham Amniotome
17-324 26 cm




Staude - Moore Tenaculum Forceps
17-325 17 cm



Barrett Tenaculum Forceps
17-326 18 cm



Braun Tenaculum Forceps Sharp
17-327 25 cm



Schroeder-Braun Tenaculum Forceps
17-328 25 cm sharp

17-329 blunt




17-330
20 cm
1 x 2th, Angled
Screw Joint

BONNEY



17-331
20 cm
1 x 2th, Curved
Screw Joint

BONNEY-BERKLEY



17-332
20 cm
Straight
Box Joint

17-333
20 cm
Curved
Box Joint


GWILLIAM



17-334
20 cm
Curved
1 tooth
Box Joint


17-335
20 cm
Curved
2 Teeth
Box Joint

HEANEY



17-336
23 cm
Curved
Box Joint

HUNTER



17-337
20 cm
Straight
Box Joint


17-338
20 cm
Curved
Box Joint

MAINGOT

PARAMETRIUM FORCEPS, HYSTERECTOMY FORCEPS

Parametriumklemmen, Hysterektomieklemmen
Pincas parametrium, Pincas à hystérectomie


Pinzas parametrium, Pinzas-clamps para hysterectomia
Pinze parametrium, Pinze per isterectomia



Parametrium
Hysterectomy
Forceps




17-339
21 cm
17-340
24 cm
17-341
30 cm
17-342
35 cm




17-343
21 cm
17-344
24 cm
17-345
30 cm
17-346
35 cm



17-347
21 cm
17-348
24 cm
17-349
30 cm
17-350
35 cm




17-351
21 cm
17-352
24 cm
17-353
30 cm
17-354
35 cm




17-355
21 cm
17-356
24 cm
17-357
30 cm
17-358
35 cm
S-curved



Rogers Parametrium
Hysterectomy Forceps




17-359
Fig. 1
22 cm




17-360
Fig. 2
21.5 cm




17-361
Fig. 3
21 cm




17-362
Fig. 4
20.5 cm




17-363
Fig. 1
25 cm



17-364
Fig. 2
25 cm




17-365
Fig. 3
24.5 cm



17-366
Fig. 4
23 cm


UTERINE DILATORS

Uterusdilatatoren | Dilataores uterinos
Dilateurs utérins | Dilatori uterini



Hank Uterine Dilators 6 Set
17-367 Complete Set
Curved, Delicate, 11-12 Fr. to 21-22 Fr.
(3.6mm to 7.33cm)

Hank Uterine Dilators, Double-ended
17-368 11-12 Fr. (3.6mm x 4.0mm), 27.5cm
17-369 13-14 Fr. (4.3mm x 4.7mm), 27.5cm
17-370 15-16 Fr. (5.0mm x 5.3mm), 27.5cm
17-371 17-18 Fr. (5.7mm x 6.0mm), 27.5cm
17-372 19-20 Fr. (6.3mm x 6.7mm), 27.5cm
17-373 21-22 Fr. (7.0mm x 7.3mm), 27.5cm



Pratt Uterine Dilators 8/Set
17-374 Complete Set
Curved, Delicate, 13-15 Fr. to
41-43 Fr., Straight

Pratt Uterine Dilators, Double-ended
17-375 13 / 15 Fr., 29.5cm
17-376 17 / 19 Fr., 29.5cm
17-377 21 / 23 Fr., 29.5cm
17-378 25 / 27 Fr., 29.5cm
17-379 29 / 31 Fr., 29.5cm
17-380 33 / 35 Fr., 29.5cm
17-381 37 / 39 Fr., 29.5cm
17-382 41 / 43 Fr., 29.5cm



Hegar uterine dilator
17-383 3 / 4 mm - 17 / 18 mm Ø
8 Pieces per set
17-384 1 / 2 mm - 25 / 26 mm Ø
13 Pieces per set

Hegar uterine dilator
17-385 1 / 2 mm Ø
17-386 3 / 4 mm Ø
17-387 5 / 6 mm Ø
17-388 7 / 8 mm Ø
17-389 9 / 10 mm Ø
17-390 11 / 12 mm Ø
17-391 13 / 14 mm Ø
17-392 15 / 16 mm Ø
17-393 17 / 18 mm Ø
17-394 19 / 20 mm Ø
17-395 21 / 22 mm Ø
17-396 23 / 24 mm Ø
17-397 25 / 26 mm Ø
double ended




Hegar Uterine Dilators 8 Set W/Roll
17-398 Complete Set
Curved, Delicate, 3mm to 18mm, Curved

Hegar Uterine Dilator, Double-ended
17-399 3mm & 4mm, (20cm)
17-400 5mm & 6mm, (20cm)
17-401 7mm & 8mm, (20cm)
17-402 9mm & 10mm, (20cm)
17-403 11mm & 12mm, (20cm)
17-404 13mm & 14mm, (20cm)
17-405 15mm & 16mm, (20cm)
17-406 17mm & 18mm, (20cm)

UTERINE DILATORS



Uterusdilatoren | Dilatadores uterinos
Dilatateurs utérins | Dilatatori uterini




17-408 18.5cm 3-4-5mm
17-409 18.5cm 4-5-6mm
17-410 18.5cm 5-6-7mm
17-411 18.5cm 6-7-8mm
17-412 18.5cm 7-8-9mm
17-413 18.5cm 8-9-10mm
17-414 18.5cm 9-10-11mm
17-415 18.5cm 10-11-12mm
17-416 18.5cm 11-12-13mm
17-417 18.5cm 12-13-14mm
17-418 18.5cm 13-14-15mm
17-419 18.5cm 14-15-16mm
17-420 18.5cm 15-16-17mm
17-421 18.5cm 16-17-18mm
17-422 18.5cm 17-18-19mm
17-423 18.5cm 18-19-20mm
Canvas Roll

Bonney-Barker Dilators
17-407 Set of 16 Pcs



17-425 27cm 9/10mm
17-426 27cm 11/12mm
17-427 27cm 13/14mm
17-428 27cm 15/16mm
17-429 27cm 17/18mm
17-430 27cm 19/20mm

Hank Double Ended
17-424 Set of 6 Pcs



17-432 17cm 3-6mm
17-433 17cm 4-7mm
17-434 17cm 5-8mm
17-435 17cm 6-9mm
17-436 17cm 7-10mm
17-437 17cm 8-11mm
17-438 17cm 9-12mm
17-439 17cm 10-13mm
17-440 17cm 11-14mm
17-441 17cm 12-15mm
17-442 17cm 13-16mm
17-443 17cm 14-17mm
17-444 17cm 15-18mm
17-445 17cm 16-19mm
17-446 17cm 17-20mm
17-447 17cm 18-21mm
Canvas Roll Only

Hawkin Ambler Dilators
17-431
Set of 16 Pcs



20cm 1/2mm
20cm 3/4mm
20cm 5/6mm
20cm 7/8mm
20cm 9/10mm
20cm 11/12mm
20cm 13/14mm
20cm 15/16mm
20cm 17/18mm
20cm 19/20mm
20cm 21/22mm
20cm 23/24mm
20cm 25/26mm
17-448 Set 8 Double Ended 3/4mm-17/18mm
17-449 Set 12 Double Ended 1/2mm-23/24mm
17-450 Set 8 Double Ended Canvas Only
17-451 Set 12 Double Ended Canvas Only

Hegar Double Ended Dilators

UTERINE DILATORS

Uterusdilatoren | Dilataores uterinos
Dilataeurs utérins | Dilatatori uterini



Hegar
with siped handle

Set	mm	pieces per set
17-452	4-17	14
17-453	1-26	26

17-454	1mm Ø	17-484	1.5mm Ø
17-455	2mm Ø	17-485	2.5mm Ø
17-456	3mm Ø	17-486	3.5mm Ø
17-457	4mm Ø	17-487	4.5mm Ø
17-458	5mm Ø	17-488	5.5mm Ø
17-459	6mm Ø	17-489	6.5mm Ø
17-460	7mm Ø	17-490	7.5mm Ø
17-461	8mm Ø	17-491	8.5mm Ø
17-462	9mm Ø	17-492	9.5mm Ø
17-463	10mm Ø	17-493	10.5mm Ø
17-464	11mm Ø	17-494	11.5mm Ø
17-465	12mm Ø	17-495	12.5mm Ø
17-466	13mm Ø	17-496	13.5mm Ø
17-467	14mm Ø	17-497	14.5mm Ø
17-468	15mm Ø	17-498	15.5mm Ø
17-469	16mm Ø	17-499	16.5mm Ø
17-470	17mm Ø	17-500	17.5mm Ø
17-471	18mm Ø	17-501	18.5mm Ø
17-472	19mm Ø	17-502	19.5mm Ø
17-473	20mm Ø	17-503	20.5mm Ø
17-474	21mm Ø	17-504	21.5mm Ø
17-475	22mm Ø	17-505	22.5mm Ø
17-476	23mm Ø	17-506	23.5mm Ø
17-477	24mm Ø		
17-478	25mm Ø		
17-479	26mm Ø		
17-480	27mm Ø		
17-481	28mm Ø		
17-482	29mm Ø		
17-483	30mm Ø		



Hegar Uterine Dilator
single ended
17-507 4 mm - 17 mm Ø
8 Pieces per set
17-508 1 mm - 26 mm Ø
24 Pieces per set



Hegar Uterine Dilator
single ended
17-509 4 mm - 17 mm Ø
24 pieces per set

RETRACTORS AND SALPINGOGRAPH



Retraktoren und Salpingograph
Écarteurs et salpingographe

Retractoires y Salpingografía
Divaricatori e salpingografi



GELPI
Vulva Self Retaining



DOYEN
Vaginal Flat Retractor
17-511
1.75mm x 10mm x 23cm



JOLL Thyroid Self Retaining
17-512 15.5cm



Forceps for I.U.D Removing
17-513
4 mm Ø x=230 mm



17-514
IUP Grasping Forceps 3 mm
22 cm



Schultze, Salpingograph
17-515 29.5 cm
complete in wooden case




CATHERES RETRACTORS & DRESSING FORCEPS

Catheres Retraktoren & Verbandzangen
Écarteurs et pinces à pansements Catheres

Retradores de catéteres y pinzas de apósito
Divaricatori Catheres e pinze per medicazione



NEEDLE EXTENSIONS



UTERINE SOUNDS



Uterus-Sonden

Histerómetros uterinos

Hystéromètres du col de l'utérus

Isterometri uterini



Sims Uterine Sound
17-538 32 cm
right, graduated in cm
17-539 32 cm
malleable graduated in cm



Sims Uterine Sound
17-540 33 cm
right, malleable



Sims Uterine Sound
17-541 32.5cm
Non-graduated



Sims Uterine Sound
17-542 32.5cm
Graduated in cm



Sims Uterine Sound
17-543 24cm
Graduated in cm



Sims-Special Uterine Sound
17-544 33 cm



Sims Uterine Sound
schwed. Modell
17-545 26 cm



Martin Uterine Sound
17-546 30 cm



Valleix Uterine Sound
17-547 25 cm



Mayo Uterine Sound
17-548 33 cm
biegsam malleable



Simpson Uterine Sound
17-549 29 cm
malleable silver-plated



Simpson Uterine Sound
17-550 32.5cm
Graduated in inches

HOOKLET, RETRACTORS

Wundhäkchen, Wundhaken
Ganchito, Separadores

Crochet, Ecarteurs
Divaricatore, Divaricatori



SMELLIE
17-551
30 cm



Braun, Midwifery
Decapitating hook
17-552
32.5 cm



Suture Retrieval
Hook
17-553
27cm



Suture Retrieval
Hook
Sims Handle
17-554
27 cm



17-555
18cm
Fine
Colour
Coated

IRIS
Skin Hooks



IRIS
Skin Hooks

17-556
25.5cm
Fine
Colour
Coated



TENACULUM EMMETT



17-557
23cm
bic



TENACULUM 3 PRONG



17-558
26cm
bic



WHITTAKER



17-559
18.5cm
bic

MYOMA INSTRUMENTS, FISTULA HOOKS



Myom-Instrumente, Fistelhaken


Instrumentos para fibroma, Ganchos para fistulas

Instruments pour myomes, Crochets pour fistules

Strumenti per miomi, Uncini per fistole



**Second
Myoma Knife
Fistula Hook
17-560**
27 cm



**Ayre
Cone Knife
Fistula Hook
17-561**
23 cm



**Doyen
Myoma Screw
And Fistula Hook
17-562 14.5cm
17-563 17 cm**



**Fig. 1
17-564**
22 cm

**Simon
Fistula Hook,
One Prong**



**Fig. 2
17-565**
22 cm

**Simon
Fistula Hook,
Two Prong**



**Emmet
Fistula Hook**



**Fig. 1
17-566**
22 cm



**Fig. 2
17-567**
22 cm



**Fig. 3
17-568**
22 cm



**Fig. 4
17-569**
22 cm



**Fig. 5
17-570**
22 cm



1:1



1:2

Uni - Randall
aspirating curette
17-571 4 mm Ø
24 cm



Novak
aspirating curette
17-576 4 mm Ø
23 cm

Mini - Novak
aspirating curette
23.5 cm
17-577 1 mm Ø
17-578 2 mm Ø
17-579 3 mm Ø



1:1

Uni - Randall
aspirating curette
17-572
30 cm
rigid, without basket

Uni - Randall
aspirating curette
17-573
30 cm
malleable,
without basket



1:1

Uni - Randall
aspirating curette
17-574
30 cm
rigid, with basket

Uni - Randall
aspirating curette
17-575
30 cm
malleable,
with basket



Specimen And
Aspirating Curette
20 cm
17-580 3 mm Ø
17-581 4 mm Ø
17-582 6 mm Ø
17-583 8 mm Ø
17-584 10 mm Ø
17-585 12 mm Ø
17-586 14 mm Ø

Coupling Piece
aspirating curette
12 cm
17-587



FLUSHING CURETTES




Küretten spülen

Curetas de lavado

Curettes de rinçage

Curette di lavaggio



17-588
Small 29cm
8mm Sharp



17-589
Medium 29cm
11mm Sharp



17-590
Large 29cm
14mm Sharp

**Rheinstadter
Flushing Curettes**



17-591
Small 29cm
8mm Blunt



17-592
Medium 29cm
11mm Blunt



17-593
Large 29cm
14mm Blunt

**Rheinstadter
Flushing Curettes**

8.0mm




17-594
12 FG 4mm
25.5cm

**Leechwilkinson
Suction Biopsy Curettes**




- 1 • **17-595**
- 2 • **17-596**
- 3 • **17-597**
- 4 • **17-598**

**Novak Suction
Biopsy Curettes**



1mm 25.5cm
2mm 25.5cm
3mm 25.5cm
4mm 25.5cm



2mm 25.5cm
3mm 25.5cm
4mm 25.5cm

- 2 • **17-599**
- 3 • **17-600**
- 4 • **17-601**

**Randall Suction
Biopsy Curettes**



**Sharman Endometrial
Biopsy Curettes**



17-602
10FG
23cm C/P

10FG 3.3mm



SCHROEDER
30 cm, 12"
sharp, rigid



17-603



17-604



17-605



17-606



17-607



17-608

KEVORKIAN



17-609
30 cm
1.2mm

Without Basket



17-610
30 cm
4.0 mm
Without Basket



17-611
30 cm
4.0 mm
With Basket



17-612
30.5cm


**Kevorkian
Endocervical
Curettes**
With Basket Round

Uteruslöffel, Absaugküretten, uterusdepressor

Curettes uterines, curettes d'aspiration, Dépresseur utérin

Cucharillas uterinas, Curetas aspiradoras, depresor uterino

Curette uterine, curette pe aspirazione, depressore uterino



17-613
30.5cm



17-614
30.5cm



Braun Uterine Depressor
17-615 28 cm



Blake Double Ended
17-616 24.5 cm
Sharp / Blunt



Zweifel Double Ended
17-617 24.5cm
Sharp / Sharp

Kevorkian
Endocervical
Curettes
With Basket Square

Kevorkian
Endocervical
Curettes
Without Basket Square



Coakley
17cm



17-618
Fig. 1



17-619
Fig. 2



17-620
Fig. 3



17-621
Fig. 4



17-622
Fig. 5



17-623
Fig. 6

UTERINE CURETTES



Uterusküretten	Curettes uterinas
Curettes utérines	Curette uterine



**Sims
Single Ended
Sharp**



17-624
Small 6mm
Sharp 25.5cm



17-625
Medium 9mm
Sharp 25.5cm



17-626
Large 12mm
Sharp 25.5cm



17-627
Small 6mm
Blunt 25.5cm




17-628
Medium 9mm
Blunt 25.5cm



17-629
Large 12mm
Blunt 25.5cm



**Sims
Single Ended
S/S**



17-630
Small 6mm X 5mm
Sharp
Sharp 25.5cm



17-631
Medium
9mm X 6mm Sharp
Sharp 25.5cm




17-632
Large
12mm X 9mm Sharp
Sharp 25.5cm



**Blake
17-633
27 cm**



Blunt



**Sims
17-634
26 cm**




**Sims
Double Ended
B/S**



17-635
Small
6mm X 5mm Sharp
Blunt 25.5mm



17-636
Medium
9mm X 6mm Sharp
Blunt 25.5mm



17-637
Large
12mm X 9mm Sharp
Blunt 25.5mm



**Sims
17-638
Double Ended
Sharp Blunt**

UTERINE CURETTES

Uterusküretten | Curettes uterinas
Curettes utérines | Curette uterine



Recamier
uterine curette
y=26 cm

Fig.	x=mm	Sharp, malleable	Blunt, malleable	Sharp, rigid	Blunt, rigid
00	5	17-639	17-653	17-667	17-681
0	6	17-640	17-654	17-668	17-682
1	7	17-641	17-655	17-669	17-683
2	8	17-642	17-656	17-670	17-684
3	9	17-643	17-657	17-671	17-685
4	11	17-644	17-658	17-672	17-686
5	12	17-645	17-659	17-673	17-687
6	14	17-646	17-660	17-674	17-688
7	15	17-647	17-661	17-675	17-689
8	16	17-648	17-662	17-676	17-690
9	19	17-649	17-663	17-677	17-691
10	20	17-650	17-664	17-678	17-692
11	21	17-651	17-665	17-679	17-693
12	23	17-652	17-666	17-680	17-694

Recamier
uterine curette
y=31 cm

Fig.	x=mm	Sharp, malleable	Blunt, malleable	Sharp, rigid	Blunt, rigid
00	5	17-695	17-709	17-723	17-737
0	6	17-696	17-710	17-724	17-738
1	7	17-697	17-711	17-725	17-739
2	8	17-698	17-712	17-726	17-740
3	9	17-699	17-713	17-727	17-741
4	11	17-700	17-714	17-728	17-742
5	12	17-701	17-715	17-729	17-743
6	14	17-702	17-716	17-730	17-744
7	15	17-703	17-717	17-731	17-745
8	16	17-704	17-718	17-732	17-746
9	19	17-705	17-719	17-733	17-747
10	20	17-706	17-720	17-734	17-748
11	21	17-707	17-721	17-735	17-749
12	23	17-708	17-722	17-736	17-750

UTERINE CURETTES



Uterusküretten | Curettes uterinas
Curettes utérines | Curette uterine



Recamier
uterine curette
y=34 cm



Fig. 13



Fig. 14



Fig. 15



Fig. 16

Fig.	x=mm	Sharp, malleable	Blunt, malleable	Sharp, rigid	Blunt, rigid
13	25	17-751	17-755	17-759	17-763
14	30	17-752	17-756	17-760	17-764
15	35	17-753	17-757	17-761	17-765
16	40	17-754	17-758	17-762	17-766



Sims
uterine curette
28 cm



Fig. 1



Fig. 2



Fig. 3



Fig. 4



Fig. 5



Fig. 6

Fig.	x=mm	Sharp, malleable	Sharp, rigid	Blunt, malleable	Blunt, rigid
1	7	17-767	17-773	17-779	17-785
2	8	17-768	17-774	17-780	17-786
3	9	17-769	17-775	17-781	17-787
4	11	17-770	17-776	17-782	17-788
5	12	17-771	17-777	17-783	17-789
6	14	17-772	17-778	17-784	17-790

UTERINE SPOONS

Scharfe Löffel
Curettes tranchantes

Cucharillas cortantes
Cucchiai tahlenti



Schröder
Uterine Spoon



Fig. 1
17-791
25 cm
Sharp, rigid



Fig. 2
17-792
25 cm
Sharp, rigid



Fig. 3
17-793
25 cm
Sharp, rigid



Fig. 4
17-794
25 cm
Sharp, rigid



Fig. 5
17-795
25 cm
Sharp, rigid



Fig. 6
17-796
25 cm
Sharp, rigid



Schröder
Uterine Spoon



Fig. 1
17-797
30 cm



Fig. 2
17-798
30 cm



Fig. 3
17-799
30 cm




Fig. 4
17-800
30 cm



Fig. 5
17-801
30 cm





Fig. 6
17-802
30 cm

OBSTETRICS

Geburtshilfe | obstetrica
obstétrique | Ostetricia




Green-Armytage
Uterine Elevating Forceps
17-803 21 cm



Smellie
Obstetrical Forceps
17-804 25 cm
Perforator



NAEGELE
17-805
25 cm



Kane
obstetrical forceps
17-806
8.5 cm



Collin
Obstetrical Forceps
17-807
9 cm



Zimmerman
Myotome,
Spoon Tip,
17-808
24 cm




Zimmerman
Myotome,
Chisel Tip,
17-809
29 cm

Abortuszangen - Abortusküretten

Pinzas para abort - Curetas para aborto

Pinces à faux germes - Curette d'avortement

Pinze per placenta ed aborto - Curette per placenta ed aborto



Winter
abortion forceps



Fig. 1
17-810
28 cm
straight

17-811
28 cm
curved



Fig. 2
17-812
28 cm
straight

17-813
28 cm
curved




Fig. 3
17-814
28 cm
straight

17-815
28 cm
curved



Cuzzi
ovum curette
blunt, back serrated



Cuzzi



14.5 mm
17-816
30 cm



18 mm
17-817
30 cm



21 mm
17-818
30 cm



Saenger
Abortion Forceps



17-819
27 cm
Straight



17-820
27 cm
Curved



Wallich
ovum curette
17-821 42 cm



1:1



Pestalozza
Ovum Curette
17-822 30 cm



1:1

PELVIMETER, OBSTETRICAL FORCEPS



Beckenzirkel, Geburtszangen	Pelvímretro, Forceps obstétricos
Pelvimètre, Forceps obstétriques	Pelvimetri, Forcipi obstetrici



Thole Pelvimeter, Obstetrical Forceps 17-823
for internal and external measurements, graduated in centimeters



De Lee Pelvimeter, Obstetrical Forceps 17-824
graduated in centimeters and inches

for internal measurements between symphysis and promontory



Collin Pelvimeter, Obstetrical Forceps 17-825
Graduated in Centimeters



Martin Pelvimeter, Obstetrical Forceps 17-826
graduated in centimeters



Collin Pelvimeter, Obstetrical Forceps 17-827
graduated in centimeters

The additional straight arm can be replaced. The pelvimeter with straight arm is used for localisation of the portio vaginalis uteri and of vaginal metastasis.



Breisky Pelvimeter, Obstetrical Forceps 17-828
graduated in centimeters

OBSTETRICAL FORCEPS

Geburtszangen
Forceps obstétriques

Forceps obstétricos
Forcipi



Simpson
Obstetrical Forceps
17-829 23 cm
17-830 30 cm
17-831 33 cm
17-832 36 cm



Simpson - Braun
Obstetrical Forceps
17-833 33 cm
17-834 36 cm



Naegle
Obstetrical Forceps
17-835 36 cm
17-836 40 cm



Piper
obstetrical forceps
17-837 44 cm



Simpson - Luikart
Obstetrical Forceps
17-838 36 cm



Wrigley
Obstetrical Forceps
17-839 25 cm
17-840 28 cm
for caesarian section



De Lee
obstetrical forceps
17-841 30 cm



De Lee
obstetrical forceps
17-842 36 cm
17-843 15 cm



Elliot
obstetrical forceps
17-844 32 cm
17-845 38 cm

OBSTETRICAL FORCEPS



Geburtszangen

Forceps obstétriques

Forceps obstétricos

Forcipi




McLean
obstetrical forceps
17-846 36 cm



Kielland
Obstetrical Forceps
17-847 41 cm



Kielland-Luikart
obstetrical forceps
17-848 40 cm



Barnes-Neville
17-849 40 cm
Simpson handles, detachable
axis traction rod and handle.



Bil Obstetrical
Holder
17-850



Barnes-Neville
Obstetrical Forceps
17-851
with traction
handle 43.5 cm



Braun
obstetrical forceps
17-852
42 cm
Cranioclast, one
jaw fenestrated



Tarnier
obstetrical forceps
17-853
46 cm
Basiotribe



Sellheim
Obstetrical Lever
17-854
28 cm

LIPOSUCTION CANNULAS AND SETS



CANNULAS FOR INILTRATION

Cannulas for Iniltration | Cánulas de iniltración
Kanülen zur Iniltration | Canules pour l'iniltration



Cannulas For Iniltration
Measurements = working length



18-001
2mmØ
25cm
10holes



18-002
3 mm Ø
25 cm
10 holes



Cannulas
18-003
9 x 70 mm Ø
18-004
1 x 120 mm Ø
Packing unit 100 pieces, steril



18-005
Luer-Lock syringes
Adaptor with Luer-Lock connection, for use of
cannulas with handle connector on



18-006
Accessories for Syringes

LIPOSUCTION CANNULAS

Liposuction Cannulas
Liposuction-Kanülen

Cánulas de liposucción
Canules pour liposuccion



Entner
Cannulas



18-007
2 mm Ø
15 cm
18 holes



18-008
3 mm Ø
15 cm
18 holes



Entner Cannulas
Measurements = working length



18-009
2 mm Ø
15 cm
18 holes



18-010
3 mm Ø
15 cm
18 holes



Entner Cannulas
Measurements = working length



18-011
3 mm Ø
20 cm
22 holes



18-012
4 mm Ø
20 cm
22 holes



Entner
Cannulas



18-013
3 mm Ø
20 cm
22 holes




18-014
4 mm Ø
20 cm
22 holes

LIPOSUCTION CANNULAS




Liposuction Cannulas
Liposuction-Kanülen

Cánulas de liposucción
Canules pour liposuction



Measurements = working length



The bayonet shaped handles allow an easier manipulation.

LIPOSUCTION CANNULAS

Liposuction Cannulas
Liposuction-Kanülen

Cánulas de liposucción
Canules pour liposuction



SET SUGGESTIONS LIPOSUCTION CANNULAS



Set Suggestions Liposuction Cannulas
Setvorschläge Liposuction-Kanülen

Proposiciones de juegos cánulas de liposucción
Propositions pour boîtes canules pour liposuction



Face Cannulas
18-031
15 cm
2 mm Ø



18-032
5 cm 1,7 mm Ø
18-033
15 cm 1,7 mm Ø



Neck Cannulas
18-034
15 cm
3 mm Ø



Neck Cannulas
18-035
20 cm
3 mm Ø



Neck Cannulas
18-036
15 cm
2 mm Ø



18-037
5 cm 2 mm Ø
18-038
5 cm 3 mm Ø
18-039
15 cm 2 mm Ø
18-040
15 cm 3 mm Ø



Thigh - Lower Leg Cannulas
18-041 20 cm 3 mm Ø
18-042 30 cm 3 mm Ø
18-043 20 cm 4 mm Ø
18-044 30 cm 4 mm Ø



Abdomen / Saddlebag Cannulas
18-045 20 cm 3 mm Ø
18-046 20 cm 4 mm Ø



Abdomen / Saddlebag Cannulas
18-047 30 cm 3 mm Ø
18-048 30 cm 4 mm Ø



18-049
with rotatable tube connection
Accessories

LIPOSUCTION CANNULAS

Liposuction Cannulas
Liposuction-Kanülen

Cánulas de liposucción
Canules pour liposuction

Fork shape, rounded tips, dissecting surface between tips / one central hole.



18-050
2 mm Ø
5 cm

18-051
2 mm Ø
15 cm



18-052
3 mm Ø
5 cm

18-053
3 mm Ø
15 cm

18-054
3 mm Ø
25 cm

18-055
3 mm Ø
30 cm



18-056
4 mm Ø
5 cm

18-057
4 mm Ø
25 cm

18-058
4 mm Ø
30 cm



18-059
5 mm Ø
30 cm

Toledo Cannulas
Measurements = working length



18-060
3 mm Ø
15 cm

18-061
3 mm Ø
25 cm



18-062
4 mm Ø
15 cm

18-063
4 mm Ø
25 cm



three lateral holes

Mercedes Cannulas
Measurements = working length



Kotzur Cannulas



Turbo

18-064
3 mm Ø
30 cm



The improved arrangement of the suction holes ensures a more effective and gentle suction in areas with much connective tissue and in case of gynecomasty.

LIPOSUCTION CANNULAS



Liposuction Cannulas
Liposuction-Kanülen

Cánulas de liposucción
Canules pour liposuction



LIPOSUCTION CANNULAS

Liposuction Cannulas
Liposuction-Kanülen

Cánulas de liposucción
Canules pour liposuction



18-075
3 mm Ø
20 cm

18-076
3 mm Ø
30 cm



18-077
4 mm Ø
30 cm

18-078
4 mm Ø
40 cm



18-079
5 mm Ø
30 cm



18-080
6 mm Ø
30 cm



18-081
4 mm Ø
30 cm



18-082
5 mm Ø
30 cm

Cannulas
three holes, diagonal
Measurements = working length

Cannulas
one central hole, two lateral holes
Measurements = working length



18-083
3 mm Ø
20 cm



18-084
3 mm Ø
30 cm

Cannulas
Point with hole, one lateral hole



18-085
30 cm

Cannulas
Measurements = working length

LIPOSUCTION CANNULAS



Liposuction Cannulas
Liposuction-Kanülen

Cánulas de liposucción
Canules pour liposuction

For facial contouring:
one lateral hole



Cannulas
with handle connector



18-086
1,7 mm Ø
5 cm



18-087
1,7 mm Ø
15 cm



18-088
2 mm Ø
5 cm



18-089
2 mm Ø
15 cm



18-090
3 mm Ø
5 cm



18-091
3 mm Ø
15 cm



18-092
1,7 mm Ø
5 cm

18-093
1,7 mm Ø
15 cm

18-094
2 mm Ø
5 cm

18-095
2 mm Ø
15 cm

18-096
3 mm Ø
5 cm

18-097
3 mm Ø
15 cm

Cannulas

with Luer-Lock connection

Measurements = working length

Handles with connection thread for cannulas and hose connection Ø 13 mm, for operation with a suction pump.

Recommended handles for the tumescence local anaesthesia technique
Hole for interrupting the suction process



+



18-098
Fix Handle



18-099
with rotatable tube connection
Handle



+



18-100
Fix Handle



18-101
with rotatable tube connection
Handle

LIPOSUCTION CANNULAS

Liposuction Cannulas

Cánulas de liposucción

Liposuction-Kanülen

Canules pour liposuction

Liposuction cannulas, stainless steel with thumb plate and hose connection Ø 13 mm, for use with a suction pump.



**Gasparotti
Cannulas**

Measurements = working length



18-102
3 mm Ø
15 cm

one central hole



18-103
4 mm Ø
25 cm

one central hole



18-104
6 mm Ø
15 cm

one central hole



18-105
6 mm Ø
25 cm

one central hole



18-106
9 mm Ø
30 cm

one central hole



**Gasparotti
Cannulas**

Measurements = working length



18-107
4 mm Ø
15 cm

two central holes



18-108
4 mm Ø
25 cm

two central holes



18-109
6 mm Ø
30 cm

two central holes



**Gasparotti
Cannulas**

Measurements = working length



18-110
3 mm Ø
15 cm

one central hole



18-111
6 mm Ø
15 cm

one central hole



18-112
6 mm Ø
25 cm

one central hole

LIPOSUCTION CANNULAS



Liposuction Cannulas
Liposuction-Kanülen

Cánulas de liposucción
Canules pour liposuction

Liposuction Cannulas Set



18-113

Liposuction Cannula Set for Abdomen and Saddlebag With Threaded Fitting
Liposuction Cannula with threaded connector, 22-Holes, diameter 3.0 mm, working length 20cm
Liposuction Cannula with threaded connector, 22-Holes, diameter 4.0 mm, working length 20cm
Liposuction Cannula with threaded connector, 30-Holes, diameter 3.0 mm, working length 30cm
Liposuction Cannula with threaded connector, 30-Holes, diameter 4.0 mm, working length 30cm
Liposuction Cannula Handle
Liposuction Cannula Handle with Interrupting Hole

18-114

Liposuction Cannula Set for Face with Luer Lock Fitting

Liposuction Cannula, with one central hole, diameter 1.7 mm, working length 5 cm
Liposuction Cannula, with one central hole, diameter 1.7 mm, working length 15 cm
Liposuction Cannula, 18-Holes, diameter 2.0 mm, working length 15 cm



18-115

Liposuction Cannula Set for Face with Threaded Fitting


Liposuction Cannula, with one central hole, diameter 1.7 mm, working length 5 cm
Liposuction Cannula, with one central hole, diameter 1.7 mm, working length 15 cm
Liposuction Cannula, 18-Holes, diameter 2.0 mm, working length 15 cm
Liposuction Cannula Handle
Liposuction Cannula Handle with Interrupting Hole

18-116

Liposuction Cannula Set for Neck

Liposuction Cannula with Luer-Lock, with one central hole, diameter 2.0 mm, working length 5cm
Liposuction Cannula with Luer-Lock, with one central hole, diameter 3.0 mm, working length 5cm
Liposuction Cannula with Luer-Lock, with one central hole, diameter 2.0 mm, working length 15cm
Liposuction Cannula with Luer-Lock, with one central hole, diameter 3.0 mm, working length 15cm
Liposuction Cannula with threaded connector, 18-Holes, diameter 3.0 mm, working length 15 cm
Liposuction Cannula with threaded connector, Angled, 22-Holes, diameter 3.0 mm, working length 20 cm
Liposuction Cannula with threaded connector, 8-Holes, diameter 2.0 mm, working length 15 cm
Liposuction Cannula Handle
Liposuction Cannula Handle with Interrupting Hole

Liposuction Cannulas Set



18-117

Liposuction Cannula Set for Thie and Lower Leg with Threaded Fitting
Liposuction Cannula with threaded connector, byonet shape, 22-Holes, diameter 3.0 mm, working length 20 cm
Liposuction Cannula with threaded connector, byonet shape, 22-Holes, diameter 4.0 mm, working length 20 cm
Liposuction Cannula with threaded connector, byonet shape, 30-Holes, diameter 3.0 mm, working length 30 cm
Liposuction Cannula with threaded connector, byonet shape, 30-Holes, diameter 4.0 mm, working length 30 cm
Liposuction Cannula Handle
Liposuction Cannula Handle with Interrupting Hole

18-118

Micro Fat Injection Cannulas Set
Open Port Injector Cannula, 1.2mm, 10cm, Luer Lock
Multi Holes Harvester Cannula Model 2, 2.5mm, 20cm, Luer Lock
Tumescent Infiltrator Cannula, 1.5mm, 12cm Luer Lock
Micro Injector Cannula, 0.8mm, 4cm, Luer Lock
Micro Injector Cannula, 1.0mm, 5cm, Luer Lock
Holding Device
Luer Lock Syringe Cap


Safe Cannula Set

18-119

1 TWO STANDARD HOLE LIPOSUCTION CANNULA, 3.5mm, 25cm, Super Handle
2 Mercedes Liposuction Cannula, 5.0mm, 30cm, Super Handle
3 Mercedes Liposuction Cannula, 4.5mm, 30cm, Super Handle
4 Mercedes Liposuction Cannula, 4.0mm, 25cm, Super Handle
5 Single Basket Liposuction Cannulas, 4.0mm, 30cm, Super Handle
6 Double Basket Liposuction Cannulas, 4.0mm, 30cm, Super Handle
7 Triple Basket Liposuction Cannulas, 4.0mm, 30cm, Super Handle
8 Spiral Three Port Harvester Cannula, 3.0mm, 20cm, Super Handle
9 Spiral Three Port Harvester Cannula, 4.0mm, 25cm, Super Handle
10 Three Port Harvester Cannula, Model 4, 4.5mm, 30cm, Super Handle
11 Three Ports Pyramid Cannula, 3.5mm, 25cm, Super Handle
12 Basket Cannula, 5.0mm, 32cm, Super Handle



Fat Injection Syringes



18-120
Fat injection Syringe
50ml



18-121
Universal Fat Injection syringe
20ml



18-122
Universal Fat Injection syringe
10ml



18-123
Fat Injection Syringes
0.1 ml per click
1 ml



18-124
Fat Injection Syringe
3 ml special
0.1 ml per click

Fat Suction Syringes



18-125
Fat Suction Syringe
50 ml



18-126
Fat Suction Syringe
50 ml



18-127
Fat Suction Syringe
20 ml



18-128
Fat Suction Syringe
10 ml



18-129
Fat Suction Syringe
5 ml



18-130
Fat Suction Syringe
3 ml

HOLLOW WARES



STERILIZING DRUMS

Sterilisieren Schlagzeug

Tambours de stérilisation

Tambores de esterilización

Tamburi di sterilizzazione



	Ø	h
	mm	mm
19-001	80	80
19-002	120	120
19-003	160	120
19-004	190	140
19-005	190	200
19-006	240	120
19-007	240	160
19-008	240	240
19-009	290	120
19-010	290	160
19-011	290	240
19-012	340	160
19-013	340	240
19-014	390	145
19-015	390	200
19-016	390	290

STERILIZING DRUMS



STERILIZING DRUM

	Ø	h
	mm	Size mm
19-017	150	150
19-018	175	175
19-019	190	190
19-020	200	140
19-021	200	200
19-022	240	120
19-023	225	225
19-024	240	240
19-025	290	160
19-026	280	240
19-027	230	150
19-028	300	250
19-029	340	250
19-030	350	250
19-031	390	200



STERILIZING DRUM

	Ø	h
	mm	Size mm
19-032	150	150
19-033	175	175
19-034	190	190
19-035	200	140
19-036	200	200
19-037	240	120
19-038	225	225
19-039	240	240
19-040	290	160
19-041	280	240
19-042	230	150
19-043	300	250
19-044	340	250
19-045	350	250
19-046	390	200



STERILIZING DRUM

	Ø	h
	mm	Size mm
19-047	150	150
19-048	175	175
19-049	190	190
19-050	200	140
19-051	200	200
19-052	240	120
19-053	225	225
19-054	240	240
19-055	290	160
19-056	280	240
19-057	230	150
19-058	300	250
19-059	340	250
19-060	350	250
19-061	390	200



STERILIZING DRUM

	Ø	h
	mm	Size mm
19-062	150	150
19-063	175	175
19-064	190	190
19-065	200	140
19-066	200	200
19-067	240	120
19-068	225	225
19-069	240	240
19-070	290	160
19-071	280	240
19-072	230	150
19-073	300	250
19-074	340	250
19-075	350	250
19-076	390	200



G LOVES BOX
19-077
300 x 150 mm



**INSTRUMENT TRAY
WITH COVER**
19-078




INSTRUMENT TRAY
19-079



CATHER TRAY

	Size
19-080	230 x 75 mm
19-081	350 x 90 mm
19-082	430 x 100 mm
19-083	200 x 100 mm
19-084	300 x 100 mm
19-085	350 x 100 mm



PERFORATED BOX

	Size
19-086	100 x 50 mm
19-087	165 x 65 mm
19-088	250 x 100 mm
19-089	250 x 150 mm
19-090	250 x 200 mm
19-091	150 x 65 mm
19-092	170 x 70 mm
19-093	200 x 100 mm
19-094	250 x 125 mm
19-095	300 x 150 mm
19-096	300 x 200 mm
19-097	300 x 130 mm
19-098	320 x 150 mm
19-099	420 x 180 mm
19-100	460 x 125 mm
19-101	355 x 180 mm
19-102	300 x 180 mm
19-103	400 x 200 mm
19-104	510 x 200 mm



BOX WITH FILTER

	Size
19-105	100 x 50 mm
19-106	165 x 65 mm
19-107	250 x 100 mm
19-108	250 x 150 mm
19-109	250 x 200 mm
19-110	150 x 65 mm
19-111	170 x 70 mm
19-112	200 x 100 mm
19-113	250 x 125 mm
19-114	300 x 150 mm
19-115	300 x 200 mm
19-116	300 x 130 mm
19-117	320 x 150 mm
19-118	420 x 180 mm
19-119	460 x 125 mm
19-120	355 x 180 mm
19-121	300 x 180 mm
19-122	400 x 200 mm
19-123	510 x 200 mm




INSTRUMENTS BOX

Without Clip	Size
19-124	100 x 50 mm
19-125	165 x 65 mm
19-126	250 x 100 mm
19-127	250 x 150 mm
19-128	250 x 200 mm
19-129	150 x 65 mm
19-130	170 x 70 mm
19-131	200 x 100 mm
19-132	250 x 125 mm
19-133	300 x 150 mm
19-134	300 x 200 mm
19-135	300 x 130 mm
19-136	320 x 150 mm
19-137	420 x 180 mm
19-138	460 x 125 mm
19-139	355 x 180 mm
19-140	300 x 180 mm
19-141	400 x 200 mm
19-142	510 x 200 mm



INSTRUMENT BOX

With Clip	Size
19-143	100 x 50 mm
19-144	165 x 65 mm
19-145	250 x 100 mm
19-146	250 x 150 mm
19-147	250 x 200 mm
19-148	150 x 65 mm
19-149	170 x 70 mm
19-150	200 x 100 mm
19-151	250 x 125 mm
19-152	300 x 150 mm
19-153	300 x 200 mm
19-154	300 x 130 mm
19-155	320 x 150 mm
19-156	420 x 180 mm
19-157	460 x 125 mm
19-158	355 x 180 mm
19-159	300 x 180 mm
19-160	400 x 200 mm
19-161	510 x 200 mm



INSTRUMENT BOILER

	Size
19-162	180 x 80 x 40 mm
19-163	200 x 100 x 50 mm
19-164	250 x 125 x 40 mm
19-165	300 x 150 x 60 mm
19-166	350 x 175 x 75 mm
19-167	400 x 200 x 100 mm



PERFORATED TRAY

	Size
19-168	180 x 80 x 40 mm
19-169	200 x 100 x 50 mm
19-170	250 x 125 x 40 mm
19-171	300 x 150 x 60 mm
19-172	350 x 175 x 75 mm
19-173	400 x 200 x 100 mm

HOLLOWARE INSTRUMENTS

Hohlware Instrumenten
Orfebrería Instrumentos

orfèverie de Instruments de
Holloware Gli strumenti



INSTRUMENT BOILER ELECTRIC

	Size
19-174	180 x 80 x 40 mm
19-175	200 x 100 x 50 mm
19-176	250 x 125 x 40 mm
19-177	300 x 150 x 60 mm
19-178	350 x 175 x 75 mm
19-179	400 x 200 x 100 mm



INSTRUMENT BOILER WITH STAND

	Size
19-180	180 x 80 x 40 mm
19-181	200 x 100 x 50 mm
19-182	250 x 125 x 40 mm
19-183	300 x 150 x 60 mm
19-184	350 x 175 x 75 mm
19-185	400 x 200 x 100 mm



NEEDLE CASE
19-186
230 x 125 mm



SCALER BOX
19-187
300 x 180 mm



SCALER TRAY
19-188
100 x 200 mm

FIMA INSTRUMENTS



INDO BOX WITH SILVER FRAME
19-189



INSTRUMENT TRAY WITHOUT COVER

Size	Size
19-191 240 x 240 mm	19-202 230 x 125 mm
19-192 355 x 400 mm	19-203 200 x 150 mm
19-193 300 x 150 mm	19-204 250 x 200 mm
19-194 300 x 180 mm	19-205 300 x 250 mm
19-195 300 x 240 mm	19-206 350 x 300 mm
19-196 340 x 320 mm	19-207 400 x 350 mm
19-197 230 x 180 mm	19-208 450 x 300 mm
19-198 250 x 130 mm	19-209 300 x 200 mm
19-199 200 x 130 mm	19-210 350 x 250 mm
19-200 180 x 130 mm	19-211 250 x 125 mm
19-201 150 x 100 mm	19-212 200 x 100 mm
	19-213 400 x 200 mm



INSTRUMENT TRAY WITHOUT COVER

Size	Size
19-214 240 x 240 mm	19-225 230 x 125 mm
19-215 355 x 400 mm	19-226 200 x 150 mm
19-216 300 x 150 mm	19-227 250 x 200 mm
19-217 300 x 180 mm	19-228 300 x 250 mm
19-218 300 x 240 mm	19-229 350 x 300 mm
19-219 340 x 320 mm	19-230 400 x 350 mm
19-220 230 x 180 mm	19-231 450 x 300 mm
19-221 250 x 130 mm	19-232 300 x 200 mm
19-222 200 x 130 mm	19-233 350 x 250 mm
19-223 180 x 130 mm	19-234 250 x 125 mm
19-224 150 x 100 mm	19-235 200 x 100 mm
	19-236 400 x 200 mm



Autoclaveable
19-190



KIDNEY TRAY WITH COVER
(English Pattern)

	Size
19-237	150 mm
19-238	200 mm
19-239	240 mm
19-240	300 mm



KIDNEY TRAY (German Pattern)
(without cover)

	Size
19-241	150 mm
19-242	200 mm
19-243	240 mm
19-244	300 mm



SPLIT LAMP
19-245
100 ml.



BOWL

	Size
19-246	80 mm
19-247	100 mm
19-248	125 mm
19-249	150 mm
19-250	200 mm
19-251	250 mm
19-252	300 mm
19-253	350 mm
19-254	400 mm



WASH BASIN

	Size
19-255	80 mm
19-256	100 mm
19-257	125 mm
19-258	150 mm
19-259	200 mm
19-260	250 mm
19-261	300 mm
19-262	350 mm
19-263	400 mm



BOWL

	Size
19-264	80 mm
19-265	100 mm
19-266	125 mm
19-267	150 mm
19-268	200 mm
19-269	250 mm
19-270	300 mm
19-271	350 mm
19-272	400 mm



DRESSING JAR WITH KNOB

	Size
19-273	80 x 80 mm
19-274	100 x 100 mm
19-275	125 x 125 mm
19-276	150 x 150 mm
19-277	200 x 200 mm



DRESSING JAR

	Size
19-278	80 x 80 mm
19-279	100 x 100 mm
19-280	125 x 125 mm
19-281	150 x 150 mm
19-282	200 x 200 mm



DRESSING JAR

	Size
19-283	80 x 80 mm
19-284	100 x 100 mm
19-285	125 x 125 mm
19-286	150 x 150 mm
19-287	200 x 200 mm



DRESSING JAR

	Size
19-288	80 x 80 mm
19-289	100 x 100 mm
19-290	125 x 125 mm
19-291	150 x 150 mm
19-292	200 x 200 mm

HOLLOWWARE INSTRUMENTS

Hohlware Instrumenten
Orfebrería Instrumentos

orfèverie de Instruments de
Holloware Gli strumenti



**FORCEPS
JAR**

	Size
19-293	55 x 100 mm
19-294	55 x 140 mm
19-295	55 x 170 mm
19-296	55 x 180 mm
19-297	80 x 100 mm



**FORCEPS
JAR**

	Size
19-298	55 x 100 mm
19-299	55 x 140 mm
19-300	55 x 170 mm
19-301	55 x 180 mm
19-302	80 x 100 mm



**FORCEPS
JAR**

	Size
19-303	55 x 100 mm
19-304	55 x 140 mm
19-305	55 x 170 mm
19-306	55 x 180 mm
19-307	80 x 100 mm



**FORCEPS
JAR**

	Size
19-308	55 x 100 mm
19-309	55 x 140 mm
19-310	55 x 170 mm
19-311	55 x 180 mm
19-312	80 x 100 mm



**FORCEPS
JAR**

	Size
19-313	55 x 100 mm
19-314	55 x 140 mm
19-315	55 x 170 mm
19-316	55 x 180 mm
19-317	80 x 100 mm



**COTTON HOLDER
19-318**



GALIPOTS

	Capacity
19-319	2 Ounces
19-320	3 Ounces
19-321	4 Ounces
19-322	5 Ounces
19-323	6 Ounces
19-324	7 Ounces
19-325	8 Ounces



GALIPOTS

	Capacity
19-326	2 Ounces
19-327	3 Ounces
19-328	4 Ounces
19-329	5 Ounces
19-330	6 Ounces
19-331	7 Ounces
19-332	8 Ounces



**SPRIT LAMP
19-333
100 ml.**



**INDO BOX ROUND
19-334**



**NEEDLE CASE
19-335**



MEASURING JUG

	Capacity
19-336	1.00 Ltr
19-337	1.50 Ltr
19-338	2.00 Ltr



MEDICINE CUP GRADUATED

	Capacity
19-339	50 ml
19-340	60 ml
19-341	100 ml



BUCKET WITH LID

	Capacity
19-342	10 Ltr
19-343	12 Ltr
19-344	15 Ltr



SPUTUM MUG
19-345
100 x 50 mm



ENEMA CANE
Capacity
19-346 0.5 Ltr
19-347 1.5 Ltr
19-348 2.0 Ltr



DUST BIN
19-349
20 Lit



ROUND BEDPAN
19-350
300 x 85 mm



BEDPAN
19-351
Size: Adult, Child



BEDPAN
19-352
Size: Adult, Child



FUNNAL
19-353
100 mm



URINAL - FEMALE
19-354



URINAL - MALE
19-355



SHAVING MUG
19-356
100 mm


HOLLOWWARE INSTRUMENTS

Hohlware Instrumenten
Orfebrería Instrumentos

orfèvrerie de Instruments de
Holloware Gli strumenti



INSTRUMENTS MANAGEMENT SYSTEM



Cassette For 4 Instruments
19-357



Cassette For 5 Instruments
19-358



Cassette For 6 Instruments
19-359



Cassette For 7 Instruments
19-360



Cassette For 10 Instruments
19-361



Cassette For 20 Instruments
19-362



Cassette For Extraction Forceps
19-363



Cassette For Root Elevator
19-364



19-365
Cassette Tray
15 Pcs



19-366
Cassette Tray
15 Pcs



Instruments Cassette
19-367 5 Pieces Cassette
19-368 7 Pieces Cassette
19-369 10 Pieces Cassette




19-370
18×11×3.5cm
Extra large cassette for 6
Instruments, with compartment
for small parts hinged lids

STERILIZING MESH TRAY

Different sizes & dimensions
are available at buyer's demand



MESH TRAYS
19-371
Different sizes & dimensions
are available at buyer's demand



MESH TRAYS
19-372
05 Pcs




MESH TRAYS
19-373
07 Pcs




MESH TRAYS
19-374
10 Pcs



MESH TRAYS
19-375
20 Pcs



19-376
Mesh Cassettes
5 Pieces, 7 Pieces, 10 Pieces,
15 Pieces & 20 Pieces



MESH TRAY
19-377 240x130x50mm
19-378 240x130x60mm
19-379 240x250x50mm
19-380 240x250x60mm
19-381 480x250x50mm
19-382 480x250x60mm
19-383 480x250x70mm



MESH TRAY
19-384 240x250x 50mm
19-385 240x250x 60mm
19-386 240x250x 70mm
19-387 480x250x 50mm
19-388 480x250x 60mm
19-389 480x250x 70mm
19-390 480x250x100mm



MESH TRAY
19-391 240x250x50mm
19-392 240x250x60mm
19-393 240x250x70mm
19-394 280x250x50mm
19-395 280x250x60mm



MESH TRAY
19-396 280x250x 70mm
19-397 280x250x100mm
19-398 280x250x 50mm
19-399 280x280x 60mm
19-400 280x280x 70mm
19-401 300x280x 60mm

DESCRIPTION	PAGE#
KALLMORGEN SPECULUM	374
KALLMORGEN SPECULUM & RETRACTOR	374
KALT NEEDLE HOLDERS	110,289
KALT NEEDLE HOLDERS T.C	120
KANADISCHES MODELL MENISCOTOMY KNIFE	205
KANE OBSTETRICAL FORCEPS	404
KAYE #3RK ROUND SCALPEL HANDLE	25
KAYE-MILLARD RETRACTOR	147
KEITH ABDOMINAL NEEDLE	125
KELLY HAEMOSTATIC FORCEPS	79
KELLY INTESTINAL RECTAL SPECULUM	351
KELLY OPERATING DISSECTING SCISSORS T/C INSERT	58
KELLY OPERATING, DISSECTING SCISSORS	39
KELLY PROCTOSCOPE INTESTINAL RECTAL SPECULUM	352
KELLY RETRACTORS	145,148,149
KELLY SPONGE & DRESSING FORCEPS	100
KERN BONE HOLDING FORCEPS	223
KERRISON BONE EAR RONGEURS	321
KERRISON CERVICAL RONGEUR	253
KERRISON RONGEUR	255,256,257,258
KEYORKIAN BIOPSY AND SPECIMEN FORCEPS	380
KEYORKIAN ENDOCERVICAL CURETTES	398,399
KEYORKIAN FLUSHING CURETTES	398
KEYORKIAN YOUNG BIOPSY FORCEPS	381
KEY HANDLE	229
KEY HERNIA DIRECTOR	367
KEY PERIOSTEAL ELEVATOR	203,204
KIDNEY STONE FORCEPS	355,356,357
KIDNEY TRAY (GERMAN PATTERN)	428
KIDNEY TRAY WITH COVER	428
KIELLAND OBSTETRICAL FORCEPS	408
KIELLAND-LUIKART OBSTETRICAL FORCEPS	408
KILIAN, SEPTUM GOUGE	187
KILIAN-CLAUS SEPTUM GOUGE	187
KILLIAN LARYNGETOMY SCISSORS	323
KILLIAN NASAL SPECULUM	304
KILLIAN NASAL SPECULUM WITH FIBER	172
KILLIAN NASAL SPECULUM WITH FIBER OPTICS ILLUMINATION	174
KILLIAN PROBES	23
KILLIAN SEPTUM ELEVATOR	165
KILLIAN SEPTUM SPECULUM	305
KILLIAN, EUSTACHIAN CATHETERS	20
KILNER HOOK	127
KILNER NEEDLE HOLDERS	105
KILNER NEEDLE HOLDERS T.C	112
KILNER-DOUGHTY MOUTH GAG	327
KILNER-DOUGHTY MOUTH GAG HANDLE	327
KILNER-DOUGHTY MOUTH GAG SET WITH 3 BLADES COMPLETE	327
KIRK MALLETS	191
KIRSCHNER EXTENSION BOWS	233
KIRSCHNER RETRACTOR	161
KIRSCHNER WIRES & STEINMANN PINS CUTTER	219
KIRSCHNER WIRES	233
KIRSCHNER-BALFOUR RETRACTOR	161
KISTLER EAR CARTILAGE RASP	183
KLAPP RETRACTOR	135
KLEINERT KUTZ ELEVATOR AND DISSECTOR	199
KLEINERT-KUTZ BONE CUTTING FORCEPS	218
KLEINERT-KUTZ BONE RASP	183
KLEINERT-KUTZ BONE RONGEUR	216
KLEINERT-KUTZ SYNOVOTOMY RONGEUR	218
KLEINERT-KUTZ TENDON RETRIEVER	206
KLINIKUM BERLIN TUBING CLAMP	94
KNIFE HANDLE (SMOOTH) FOR BLADES NOS. 60 AND 70	24
KNIGHT NASAL POLYPUS AND SEPTUM FORCEPS	317
KNOWLES PLASTER SCISSORS	52
KNUTSON URETHROGRAPHIE-INSTRUMENTS	364
KOCHER GOITRE ENUCLEATOR	330
KOCHER HAEMOSTATIC FORCEPS	83
KOCHER HAEMOSTATIC FORCEPS STRONG PATTERN	82
KOCHER HOOKLET, RETRACTORS	130
KOCHER INTESTINAL FORCEPS CLAMPS	343,344
KOCHER RETRACTOR	133,134,136
KOCHER-LANGENBECK RETRACTOR	137
KOCHER-OCHSNER HAEMOSTATIC FORCEPS	84
KOCHER'S BLADDER RETRACTOR	367
KOCHER-WAGNER RETRACTOR	136
KOERTE RETRACTOR	134,135
KOGAN - MODIF. SPECULUM	372
KOGAN ENDOCERVICAL SPATULA	372
KOGAN SPECULUM	372
KOLBEL SOFT TISSUE RETRACTOR	155
KOMPLETT	01
KONIG ELEVATOR	202
KONIG RETRACTOR	132
KONIG REVERDIN NEEDLES	122
KORTE-WAGNER RETRACTOR	134
KOTZUR CANNULAS	414
KRATZ - BARRAQUER WIRE SPECULUM	294
KRAUS-VOSSE NASAL POLYPOUS SNARES	301
KRAUS-WILDE NASAL POLYPOUS SNARES	301
KRAUSE EAR POLYPOUS SNARES	301
KRAUSE NASAL POLYPOUS SNARES	301
KRAUSE, TROCER	15
KRISTELLER SPECULUM	373
KRISTELLER RETRACTOR	373
KRISTELLER SPECULUM & RETRACTOR	373
KRISTELLER SPECULUM TAPERED, FOR VIRGINS	373

DESCRIPTION	PAGE#
KRONECKER NEEDLE	123
KUEHNE CILIA FORCEPS	74
K-WIRE EXTRACTION PLIERS	210
L	
LA FORCE ADENOTOME PERFORATED	316
LA FORCE ADENOTOME UNPERFORATED	316
LABORDE TRACHEA DILATOR 3 PRONG	331
LACHRYMAL DILATOR	292
LACHRYMAL PROBES	292
LAHEY (SWEET) HAEMOSTATIC FORCEPS	87
LAHEY GOITER RETRACTOR	135
LAHEY GOITRE GRASPING FORCEPS	330
LAHEY NERVE DURA HOOK RETRACTOR	243
LAHEY RETRACTOR	135
LAMBERT-LOWMAN BONE CLAMP	221
LAMBOTTE BONE HOLDING FORCEPS	224
LAMBOTTE ELEVATOR	201
LAMBOTTE OSTEOTOME	187
LAMBOTTE RETRACTORS	131
LAMINECTOMY PUNCHES	261
247,248,249,250,251,252,253,254,255,256,257,258,259,260,	261
LAMINECTOMY RETRACTORS	152
LANCASTER EYE SPECULUM	151
LANCET TIP	20
LANDAU SPECULUM	375
LANDAU, TROCER	15
LANDOLT DISSECTING SUCTION TUBES	17
LANDOLT FORCEPS	69
LANDOLT OPERATING	35
LANDON VAGINAL RETRACTOR	374,375
LANE BONE ELEVATOR	204
LANE BONE HOLDING FORCEPS	224
LANE BONE LEVER	207
LANE BONE PLATE HOLDING FORCEPS	225
LANE STOMACH INTESTINAL CLAMP FORCEPS	346
LANE TISSUE FORCEPS	64
LANE TISSUE GRASPING FORCEPS	342
LANGE BANDAGE SCISSORS	52
LANGE TISSUE FORCEPS	75
LANGENBECK AMPUTATION	30
LANGENBECK BONE ELEVATOR	164
LANGENBECK BONE HOLDING FORCEPS	223
LANGENBECK DISSECTING KNIVES	30
LANGENBECK EAR CURETTE	308
LANGENBECK METACARPAL SAW	182
LANGENBECK NEEDLE HOLDERS T.C	121
LANGENBECK PERIOSTEAL ELEVATOR	202
LANGENBECK PERIOSTEAL ELEVATOR	203,204
LANGENBECK POLYPOUS GRASPING FORCEPS	383
LANGENBECK RETRACTOR	137,138
LARGE PIN CUTTER	220
LARRY PROBE AND GROOVED DIRECTOR	21
LARYNGEAL MIRRORS	330
LARYNGEAL MIRRORS HANDLE	330
LARYNGEAL MIRRORS WITH HANDLE	330
LARYNGEAL POLYPUD FORCEPS	337
LARYNGOSCOPE BATTERY HANDLE FOR	05
LARYNGOSCOPE BATTERY HANDLE FOR ADULT	03,04
LARYNGOSCOPE BATTERY HANDLE FOR ADULT AA BATTERIES	06
LARYNGOSCOPE BATTERY HANDLE FOR ADULT C TYPE	05
LARYNGOSCOPE FLEXIBLE (MC-COY) CONVENTIONAL SET	06
LARYNGOSCOPE SET OF CONVENTIONAL MILLER	05
LARYNGOSCOPES	03,04,05,06
LASTER-BURCH SPECULUM	294
LATERAL VAGINAL RETRACTOR	371
LATERAL VAGINAL RETRACTOR BLUE COATED LEEP / LIETZ	371
LATERALWALL RETRACTOR	371
LATHBURY-FARREL TRIANGULAR SERRATED END	23
LEAD HAND	209
LEBSCHKE STERNUM CHISEL	283
LEECHWILKINSON BIOPSY FORCEPS	381
LEECHWILKINSON SUCTION BIOPSY CURETTES	398
LEEDHAM TISSUE FORCEPS	63
LEGUEU BLADDER RETRACTOR	367
LEGUEU BLADDER SPATULA	366
LEGUEU RETRACTOR	160
LEKSELL BONE RONGEUR	214,216
LEKSELL RONGEUR	214
LEKSELL-STILLE BONE RONGEUR	215
LELAND-JONES PERIPHERAL VASCULAR CLAMPS	269
LEMPERT BONE CURETTE	197
LEMPERT BONE RONGEUR	212
LEMPERT SEPTUM ELEVATORS RASPARTORY	311
LENS MANIPULATING HOOK	293
LENSE EXPRESSOR	292
LERCHE TISSUE FORCEPS	62
LESTER RETRACTOR	176
LEUR NEEDLE EXTENSIONS	393
LEVASSEUR MERRILL FIBER	175
LEWIN PERFORATING BONE HOLDING FORCEPS	222
LEWIS TONSIL SCREW	313
LEWKOWITZ LITHOTOMY FORCEPS	365
LEXER CHISEL	185
LEXER GOUGE	186,185
LEXER GOUGE TUFFNEL HANDLE	186
LEXER MODIF. CHISEL	185
LEXER MODIF. OSTEOTOME	185

DESCRIPTION	PAGE#	DESCRIPTION	PAGE#
LEXER OPERATING SUPER CUT SCISSORS	55	MANDIBULAR CHANNEL RETRACTOR	174
LEXER OSTEOTOME	185	MANNERFELT HOOK	128
LEXER OSTEOTOME FLAT	186	MANNERFELT HOOKLET RETRACTORS	133
LEYLA SELF-RETAINING BRAIN RETRACTOR	239,262	MANNERFELT MODIF. BALLENGER	187
L-G INTERVERTEBRAL DISC RONGEURS	248	MANNERFELT-LANGENBECK HOOK	128
LICHTENBERG NEEDLE HOLDER	118	MARBACH PLASTER SCISSORS	53
LICHTENBERG NEEDLE HOLDERS T.C	121	MARKHAM/MEYERDING RETRACTOR	152
LICHTWITZ ANTRUM TROCAR	15	MARKWALDER BONE CUTTING RONGEUR FORCEPS	217
LIEBERMAN SPECULUM	294	MARQUARDT BONE RONGEUR	213
LIGATOR RUBBER RINGS IN PACKAGE OF 100 PIECES	353	MARSHNIK TONSIL SEIZING FORCEPS	329
LIGATURE NEEDLE	123	MARTEL HAEMOSTATIC FORCEPS	95
LIGATURE REVERDIN NEEDLES	122	MARTIN CARTILAGE BONE HOLDING CLAMP FORCEPS	225
LIGATUTR NEEDLES	124	MARTIN NEEDLE HOLDER	118
LINDERMAN AMPUTATION SAW	182	MARTIN NERVE ROOT RETRACTOR	242
LINK RONGEUR	259	MARTIN PELVIMETER, OBSTETRICAL FORCEPS	406
LINK SIMPLECLEAN KERRISON RONGEUR	260	MARTIN TONGUE SPATULA	325
LINK-BEATTY RONGEUR	259	MARTIN UTERINE SOUND	394
LIPOSUCTION CANNULA SET	419,420,421	MARTIN VAGINAL SPECULA	378
LIPOSUCTION CANNULA SET FOR NECK	420	MARTINI CURETTE	198
LIPOSUCTION CANNULAS	410,411,412,414,415,416,417,418,419,40,421,422	MARX BREAST RETRACTOR	166
LISTER BANDAGE PLASTER SCISSORS	51	MASING PERIOSTEAL SEPTUM ELEVATORS	312
LISTER BANDAGE PLASTER SCISSORS, T/C INSERT	59	MASING SEPTUM ELEVATOR	165
LISTER BOUGIES	362	MASS. GENERAL HOSPITAL NEEDLE HOLDERS	103
LISTER SINUS FORCEPS	95	MASSON NEEDLE HOLDER	106
LISTER-EXCENTRIC PLASTER SCISSORS	51	MASSON NEEDLE HOLDERS T.C	118
LISTON BONE CUTTING RONGEUR FORCEPS	216,217,218	MASSON-LUETHY NEEDLE HOLDERS	104
LISTON DISSECTING KNIVES	30	MASTIN GOITRE GRASPING FORCEPS	330
LISTON KEY HORSLEY BONE CUTTING RONGEUR FORCEPS	217	MATHIEU FOR CHILDREN AND VIRGINS	375
LISTON-STILLE BONE CUTTING RONGEUR FORCEPS	218	MATHIEU FOR CHILDREN AND VIRGINS FITTING	375
LITTAUER / LISTON BONE CUTTING RONGEUR FORCEPS	216	MATHIEU NEEDLE HOLDER T.C	120
LITTAUER CILIA FORCEPS	74	MATHIEU NEEDLE HOLDERS	108
LITTAUER EAR POLYPUS FORCEPS	318	MATHIEU RECTAL SPECULUM	350
LITTAUER NASAL POLYPUS AND SEPTUM FORCEPS	317	MATHIEU RETRACTOR	139
LITTAUER STITCH SCISSORS	43	MATHIEU VAGINAL RETRACTOR	374
LITTLE RETRACTOR	135	MATHIEU VAGINAL SPECULUM	374
LITTLE WOOD TISSUE GRASPING FORCEPS	342	MATHIEU-RYDER NEEDLE HOLDERS T.C	120
LOCKE PHALANGEAL FORCEPS	222	MATHIEU-RYDER VASCULAR NEEDLE HOLDERS T.C	120
LOCKHART-MUMMERY DIRECTORS	21	MATSON MENISCOTOMY KNIFE	205
LOCKING PLIER	211	MAUNDER MOUTH GAG	326
LOCKLIN DISSECTING SCISSORS	42	MAURER-COLLIN RIB SHEARS	284
LOCKWOOD GRASP FORCEPS	341	MAYFIELD BONE RONGEUR	214
LONG ANTERIOR KERRISON LAMINECTOMY PUNCH	261	MAYO DISSECTING & GYNECOLOGICAL SCISSORS	36
LONG CHROME HANDLE INTESTINAL RECTAL SPECULUM	352	MAYO DISSECTING & GYNECOLOGICAL SCISSORS T/C INSERT	58
LONG HAEMOSTATIC TISSUE FORCEPS	89	MAYO ROCHESTER NEEDLE HOLDERS	108
LONGE OPERATING DISSECTING SCISSORS	35	MAYO SCISSORS, SUPERCUT PLUS T.C.	170
LOPEZ-REINKE SUCTION TUBES	17	MAYO UTERINE SOUND	394
LORENZ PLASTER SCISSORS	52	MAYO VESSEL CLAMP	343
LORNA TOWEL FORCEPS	93	MAYO-BLAKE GALL STONE FORCEPS	358
LOVE NERVE DURA HOOK RETRACTOR	243	MAYO-COLLIN RETRACTOR	141
LOVE-GRUENWALD INTERVERTEBRAL DISC RONGEURS	247,248	MAYO-HARRINGTON DISSECTING & GYNECOLOGICAL SCISSORS	36
LOVELACE LUNG TISSUE FORCEPS	341	MAYO-HEGAR NEEDLE HOLDERS	101
LOWER HAEMOSTATIC FORCEPS	87	MAYO-HEGAR NEEDLE HOLDERS T.C	113
LOWMANN BONE CLAMP	221	MAYO-NOBLE DISSECTING & GYNECOLOGICAL SCISSORS	36
LUBET-BARBON NASAL POLYPUS AND SEPTUM FORCEPS	317	MAYO-ROBSON INTESTINAL CLAMP	344
LUC NASAL CUTTING FORCEPS	319	MAYO-ROBSON INTESTINAL FORCEPS CLAMPS	343
LUCAE / ENT DRESSING FORCEPS	67	MAYO-STILLE DISSECTING & GYNECOLOGICAL SCISSORS	36
LUCAE DRESSING FORCEPS	318	MAYO-STILLE DISSECTING & GYNECOLOGICAL SCISSORS T/C INSERT	58
LUCAE EAR LOOPS - HOOKS	308	MC DONALD ELEVATORS	200
LUCAE HOOKLET & SKIN HOOKS	128	MC INDOE BONE CUTTING FORCEPS	218
LUCAE MALLETS	191	MC INDOE BONE CUTTING RONGEUR FORCEPS	218
LUCAE TISSUE & DRESSING FORCEPS	67	MC INDOE DRESSING FORCEPS	65
LUCAE TUNING FORK	299	MC INDOE DRESSING FORCEPS T.C	78
LUCAS CHISEL	184	MCGEE LATERAL VAGINAL RETRACTOR	371
LUCAS CURETTE	199	MCGIVNEY HEMORRHOID GRASPING FORCEPS	353
LUCAS GOUGE	184	MCGIVNEY HEMORRHOIDAL LIGATOR	353
LUER BONE RONGEUR	211	MCGIVNEY HEMORRHOIDAL LIGATOR FORCEPS	354
LUER LOCK ADAPTER	14	MCINNIS POLYURETTE BREAST INSTRUMENT	163
LUER RETRACTOR	139	MC-INTOSH FIBER OPTIC LARYNGOSCOPE BLADE	03
LUER SUCTION TUBE WORKING	19	MC-INTOSH FIBER OPTIC LARYNGOSCOPE SET	03
LUER, TROCAR	14	MC-INTOSH FIBER OPTIC LARYNGOSCOPE SET HALOGEN ILLUMINATION	03
LUER-KOERTE MODIF WITH MALLEABLE SHAFT	361	MC-INTOSH LARYNGOSCOPE	03,04,05
LUER-LOCK SUCTION & IRRIGATION CANNULA	19	MC-INTOSH LARYNGOSCOPE CONVENTIONAL	05
LUER-STANDARD TRACHEA TUBES	331	MC-INTOSH LARYNGOSCOPE SET	03,04
LUER-STILLE BONE RONGEUR	215	MCKENZIE BRAIN CLIP RACK	240
LUETHY-BECK NEEDLE HOLDERS T.C	121	MCKENZIE CLIP APPLYING FORCEPS	240
LUNG SPATULAS, AORTIC DILATORS	278	MCKENZIE HAND DRILL BIT	235
LUNIATSCHEK TONSIL TAMPON PLUGGER	313	MCKENZIE INSTRUMENTS	240
LURE POLYOUS GRASPING FORCEPS	383	MCKEAN OBSTETRICAL FORCEPS	408
LURZ HEMORRHOIDAL SUCTION LIGATOR	353	MCLVOR MOUTH GAG COMPLETE SET	328
		MCLVOR MOUTH GAG HANDLE	328
		MCPHERSON - VANNAS	289
M		MCPHERSON FORCEPS	285,287
MACCOLLUM DINGMAN BREAST DISSECTOR	162,163	MCPHERSON NEEDLE HOLDER	289
MACCOLLUM DINGMAN SUBMAMMARY DISSECTOR	162	MEAD BONE RONGEUR	213
MACPHAIL NEEDLE HOLDERS	108	MEAD MALLET	192
MAGILL CATHETER INTUBATING FORCEP GOLD	01	MEASURING GAUGE	234
MAGILL CATHETER INTUBATING FORCEPS	01	MEASURING INSTRUMENTS	13
MAGILL SUCTION TUBE	20	MEASURING JUG	429
MAINGOT UTERUS HOLDING FORCEPS	387	MEDIAL RAMUS RETRACTOR	174
MALE CATHETERS CHARR., METAL	362	MEDICINE CUP GRADUATED	429
MALINIAC NASAL RETRACTOR WITH FIBER OPTIC	177	MEMORAL HEAD EXTRATOR	229
MALLEABLE RETRACTOR	150	MENISCOTOMY KNIFE	205
MALLEOLLAR BONE HOLDING FORCEPS	227	MENISCUS CLAMP	222
MALLET #49	192	MERCEDES CANNULAS	414
MALLETS	191,192,193	MESH BASKET	433
MALLET, LIGHTWEIGHT	192	MESH CASSETTES 5 PIECES, 7 PIECES, 10 PIECES, 15PIECES & 20 PIECES	432
MALONEY INTRA-OPERATIVE KERATOMETER TITANIUM	294	MESH TRAYS 05 PCS	432
MALTZ NASAL SAW	314	MESH TRAYS 07 PCS	432
MALTZ RASP	183,313		

DESCRIPTION	PAGE#
MESH TRAYS 10 PCS	432
MESH TRAYS 20 PCS	432
MESH TRAYS	432
METZENBAUM NEEDLE HOLDERS	106
METZENBAUM OPERATING SCISSORS T/C INSERT	57
METZENBAUM OPERATING, DISSECTING SCISSORS	37
METZENBAUM-FINO OPERATING SCISSORS T/C INSERT	57
METZENBAUM-FINO SCISSORS	38
METZENBAUM-NELSON OPERATING SCISSORS T/C INSERT	57
METZENBAUM-NELSON SCISSORS	38
MEYERDING FINGER RETRACTOR	138
MEYERDING RETRACTOR	146,155,208
MICHEL CLIPS	76
MICHEL FORCEPS	76
MICHEL MOVEABLE MIRROR	330
MICHEL REVERDIN NEEDLES	122
MICHEL TISSUE FORCEPS	75
MICHELE TREPHINE	238
MICRO NEEDLE HOLDER	289
MICRO ADSON TISSUE & DRESSING FORCEPS T.C	77
MICRO ADSON TISSUE FORCEPS T.C	77
MICRO ALLIGATOR	334
MICRO ALLIGATOR SCISSORS	334
MICRO CUP FORCEPS	249
MICRO EJECTOR KERRISON LAMINECTOMY RONGEURS	251,253
MICRO FAT INJECTION CANNULAS SET	421
MICRO INCISION FORCEPS	290
MICRO INCISION SCISSORS	290
MICRO INSTRUMENTS FOR LARYNGO BRONCHO OESOPHAGOSCOPY	334,335,336
MICRO NEEDLE HOLDERS	110,119,288
MICRO NEEDLE TISSUE FORCEPS T.C	78
MICRO RASP	183
MICRO RYDER NEEDLE HOLDERS	107
MICRO RYDER NEEDLE HOLDERS T.C	111
MICRO SCISSORS	48,49,50,288
MICRO TYING FORCEPS	67
MICRO-ADSON TISSUE & DRESSING FORCEPS	70
MICRO-ADSON-BIEMER TISSUE & DRESSING FORCEPS	70
MICROVASCULAR NEEDLE HOLDERS	106
MICROVASCULAR NEEDLE HOLDERS T.C	111
MIDDLEDORFF RETRACTOR	135
MIDDLETON-JANSEN SEPTUM FORCEPS	323
MIDFACE FACIAL DISSECTOR	163
MIKULICZ HAEMOSTATIC TISSUE FORCEPS	84,90
MIKULICZ RETRACTOR	145,146
MILLER ALVEOTOMY FILE	183
MILLER RASP	184
MILLER RETRACTOR	129
MILLER-COLBURN BONE FILE	163
MILLIGAN DURA DISSECTOR	242
MILLIGAN RETRACTOR	152
MILLIN BLADDER NECK SPREADER	365
MILLIN BLADDER RETRACTOR	366
MILLIN LIGATURE CARRYING FORCEPS ANGULAR	359
MILLIN PROSTECTOMY FORCEPS	359
MILLIN SUCTION TUBE	367
MILLIN VULSELLUM FORCEPS	359
MINI - NOVAK ASPIRATING CURETTE	397
MINI - TISCHLER BIOPSY AND SPECIMEN FORCEPS	380
MINI DRILL SLEEVE	230
MINI LEXER CHISEL	186
MINI LEXER OSTEOTOME	186
MINI NON ROTATING	241
MINI OPHTHALMOSCOPE	10
MINI ROTATING	241
MINI-HOHMANN BONE ELEVATOR	207
MINI-LAMBOTTE OSTEOTOME	187
MINI-ROUX RETRACTOR	140
MIXTER FORCEPS, LEE MODIF	88
MIXTER GALL STONE FORCEPS	358
MIXTER HAEMOSTATIC GALL DUCT FORCEPS	87
MIXTER-O'SHAUGNESSY HAEMOSTATIC FORCEPS	85
MOD. BONN SCISSORS	45
MOD. MARTIN DRESSING FORCEPS	61
MOD. MICHEL FORCEPS	76
MOD. USA, DRESSING FORCEPS	61
MOD. USA, TISSUE FORCEPS	61
MOD. WIEN NASAL SPECULUM	305
MODEL USA CAST SPREADER PLASTER	54
MODELL PICCOLO ROUND SERRATED	23
MODELL PICCOLO TRIANGULAR SERRATED END	23
MODIFIED-USA PLASTER SCISSORS	53
MÖELTGEN FLEXOMETER	
MOGEN CIRCUMCISION CLAMP	364
MOLLISON RETRACTOR	152
MOLT MOUTH GAG	326
MOORE DRILL AND BIT SET COMPLETE	228
MOREFIELD FORCEPS	285
MORGEN PROCTOSCOPE INTESTINAL RECTAL SPECULUM	352
MORRIS RETRACTOR	149
MORSE STERNAL RETRACTOR	279
MORSE-FAVALORO STERNAL RETRACTOR	280
MORTISING CHISEL	230
MOSQUITO HAEMOSTATIC FORCEPS	79
MOUTH GAGS	326,327,328
MOUTH MIRRORS HEAD	06
MOYNIHAN GALL DUCT DILATOR	361
MOYNIHAN HAEMOSTATIC FORCEPS	83

DESCRIPTION	PAGE#
MULTI PURPOSE VASCULAR CLAMPS	272,273,274,276
MURPHY RETRACTORS	131
MUSEUX TENACULUM FORCEPS	385
MUSEUX TONSIL HOLDING FORCEPS	329
MYOMA INSTRUMENTS, FISTULA HOOKS	396
MYRTLE LEAF PROBE	23
N	
NABATOFF CABLE	263
NABATOFF HANDLE	263
NABATOFF METAL CASE	263
NABATOFF METAL OLIVE	263
NABATOFF METAL PROBE	263
NABATOFF PROBE TIP IN PLASTIC	263
NABATOFF VARICOSE VEIN STRIPPER	263
NAEGELE OBSTETRICAL FORCEPS	407
NAEGELE OBSTETRICALS	404
NAGAHARA CHOPPER	292
NAGATA AURICULAR SEPTUM ELEVATOR	165
NAIL EXTRACTING FORCEPS	208
NAIL EXTRACTOR	
NAIL NIPPERS	47
NAIL RASP	183
NAIL SCISSORS	47
NAIL SPLITTING FORCEPS	209
NASAL CUTTING FORCEPS	320
NASAL KNIVES	309
NASAL POLYPUS & NASAL DRESSING FORCEPS	319
NASAL POLYPUS AND SEPTUM FORCEPS	317
NASAL POLYPUS SNARES	301
NASAL SAWS	313,314
NASAL SPECULA	302,303,304,305,306
NAUNTON MORGAN NEEDLE HOLDERS	103
NECK CANNULAS	413
NEEDLE HOLDER	287
NEEDLE CASE	427,429
NEEDLE HOLDERS (SCHWED. MODEL)	104
NEEDLE HOLDERS	101,102,104-121
NEEDLE NOSE PLIERS	209,210
NEEDLE-NOSE LOCKING PLIERS	211
NEFF MENISCOTOMY KNIFE	205
NEGUS HAEMOSTATIC FORCEPS	83
NEGUS MOUTH GAG HANDLE	327
NEGUS MOUTH GAG COMPLETE SET	327
NEIVERT NEEDLE HOLDERS	105
NEIVERT NEEDLE HOLDERS T.C	117
NEIVERT-ANDERSON OSTEOTOME	189
NELATON DIRECTORS	21
NELATON PROBES & DIRECTORS	21
NELSON TISSUE FORCEPS	64
NELSON TROCER THORAX	14
NEONATAL VASCULAR CLAMPS	276
NERVE HOOKS - DURA HOOKS	243
NEUFELD HAMMER	191
NEURO-DUOFLEX HAMMER	02
NEURO-MICRO DISSECTORS	245,246
NEUTRAL & LOAD DRILL GUIDE	231
NEW INSTRUTOR'S DISSECTING SET	11
NEW ORLEANS NEEDLE HOLDERS	103
NEW ORLEANS NEEDLE HOLDERS T.C	111
NEW PORT RETRACTORS	371
NICHOLS WEIGHTED SPECULUM	378
NICHOLS-VAGINALWEIGHTED SPECULUM	378
NO SCALPEL VASECTOMY HEMOSTAT FORCEPS	367
NO SCALPEL VASECTOMY RINGED FORCEPS CLAMPS	367
NOISE DRUM	299
NON-PERFORATING TOWEL FORCEPS	93
NORFOLK & NORWICH RETRACTOR	156
NORTHBENT STITCHES SCISSORS	43
NOSE PLIERS	210
NOSE SPECULUM	304
NOTO SPONGE & DRESSING FORCEPS	100
NOTT VIRGIN SPECULAUMS	370
NOVAK ASPIRATING CURETTE	397
NOVAK SUCTION	398
NOYES MICRO SCISSORS	50,287
NOYES NASAL DRESSING FORCEPS	319
NOYES SCISSORS	49
NUSSBAUM INTESTINAL FORCEPS CLAMPS	345
NUSSBAUM STOMACH INTESTINAL TWIN CLAMP FORCEPS	345
NYLON CAP	192
NYLON CAPS LEFT & RIGHT	192
NYLON HEAD Mallet	191
O	
O'BRIEN NEEDLE HOLDERS T.C	118
O'BRIEN VASCULAR SCISSORS	42
O'SHAUGNESSY HAEMOSTATIC FORCEPS	82
O'SULLIVAN - O'CONNOR VAGINAL RETRACTOR	379
O'SULLIVAN/O'CONNOR RETRACTOR	159
OBSTETRICAL FORCEPS	407,408
OBSTETRICALS	404
OBWEGESER - FREER	202
OBWEGESER ALVEOLAR (WISDOM TOOTH) RETRACTOR	173
OBWEGESER BLADE OSTEOTOME	190
OBWEGESER CHANNEL RETRACTOR	325
OBWEGESER CHEEK RETRACTOR	324
OBWEGESER CHIN RETRACTOR	174,325

DESCRIPTION	PAGE#	DESCRIPTION	PAGE#
OBWEGESER ELEVATOR	199,202	PERIOSTEAL ELEVATORS, SEPTUM ELEVATORS	312
OBWEGESER ILIAC CREST RETRACTOR	173	PERIOSTEAL ELEVATORS, TAMPERS	194
OBWEGESER MANDIBULAR CHANNEL RETRACTOR	174	PERIOSTEAL RASPATORIES	203
OBWEGESER NASAL SPINE RETRACTOR	173,325	PERIOSTEAL SPREADER	164
OBWEGESER NEEDLE	200	PERIPHERAL VASCULAR CLAMPS	269,270
OBWEGESER OSTEOTOME	190	PERKINS RETRACTOR	154
OBWEGESER RAMUS RETRACTOR	325	PERTHES	230
OBWEGESER RIB RASPATORY	165	PESTALOZZA OVUM CURETTE	405
OBWEGESER SOFT TISSUE RETRACTOR	325	PETRI PTERYGOID RETRACTOR WITH FIBER OPTIC	177
OBWEGESER-FREER ELEVATOR	199	PHANEUF HAEMOSTATIC TISSUE FORCEPS	89
OCHSNER DRESSING FORCEPS	65	PHENOLIC HANDLE Mallet	192
OCHSNER FORCEPS	75	PIERCE ANTRUM TROCER	15
OCHSNER TROCER THORAX	15	PILE FORCEPS	383
OCHSNER VASCULAR FORCEPS	266	PIN CUTTER MEDIUM	220
OLDBERG INTERVERTEBRAL DISC RONGEURS	248	PIN CUTTER SMALL	220
OLIVECRONA BRAIN SPATULA	241	PIN PULLER W/ GROOVE TC	219
OLIVECRONA DURA DISSECTOR	242	PINARD ALUMINIUM STETHOSCOPE	01
OLIVECRONA TOENNIS (SCOVILLE), CLIP FORCEPS	240	PINARD WOODEN STETHOSCOPE	01
OLIVECRONA TRIGEMINAL KNIFE	31	PIPER OBSTETRICAL FORCEPS	407
OLIVECRONA VASCULAR SCISSORS	44	PLAIN FORCEPS	68
OLIVECRONA WIRE SAW	238	PLASTER SAW ELECTRONIC POWER	228
OLLIER RETRACTOR	135	PLASTER SCISSORS	51-53
OLSEN-HEGAR NEEDLE HOLDERS	104	PLASTER SHEARS	54
OLSEN-HEGAR NEEDLE HOLDERS T.C	112	PLASTIC EAR CANULAS	06
OMBREDANNE MALLETS	192	PLATE BENDER	232
OPERATING & DISSECTING SCISSORS	35,37,38,39	PLATE BENDING PLIER	211
OPERATING SCISSORS (FLAT PATTERN)	33	PLESTER DAUMENPLATT	18
OPERATING SCISSORS	33,34,41	PLESTER RETRACTOR	154
OPERATING SCISSORS T.C.	56-59	PLESTER SUCTION TUBES 17 CM	18
OPERATING SCISSORS, B/B CVD. T/C INSERT	56	POLITZER EAR SPECULUM	295
OPERATING SCISSORS, B/B STR. T/C INSERT	56	POLYOUS GRASPING FORCEPS, ANASTOMOSIS	345,347,348,383
OPERATING SCISSORS, S/B STANDARD CVD. T/C INSERT	56	POOL SUCTION CANNULA	16
OPERATING SCISSORS, S/B STANDARD STR. T/C INSERT	56	POTTS COARCTATION CLAMPS	267
OPERATING SCISSORS, S/S CVD. T/C INSERT	56	POTTS LIGATURE FORCEPS	268
OPERATING SCISSORS, S/S STR. T/C INSERT	56	POTTS MULTI PURPOSE VASCULAR CLAMPS	272
OPHTHALMOSCOPE	09	POTTS MULTI PURPOSE VASCULAR CLAMPS	272
OPHTHALMOSCOPE HEAD	07,08	POTTS-SMITH DRESSING FORCEPS	61
OPTIC RETRACTOR FOR VERTICAL OSTEOTOMIES	175	POTTS-SMITH DRESSING FORCEPS T.C	78
ORBITAL RIM DISSECTOR	164	POTTS-SMITH T.C DRESSING FORCEPS	78
ORINGER LIP AND CHEEK RETRACTOR	324	POTTS-SMITH TISSUE FORCEPS	62
ORTHO MALLETS	192	POTTS-SMITH VASCULAR SCISSORS	41
ORTHOPAEDICS, TENDON STRIPPERS	209	POZZI TENACULUM FORCEPS	386
ORTHOPEDIC Mallet	192	PRATT ANTRUM CURETTE	310
OSTEOTOME BLADES, INTERCHANGEABLE	184	PRATT HAEMOSTATIC FORCEPS	95
OSTEOTOME ELEVATOR	202	PRATT RECTAL SPECULUM	350
OSTEOTOMES	184,187	PRATT TENACULUM FORCEPS	385
OTOSCOPE + OPHTHALMOSCOPE SET	10	PRATT UTERINE DILATORS 8/SET	389
OTOSCOPE HEAD	06,08	PRATT UTERINE DILATORS, DOUBLE-ENDED	389
P		PRESBYTERIAN TUBING CLAMP	94
PACKER MOSQUITO FORCEPS	81	PRESTON LOGAMENTUM FLAVUM FORCEPS	222
PANTINGS BLADDER RETRACTOR	366	PRINCE SCISSORS	45
PAPER DRESSING FORCEPS	67	PROBE MYRTLE-SHAPED	19
PARALLEL TARGET DEVICE	230	PROBE WITH EYE	19,23
PARAMETRIUM FORCEPS, HYSTERECTOMY FORCEPS	388	PROBE, DOUBLE ENDED	19
PARAMETRIUM HAEMOSTATIC FORCEPS	89	PROBES AND DIRECTORS	21-23
PARAMETRIUM HYSTERECTOMY FORCEPS	388	PROFESSIONAL OTOSCOPE & OPHTHALMOLOGY SET	07
PARHAM MARTIN BANDS	221	PROSTATIC RETRACTOR	360
PARHAM-MARTIN BONE CLAMP	221	PUNCH BRONCHOS OESOPHAGOSCOPY FORCEPS	335
PARKER RETRACTOR	141	PUNCHES SPHENOID BONE	321
PARKER-LANGENBECK RETRACTOR	141	PUTTI RASP	182,183
PARKER-MOTT RETRACTOR	141	PUTTI BONE FILE	182
PARTSCH CHISEL	184	PUTTI BONE RASP	182
PARTSCH GOUGE	185	Q	
PARTSCH MALLETS	191	QUERVAIN ELEVATOR	201
PASSOW BONE TAMPER	184	QUIRE EAR POLYPUS FORCEPS	318
PATELLA BONE SAW	183	R	
PATELLA MALLEOLAR FORCEPS	227	RABINER HAMMER	02
PAUL HOOKLET RETRACTORS	129	RAGNELL-DAVIS RETRACTOR	132
PAYR CRUSHING INTESTINAL CLAMP	349	RALK DRILL	228
PAYR REVERDIN NEEDLES	122	RAMATHIBODI TUBAL HOOK	372
PAYR-BABY CRUSHING INTESTINAL CLAMP	349	RAMATHIBODI UTERINE MANIPULATOR	372
PEAN HAEMOSTATIC FORCEPS SLENDER TYPE	82	RAMSAY DISSECTING FORCEPS	64
PEAPOD INTERVERTEBRAL DISC RONGEURS	254	RAMSAY TISSUE FORCEPS	64
PECK RAKE RETRACTOR	130	RAMUS RETRACTOR	174
PEDERSON VAGINAL SPECULUM	369	RANDALL KIDNEY STONE CLAMPS FORCEPS	355
PEDERSON VAGINAL SPECULUM LEEP / LIETZ	370	RANDALL KIDNEY STONE FORCEPS	356
PEDERSON VIRGIN SPECULAUMS	370	RANDALL SUCTION BIOPSY CURETTES	398
PEDIATRIC VASCULAR CLAMPS	275	RANEY APPLYIG & REMOVING FORCEPS	240
PEERS-BERTMAN TOWEL FORCEPS	93	RASPATORIES	186
PEET SPLINTER FORCEPS	76	RASPATORIES	186
PELKMANN SPONGE FORCEPS	98,99	RASPS WITH TUNGSTEN CARBIDE	313
PELTESON COTTON APPLICATOR	22	RASPS, FILES	184
PELVIC REDUCTION FORCEPS ASYMMETRIC	227	RAY BONE CURETTE	198
PELVIC REPOSITION FORCEPS	227	RAZOR BLADE BREAKER & HOLDER, STRAIGHT	26
PELVIMETER, OBSTETRICAL FORCEPS	406	RAZOR BLADE HOLDER & BREAKERS	26
PENFIELD DURA DISSECTOR	245	RAZOR HOLLOW	32
PENFIELD ELEVATORS	200	RECAMIER UTERINE CURETTE	401,402
PENN BREAST DISSECTOR	163	RECTAL BIOPSY FORCEPS	347
PENNINGTON / SIEZING HAEMOSTATIC FORCEPS	95	RECTUM VASECTOMY SET	368
PENNINGTON SEPTUM ELEVATORS RASPATORY	311	RECTUM VASECTOMY SET OF 2 PCS	368
PENNYBACKER PERIOSTEAL ELEVATOR	203	RECTUM VASECTOMY SET OF 3 PCS	368
PERCUSSION HAMMERS	02	RECTUM VASECTOMY SET OF 4 PCS	368
PERCY AMPUTATION SHIELD	182	REDUCTION FORCEPS	227
PERFORATED BASKETS	433	REICHERT EPIGLOTIS RETRACTOR	323
PERFORATED TRAY	426	REIL BONE HOLDING FORCEPS	227
PERIOSTEAL ELEVATORS	164,199-204,207	REILL BONE REDUCTION FORCEPS	226
PERIOSTEAL ELEVATORS, MENISCOTOMY KNIVES	205		

DESCRIPTION	PAGE#
REINER EAR SYRINGE	298
REINER EAR SYRINGE SHIELD	298
REPOSITION FORCEPS	227
REPROCUSSION FREE Mallet	193
RESANO DRESSING FORCEPS T.C	78
RESANO INTESTINAL POLYOUS GRASPING FORCEPS	345
RESECTION KNIFE	30
RETRACTOR, INSULATED BLADE	169
RETRACTORS & DRESSING FORCEPS	393
RETRACTORS	133-161,174-177,237,325,262,332,395
RETRACTORS INSTRUMENTS	173
RETRACTORS, VEIN HOOK	147
REUTHER LLIAC CREST RETRACTOR	176
REVERDIN NEEDLES	122,123,125,126
REVERDIN RETRACTOR	150
REYNOLDS OPERATING SCISSORS T/C INSERT	59
REYNOLDS TRANSAUXILLARY BREAST DISSECTOR, FENESTRATED,	162
RHEINSTADTER FLUSHING CURETTES	398
RHINOPLASTY INSTRUMENTS	178-181
RHINOPLASTY RETRACTORS	172
RHINOPLASTY SET	178-181
RIB CONTRACTORS	277
RIB RETRACTORS	279-282
RIB SHEAR	283
RIB SHEARS, STERNUM CHISEL	283,284
RIBBON RETRACTOR	150
RICHARDSON EASTMAN RETRACTOR	139
RICHARDSON RETRACTOR	148,149,176
RICHARDSON/KELLY RETRACTOR	147
RICHARDSON-EASTMAN RETRACTOR	149
RICHNAU HOLMGREN EAR SPECULUM	298
RICHTER OPERATING, DISSECTING SCISSORS	39
RICHTER-HEATH FORCEPS	76
RIECHERT SELF RETAINING BRAIN RETRACTOR	239
RIGBY RETRACTOR	154
RING CUTTER	209
RINGMESSER	315,316
RISKIN BONE RONGEUR	214
RITTER-HALLE FORNTAL SINUS BOUGIES	322
ROBERTS HAEMOSTATIC FORCEPS	86
ROBINSON'S KNIFE	205
ROCHARD RETRACTOR	160
ROCHESTER NEEDLE HOLDERS T.C	114
ROCHESTER-CARMALT HAEMOSTATIC FORCEPS	83,84
ROCHESTER-CHSNER HAEMOSTATIC FORCEPS	84
ROCHESTER-PEAN HAEMOSTATIC FORCEPS CURVED	84
ROCHESTER-PEAN HAEMOSTATIC FORCEPS STRAIGHT	84
ROD BENDER	232
ROEDER TOWEL FORCEPS	93
ROGERS PARAMETRIUM HYSTERECTOMY FORCEPS	388
ROGGE STERILIZING FORCEPS	434
ROLLET (BLAIR) HOOKLET RETRACTORS	129
ROLLET HOOK	127
ROLLET HOOKLET RETRACTORS	129
ROOSEVELT STOMACH INTESTINAL TWIN CLAMP FORCEPS	346
ROSE RETRACTOR	139
ROSEN EAR SPECULUM	296
ROSEN PHCO SPLITTER	292
ROSEN SUCTION TUBE	20
ROSER-KONIG MOUTH GAG	326,328
ROSS RETRACTOR	147
ROTTGEN-RUSKIN BONE RONGEUR	215
ROTUNDA CATHERERS	393
ROUND BEDPAN	430
ROUND KNURLED SCALPEL HANDLE	25
ROUND SCALPEL HANDLE	26
ROUX RETRACTOR	140
ROUX-MAGNUS RETRACTOR	140
RUBIN OSTEOTOME	189
RUBIN SEPTAL MORSELIZER FORCEPS	170
RUBIN SEPTAL MORSELIZER FORCEPS T.C.	170
RUBIN SEPTUM FORCEPS	323
RUBIO WIRE TWISTER NEEDLE HOLDERS T.C	117
RUSKIN BONE CUTTING RONGEUR FORCEPS	217
RUSKIN RONGEUR	214
RUSKIN-LISTON BONE CUTTING RONGEUR FORCEPS	218
RUSKIN-MINI BONE RONGEUR	216
RUSSIAN-MODELL TISSUE FORCEPS	65
RYDEL-SEIFFER TUNING FORK	299
RYDER FRENCH-EYE/STANDARD/MICRO NEEDLE HOLDERS T.C	115
RYDER NEEDLE HOLDERS	107
RYDER NEEDLE HOLDERS T.C	115

S

SACHS DURA DISSECTOR	244
SACHS NERVE ELEVATOR	199
SAENGER ABORTION FORCEPS	405
SAFE CANNULA SET	422
SAFETY SUPER CUT SCISSORS	55
SAILER SCISSORS	44
SALDANHA SOLZ ABDOMINAL RETRACTOR	174
SALPINGOGRAPH	392
SAMII FORCEPS	67
SAROOT HAEMOSTATIC FORCEPS	81
SAROT NEEDLE HOLDERS	102,107
SAROT NEEDLE HOLDERS T.C	113
SATINSKY TANGENTIAL OCCLUSION FORCEPS	271
SATTERLEE AMPUTATION & RESECTION	32

DESCRIPTION	PAGE#
SATTERLEE SAW	182
SAUERBRUCH BONE RONGEUR	215
SAUERBRUCH RETRACTOR	142
SAUERBRUCH RIB SHEARS	284
SAUERBRUCH-FREY RIB SHEARS	283
SAW BLADE ONLY	209
SAWYER RETRACTORS	148
SAYRE PERIOSTEAL ELEVATOR	204
SCALER BOX	427
SCALER TRAY	427
SCALP ELEVATOR	164
SCALPEL #3P HANDLE	25
SCALPEL #3RC HANDLE	25
SCALPEL BLADES	27
SCALPEL GRIP	25
SCALPEL HANDLE #4 WITH ADJUSTABLE	25
SCALPEL HANDLE #4 WITH SCALE	24
SCALPEL HANDLE #7	24
SCALPEL HANDLE / LONG ANGLED #3LA	24
SCALPEL HANDLE 3L LONG	24
SCALPEL HANDLE 4L LONG BLADE	24
SCALPEL HANDLE 4LA	24
SCALPEL HANDLE DOUBLE ENDED #3, #4	24
SCALPEL HANDLE/ #3	24
SCALPEL HANDLE/ #4	24
SCALPEL HANDLES	24-26
SCALPEL KNIVES	28-29
SCALPEL PLASTIC HANDLE # 5	26
SCALPEL PLASTIC HANDLE # 6	26
SCALPEL PLASTIC HANDLE #8 BLADE	26
SCALPEL PLASTIC HANDLE	26
SCHAEDEL TOWEL FORCEPS	93
SCHAEDE BONE CURETTE	195
SCHAEDE ORGAN GRASPING FORCEPS	341
SCHERBACK SPECULUM	376
SCHINDLER HAEMOSTATIC TISSUE FORCEPS	88
SCHINK DERMATOMES	32
SCHLESINGER CERVICAL RONGEUR	253
SCHLESINGER INTERVERTEBRAL DISC RONGEURS	248
SCHMIEDEN NEEDLE	123
SCHMIEDEN-DICK NEEDLE	123
SCHOENBORE RETRACTOR	132
SCHROEDER UTERINE SPOON	403
SCHROEDER BONE CURETTE	196
SCHROEDER FLUSHING CURETTES	398
SCHROEDER TENACULUM FORCEPS	386
SCHROEDER VULSELLUM TENACULUM FORCEPS	385
SCHROEDER-BRAUN TENACULUM FORCEPS	386
SCHUARDT RIB SHEARS	284
SCHUBERT BIOPSY FORCEPS	384
SCHUHKNECHT RETRACTOR	153
SCHULTZE, SALPINGOGRAPH	392
SCHULZE-BERGMANN INTESTINAL RECTAL SPECULUM	350
SCHUMACHER BIOPSY AND SPECIMEN FORCEPS	380
SCHUMACHER BIOPSY FORCEPS	380,384
SCHUMACHER PLASTER SCISSORS	53
SCHUMACHER RIB SHEARS	283
SCHUMACHER STERNUM SHEARS	220
SCHWARTZ CLIP APPLYING FORCEPS	240
SCHWARTZ VESSEL CLIP	240
SCHWARZ NASAL SPECULUM	304
SCHWENZER CHEEK RETRACTOR	176
SCISSORS	35,290
SCLERAL PLUG FORCEPS	286
SCOVILLE BONE CURETTE	197
SCOVILLE NERVE RETRACTOR	208
SCOVILLE, NERVE ROOT RETRACTOR	242
SCOVILLE-RICHTER RETRACTOR	159
SCREW DEPTH GAUGES	234
SCREW DRIVERS	230
SCUDDER INTESTINAL FORCEPS CLAMPS	345
SCULPO BREAST RETRACTORS	169
SEBILLEAU PERIOSTEAL ELEVATOR	194
SEDILOT ELEVATOR	200
SEGOND MYOMA KNIFE	396
SEGOND MYOMATOME	31
SEGOND-LANDAU HAEMOSTATIC FORCEPS	88
SEIFFERT PROBES	23
SELETZ-GELPI RETRACTOR	153
SELF RETAINING BRAIN RETRACTOR	239
SELF-RETAINING RETRACTORS	286
SELLHEIM OBSTETRICAL LEVER	408
SELLORS RIB CONTRACTOR	277
SILVERSTONE INTERVERTEBRAL DISC RONGEURS	248
SEMB BONE CUTTING RONGEUR FORCEPS	217
SEMB BONE HOLDING FORCEPS	223,225
SEMB BONE RONGEUR	215
SEMB PERIOSTEAL ELEVATOR	202
SEMB RIB SHEARS	283
SEMKEN DRESSING FORCEPS	61
SEMKEN TISSUE FORCEPS	76
SEMM ABDOMINAL RETRACTORS	161
SEMM VAGINAL SPECULUM	369
SENN RETRACTOR	131
SENN-GREEN RETRACTOR	132
SENNING NEEDLE HOLDERS	104
SENNING NEEDLE HOLDERS T.C (SCHWED. MODEL)	111
SENN-MUELLER RETRACTOR	131

DESCRIPTION

PAGE#

SEPTUM ELEVATORS - RASPATORIES	311,324,325
SEPTUM FORCEPS	323
SEPTUM INSTRUMENTS	322
SEPTUM NEEDLE SUCTION & IRRIGATION CANNULA	19
SEQUESTUM (NIRCOSIS) BONE HOLDING FORCEPS	223
SERRA GIGLI	182
SET SUGGESTIONS LIPOSUCTION CANNULAS	413
SEUTIN PLASTER SCISSORS	53
SEYFFERT VIRGIN SPECULUMS	370
SHARMAN ENDOMETRIAL BIOPSY CURETTES	398
SHAVING MUG	430
SHEA WITH FIBER OPTIC	175
SHEEHAN OSTEOHOME	190
SHORT REVERDIN NEEDLES	122
SHORTBENT STITCH SCISSORS	43
SHOULDER PENETRATING GOUGE	187
SHOULDER PERCUSSION AWL	222
SIDE CUTTER	220
SIDE PIN CUTTER	219
SIDE WALL RETRACTORS	371
SIEGEL ROUND	25
SIEGEL ROUND SCALPEL HANDLE #3 WITH ADJUSTABLE	25
SIEGLE MAGNIFIER EAR SPECULUM COMPLETE SET	297
SIEGLE MAGNIFIER OTOSCOPE	297
SIGMOID NPTCH RETRACTOR WITH FIBER OPTIC	177
SILVER DERMATOMES	32
SILVER RHINOPLASTIC CHISEL	189
SILVERSTEIN BREAST DISSECTOR	163
SIMON BONE CURETTE	196
SIMON DISSECTING KNIVES	30
SIMON FISTULA HOOK	396
SIMON VAGINAL RETRACTOR	374
SIMPSON - BRAUN OBSTETRICAL FORCEPS	407
SIMPSON - LUKART OBSTETRICAL FORCEPS	407
SIMPSON OBSTETRICAL FORCEPS	407
SIMPSON UTERINE SOUND	394
SIMS DISSECTING & GYNECOLOGICAL SCISSORS	40
SIMS HANDLE	395
SIMS RECTAL SPECULUM	350
SIMS SPECULUM & RETRACTOR	373
SIMS UTERINE CURETTES	400,402
SIMS UTERINE SOUND	394
SIMS-MAIER SPONGE & DRESSING FORCEPS	98
SIMS-SPECIAL UTERINE SOUND	394
SINGLEY-TUTTLE FORCEPS	65
SINSKEY HOOK	292
SINSKEY II LENS MANIPULATING HOOK	293
SINUSCOPY TROCAR	14
SKENE TENACULUM FORCEPS	386
SKIN FIXING PLATE	32
SLOTTED HAMMER	232
SLUDER-BALLENGER TONSIL GUILLOTINES & PUNCH FORCEPS	316
SLUDER-BALLENGER TONSIL GUILLOTINES & PUNCH FORCEPS HANDLE	316
SMELLIE OBSTETRICAL FORCEPS	404
SMELLIE RETRACTORS	395
SMILLIE MENISCOTOMY KNIFE	205
SMITH (MOD. USA) PLASTER SCISSORS	52
SMITH PETERSON	185
SMITH PETERSON OSTEOHOME	185
SMITH-BLUE INTESTINAL RECTAL SPECULUM	350
SMITH-PETERSEN	188
SMITH-PETERSEN GOUGE	188
SMITH-PETERSEN OSTEOHOME	188
SMITH-PETERSON LAMINECTOMY RONGEUR	213
SMITHWICK DURA DISSECTOR	244
SMITHWICK CLIP APPLYING FORCEPS	240
SMITHWICK NERVE DURA HOOK	242
SOFIELD RETRACTOR	139
SOLID COPPER HEAD Mallet	191
SOLZ AREOLA MARKER	162
SOLZ ATRAUMATIC BREAST DISSECTOR	162
SOLZ BREAST HOOK DISSECTOR	162
SOMER KIDNEY ELEVATING STONE CLAMPS FORCEPS	357
SPARE BLADE	32
SPATULA	23
SPATULA BREAST RETRACTOR DISSECTOR	163
SPECIMEN AND ASPIRATING CURETTE	397
SPECULUMS & CALIPERS	294
SPENCER CORONARY NEEDLE HOLDERS T.C	119
SPENCER STITCH SCISSORS	43
SPENCER-WELLS HAEMOSTATIC FORCEPS SLENDER TYPE	83
SPINAL DISTRACTOR	236
SPINAL RASPATORY	204
SPLINTER FORCEPS	74,75,76
SPLIT LAMP	428
SPRATT BONE CURETTE / HOLLOW HANDLE	195
SPREADER PLASTER	54
SPRIT LAMP	429
SPURLING KERRISON LAMINECTOMY RONGEURS	251
SPURLING INTERVERTEBRAL DISC RONGEURS	247,248,253
SPURLING LAMINECTOMY RONGEURS	252,253
SPUTUM MUG	430
S-SHAPED BLUNT CANNULA	16
S-SHAPED RETRACTOR	144
ST. CLAIR-THOMSON ADENOID CURETTE	315
ST. THOMAS KNIFE	205
STACKE PROBES	23
STAGBEETLE BONE REDUCTION FORCEPS	222

DESCRIPTION

PAGE#

STANDARD DRESSING FORCEPS	61
STANDARD DRESSING FORCEPS NARROW POINT	61
STANDARD DRESSING FORCEPS SERRATED	61
STANDARD DRESSING FORCEPS T.C	78
STANDARD NON ROTATING	241
STANDARD OPERATING SUPER CUT SCISSORS	55
STANDARD OPHTHALMOSCOPE	10
STANDARD OPHTHALMOSCOPE SET	09
STANDARD OTOSCOPE + OPHTHALMOSCOPE SET	09
STANDARD OTOSCOPE SET	09
STANDARD ROTATING	241
STANDARD TISSUE FORCEPS	61
STAUE - MOORE TENACULUM FORCEPS	386
ST-CLAIR THOMSON NASAL SPECULUM	302
STELLBRINK BONE RONGEUR	212
STERILE SCALPEL BLADE SIZE # 10	27
STERILE SCALPEL BLADE SIZE # 10A	27
STERILE SCALPEL BLADE SIZE # 11	27
STERILE SCALPEL BLADE SIZE # 11P	27
STERILE SCALPEL BLADE SIZE # 12	27
STERILE SCALPEL BLADE SIZE # 13	27
STERILE SCALPEL BLADE SIZE # 15	27
STERILE SCALPEL BLADE SIZE # 15A	27
STERILE SCALPEL BLADE SIZE # 17	27
STERILE SCALPEL BLADE SIZE # 18	27
STERILE SCALPEL BLADE SIZE # 20	27
STERILE SCALPEL BLADE SIZE # 21	27
STERILE SCALPEL BLADE SIZE # 22	27
STERILE SCALPEL BLADE SIZE # 23	27
STERILE SCALPEL BLADE SIZE # 24	27
STERILE SCALPEL BLADE SIZE # 25	27
STERILE SCALPEL BLADE SIZE # 36	27
STERILE SCALPEL BLADE SIZE # 40	27
STERILE SCALPEL BLADE SIZE # 6	27
STERILE SCALPEL BLADE SIZE # 60	27
STERILE SCALPEL BLADE SIZE # 70	27
STERILE SCALPEL BLADE SIZE # E11	27
STERILIZING BOXES (WITH SILICON PADS)	291
STERILIZING DRUMS	424
STERILIZING FORCEPS	434
STERNBERG LIP AND CHEEK RETRACTOR	324
STERNBERG RETRACTOR	132
STERNUM NEEDLE HOLDERS T.C	116
STETHOSCOPES	01
STEVENS OPERATING SCISSORS T/C INSERT	59
STEVENS SCISSORS	47
STEVENS TENOTOMY SUPER CUT SCISSORS	55
STEVENS TISSUE FORCEPS	72
STIEGLITZ FORCEPS WITH MEGNIFYING MIRRORS	75
STIEGLITZ SPLINTER FORCEPS	76
STILLE	185,188
STILLE BONE RONGEUR	214
STILLE CHISEL	185,188
STILLE GOUGE	185,188
STILLE HAEMOSTATIC FORCEPS	87
STILLE OSTEOHOME	185,188
STILLE PATTERN DRILL SET	228
STILLE PLASTER SHEARS	54,220
STILLE-BARRAYA TISSUE FORCEPS	65
STILLE-HORSLEY BONE CUTTING FORCEPS	218
STILLE-LEUR RONGEUR	215
STILLE-LISTON BONE CUTTING RONGEUR FORCEPS	217,218
STILLE-LUER ECHLIN DUCKBILL RONGEUR	215
STILLE-LUER RONGEUR FORCEPS	214
STILLE-RUSKIN BONE RONGEUR	214
STITCH SCISSORS	43
STOCKMANN PENILE CLAMP	364
STONE TISSUE FORCEPS	62
STONE-WATT CLOSING FORCEPS	342
STONE-WATT INTESTINAL ANASTOMOSIS CLAMP	348
STRABISMUS HOOK	292
STRANDELL-STILLE RETRACTOR	132
STRASSMANN KIDNEY ELEVATING STONE CLAMPS FORCEPS	357
STRATTE NEEDLE HOLDERS	102,359
STRATTE NEEDLE HOLDERS T.C	114
STRAUSS PENILE CLAMP	365
STRUEMPFLER EAR FORCEPS	333
SUBACROMION SPREADER	224
SUCTION AND IRRIGATION CANNULAS	19
SUCTION CANNULA FOR TONSILLECTOMY, ORAL	16
SUCTION ELEVATOR	164
SUCTION TUBE, GUIDE NEEDLE	20
SUCTION TUBES	16-18
SULLIVAN SCISSORS	45
SUPER CUT SCISSORS	55
SUPRA CUT DISTRACTION SCREW	237
SUPRA CUT DISTRACTION SCREW DRIVER	237
SURGEONS NEEDLES	125,126
SUTURE RETRIEVAL HOOK	395
SUTURE TYING FORCEPS	285,287
SWEDISH PATTERN DRESSING FORCEPS	65
SYPERT RONGEUR	216
SYSTRUNK SCISSORS	46
T.C OEHLER TISSUE FORCEPS	76
TAKAHASHI NASAL CUTTING FORCEPS	19
TANGENTIAL OCCLUSION FORCEPS	271

DESCRIPTION	PAGE#
TAP SLEEVE	231
TARNIER OBSTETRICAL FORCEPS	408
TAY-BRONSON BONE CURETTE	197
TAYLOR HAMMER	02
TAYLOR RETRACTOR	208
TAYLOR SPINAL RETRACTORS	139
TEALE HERNIA DIRECTOR	367
TEALE TENACULUM FORCEPS	385
TEBBETTS BREAST AUGMENTATION FIBER OPTIC	169
TEBBETTS BREAST AUGMENTATION FIBER OPTIC RETRACTOR	168
TEBBETTS FIBER OPTIC RETRACTOR	168,169
TEBBETTS RETRACTOR	168
TECHNICAL FORCEPS	211
TEMPORAL T DISSECTOR	163
TENACULUM EMETT	395
TENACULUM FORCEPS	385,386
TENDON HOOKS RETRACTOR	135
TENDON SEIZING FORCEPS	99
TENDON STRIPPERS, BONE LEVERS	206
TENDON STRIPPERS, TECHNICAL FORCEPS	210
TENDON TUNELLING FORCEPS	209
TENNANT TYING FORCEPS	286
TENOTOMY, CAPSULOTOMY VANNAS SCISSORS	289
TESSIER RIB GRAFT RETRACTOR ADULT	176
T-HANDLE AO COUPLING	229
T-HANDLE FOR AO COUPLING	229
THIERSCH DERMATOMES	32
THIGH - LOWER LEG CANNULAS	413
THIN FOOTPLATE KERRISON LAMINECTOMY PUNCH	261
THOLE PELVIMETER, OBSTETRICAL FORCEPS	406
THOMS-ALLIS INTESTINAL AND TISSUE GRASPING FORCEPS	339
THOMS-GAYLOR BIOPSY FORCEPS	384
THOMSON WALKER BLADDER RETRACTOR	366
THOMSON WALKER NEEDLE HOLDERS	103
THOMSON WALKER PENILE CLAMP	365
THOMSON-WALKER NEEDLE HOLDERS T.C	118
THOMSON-WALKER'S PENILE CLAMPS	93
THORACIC & I.M.A RIB RETRACTORS	281
THORACIC RIB RETRACTORS	281
THORPE TISSUE & DRESSING FORCEPS	70
THUDICHUM (GOLDSMITH) NASAL SPECULUM	302,303
TIECK-HALLE NASAL SPECULUM	304
TILLEY NASAL POLYPUS AND SEPTUM FORCEPS	317
TISCHLER - MORGAN BIOPSY AND SPECIMEN FORCEPS	380
TISCHLER BABY BIOPSY FORCEPS	380
TISCHLER BIOPSY FORCEPS	381
TISSUE & DRESSING FORCEPS	61-72
TISSUE & DRESSING FORCEPS T.C	77
TISSUE FORCEPS	63,68,71
TOBOLD TONGUE DEPRESSOR	324
TOENNIS DURA DISSECTOR	244
TOENNIS HOOK & SKIN HOOK	127
TOENNIS NEEDLE HOLDERS	109
TOENNIS NEEDLE HOLDERS T.C	119
TOENNIS NERVE DURA HOOK	242
TOHOKU TOWEL FORCEPS	93
TOLEDO CANNULAS	414
TONGUE SPATULA	325
TONOMETER SET	302,318,329
TONSIL INSTRUMENTS	308
TONSIL KNIFE	308
TONSIL KNIVES	308
TONSIL NEEDLE SUCTION & IRRIGATION CANNULA	19
TONSIL, NASAL SCISSORS	44
TOWEL FORCEPS	93
TOWNLEY CALIPER	13
TOYNBEE EAR SPECULUM	295
TRACHEA DILATORS, TRACHEA TUBES	331
TRACHEAL RETRACTOR	129
TRACHEOTOMY BRONCHO-OESOPHAGOSCOPY FORCEPS SET	332
TRACHEOTOMY SET	332
TRANSORAL DISSECTOR	163
TRAUBE HAMMER	02
TRAVERS RETRACTOR	153
TRELAT ELEVATOR	201
TRELAT VAGINAL SPECULUM	369
TRETHOWAN'S BONE ELEVATOR	204
TREVES TISSUE FORCEPS	62
TRIGEMINAL KNIFE	31
TROCAR	15
TROCAR BREAST TUBE	163
TROCARS	14-16
TROELTSCH DRESSING FORCEPS	318
TROELTSCH FORCEPS	69
TROEMNER HAMMER	02
TROEMNER-SIMPLEX HAMMER	02
TROELTSCH EUSTACHIAN CATHETERS	20
TROUSSEAU TRACHEA DILATOR 2 PRONG	331
TROUTMAN MICRO SCISSORS	49
TROUTMAN NEEDLE HOLDERS	109
TROUTMAN-CASTROVIEJO SCISSORS	290
TROUTMAN-CHRIS BLADE BREAKER & HOLDER ANGLED, WITH PUSH-PULL LOCK	26
TROUTMAN-CHRIS BLADE BREAKER & HOLDER WITH PUSH-PULL LOCK	26
TUBBS MITRAL VALVE AORTIC DILATOR	278
TUBING CLAMP	94
TUBING CLUMP WITH PROTECTOR	94
TUFFIER / HABERER RETRACTOR	150

DESCRIPTION	PAGE#
TUFFIER HAEMOSTATIC FORCEPS	95
TUFFIER RETRACTOR	145,146
TUFFIER RIB RETRACTORS	280,281
TUFFIER TISSUE FORCEPS	339
TUMARKIN EAR SPECULUM	296
TUNING FORKS	299,300
TUNING FORKS NEUROL WITH HAMMER	299
TUNING FORKS SET	300
TUPPERS HAND RETRACTOR SET	173
TURNBULL NIPPERS	47
TURNER-DOYEN RETRACTORS	148
TYDING TONSIL SNARE	302
U	
U.S. ARMY PATTERN	184
UCKERMANN TRACHEA TUBES	331
UCKERMANN TROCER	14
UFFENORDE BONE CURETTE LEFT	195
UFFENORDE BONE CURETTE RIGHT	195
ULRICH BONE HOLDING FORCEPS	224
ULRICH BONE HOLDING SELF RETAINING FORCEPS	224
ULRICH SPONGE FORCEPS	98
UNI - BURKE BIOPSY AND SPECIMEN FORCEPS	382
UNI - KEVORKIAN BIOPSY AND SPECIMEN FORCEPS	382
UNI - OVAL BIOPSY AND SPECIMEN FORCEPS	382
UNI - RANDALL ASPIRATING CURETTE	397
UNI - RANDALL ASPIRATING CURETTE	397
UNI - TISCHLER BIOPSY AND SPECIMEN FORCEPS	382
UNI - TOWNSEND BIOPSY AND SPECIMEN FORCEPS	382
UNIVERSAL ABDOMINAL TROCER SET	14
UNIVERSAL BONE RONGEUR	213
UNIVERSAL FAT INJECTION SYRINGE	423
UNIVERSAL HAMMER	192
UNIVERSAL OTOSCOPE & OPHTHALMOLOGY SET	08
UNIVERSAL PARAMEDIC SHEAR	60
UNIVERSAL ROD BENDER	232
UNIVERSAL SCISSORS	60
UNIVERSAL WIRE CUTTING SCISSORS	46
URETHRAL CATHETERS	362,363
URINAL - FEMALE	430
URINAL - MALE	430
URO-TANGENTIAL GALL STONE FORCEPS	358
USA MODEL SPECIMEN FORCEPS	380
US-ARMY RETRACTOR	139
UTERINE CURETTES	400-402
UTERINE DILATORS	389
UTERINE SCOOPS, SUCTION CURETTES & UTERINE DEPRESSOR	399
UTERINE SOUNDS	390,391,394
UTERINE SPOONS	403
UTERUS HOLDING FORCEPS	387
UTILITY STERILIZING FORCEPS	434
UTILITY TECHNICAL FORCEPS PLIERS	210
UTRATA CAPSULORHEXIS FORCEPS	285-287
V	
VACHER NASAL SPECULUM	298
VAGINAL DRESSING FORCEPS	393
VAGINAL SPECULA	369-379
VAGINAL SPECULA, VAGINAL RETRACTORS	374
VALLIEUX UTERINE SOUND	394
VAN BUREN SEQUESTERUM FORCEPS	223
VAN BUREN URETHRAL SOUND	362
VAN DOREN BIOPSY FORCEPS	384
VANNAS SCISSORS	288,289
VARICOSE VEIN PROBE SET	263
VASCULAR FORCEPS	265,266
VASCULAR SCISSORS	41,42
VEBRUGGE SELF-CENTERING BONE HOLDING FORCEPS	225
VEBRUGGE BONE HOLDING FORCEPS	225,226
VERNIER CALIPER	13
VERNON-DAVID INTESTINAL RECTAL SPECULUM	351
VERONIKIS VAGINAL RETRACTORS	375
VESTIBULUM RETRACTOR	324
VICKERS MALLETS	193
VIENNA SPECULUM	304
VIRCHOW DISSECTING KNIVES	31
VITAL HEAVY DUTY PLIERS T.C	219
VITAL WIRE CUTTER	219
VITAL WIRE EXTRACTION PLIERS T.C	219
VITRECTOMY SCISSORS	290
VITREORETINAL FORCEPS	290
VOLKMAN RETRACTOR	131
VOLKMANN BONE CURETTE	195
VOLKMANN CURETTE	198
VOLKMANN RETRACTORS	131,134
VOLTOLINI NASAL SPECULUM	298
VON GRAEFE KNIFE	293
VON GRAEFE STRABISMUS & SKIN HOOKS	128
W	
WACHENFELDT TISSUE FORCEPS	75,76
WAGENER EAR LOOPS - HOOKS	308
WAGENER HOOKLET & SKIN HOOKS	128
WALB DISSECTING KNIVES	31
WALDMANN PLASTER SCISSORS	52
WALDMANN SCISSORS	53
WALLICH OVUM CURETTE	405
WALTER ALAR AND NASAL RETRACTOR	171
WALTER GOUGE	189

DESCRIPTION PAGE#

WALTON-LISTON BONE CUTTING RONGEUR FORCEPS	218
WANGENSTEEN DRESSING FORCEPS	62
WANGENSTEEN DRESSING FORCEPS T.C	78
WANGENSTEEN NEEDLE HOLDERS	106
WANGENSTEEN NEEDLE HOLDERS T.C	113
WARTENBERG PINWHEEL	02
WASH BASIN	428
WASSMUND RETRACTOR	135
WATSON DERMATOMES	32
WATSON-CHEYNE ELEVATORS	200
WATSON-CHEYNE RASPATORY	312
WATSON-WILLIAMS NASAL CUTTING FORCEPS	319
WAUGH DRESSING FORCEPS	63
WAUGH TISSUE FORCEPS	63
WEBER AESTHESIOMETER	01
WEBSTER NEEDLE HOLDER	101
WEBSTER NEEDLE HOLDERS T.C	117
WECKER SCISSORS	49,290
WEIL NASAL CUTTING FORCEPS	319
WEIL-BLAKESLEY NASAL CUTTING FORCEPS	320
WEINBERG VAGOTOMY RETRACTOR (JOE'S HOE)	147
WEINGÄRTNER EAR FORCEPS	333
WEISSBARTH SPECULUM & RETRACTOR	377
WEITLANER RETRACTOR	154,155
WELL BLAKESLEY LAMINECTOMY RONGEURS	249
WELLENSCHLIFF NASAL SAW	182
WERTHEIM HAEMOSTATIC FORCEPS	92
WERTHEIM HAEMOSTATIC TISSUE FORCEPS	91
WERTHEIM HYSTERECTOMY FORCEPS	91
WERTHEIM NEEDLE HOLDERS	104
WERTHEIM NEEDLE HOLDERS T.C	112
WERTHEIM PARAMETRIUM CLAMPS	91
WERTHEIM RETRACTOR	375
WERTHEIM-CULLEN HAEMOSTATIC FORCEPS	92
WESTCOTT MICRO SCISSORS	49
WESTCOTT SCISSORS	288
WESTERMARK HAEMOSTATIC FORCEPS	89
WESTMACOTT NASAL POLYPUS AND SEPTUM FORCEPS	317
WHITE TONSIL SEIZING FORCEPS	329
WHITEHEAD MOUTH GAG	326
WHITTAKER TENACULUM 3 PRONG	395
WIDE MOUTH EJECTOR KERRISON LAMINECTOMY RONGEURS	251,253
WIEDER BREAST AND FACIAL RETRACTOR / TONGUE DEPRESSOR, SERRATED,	166
WIEDER BREAST AND FACIAL RETRACTOR / TONGUE DEPRESSOR, SERRATED,	166
WIEDER TONGUE DEPRESSOR	316
WIENER HOOK & SKIN HOOK	127
WIENER NASAL SAW	314
WILDE FORCEPS	69
WILDE INTERVERTEBRAL DISC RONGEURS	248
WILDE NASAL DRESSING FORCEPS	318
WILDER LACRIMAL DILATOR	23
WILLET-MARTEL HAEMOSTATIC FORCEPS	95
WILLIAM LACRIMAL PROBES DOUBLE ENDED, STERLING SILVER	22
WILLIAM PROBE	22
WILLIAMS GRASP FORCEPS	338
WILLIAMS STYLE MICRO DISC INTERVERTEBRAL DISC RONGEURS	249
WILLIGER CURETTE	198,199
WILLIGER ELEVATOR	201
WILLIGER HAMMER	191
WILLIGER MALLETS	192
WINKELMANN CIRCUMCISION CLAMPS	364
WINTER ABORTION FORCEPS	405
WINTER TENACULUM FORCEPS	385
WINTERTON VAGINAL SPECULUM	369
WIRE BENDING FORCEPS	210
WIRE CUTTING INSTRUMENTS	219,220
WIRE CUTTING SCISSORS	46
WIRE CUTTING SCISSORS T/C INSERT	58
WIRE FORCEPS	220
WIRE GUIDE	230
WIRE INSTRUMENTS	232
WIRE RETAINING RETRACTOR	151
WIRE SPECULUMS	294
WIRE TIGHTENING NEEDLE HOLDERS T.C	117
WIRE TWISTERS NEEDLE HOLDERS T.C	116
WIRE VECTUS	292
WITH FIBER OPTIC ILLUMINATION	171
WITTNER BIOPSY AND SPECIMEN FORCEPS	380
WITTNER BIOPSY FORCEPS	380
WOLFF SPREADER PLASTER	54
WOODSON DURAL SEPARATOR/PACKER	200
WOODSON ELEVATOR/SPATULA	200
WRENCH KIRSCHNER EXTENSION BOWS	233
WRIGLEY OBSTETRICAL FORCEPS	407
WULLSTEIN RETRACTOR	153
WULLSTEIN SUCTION TUBE COMPLETE	17
Y	
YANCOSKIE RETRACTOR	177
YANKAUER NASAL SPECULUM	298
YANKAUER SUCTION TUBE	17
YASARGIL FORCEPS	68,69
YASARGIL NERVE DURA HOOK RETRACTOR	243
YASARGIL NEURO-MICRO DISSECTOR	246
YASARGILL SCISSORS	50
YEOMAN BIOPSY FORCEPS	335
YEOMANN INTESTINAL POLYOUS GRASPING FORCEPS	347

DESCRIPTION PAGE#

YEOMANN POLYOUS GRASPING FORCEPS	383
YOUNG LOBE HOLDING FORCEPS	365
YOUNG NEEDLE	123
YOUNG PROSTATIC RETRACTOR	360
Z	
ZAUFAL-JANSEN BONE RONGEUR	213
ZIEGLER CILIA FORCEPS	74
ZIEGLER PROTECTION CAP	01
ZIESER MEAEUS CLAMP	365
ZIMMERMAN MYOTOME, CHISEL TIP	404
ZIMMERMAN MYOTOME, SPOON TIP	404
ZOLLNER EAR SPECULUM	296
ZOLLNER, SUCTION TUBE	20
ZWEIFEL DOUBLE ENDED	399
ZWEIFEL NEEDLE HOLDERS T.C	121
ZYGOMATIC ARCH ELEVATOR	164



FIMA INSTRUMENTS

Manufacturer and Exporter of Surgical / Dental/ Manicure and
Veterinary Instruments.

H No. 8/121 Pura Hairan Rang Pura Road, Postal code : 51310 Sialkot Pakistan.
Cell & WhatsApp No. +92-3006106380 Cell No. +92-3349791619
Phone No. +92-52 4568278
E-mail: address:fima@skt.comsats.net.pk fimastruments@gmail.com
www.fimastruments.com